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1. Lawei Yang,Gang Liu,Ziyang Lin,Yahong Wang,Huijuan He,Tie Liu,David W Kamp. Pro-inflammatory response and oxidative stress induced by specific components in ambient particulate matter in human bronchial epithelial cells. *Environ Toxicol*. 2016 Aug;31(8):923-36.;doi: 10.1002/tox.22102
2. Wen D,Cui C,Duan W,Wang W,Wang Y,Liu Y,Li Z,Li C.The role of insulin-like growth factor 1 in ALS cell and mouse models: A mitochondrial protector. *Brain Res Bull*. 2019 Jan;144:1-13
3. Gao X,Zhou Y,Sun H,Liu D,Zhang J,Zhang J,Liu W,Pan X.Effects of a spiroketal compound Peniciketol A and its molecular mechanisms on growth inhibition in human leukemia. *TOXICOL APPL PHARM* . 2019 Mar 1;366:1-9
4. Sun L,Chen Y,Luo H,Xu M,Meng G,Zhang W.Ca2+/calmodulin-dependent protein kinase II regulation by inhibitor 1 of protein phosphatase 1 alleviates necroptosis in high glucose-induced cardiomyocytes injury. *Biochem Pharmacol*. 2019 May;163:194-205
5. Deng Z,Chai J,Zeng Q,Zhang B,Ye T,Chen X,Xu X.The anticancer properties and mechanism of action of tablysin-15, the RGD-containing disintegrin, in breast cancer cells. *Int J Biol Macromol*. 2019 May 15;129:1155-1167
6. Jia J,Yuan X,Peng X,Yan B.Cr(VI)/Pb2+ are responsible for PM2.5-induced cytotoxicity in A549 cells while pulmonary surfactant alleviates such toxicity. *ECOTOX ENVIRON SAFE*. 2019 May 15;172:152-158

7. Xiao X, Li S, Zhang X, Lu J, Wang W, Zhou S, Zhang J, Wang R, Li A. HHQ-4, a quinoline derivate, preferentially inhibits proliferation of glucose-deprived breast cancer cells as a GRP78 down-regulator. *TOXICOL APPL PHARM*. 2019 Jun 15;373:10-25
8. Cao YQ, Yuan L, Zhao Q, Yuan JL, Miao C, Chang YF, Wen XT, Wu R, Huang XB, Wen YP, Yan QG, Huang Y, Han XF, Ma XP, Cao SJ. Hsp40 Protein DNAJB6 Interacts with Viral NS3 and Inhibits the Replication of the Japanese Encephalitis Virus. *Int J Mol Sci*. 2019 Nov 14;20(22). pii: E5719
9. Ke Song, Weicheng Zhang, Caiyun Sun, Xiaomin Hu, Jiangzheng Wang, Lunguang Yao. Dynamic cytotoxicity of ZnO nanoparticles and bulk particles to *Escherichia coli*: A view from unfixed ZnO particle:Zn<sup>2+</sup> ratio. *Aquat Toxicol*. 2020 Mar;220:105407.;doi: 10.1016/j.aquatox.2020.105407
10. Xiaodong Sun, Guohui Sun, Yaxin Huang, Shufen Zhang, Xiaoyu Tang, Na Zhang, Lijiao Zhao, Rugang Zhong, Yongzhen Peng. Glycolytic inhibition by 3-bromopyruvate increases the cytotoxic effects of chloroethylnitrosourea to human glioma cells and the DNA interstrand cross-links formation. *Toxicology*. 2020 Apr 15;435:152413.;doi: 10.1016/j.tox.2020.152413
11. Jingyi Wang, Minghui Li, Xu Han, Hui Wang, Xinyang Wang, Ge Ma, Tiansong Xia, Shui Wang. MiR-1976 knockdown promotes epithelial-mesenchymal transition and cancer stem cell properties inducing triple-negative breast cancer metastasis. *Cell Death Dis*. 2020 Jul 3;11(7):500.;doi: 10.1038/s41419-020-2711-x
12. Jingyi Liu, Huimin Li, Tianqi Xia, Pengfei Du, Bikash Giri, Xue Li, Xuxin Li, Guofeng Cheng. Identification of *Schistosoma japonicum* GSK3 $\beta$  interacting partners by yeast two-hybrid screening and its role in parasite survival. *Parasitol Res*. 2020 Jul;119(7):2217-2226.;doi: 10.1007/s00436-020-06731-2
13. Zhenhua Guan, Jing Chen, Xueliang Li, Na Dong. Tanshinone IIA induces ferroptosis in gastric cancer cells through p53-mediated SLC7A11 down-regulation. *BIOSCIENCE REP*. 2020 Aug 28;40(8):BSR20201807.;doi: 10.1042/BSR20201807
14. Caiyun Sun, Weicheng Zhang, Runrun Ding, Jianzheng Wang, Lunguang Yao. Mechanism of low concentrations of polystyrene microplastics influence the cytotoxicity of Ag ions to *Escherichia coli*. *Chemosphere*. 2020 Aug;253:126705.;doi: 10.1016/j.chemosphere.2020.126705
15. Haoyuan Luo, Wenjun Zhu, Weiyan Mo, Min Liang. High-glucose concentration aggravates TNF-alpha-induced cell viability reduction in human CD146-positive periodontal ligament cells via TNFR-1 gene demethylation. *Cell Biol Int*. 2020 Dec;44(12):2383-2394.;doi: 10.1002/cbin.11445
16. Cheng Zhang, Pengqing Nie, Chunliu Zhou, Yue Hu, Suling Duan, Meijia Gu, Dongxu Jiang, Yunfu Wang, Zixin Deng, Jincao Chen, Shi Chen, Lianrong Wang. Oxidative stress-induced mitophagy is suppressed by the miR-106b-93-25 cluster in a protective manner. *Cell Death Dis*. 2021 Feb 24;12(2):209.;doi: 10.1038/s41419-021-03484-3.
17. Jin Hou, Ling Zhao, Huaqiao Tang, Xiaoli He, Gang Ye, Fei Shi, Min Kang, Helin Chen, Yinglun Li. Silver Nanoparticles Induced Oxidative Stress and Mitochondrial Injuries Mediated Autophagy in HC11 Cells Through Akt/AMPK/mTOR Pathway. *Biol Trace Elem Res*. 2021 Mar;199(3):1062-1073.;doi: 10.1007/s12011-020-02212-w.
18. Minmin Zhou, Shaobo Wang, Jiao Guo, Yang Liu, Junyuan Cao, Xiaohao Lan, Xiaoying Jia, Bo Zhang, Gengfu Xiao, Wei Wang. RNA Interference Screening Reveals Requirement for Platelet-Derived Growth Factor Receptor Beta in Japanese Encephalitis Virus Infection. *Antimicrob Agents Chemother*. 2021 May 18;65(6):e00113-21.;doi: 10.1128/AAC.00113-21.
19. Shunjie Zhang, Xinyu Wang, Ru Zhang, Yujun Cui, Heming Zhang, Weijie Song, Xiaohan Hou, Shibo Fu, Qingzhi Gao, Shengnan Liu. A GLUT1 inhibitor-based fluorogenic probe for Warburg effect-targeted drug screening and diagnostic imaging of hyperglycolytic cancers. *Anal Chim Acta*. 2021 Jul 4;1167:338593.;doi: 10.1016/j.aca.2021.338593.
20. Yi Han, Yanchao Liu, Xiaoyi Ma, Aijun Shen, Yiqiong Liu, Nopphon Weeranoppanant, Haiqing Dong, Yan Li, Tianbin Ren, Le Kuai, Bin Li, Maomao An, Yongyong Li. Antibiotics armed neutrophils as a potential therapy for brain fungal infection caused by chemotherapy-induced neutropenia. *Biomaterials*. 2021 Jul;274:120849.;doi: 10.1016/j.biomaterials.2021.120849.
21. Long Kong, Yanxin Wu, Cong Li, Jian Liu, Jianbo Jia, Hongyu Zhou, Bing Yan. Nano-cell and nano-pollutant interactions constitute key elements in nanoparticle-pollutant combined cytotoxicity. *J Hazard Mater*. 2021 Sep 15;418:126259.;doi: 10.1016/j.jhazmat.2021.126259.
22. Huihui Yang, Jun Chen, Yuling Liang, Yanfei Zhang, Wen Yin, Yuzhi Xu, Si-Yang Liu, Zong Dai, Xiaoyong Zou. A MOF-Shell-Confined I-Motif-Based pH Probe (MOFC-i) Strategy for Sensitive and Dynamic Imaging of Cell Surface pH. *ACS Appl Mater Interfaces*. 2021 Sep 29;13(38):45291-45299.;doi: 10.1021/acsami.1c13720.
23. Ya-Jing Zhang, Zhen-Lin Mu, Ping Deng, Yi-Dan Liang, Li-Chuan Wu, Ling-Ling Yang, Zhou Zhou, Zheng-Ping Yu. 8-Formylphlopiogonanone B induces ROS-mediated apoptosis in nasopharyngeal carcinoma CNE-1 cells. *Toxicol Res (Camb)*. 2021 Sep 6;10(5):1052-1063.;doi: 10.1093/toxres/tfab087.
24. Yawen Bu, Qingyuan Teng, Delan Feng, Lu Sun, Jia Xue, Guozhong Zhang. YLMY Tyrosine Residue within the Cytoplasmic Tail of Newcastle Disease Virus Fusion Protein Regulates Its Surface Expression to Modulate Viral Budding and Pathogenicity. *Microbiol Spectr*. 2021 Dec 22;9(3):e0217321.;doi: 10.1128/spectrum.02173-21.
25. Fenglian Yan, Dalei Cheng, Haiyan Wang, Min Gao, Junfeng Zhang, Hongyan Cheng, Changying Wang, Hui Zhang, Huabao Xiong. Corilagin Ameliorates Con A-Induced Hepatic Injury by Restricting M1 Macrophage Polarization. *Front Immunol*. 2022 Jan 13;12:807509.;doi: 10.3389/fimmu.2021.807509.
26. Ruitu Lyu, Xuguo Zhu, Yinghui Shen, Lijun Xiong, Lu Liu, Hang Liu, Feizhen Wu, Christian Argueta, Li Tan. Tumour suppressor TET2 safeguards enhancers from aberrant DNA methylation and epigenetic reprogramming in ER $\alpha$ -positive breast cancer cells. *Epigenetics*. 2022 Oct;17(10):1180-1194.;doi: 10.1080/15592294.2021.1997405.

Version 2023.03.02