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1. Wang J,Chen Y,Mo PL,Wei YJ,Liu KC,Zhang ZG,Zhang ZW,Chen XP,Zhang L. 1 α ,25-Dihydroxyvitamin D3 inhibits aflatoxin B1-induced proliferation and dedifferentiation of hepatic progenitor cells by regulating PI3K/Akt and Hippo pathways. J STEROID BIOCHEM. 2018 Oct;183:228-237.
2. Tang Y,Chen M,Zhou L,Ma J,Li Y,Zhang H,Shi Z,Xu Q,Zhang X,Gao Z,Zhao Y,Cheng Y,Jiao S,Zhou Z. Architecture, substructures, and dynamic assembly of STRIPAK complexes in Hippo signaling. Cell Discov. 2019 Jan 8;5:3.
3. Li Y,Ruan X,Xu X,Li C,Qiang T,Zhou H,Gao J,Wang X. Shengmai Injection Suppresses Angiotensin II-Induced Cardiomyocyte Hypertrophy and Apoptosis via Activation of the AMPK Signaling Pathway Through Energy-Dependent Mechanisms. Front Pharmacol. 2019 Sep 20;10:1095.
4. Naiqiang Zhu,Jingyi Hou. Molecular mechanism of the anti-inflammatory effects of Sophorae Flavescentis Aiton identified by network pharmacology. SCI REP-UK. 2021 Jan 13;11(1):1005.;doi: 10.1038/s41598-020-80297-y.

5. Siqing Jiang,Hao Zhang,Xin Li,Bin Yi,Lihua Huang,Zhaoxin Hu,Aimei Li,Jie Du,Yanchun Li,Wei Zhang. Vitamin D/VDR attenuate cisplatin-induced AKI by down-regulating NLRP3/Caspase-1/GSDMD pyroptosis pathway. *J Steroid Biochem Mol Biol.* 2021 Feb;206:105789.;doi: 10.1016/j.jsbmb.2020.105789.
6. Weibin Zhang,Chunhua Zhu,Fangnan Xiao,Xiaodong Liu,Anhua Xie,Fangman Chen,Panpan Dong,Pingdong Lin,Chenyang Zheng,Hong Zhang,Hui Gong,Yunkun Wu. pH-Controlled Release of Antigens Using Mesoporous Silica Nanoparticles Delivery System for Developing a Fish Oral Vaccine. *Front Immunol.* 2021 Apr 19;12:644396.;doi: 10.3389/fimmu.2021.644396.
7. Hao Liang,Rongchuan Yue,Chao Zhou,Mengyu Liu,Xi Yu,Shengzhong Lu,Jing Zeng,Zhengping Yu,Zhou Zhou,Houxiang Hu. Cadmium exposure induces endothelial dysfunction via disturbing lipid metabolism in human microvascular endothelial cells. *J Appl Toxicol.* 2021 May;41(5):775-788.;doi: 10.1002/jat.4115.
8. Xiaofan Jia,Yan Zhang,Ting Wang,Yuan Fu. Highly Efficient Method for Intracellular Delivery of Proteins Mediated by Cholera Toxin-Induced Protein Internalization. *Mol Pharm.* 2021 Nov 1;18(11):4067-4078.;doi: 10.1021/acs.molpharmaceut.1c00479.
9. Jian Luo,Qiang He,Jin-Zhi Xu,Chen Xu,Yin-Ze Han,Hai-Long Gao,Xian-Zhi Meng,Guo-Qing Pan,Tian Li,Ze-Yang Zhou. Microsporidia infection upregulates host energy metabolism but maintains ATP homeostasis. *J Invertebr Pathol.* 2021 Nov;186:107596.;doi: 10.1016/j.jip.2021.107596.
10. Hong Wen,Yang Fei,Ruisi Cai,Xuemei Yao,Yanan Li,Xuan Wang,Chencheng Xue,Yan Hu,Menghuan Li,Zhong Luo. Tumor-activatable biomaterialized nanotherapeutics for integrative glucose starvation and sensitized metformin therapy. *Biomaterials.* 2021 Nov;278:121165.;doi: 10.1016/j.biomaterials.2021.121165.

Version 2023.03.02