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1. Wang R, Sun Y, Yu W, Yan Y, Qiao M, Jiang R, Guan W, Wang L. Downregulation of miRNA-214 in cancer-associated fibroblasts contributes to migration and invasion of gastric cancer cells through targeting FGF9 and inducing EMT. *J EXP CLIN CANC RES*. 2019 Jan 15;38(1):20.
2. Ma J, Hu X, Liao C, Xiao H, Zhu Q, Li Y, Liu Z, Tao A, He Z, Xu C, Zheng K, Gypenoside L Inhibits Proliferation of Liver and Esophageal Cancer Cells by Inducing Senescence. *Molecules*. 2019 Mar 18;24(6). pii: E1054
3. hao YF, Zhao JY, Hu WZ, Ma K, Chao Y, Sun PJ, Fu XB, Zhang H. Synthetic poly(vinyl alcohol)-chitosan as a new type of highly efficient hemostatic sponge with blood-triggered swelling and high biocompatibility. *J Mater Chem B*. 2019 Mar 21;7(11):1855-1866
4. Lu W, Gao B, Fan J, Cheng P, Hu Y, Jie Q, Luo Z, Yang L. Mesenchymal Progenitors Derived from Different Locations in Long Bones Display Diverse Characteristics. *Stem Cells Int*. 2019 Apr 4;2019:5037578
5. Guo Z, Sui L, Qi J, Sun Q, Xu Y, Zou N, Xie Y, Kong Y. miR-196b inhibits cell migration and invasion through targeting MAP3K1 in hydatidiform mole. *Biomed Pharmacother*. 2019 May;113:108760
6. Sun B, Wang X, Pan Y, Jiao Y, Qi Y, Gong H, Jiang D. Antitumor effects of conditioned media of human fetal dermal mesenchymal stem cells on melanoma cells. *ONCOTARGETS THER*. 2019 May 28;12:4033-4046
7. Liu C, Jiang YH, Zhao ZL, Wu HW, Zhang L, Yang Z, Hoffman RM, Zheng JW. Knockdown of Histone Methyltransferase WHSC1 Induces Apoptosis and Inhibits Cell Proliferation and Tumorigenesis in Salivary Adenoid Cystic Carcinoma. *Anticancer Res*. 2019 Jun;39(6):2729-2737
8. Chen X, Tan W, Li W, Li W, Zhu S, Zhong J, Shang C, Chen Y. miR-1226-3p Promotes Sorafenib Sensitivity of Hepatocellular Carcinoma via Downregulation of DUSP4 Expression. *J Cancer*. 2019 Jun 2;10(12):2745-2753
9. Wang K, Chen B, Yin T, Zhan Y, Lu Y, Zhang Y, Chen J, Wu W, Zhou S, Mao W, Tan Y, Du B, Liu X, Ho HI, Xiao J. N-Methylparoxetine Blocked Autophagic Flux and Induced Apoptosis by Activating ROS-MAPK Pathway in Non-Small Cell Lung Cancer Cells. *Int J Mol Sci*. 2019 Jul 11;20(14). pii: E3415
10. Li P, Liu Y, Liu W, Li G, Tang Q, Zhang Q, Leng F, Sheng F, Hu C, Lai W, Liu Y, Zhou M, Huang J, Zhou H, Zhang R, Zhao Y. IR-783 inhibits breast cancer cell proliferation and migration by inducing mitochondrial fission. *Int J Oncol*. 2019 Aug;55(2):415-424
11. Shuiyuan Liu, Yingqiang Du, Kexin Shi, Yaqing Yang, Zhijian Yang. Resveratrol improves cardiac function by promoting M2-like polarization of macrophages in mice with myocardial infarction. *Am J Transl Res*. 2019 Aug 15;11(8):5212-5226
12. Zeng RJ, Zheng CW, Gu JE, Zhang HX, Xie L, Xu LY, Li EM. RAC1 inhibition reverses cisplatin resistance in esophageal squamous cell carcinoma and induces downregulation of glycolytic enzymes. *Mol Oncol*. 2019 Sep;13(9):2010-2030
13. Yang L, Xiao B, Hou L, Zhou G, Mo B, Yao D. Bioactive molecule Icaritin inhibited proliferation while enhanced apoptosis and autophagy of rat airway smooth muscle cells in vitro. *Cytotechnology*. 2019 Dec;71(6):1109-1120
14. Chu Q, Yu X, Jia R, Wang Y, Zhang Y, Zhang S, Liu Y, Li Y, Chen W, Ye X, Zheng X. Flavonoids from *Apios americana* Medikus Leaves Protect RAW264.7 Cells against Inflammation via Inhibition of MAPKs, Akt-mTOR Pathways, and Nfr2 Activation. *Oxid Med Cell Longev*. 2019 Dec 13;2019:1563024
15. Liu B, Liu J, Yu H, Wang C, Kong C. Transcription factor RUNX2 regulates epithelial-mesenchymal transition and progression in renal cell carcinomas. *Oncol Rep*. 2020 Feb;43(2):609-616
16. Huang ZX, Mao XM, Lin DC, Hong YH, Liang GS, Chen QX, Chen QH. Establishment and characterization of immortalized human eutopic endometrial stromal cells. *Am J Reprod Immunol*. 2020 Mar;83(3):e13213
17. Yumei Li, Xinghe Chen, Xixi Zeng, Shaokun Chen, Xi Yang, Li Zhang. Galectin-3 mediates pulmonary vascular endothelial cell dynamics via TRPC1/4 under acute hypoxia. *J BIOCHEM MOL TOXIC*. 2020 May;34(5):e22463.;doi: 10.1002/jbt.22463
18. Wenlong Yang, Shuning Zhang, Tiantong Ou, Hao Jiang, Daile Jia, Zhiyong Qi, Yunzeng Zou, Juying Qian, Aijun Sun, Junbo Ge. Interleukin-11 regulates the fate of adipose-derived mesenchymal stem cells via STAT3 signalling pathways. *CELL PROLIFERAT*. 2020 May;53(5):e12771.;doi: 10.1111/cpr.12771
19. Dan-Wen Wang, Fei Su, Tao Zhang, Tie-Cheng Yang, Hua-Qiao Wang, Li-Jie Yang, Fen-Fang Zhou, Mao-Hui Feng. The miR-370/UQCRC2 axis facilitates tumorigenesis by regulating epithelial-mesenchymal transition in Gastric Cancer. *J Cancer*. 2020 Jun 23;11(17):5042-5055.;doi: 10.7150/jca.45553
20. Jun Xu, Hua Yao, Shichen Wang, Huanrong Li, Xiaolin Hou. Mangiferin Inhibits Apoptosis and Autophagy Induced by *Staphylococcus aureus* in RAW264.7 Cells. *INFLAMM RES*. 2020 Nov 3;13:847-857.;doi: 10.2147/JIR.S280091

Version 2021.09.01