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- Liu H, Zhang C, Liang T, Song J, Hao J, Hou G. Inhibition of Pim2-prolonged skin allograft survival through the apoptosis regulation pathway. *Cell Mol Immunol.* 2012 Nov;9(6):503-10.
- Luo J, Zhao X, Tan Z, Su Z, Meng F, Zhang M. Mesenchymal-like progenitors derived from human embryonic stem cells promote recovery from acute kidney injury via paracrine actions. *Cytotherapy.* 2013 Jun;15(6):649-62.
- Shi H, Tan B, Ji G, Lu L, Cao A, Shi S, Xie J. Zedoary oil (Ezhu You) inhibits proliferation of AGS cells. *Chin Med.* 2013 Jun 28;8(1):13.
- Xi S, Yue L, Shi M, Peng Y, Xu Y, Wang X, Li Q, Kang Z, Li H, Wang Y. The Effects of Taoren-Honghua Herb Pair on Pathological Microvessel and Angiogenesis-Associated Signaling Pathway in Mice Model of CCl4-Induced Chronic Liver Disease. *Evid Based Complement Alternat Med.* 2016;2016:2974256.
- Zhou X, da Li W, Yuan BL, Niu LJ, Yang XY, Zhou ZB, Chen XH, Feng X. Lithium Treatment Prevents Apoptosis in Neonatal Rat Hippocampus Resulting from Sevoflurane Exposure. *Neurochem Res.* 2016 Aug;41(8):1993-2005.
- Wang B, Zhang Y, Cao W, Wei X, Chen J, Ying W. SIRT2 Plays Significant Roles in Lipopolysaccharides-Induced Neuroinflammation and Brain Injury in Mice. *Neurochem Res.* 2016 Sep;41(9):2490-500.
- Zhang JH, Jiao LY, Li TJ, Zhu YY, Zhou JW, Tian J. GSK-3 $\beta$  suppresses HCC cell dissociation in vitro by upregulating epithelial junction proteins and inhibiting Wnt/ $\beta$ -catenin signaling pathway. *J Cancer.* 2017 Jun 3;8(9):1598-1608.
- Chen YY, Sun ZW, Jiang JP, Kang XD, Wang LL, Shen YL, Xie XD, Zheng LR.  $\alpha$ -adrenoceptor-mediated enhanced inducibility of atrial fibrillation in a canine system inflammation model. *Mol Med Rep.* 2017 Jun;15(6):3767-3774.
- He Z, Yu Y, Nong Y, Du L, Liu C, Cao Y, Bai L, Tang H. Hepatitis B virus X protein promotes hepatocellular carcinoma invasion and metastasis via upregulating thioredoxin interacting protein. *Oncol Lett.* 2017 Aug;14(2):1323-1332.
- Wang X, Zhang W, Sun X, Lin Y, Chen W. Cancer-associated fibroblasts induce epithelial-mesenchymal transition through secreted cytokines in endometrial cancer cells. *Oncol Lett.* 2018 Apr;15(4):5694-5702.
- Xia HB, Wang HJ, Fu LQ, Wang SB, Li L, Ru GQ, He XL, Tong XM, Mou XZ, Huang DS. Decreased CRHBP expression is predictive of poor prognosis in patients with hepatocellular carcinoma. *Oncol Lett.* 2018 Sep;16(3):3681-3689.
- Han B, Feng D, Yu X, Liu Y, Yang M, Luo F, Zhou L, Liu F. MicroRNA-144 mediates chronic inflammation and tumorigenesis in colorectal cancer progression via regulating C-X-C motif chemokine ligand 11. *Exp Ther Med.* 2018 Sep;16(3):1935-1943.
- Yin XF, Zhang Q, Chen ZY, Wang HF, Li X, Wang HX, Li HX, Kang CM, Chu S, Li KF, Li Y, Qiu YR. NLRP3 in human glioma is correlated with increased WHO grade, and regulates cellular proliferation, apoptosis and metastasis via epithelial-mesenchymal transition and the PTEN/AKT signaling pathway. *Int J Oncol.* 2018 Sep;53(3):973-986.
- Peng H, Zhang K, Liu Z, Xu Q, You B, Li C, Cao J, Zhou H, Li X, Chen J, Cheng G, Shi R, Zhang G. VPO1 Modulates Vascular Smooth Muscle Cell Phenotypic Switch by Activating Extracellular Signal-regulated Kinase 1/2 (ERK 1/2) in Abdominal Aortic Aneurysms. *J Am Heart Assoc.* 2018 Sep 4;7(17):e010069.
- Zhang C, Hao Y, Wu L, Dong X, Jiang N, Cong B, Liu J, Zhang W, Tang D, De Perrot M, Zhao X. Curcumin induces apoptosis and inhibits angiogenesis in murine malignant mesothelioma. *Int J Oncol.* 2018 Dec;53(6):2531-2541.
- Xia YF, Chang J, Yang JF, Ouyang W, Pitt B, Billiar T, Zhang LM. Non-canonical Wnt signaling contributes to ventilator-induced lung injury through upregulation of WISP1 expression. *Int J Mol Med.* 2019 Mar;43(3):1217-1228.
- Zhu Z, Yu J, Lin W, Tang H, Zhang W, Lu B. Molecular hydrogen accelerates the reversal of acute obstructive cholangitis-induced liver dysfunction by restoring gap and tight junctions. *Mol Med Rep.* 2019 Jun;19(6):5177-5184.
- Chen L, Wu M, Zhang S, Tan W, Guan M, Feng L, Chen C, Tao J, Chen L, Qu L. Estrogen-related receptor  $\gamma$  regulates hepatic triglyceride metabolism through phospholipase A2 G12B. *FASEB J.* 2019 Jul;33(7):7942-7952.
- Diao H, Cheng N, Zhao Y, Xu H, Dong H, Thamm DH, Zhang D, Lin D. Ivermectin inhibits canine mammary tumor growth by regulating cell cycle progression and WNT signaling. *BMC Vet Res.* 2019 Aug 2;15(1):276.
- Wang Z, Zhang M, Shan R, Wang YJ, Chen J, Huang J, Sun LQ, Zhou WB. MTMR3 is upregulated in patients with breast cancer and regulates proliferation, cell cycle progression and autophagy in breast cancer cells. *Oncol Rep.* 2019 Nov;42(5):1915-1923.
- Hu Z, Niu G, Ren J, Wang X, Chen L, Hong R, Ke C. TFA5 promotes proliferation and migration in gastric cancer. *Mol Med Rep.* 2019 Nov;20(5):4477-4488.
- Li J, Xiang X, Xu H, Shi Y. Cilostazol Promotes Angiogenesis and Increases Cell Proliferation After Myocardial Ischemia-Reperfusion Injury Through a cAMP-Dependent Mechanism. *Cardiovasc Eng Technol.* 2019 Dec;10(4):638-647.
- Hai-Bing Xia, Hui-Ju Wang, Shu-Shu Song, Jun-Gang Zhang, Xiang-Lei He, Zhi-Ming Hu, Cheng-Wu Zhang, Dong-Sheng Huang, Xiao-Zhou Mou. Decreased DC-SIGNR expression in hepatocellular carcinoma predicts poor patient prognosis. *Oncol Lett.* 2020 Jan;19(1):69-76.;doi:10.3892/ol.2019.11074
- Yan Sun, Chuang Guo, Lin Yuan, Wen Li, Zhan-You Wang, Feng Yue, Jia-Yi Li. Cynomolgus Monkeys With Spontaneous Type-2-Diabetes-Mellitus-Like Pathology Develop Alpha-Synuclein Alterations Reminiscent of Prodromal Parkinson's Disease and Related Diseases. *FRONT NEUROSCI-SWITZ.* 2020 Feb 7;14:63.;doi:10.3389/fnins.2020.00063
- Ning-Wei Zhu, Xiao-Lan Yin, Ren Lin, Xiao-Lan Fan, Shi-Jie Chen, Yuan-Ming Zhu, Xiao-Zhen Zhao. Possible mechanisms of lycopene amelioration of learning and memory impairment in rats with vascular dementia. *Neural Regen Res.* 2020 Feb;15(2):332-341.;doi:10.4103/1673-5374.265565
- Wang J, Zuo Y, Zhuang K, Luo K, Yan X, Li J, Zhang JH, Liu F. Recombinant Human Milk Fat Globule-Epidermal Growth Factor 8 Attenuates Cerebrothrombosis after Subarachnoid Hemorrhage in Rats. *J STROKE CEREBROVASC.* 2020 Mar;29(3):104536
- Yuling Lian, Fangui Zhao, Wenjun Wang. Use of Bao Gui capsule in treatment of a polycystic ovary syndrome rat model. *Mol Med Rep.* 2020 Mar;21(3):1461-1470.;doi:10.3892/mmr.2020.10953
- Z-Z Wang, S-C Chen, X-B Zou, L-L Tian, S-H Sui, N-Z Liu. Nesfatin-1 alleviates acute lung injury through reducing inflammation and oxidative stress via the regulation of HMGB1. *EUR REV MED PHARMACO.* 2020 May;24(9):5071-5081.;doi:10.26355/eurrev\_202005\_21200
- Lingqiao Wang, Lixiong He, Hui Zeng, Wenjuan Fu, Jia Wang, Yao

- Tan,Chuanfen Zheng,Zhiqun Qiu,Jiaohua Luo,Chen Lv,Yujing Huang,Wei qun Shu.Low-dose microcystin-LR antagonizes aflatoxin B1 induced hepatocarcinogenesis through decreasing cytochrome P450 1A2 expression and aflatoxin B1-DNA adduct generation *Chemosphere*. 2020 Jun;248:126036.;doi: 10.1016/j.chemosphere.2020.126036
30. Xidan Zhu,Jia Feng,Wenguang Fu,Xiaoqia Shu,Xue Wan,Jinbo Liu.Effects of cisplatin on the proliferation, invasion and apoptosis of breast cancer cells following  $\beta$ -catenin silencing *Int J Mol Med*. 2020 Jun;45(6):1838-1850.;doi: 10.3892/ijmm.2020.4543
  31. Xiuming Zhang,Zheng Kang,Xiaodong Xie,Wei Qiao,Lei.Zhang,Zhen Gong,Yan Chen,Wenrong Shen Silencing of HIF-1  $\alpha$  inhibited the expression of lncRNA NEAT1 to suppress development of hepatocellular carcinoma under hypoxia *Am J Transl Res*. 2020 Jul 15;12(7):3871-3883
  32. Yiling Yang,Qianye Chen,Siru Zhou,Xinyi Gong,Hongyuan Xu,Yueyang Hong,Qinggong Dai,Lingyong Jiang.Skeletal Phenotype Analysis of a Conditional Stat3 Deletion Mouse Model *JOVE-J VIS EXP*. 2020 Jul 3;(161).;doi: 10.3791/61390
  33. Shubo Zhai,Baichao Sun,Yan Zhang,Lengyue Zhao,Li Zhang.IL-17 aggravates renal injury by promoting podocyte injury in children with primary nephrotic syndrome *Exp Ther Med*. 2020 Jul;20(1):409-417.;doi: 10.3892/etm.2020.8698
  34. Chengjie Ma,Peipei Peng,Yan Zhou,Tianya Liu,Lijuan Wang,Chen Lu.MicroRNA-93 promotes angiogenesis and attenuates remodeling via inactivation of the Hippo/Yap pathway by targeting Lats2 after myocardial infarction  $\omega$ . *Mol Med Rep*. 2020 Jul;22(1):483-493.;doi: 10.3892/mmr.2020.11085
  35. Zhen Wang,Han Jiang,Lu-Yao Cai,Ning Ji,Xin Zeng,Yu Zhou,Ying-Qiang Shen,Qian-Ming Chen.Repurposing disulfiram to induce OSCC cell death by cristae dysfunction promoted autophagy *Oral Dis*. 2020 Sep 28.;doi: 10.1111/odi.13652
  36. Xiuhui Wu,Leng Liu,Hongxia Zhang.miR-802 inhibits the epithelial-mesenchymal transition, migration and invasion of cervical cancer by regulating BTF3 *Mol Med Rep*. 2020 Sep;22(3):1883-1891.;doi: 10.3892/mmr.2020.11267
  37. Bin Zhan,Linjin Huang,Yachun Chen,Wen Ye,Jingkun Li,Jianhui Chen,Sheng Yang,Wei Jiang.miR-196a-mediated downregulation of p27 kipl protein promotes prostate cancer proliferation and relates to biochemical recurrence after radical prostatectomy *Prostate*. 2020 Sep;80(12):1024-1037.;doi: 10.1002/pros.24036.
  38. Wenjing Wang,Wenxiao Ding,Hanpeng Huang,Yina Zhu,Ning Ding,Ganzhu Feng,Xilong Zhang.The role of mitophagy in the mechanism of genioglossal dysfunction caused by chronic intermittent hypoxia and the protective effect of adiponectin.*Sleep Breath*. 2020 Oct 8.;doi: 10.1007/s11325-020-02211-0
  39. Juan Tan,Ling Liu,Zhihua Zuo,Bin Song,Tingting Cai,Dafa Ding,Yibing Lu,Xiaolong Ye.Overexpression of novel long intergenic non-coding RNA LINC02454 is associated with a poor prognosis in papillary thyroid cancer *Oncol Rep*. 2020 Oct;44(4):1489-1501.;doi: 10.3892/or.2020.7712
  40. Yun Chen,Chaojin Hong,Xiaochen Chen,Zhiquan Qin.Demethoxycurcumin increases the sensitivity of cisplatin-resistant non-small lung cancer cells to cisplatin and induces apoptosis by activating the caspase signaling pathway *Oncol Lett*. 2020 Nov;20(5):209.;doi: 10.3892/ol.2020.12072
  41. Hangshun Li,Haina Yang,Dandan Wang,Lixin Zhang,Tao Ma.Peroxiredoxin2 (Prdx2) Reduces Oxidative Stress and Apoptosis of Myocardial Cells Induced by Acute Myocardial Infarction by Inhibiting the TLR4/Nuclear Factor kappa B (NF- $\kappa$ B) Signaling Pathway *MED SCI MONITOR*. 2020 Dec 3;26:e926281.;doi: 10.12659/MSM.926281
  42. Zhi Yang,Weigang Wu,Pengcheng Ou,Minna Wu,Furong Zeng,Boping Zhou,Shipin Wu.MiR-122-5p knockdown protects against APAP-mediated liver injury through up-regulating NDRG3 *Mol Cell Biochem*. 2021 Feb;476(2):1257-1267.;doi: 10.1007/s11010-020-03988-0

Version 2021.09.01