碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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1. Zhang J, Yan H, Wu YP, Li C, Zhang GY. Activation of GluR6-containing kainate receptors induces ubiquitin-dependent Bcl-2 degradation via denitrosylation in the rat hippocampus after kainate treatment. *J Biol Chem*. 2011 Mar 4;286(9):7669-80.
2. Liu H, Zhang S, Lin H, Jia R, Chen Z. Identification of microRNA-RNA interactions using tethered RNAs and streptavidin aptamers. *Biochem Biophys Res Commun*. 2012 Jun 8;422(3):405-10.
3. Sun N, Hao JR, Li XY, Yin XH, Zong YY, Zhang GY, Gao C. GluR6-FasL-Trx2 mediates denitrosylation and activation of procaspase-3 in cerebral ischemia/reperfusion in rats. *Cell Death Dis*. 2013 Aug 15;4:e771.
4. Yang H, Zhao N, Lv L, Yan X, Hu S, Xu T. Functional research and molecular mechanism of Kainic acid-induced denitrosylation of thioredoxin-1 in rat hippocampus. *Neurochem Int*. 2017 Sep;108:448-456.
5. Zhang J, Hutchings P. Taxonomy and distribution of Terebellidae (Polychaeta: Trichobranchidae) in the northern South China Sea, with description of three new species. *Zootaxa*. 2018 Feb 2;4377(3):387-411.
6. Ping Yang, Li-Lei Zhang, Zi-Zhou Wang, Dan-Dan Zhang, Ya-Min Liu, Qing-Shan Shi, Xiao-Bao Xie. Nickel complexes of aroylhydrazone ligand: synthesis, crystal structure and DNA binding properties. *J Inorg Biochem*. 2020 Feb;203:110919.;doi: 10.1016/j.jinorgbio.2019.110919
7. Ping Yang, Hong Chen, Zi-Zhou Wang, Li-Lei Zhang, Dan-Dan Zhang, Qing-Shan Shi, Xiao-Bao Xie. Crystal structures and biological properties of aroylhydrazone Ni(II) complexes. *J Inorg Biochem*. 2020 Dec;213:111248.;doi: 10.1016/j.jinorgbio.2020.111248

Version 2021.09.01