碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

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1. Zheng LN, Wang M, Wang B, Chen HQ, Ouyang H, Zhao YL, Chai ZF, Feng WY. Determination of quantum dots in single cells by inductively coupled plasma mass spectrometry. *Talanta*. 2013 Nov 15;116:782-7.
2. Sun SN, Jia WD, Chen H, Ma JL, Ge YS, Yu JH, Li JS. Docosahexaenoic acid (DHA) induces apoptosis in human hepatocellular carcinoma cells. *Int J Clin Exp Pathol*. 2013;6(2):281-9.
3. Liu J, Mao Z, Huang J, Xie S, Liu T, Mao Z. Blocking the NOTCH pathway can inhibit the growth of CD133-positive A549 cells and sensitize to chemotherapy. *Biochem Biophys Res Commun*. 2014 Feb 21;444(4):670-5.
4. Wang HW, Wang JQ, Zheng BQ, Li SL, Zhang YD, Li FD, Zheng N. Cytotoxicity induced by ochratoxin A, zearalenone, and α -zearalenol: effects of individual and combined treatment. *Food Chem Toxicol*. 2014 Sep;71:217-24.
5. Sun L, Li H, Qu L, Zhu R, Fan X, Xue Y, Xie Z, Fan H. Immobilized lentivirus vector on chondroitin sulfate-hyaluronate acid-silk fibroin hybrid scaffold for tissue-engineered ligament-bone junction. *Biomed Res Int*. 2014;2014:816979.
6. Min Chen, Hongwei Lin, Yanjun Gao, Zaiqiang Wang, Yujuan Li, Faguang Jin. Ghrelin attenuates drowning injury via dual effects on damage protection and immune repression. *Ann Transl Med*. 2021 Jun;9(11):920.
7. Xuancheng Ou, Tianyong Wen, Jinwei Ying, Qing He, Anwu Xuan, Dike Ruan. MCP-1/CCR2 axis inhibits the chondrogenic differentiation of human nucleus pulposus mesenchymal stem cells. *Mol Med Rep*. 2022 Sep;26(3):277.

Version 2024.04.29