碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

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1. Duan C,Liu Y,Li Y,Chen H,Liu X,Chen X,Yue J,Zhou X,Yang J.Sulfasalazine alters microglia phenotype by competing endogenous RNA effect of miR-136-5p and long non-coding RNA HOTAIR in cuprizone-induced demyelination.Biochem Pharmacol. 2018 Sep;155:110-123.
2. Fenghua Lv,Zhuo Wang,Yanli Huang,Aoyang Si,Yulei Chen. CLEC3B protects H9c2 cardiomyocytes from apoptosis caused by hypoxia via the PI3K/Akt pathway Braz J Med Biol Res. 2020;53(9):e9693.;doi: 10.1590/1414-431x20209693.
3. Yang Gu,Shulan Zhang.High-throughput sequencing identification of differentially expressed microRNAs in metastatic ovarian cancer with experimental validations Cancer Cell Int. 2020 Oct 21;20:517.;doi: 10.1186/s12935-020-01601-4.
4. Jiayin Lu, Jiaqiang Huang, Shisu Zhao, Wenjiao Xu, Yaoxing Chen, Yuanyuan Li, Zixu Wang, Yanjun Dong, Renrong You, Jing Cao, Yulan Dong. FOXO1 Is a Critical Switch Molecule for Autophagy and Apoptosis of Sow Endometrial Epithelial Cells Caused by Oxidative Stress. Oxid Med Cell Longev. 2021 Dec 21;2021:1172273.
5. Fengchan Li, Kunmin Yan, Lili Wu, Zhong Zheng, Yun Du, Ziting Liu, Luyao Zhao, Wei Li, Yulan Sheng, Lijie Ren, Chaojun Tang, Li Zhu. Single-cell RNA-seq reveals cellular heterogeneity of mouse carotid artery under disturbed flow. Cell Death Discov. 2021 Jul 16;7(1):180.
6. Xuan Ma, Yang Liu, Hao Tian, Bo Zhang, Meiling Wang, Xia Gao. LINC01272 Suppressed Cell Multiplication and Induced Apoptosis Via Regulating MiR-7-5p/CRLS1 Axis in Lung Cancer. J Microbiol Biotechnol. 2021 Jul 28;31(7):921-932.
7. Peng Li, Zhiping Lin, Qianzheng Liu, Siyuan Chen, Xiang Gao, Weixiong Guo, Fan Gong, Jinsong Wei, Hao Lin. Enhancer RNA SLIT2 Inhibits Bone Metastasis of Breast Cancer Through Regulating P38 MAPK/c-Fos Signaling Pathway. Front Oncol. 2021 Oct 15;11:743840.
8. Tingting Zhang, Hui Liu, Min Liu, Chunhong Wang. Farrerol suppresses the progression of laryngeal squamous cell carcinoma via the mitochondria-mediated pathway. Eur J Pharmacol. 2021 Dec 15;913:174636.
9. Zongting Shi, Jun He, Jian He, Yuan Xu. High hydrostatic pressure (30 atm) enhances the apoptosis and inhibits the proteoglycan synthesis and extracellular matrix level of human nucleus pulposus cells via promoting the Wnt/ β -catenin pathway. Bioengineered. 2022 Feb;13(2):3070-3081.
10. Chengzhi Zhang, Yanmei Yang, Yanhui Gao, Dianjun Sun. NaF-induced neurotoxicity via activation of the IL-1 β /JNK signaling pathway. Toxicology. 2022 Mar 15;469:153132.
11. Xiaojin Liu, Hao Lv, Haibin Shen. Vitamin D enhances the sensitivity of breast cancer cells to the combination therapy of photodynamic therapy and paclitaxel. Tissue Cell. 2022 Aug;77:101815.
12. Cheng Yu, Danna Wu, Chong Zhao, Chaoguang Wu. Overexpressed Thrombospondin 2 Induced Osteogenic Differentiation of Valve Interstitial Cells via Inhibition of Akt/NF- κ B Signaling Pathway to Promote Calcific Aortic Valve Disease Development. Dis Markers. 2022 Sep 8;2022:2022958.
13. Juan Zhao, Yong Tan, Zhe Feng, Yahong Zhou, Feihong Wang, Ge Zhou, Jing Yan, Xiaowei Nie. Catalpol attenuates polycystic ovarian syndrome by regulating sirtuin 1 mediated NF- κ B signaling pathway. Reprod Biol. 2022 Sep;22(3):100671.
14. Yuqi Zhang, Enze Yao, Yijing Liu, Yining Zhang, Mengyang Ding, Jingyu Liu, Xiaoming Chen, Sairong Fan. FUT2 Facilitates Autophagy and Suppresses Apoptosis via p53 and JNK Signaling in Lung Adenocarcinoma Cells. Cells. 2022 Dec 13;11(24):4031.
15. Xu Han, Minghui Li, Jin Xu, Jingyue Fu, Xinyang Wang, Jingyi Wang, Tiansong Xia, Shui Wang, Ge Ma. miR-1275 targets MDK/AKT signaling to inhibit breast cancer chemoresistance by lessening the properties of cancer stem cells. Int J Biol Sci. 2023 Jan 1;19(1):89-103.
16. Xudong Dai, Yong Chen, Ning Chen, Jin Dou, Haiwen Zhuang, Jian Wang, Xin Zhao, Xiaoyu Zhang, Haijian Zhao. KLF5-mediated aquaporin 3 activated autophagy to facilitate cisplatin resistance of gastric cancer. Immunopharmacol Immunotoxicol. 2023 Apr;45(2):140-152.

Version 2024.04.29