碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

BeyoEmbryo™ Human Vitronectin (胚胎细胞培养级)

产品编号	产品名称	包装
C0318-500µg	BeyoEmbryo™ Human Vitronectin (胚胎细胞培养级)	500µl
C0318-2mg	BeyoEmbryo™ Human Vitronectin (胚胎细胞培养级)	2ml
C0318-10mg	BeyoEmbryo™ Human Vitronectin (胚胎细胞培养级)	10ml

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Species	Gene ID	Accession	Source	Length	MW	Tag
Human	7448	P04004	<i>E. coli</i>	~400aa	46kDa	6X His

About this protein	
Name	BeyoEmbryo™ Human Vitronectin (胚胎细胞培养级); BeyoEmbryo™ Human Vitronectin (Embryo cell culture grade); BeyoEmbryo™ Human Vitronectin (Embryo cell culture tested); Recombinant Human Vitronectin; 重组人玻连蛋白
Synonyms	VN; V75; VNT; VTN-N; rh-VTN; rhVTN; Vitronectin; S-Protein; Serum-Spreading Factor; Epibolin
Source	Recombinant protein containing the truncated human vitronectin (62-478aa) is expressed in <i>E. coli</i> .
Concentration	1mg/ml
Purity	Sterile, ≥ 95% by SDS-PAGE
Biological Activity	Used as a substrate for the feeder-free culture of human pluripotent stem cells (PSCs) and has been proven to maintain pluripotency and normal growth characteristics in multiple PSC lines such as H9 human embryonic stem cells. Measured by the ability of the immobilized protein to support the adhesion of DU145 human prostate carcinoma cells. When cells are added to Vitronectin-coated plates (10µg/ml and 100µl/well), approximately >40% cells will adhere specifically in 30 minutes at 37°C.
Physical Appearance	Liquid
Formulation	Liquid in 10mM Sodium Phosphate, 150mM NaCl, 0.1mM DTT, 0.1mM PMSF and 50% Glycerol, pH7.1
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in ultrapure water or aqueous buffer containing 0.1% BSA to a concentration of 0.1-1.0mg/ml. Stock solutions should be apportioned into working aliquots and stored at ≤ -20°C. Further dilutions should be made in appropriate buffered solutions.
Category	Extracellular Matrix (ECM), Cell Culture, Embryo cell culture
Background	Vitronectin has well-defined activity in cell adhesion, thrombosis, complement activation, fibrinolysis, inflammation, and platelet adhesion. Vitronectin binds to a variety of macromolecules and participates in regulatory circuits to control these processes. Some target macromolecules that interact with vitronectin are heparin, PAI-1, proteases such as thrombin and urokinase-type plasminogen activator (uPA), serine protease inhibitor-protease complexes, the uPA receptor, and a sub-class of integrin receptors on the surface of cells. Vitronectin is also applied as a thin coating on tissue-culture surfaces, which is useful to promote cell attachment, spreading, proliferation, and differentiation of many normal and neoplastic cells, and to study cell migration [1,2].

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