

## 细胞衰老β-半乳糖苷酶染色试剂盒

产品编号	产品名称	包装
C0602	细胞衰老β-半乳糖苷酶染色试剂盒	>100次

### 产品简介:

- 细胞衰老β-半乳糖苷酶染色试剂盒(Senescence β-Galactosidase Staining Kit)是一种基于衰老时衰老相关β-半乳糖苷酶(senescence-associated β-galactosidase; SA-β-Gal)活性水平上调而对衰老细胞或组织进行染色检测的试剂盒。在普通的光学显微镜下就可以观测到细胞或组织的衰老情况。本试剂盒可以用于培养细胞的衰老检测,也可以用于冰冻切片的衰老检测。
- 绝大多数正常细胞被认为仅有有限的分裂能力,在不能分裂后就进入衰老(senescence)状态。此时细胞仍然是存活的,但细胞的基因和蛋白的表达谱发生了很大改变。衰老的细胞不能在一些常规的刺激下再诱导细胞分裂,并且衰老细胞的细胞周期分布也比较特殊,不同于一些损伤诱导的细胞休眠,也不同于细胞生长接触抑制的情况。衰老细胞通常体积变大,并表达在pH6.0时有高酶活性的β-半乳糖苷酶(β-galactosidase)。细胞衰老也被认为是生物体抑制肿瘤的一种方式,同时也是生物体老化(aging)的一种潜在原因。
- 碧云天生产的细胞衰老β-半乳糖苷酶染色试剂盒,以X-Gal为底物,在衰老特异性的β-半乳糖苷酶催化下会生成深蓝色产物。从而在光学显微镜下很容易观察到变成蓝色的表达β-半乳糖苷酶的细胞或组织。本产品染色的染色效果请参考图1。




图1. 细胞衰老β-半乳糖苷酶染色试剂盒用于MCF-7细胞的染色效果图。图A是正常的MCF-7细胞,图B是MCF-7细胞用Etoposide(SC0173)处理诱导的衰老组。实际染色效果会因样品及实验条件的不同而存在差异,本图仅供参考。

- 本试剂盒仅染色衰老细胞,不会染色衰老前的细胞(presenescent cells)、静止期细胞(quiescent cells)、永生细胞(immortal cells)或肿瘤细胞。
- **主要特点:** 本试剂盒经过多方面的优化,是同类产品中首创的能兼容普通的细胞培养用多孔板、移液管等聚苯乙烯类材质耗材或容器的试剂盒。本试剂盒可以有效避免由于和多孔板、移液管等的不兼容导致的染色偏弱、染色效果不稳定等情况。通常同类试剂盒要求使用可高温高压灭菌的聚丙烯(polypropylene)材质的耗材、容器或玻璃容器进行溶液的配制,而不能使用普通的多孔板、移液管等聚苯乙烯(polystyrene)类材质的容器或耗材,否则可能会出现絮状沉淀,影响实验观察。即使严格按照要求操作,也会在染色时间比较长的情况下,容易出现絮状物沉淀(参考图2)。本试剂盒经过多方面的优化,对耗材或容器的材质无特殊要求,可以兼容普通的多孔板和移液管等常用耗材和容器。而且配制的工作液不会产生沉淀或不溶物,使用更加便捷。




图2. 本试剂盒优化前后的对比图。优化前, X-Gal溶液直接接触聚苯乙烯类材质的材料如移液管、多孔板等会产生明显的腐蚀(A图), 使用聚苯乙烯容器配制染色工作液后, 在显微镜下观察有异常的絮状不溶物(C图); 优化后, X-Gal溶液直接接触聚苯乙烯类材质的材料观察不到有任何异常情况(B图), 使用聚苯乙烯容器配制染色工作液后, 在显微镜下观察也没有任何异常情况(D图)。

- 如果使用6孔板检测, 足够测定100个样品; 使用24孔板测定, 足够测定400个样品; 使用96孔板测定, 足够测定1000个样品。对于组织切片或组织块, 可以检测的样品数量视样品的大小而定。对于普通切片的滴染足够检测100个样品。

### 包装清单:

产品编号	产品名称	包装
C0602-1	β-半乳糖苷酶染色固定液	100ml
C0602-2	X-Gal溶液	5ml
C0602-3	β-半乳糖苷酶染色液A	1ml
C0602-4	β-半乳糖苷酶染色液B	1ml
C0602-5	β-半乳糖苷酶染色液C	100ml
—	说明书	1份

### 保存条件:

-20°C保存, 一年有效。其中X-Gal溶液需避光保存。

### 注意事项:

- 加入染色工作液后, 由于溶液蒸发或其它未知原因等因素, 可能会有结晶形成而影响观察和拍摄照片, 此时建议吸除染色工作液, 加入适量70%乙醇进行洗涤, 70%乙醇可在短时间内溶解结晶, 待结晶溶解消失后再更换成PBS或生理盐水。70%乙醇的洗涤对染色效果没有任何影响。
- 在石蜡包埋的过程中, 温度、固定液等因素可能会导致β-半乳糖苷酶失活, 从而造成染色失败, 因此本试剂盒不建议用于石蜡切片的衰老检测。如果一定要用于石蜡切片的检测, 建议自行对实验条件进行一定的优化。
- β-半乳糖苷酶染色固定液对人体有毒、有腐蚀性, 操作时请特别小心, 并注意有效防护以避免直接接触人体或吸入体内。
- X-Gal溶液在-20°C或4°C保存会冻结, 室温或37°C水浴2-5分钟并适当摇动即可完全融解。
- 细胞衰老β-半乳糖苷酶染色反应依赖于特定的pH条件, 不能在二氧化碳培养箱中进行染色反应。用于细胞培养的二氧化碳培养箱中较高浓度的二氧化碳会影响染色工作液的pH值, 而导致染色失败。
- 试剂解冻后或使用前如果有沉淀, 必须在使用前确保沉淀全部溶解。β-半乳糖苷酶染色液B刚从试剂盒中取出时, 管底可能存在少量沉淀, 属正常现象, 充分混匀或Vortex后, 沉淀会全部溶解, 并须确保在全部溶解后使用。配制染色工作液时, 也可能有少量絮状沉淀出现, 震荡混匀后就会完全溶解, 且须确保全部溶解后才能使用。
- 使用96孔板等多孔板进行检测时, 如果孵育过夜容易产生所谓的‘边缘效应’(edge effect), 即多孔板四周的孔由于和外界最直接接触, 易受外界环境影响, 其中最明显的是四周细胞培养孔的蒸发效应。边缘效应会导致细胞生长不均匀、细胞分布不均一、培养液体积不一致、培养液中相关成分的浓度、pH值不一致。建议采取以下方法避免96孔板等多孔板的边缘效应: 避免孵育过长时间, 以避免蒸发等带来的边缘效应; 弃用边缘孔并在弃用的边缘孔中加入等量的水、PBS或其他适当溶液; 在多孔板非孔的凹陷处加入适量的水或其他适当溶液; 将整块板放在湿盒中; 使用防挥发盖; 在实验设计时, 实验样品最好进行随机分配, 不要将某一组样品固定放在某个位置而引入可能的系统性误差。
- 需自备PBS或HBSS(Hanks Balanced Salt Solution)。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

### 使用说明:

#### 1. 对于贴壁细胞:

- 对于6孔板中培养的细胞, 吸除细胞培养液, 用PBS或HBSS洗涤1次, 加入1毫升β-半乳糖苷酶染色固定液, 室温固定15分钟。对于其它类型的培养板, 固定液及后续溶液的用量参照此比例进行操作。
- 吸除细胞固定液, 用PBS或HBSS洗涤细胞3次, 每次3分钟。
- 吸除PBS或HBSS, 每孔加入1毫升染色工作液。染色工作液的配制方法参考表1。

表1. 染色工作液的配制方法。

β-半乳糖苷酶染色液A	10μl
β-半乳糖苷酶染色液B	10μl
β-半乳糖苷酶染色液C	930μl
X-Gal溶液	50μl

- 37°C孵育过夜, 可以用parafilm或保鲜膜封住6孔板防止蒸发。注意: 37°C孵育不能在二氧化碳培养箱中进行。
- 普通光学显微镜下观察。如不能及时观察计数, 可以去除染色工作液, 加入2毫升PBS, 4°C可以保存数天; 或者加上封片液封片后, 4°C可以保存较长时间。注: 如果有结晶形成, 请参考注意事项中的建议使用70%乙醇进行洗涤处理。

#### 2. 对于悬浮细胞:

- 离心收集细胞至1.5ml离心管内, 用PBS或HBSS洗涤1次, 加入1毫升β-半乳糖苷酶染色固定液, 室温固定15分钟。固定时可

以在摇床上缓慢摇动，以避免细胞结成团块。

- b. 离心，吸除细胞固定液，用PBS或HBSS洗涤细胞3次，每次3分钟。
- c. 离心，吸除PBS或HBSS，每管加入0.5-1毫升染色工作液。染色工作液的配制方法参考表1。
- d. 37°C孵育过夜。注意：37°C孵育不能在二氧化碳培养箱中进行。
- e. 取部分染色后的细胞，滴加到载玻片上或6孔板内，普通光学显微镜下观察。如不能及时观察计数，可以离心，去除染色工作液，然后加入1毫升PBS，4°C可以保存数天。如果离心，取细胞用于涂片，加上封片液封片后，4°C可以保存较长时间。注：如果有结晶形成，请参考注意事项中的建议使用70%乙醇进行洗涤处理。

### 3. 对于冰冻切片：

- a. 冰冻切片先进行复温，用PBS浸泡洗涤组织3次，每次不少于5分钟。
- b. 加入适当体积的β-半乳糖苷酶染色固定液，以充分盖住组织为宜，室温固定不少于15分钟。
- c. 用PBS浸泡洗涤组织3次，每次不少于5分钟。
- d. 吸除PBS，加入适当量的染色工作液。染色工作液的配制方法参考表1。
- e. 37°C孵育过夜，可以用parafilm或保鲜膜封住防止蒸发。最好把整个切片浸泡在染色工作液中。注意：37°C孵育不能在二氧化碳培养箱中进行。
- f. 普通光学显微镜下观察。如不能及时观察，加上封片液封片后4°C可以保存较长时间。注：如果有结晶形成，请参考注意事项中的建议使用70%乙醇进行洗涤处理。

### 使用本产品的文献：

1. Tao Q, Lv B, Qiao B, Zheng CQ, Chen ZF. Immortalization of ameloblastoma cells via reactivation of telomerase function: Phenotypic and molecular characteristics. *Oral Oncol*. 2009 Dec;45(12):e239-44.
2. Zhu C, Chen Q, Xie Z, Ai J, Tong L, Ding J, Geng M. The role of histone deacetylase 7 (HDAC7) in cancer cell proliferation: regulation on c-Myc. *J MOL MED*. 2011 Mar;89(3):279-89.
3. Liu XR, Cai Y, Cao X, Wei RC, Li HL, Zhou XM, Zhang KJ, Wu S, Qian QJ, Cheng B, Huang K, Liu XY. A new oncolytic adenoviral vector carrying dual tumour suppressor genes shows potent anti-tumour effect. *J Cell Mol Med*. 2012 Jun;16(6):1298-309.
4. Tang Y, Xu J, Qu W, Peng X, Xin P, Yang X, Ying C, Sun X, Hao L. Resveratrol reduces vascular cell senescence through attenuation of oxidative stress by SIRT1/NADPH oxidase-dependent mechanisms. *J Nutr Biochem*. 2012 Nov;23(11):1410-6.
5. Guo D, Zhang X, Huang Z, Zhou X, Zhu L, Zhao Y, Gu N. Comparison of cellular responses across multiple passage numbers in Ba/F3-BCR-ABL cells induced by silver nanoparticles. *Sci China Life Sci*. 2012 Oct;55(10):898-905.
6. Wang D, Hu Z, Hao J, He B, Gan Q, Zhong X, Zhang X, Shen J, Fang J, Jiang W. SIRT1 inhibits apoptosis of degenerative human disc nucleus pulposus cells through activation of Akt pathway. *AGE*. 2013 Oct;35(5):1741-53.
7. Ruan Q, Liu F, Gao Z, Kong D, Hu X, Shi D, Bao Z, Yu Z. The anti-inflammatory and hepatoprotective effects of huperzine A in D-galactose-treated rats. *Mech Ageing Dev*. 2013 Mar;134(3-4):89-97.
8. Zhang G, Morin C, Zhu X, Bao Huynh M, Ouidir Ouidja M, Sepulveda-Diaz JE, Raisman-Vozari R, Li P, Papy-Garcia D. Self-evolving oxidative stress with identifiable pre- and postmitochondrial phases in PC12 cells. *J Neurosci Res*. 2013 Feb;91(2):273-84.
9. Ye C, Zhang X, Wan J, Chang L, Hu W, Bing Z, Zhang S, Li J, He J, Wang J, Zhou G. Radiation-induced cellular senescence results from a slippage of long-term G2 arrested cells into G1 phase. *Cell Cycle*. 2013 May 1;12(9):1424-32.
10. Zhu Z, Li K, Xu D, Liu Y, Tang H, Xie Q, Xie L, Liu J, Wang H, Gong Y, Hu Z, Zheng J. Erratum to: ZFX regulates glioma cell proliferation and survival in vitro and in vivo. *NEURO-ONCOLOGY*. 2013 Jul;113(3):533.
11. Yu A, Zheng Y, Zhang R, Huang J, Zhu Z, Zhou R, Jin D, Yang Z. Resistin impairs SIRT1 function and induces senescence-associated phenotype in hepatocytes. *Mol Cell Endocrinol*. 2013 Sep 5;377(1-2):23-32.
12. Zhang JW, Zhang SS, Song JR, Sun K, Zong C, Zhao QD, Liu WT, Li R, Wu MC, Wei LX. Autophagy inhibition switches low-dose camptothecin-induced premature senescence to apoptosis in human colorectal cancer cells. *Biochem Pharmacol*. 2014 Aug 1;90(3):265-75.
13. Chen B, Zhang W, Gao J, Chen H, Jiang L, Liu D, Cao Y, Zhao S, Qiu Z, Zeng J, Zhang S, Li W. Downregulation of ribosomal protein S6 inhibits the growth of non-small cell lung cancer by inducing cell cycle arrest, rather than apoptosis. *Cancer Lett*. 2014 Nov 28;354(2):378-89.
14. Chen H, Liu X, Zhu W, Chen H, Hu X, Jiang Z, Xu Y, Wang L, Zhou Y, Chen P, Zhang N, Hu D, Zhang L, Wang Y, Xu Q, Wu R, Yu H, Wang J. SIRT1 ameliorates age-related senescence of mesenchymal stem cells via modulating telomere shelterin. *Front Aging Neurosci*. 2014 Jun 3;6:103.
15. Gong L, Gong HJ, Pan X, et al. p53 isoform Δ113p53/Δ133p53 promotes DNA double-strand break repair to protect cell from death and senescence in response to DNA damage. *Cell Research*. 2015 Feb(25):351-369.
16. Zhou XX, Yang Q, Xie YH, Sun JY, Hu J, Qiu PC, Cao W, Wang SW. Tetrahydroxy stilbene glucoside extends mouse life span via upregulating neural klotho and downregulating neural insulin or insulin-like growth factor 1. *Neurobiol Aging*. 2015 Mar;36(3):1462-1470.
17. Qi XF, Chen ZY, Xia JB, Zheng L, Zhao H, Pi LQ, Park KS, Kim SK, Lee KJ, Cai DQ. FoxO3a suppresses the senescence of cardiac microvascular endothelial cells by regulating the ROS-mediated cell cycle. *J Mol Cell Cardiol*. 2015 Apr;81:114-26.
18. Li Z, Han S, Wang X, Han F, Zhu X, Zheng Z, Wang H, Zhou Q, Wang Y, Su L, Shi J, Tang C, Hu D. Rhokinase inhibitor Y-27632 promotes the differentiation of human bone marrow mesenchymal stem cells into keratinocyte-like cells in xeno-free conditioned medium. *Stem Cell Res Ther*. 2015 Mar 11;6:17.
19. Wu RX, Bi CS, Yu Y, Zhang LL, Chen FM. Age-related decline in the matrix contents and functional properties of human periodontal ligament stem cell sheets. *Acta Biomater*. 2015 Aug;22:70-82.
20. Xia N, Cai J, Wang F, Dong B, Liu S, Chen F, Cheng J, Zuo Y. SENP1 is a crucial regulator for cell senescence through deSUMOylation of Bmi1. *SCI REP-UK*. 2016 Sep 23;6:34099.
21. Zhao N, Han D, Liu H, Li Y, Wong SW, Cao Z, Xu J, Zhang X, Cai T, Wang Y, Feng H. Senescence: novel insight into DLX3 mutations leading to enhanced bone formation in Tricho-Dento-Osseous syndrome. *SCI REP-UK*. 2016 Dec 7;6:38680.
22. Niu T, Tian Y, Mei Z, Guo G. Inhibition of autophagy enhances curcumin United light irradiation-induced oxidative stress and tumor growth suppression in human melanoma cells. *SCI REP-UK*. 2016 Aug 9;6:31383.
23. Gao Q, Chen K, Gao L, Zheng Y, Yang YG. Thrombospondin-1 signaling through CD47 inhibits cell cycle progression and induces senescence in endothelial cells. *Cell Death Dis*. 2016 Sep 8;7(9):e2368.
24. Meng J, Lv Z, Qiao X, Li X, Li Y, Zhang Y, Chen C. The decay of redox-stress response capacity is a substantive characteristic of aging: revising the redox theory of aging. *Redox Biol*. 2016 Dec 28;11:365-374.
25. Wang C, Chen H, Zhang M, Zhang J, Wei X, Ying W. Malate-aspartate shuttle inhibitor aminooxyacetic acid leads to decreased intracellular ATP levels

- unaltered cell cycle of C6 glioma cells by inhibiting glycolysis. *Cancer Lett.* 2016 Aug 1;378(1):1-7.
26. Chen Y, Pan K, Wang P, Cao Z, Wang W, Wang S, Hu N, Xue J, Li H, Jiang W, Li G, Zhang X. HBP1-mediated Regulation of p21 Protein through the Mdm2/p53 and TCF4/EZH2 Pathways and Its Impact on Cell Senescence and Tumorigenesis. *J Biol Chem.* 2016 Jun 10;291(24):12688-705.
  27. Sui B, Hu C, Liao L, Chen Y, Zhang X, Fu X, Zheng C, Li M, Wu L, Zhao X, Jin Y. Mesenchymal progenitors in osteopenias of diverse pathologies: differential characteristics in the common shift from osteoblastogenesis to adipogenesis. *SCI REP-UK.* 2016 Jul 22;6:30186.
  28. Yang X, Liu S, Huang C, Wang H, Luo Y, Xu W, Huang K. Ochratoxin A induced premature senescence in human renal proximal tubular cells. *Toxicology.* 2017 May 1;382:75-83.
  29. Li P, Hou G, Zhang R, Gan Y, Xu Y, Song L, Zhou Q. High-magnitude compression accelerates the premature senescence of nucleus pulposus cells via the p38 MAPK-ROS pathway. *Arthritis Res Ther.* 2017 Sep 18;19(1):209.
  30. Li P, Gan Y, Xu Y, Song L, Wang L, Ouyang B, Zhang C, Zhou Q. The inflammatory cytokine TNF- $\alpha$  promotes the premature senescence of rat nucleus pulposus cells via the PI3K/Akt signaling pathway. *SCI REP-UK.* 2017 Feb 17;7:42938.
  31. Hu Y, Zhou KY, Wang ZJ, Lu Y, Yin M. N-stearoyl-l-Tyrosine inhibits the cell senescence and apoptosis induced by H<sub>2</sub>O<sub>2</sub> in HEK293/Tau cells via the CB2 receptor. *CHEM-BIOL INTERACT.* 2017 Jun 25;272:135-144.
  32. Tao YF, Wang NN, Xu LX, Li ZH, Li XL, Xu YY, Fang F, Li M, Qian GH, Li YH, Li YP, Wu Y, Ren JL, Du WW, Lu J, Feng X, Wang J, He WQ, Hu SY, Pan J. Molecular mechanism of G1 arrest and cellular senescence induced by LEE011, a novel CDK4/CDK6 inhibitor, in leukemia cells. *Cancer Cell Int.* 2017 Mar 6;17:35.
  33. Li P, Gan Y, Xu Y, Wang L, Ouyang B, Zhang C, Luo L, Zhao C, Zhou Q. 17 $\beta$ -estradiol Attenuates TNF- $\alpha$ -Induced Premature Senescence of Nucleus Pulposus Cells through Regulating the ROS/NF- $\kappa$ B Pathway. *Int J Biol Sci.* 2017 Jan 15;13(2):145-156.
  34. You J, Sun J, Ma T, Yang Z, Wang X, Zhang Z, Li J, Wang L, Li M, Yang J, Shen Z. Curcumin induces therapeutic angiogenesis in a diabetic mouse hindlimb ischemia model via modulating the function of endothelial progenitor cells. *Stem Cell Res Ther.* 2017 Aug 3;8(1):182.
  35. Liu C, Zhang S, Wang Q, Zhang X. Tumor suppressor miR-1 inhibits tumor growth and metastasis by simultaneously targeting multiple genes. *ONCOTARGET.* 2017 Jun 27;8(26):42043-42060.
  36. Wu W, Hou CL, Mu XP, Sun C, Zhu YC, Wang MJ, Lv QZ. H<sub>2</sub>S Donor NaHS Changes the Production of Endogenous H<sub>2</sub>S and NO in D-Galactose-Induced Accelerated Ageing. *Oxid Med Cell Longev.* 2017;2017:5707830.
  37. Zhao Y, Jia Z, Huang S, Wu Y, Liu L, Lin L, Wang D, He Q, Ruan D. Age-Related Changes in Nucleus Pulposus Mesenchymal Stem Cells: An In Vitro Study in Rats. *Stem Cells Int.* 2017;2017:6761572.
  38. Yan Y, Du Y, Zheng H, Wang G, Li R, Chen J, Li K. NS1 of H7N9 Influenza A Virus Induces NO-Mediated Cellular Senescence in Neuro2a Cells. *CELL PHYSIOL BIOCHEM.* 2017;43(4):1369-1380.
  39. Guo J, Wu Y, Yang L, Du J, Gong K, Chen W, Dai J, Li X, Xi S. Repression of YAP by NCTD disrupts NSCLC progression. *ONCOTARGET.* 2017 Jan 10;8(2):2307-2319.
  40. Li X, Song Y, Liu D, Zhao J, Xu J, Ren J, Hu Y, Wang Z, Hou Y, Zhao G. MiR-495 Promotes Senescence of Mesenchymal Stem Cells by Targeting Bmi-1. *CELL PHYSIOL BIOCHEM.* 2017;42(2):780-796.
  41. Xu W, Li Y, Yuan WW, Yin Y, Song WW, Wang Y, Huang QQ, Zhao WH, Wu JQ. Membrane-Bound CD40L Promotes Senescence and Initiates Senescence-Associated Secretory Phenotype via NF- $\kappa$ B Activation in Lung Adenocarcinoma. *CELL PHYSIOL BIOCHEM.* 2018;48(4):1793-1803.
  42. Zhao D, Sun Y, Tan Y, Zhang Z, Hou Z, Gao C, Feng P, Zhang X, Yi W, Gao F. Short-Duration Swimming Exercise after Myocardial Infarction Attenuates Cardiac Dysfunction and Regulates Mitochondrial Quality Control in Aged Mice. *Oxid Med Cell Longev.* 2018 Apr 11;2018:4079041.
  43. Sun J, Zhang Z, Ma T, Yang Z, Zhang J, Liu X, Lu D, Shen Z, Yang J, Meng Q. Endothelial progenitor cell-derived exosomes, loaded with miR-126, promoted deep vein thrombosis resolution and recanalization. *Stem Cell Res Ther.* 2018 Aug 23;9(1):223.
  44. Quan XJ, Liang CL, Sun MZ, Zhang L, Li XL. Overexpression of steroid receptor coactivators alleviates hyperglycemia-induced endothelial cell injury in rats through activating the PI3K/Akt pathway. *Acta Pharmacol Sin.* 2018 Aug 8.
  45. Chen Y, Li Y, Peng Y, Zheng X, Fan S, Yi Y, Zeng P, Chen H, Kang H, Zhang Y, Xiao ZX, Li C.  $\Delta$ Np63 $\alpha$  down-regulates c-Myc modulator MM1 via E3 ligase HERC3 in the regulation of cell senescence. *Cell Death Differ.* 2018 Dec;25(12):2118-2129.
  46. Jiang W, Wei K, Pan C, Li H, Cao J, Han X, Tang Y, Zhu S, Yuan W, He Y, Xia Y, Chen L, Chen Y. MicroRNA-1258 suppresses tumour progression via GRB2/Ras/Erk pathway in non-small-cell lung cancer. *CELL PROLIFERAT.* 2018 Dec;51(6):e12502.
  47. Li R, Dong J, Bu XQ, Huang Y, Yang JY, Dong X, Liu J. Interleukin-6 promotes the migration and cellular senescence and inhibits apoptosis of human intrahepatic biliary epithelial cells. *J Cell Biochem.* 2018 Feb;119(2):2135-2143.
  48. Li X, Zhang Z, Zhang Y, Cao Y, Wei H, Wu Z. Upregulation of lactate-inducible snail protein suppresses oncogene-mediated senescence through p16INK4a inactivation. *J EXP CLIN CANC RES.* 2018 Feb 26;37(1):39.
  49. Deng X, Zhang X, Li W, Feng RX, Li L, Yi GR, Zhang XN, Yin C, Yu HY, Zhang JP, Lu B, Hui L, Xie WF. Chronic Liver Injury Induces Conversion of Biliary Epithelial Cells into Hepatocytes. *Cell Stem Cell.* 2018 Jul 5;23(1):114-122.e3.
  50. Daguia Zambé JC, Zhai Y, Zhou Z, Du X, Wei Y, Ma F, Hua J. miR-19b-3p induces cell proliferation and reduces heterochromatin-mediated senescence through PLZF in goat male germline stem cells. *J Cell Physiol.* 2018 Jun;233(6):4652-4665.
  51. Li XC, Wang MS, Liu W, Zhong CF, Deng GB, Luo SJ, Huang CM. Co-culturing nucleus pulposus mesenchymal stem cells with notochordal cell-rich nucleus pulposus explants attenuates tumor necrosis factor- $\alpha$ -induced senescence. *Stem Cell Res Ther.* 2018 Jun 26;9(1):171.
  52. Xie J, Li B, Zhang P, Wang L, Lu H, Song X. Osteogenic protein-1 attenuates the inflammatory cytokine-induced NP cell senescence through regulating the ROS/NF- $\kappa$ B pathway. *Biomed Pharmacother.* 2018 Mar;99:431-437.
  53. Han X, Chen H, Zhou J, Tai H, Gong H, Wang X, Huang N, Qin J, Fang T, Wang F, Xiao H. The inhibitory effect of Fraxinellone on oxidative stress-induced senescence correlates with AMP-activated protein kinase-dependent autophagy restoration. *J Cell Physiol.* 2018 May;233(5):3945-3954.
  54. Chen C, Tang Q, Zhang Y, Yu M, Jing W, Tian W. Physioxia: a more effective approach for culturing human adipose-derived stem cells for cell transplantation. *Stem Cell Res Ther.* 2018 May 24;9(1):148.
  55. Jin D, Guo J, Wang D, Wu Y, Wang X, Gao Y, Shao C, Xu X, Tan S. The antineoplastic drug metformin downregulates YAP by interfering with IRF-1 binding to the YAP promoter in NSCLC. *EBioMedicine.* 2018 Nov;37:188-204.
  56. Sun K, Yang P, Zhao R, Bai Y, Guo Z. Matrine Attenuates D-Galactose-Induced Aging-Related Behavior in Mice via Inhibition of Cellular Senescence and Oxidative Stress. *Oxid Med Cell Longev.* 2018 Nov 27;2018:7108604.
  57. Guo J, Jin D, Wu Y, Yang L, Du J, Gong K, Chen W, Dai J, Miao S, Xi S. The miR 495-UBE2C-ABCG2/ERCC1 axis reverses cisplatin resistance by downregulating drug resistance genes in cisplatin-resistant non-small cell lung cancer cells. *EBioMedicine.* 2018 Sep;35:204-221.
  58. Ding J, Yu C, Sui Y, Wang L, Yang Y, Wang F, Yao H, Xing F, Liu H, Li Y, Shah JA, Cai Y, Jin J. The chromatin remodeling protein INO80 contributes to the removal of H2A.Z at the p53-binding site of the p21 gene in response to doxorubicin. *FEBS J.* 2018 Sep;285(17):3270-3285.
  59. Ge H, Ke J, Xu N, Li H, Gong J, Li X, Song Y, Zhu H, Bai C. Dexamethasone alleviates pemetrexed-induced senescence in Non-Small-Cell Lung Cancer. *Food Chem Toxicol.* 2018 Sep;119:86-97.
  60. Liu H, Kang H, Song C, Lei Z, Li L, Guo J, Xu Y, Guan H, Fang Z, Li F. Urolithin A Inhibits the Catabolic Effect of TNF $\alpha$  on Nucleus Pulposus Cell and Alleviates Intervertebral Disc Degeneration in vivo. *Front Pharmacol.* 2018 Sep 18;9:1043.
  61. He Y, Xu H, Xiang Z, Yu H, Xu L, Guo Y, Tian Y, Shu R, Yang X, Xue C, Zhao M, He Y, Han X, Bai D. YAP regulates periodontal ligament cell differentiation into

- myofibroblast interacted with RhoA/ROCK pathway. *J Cell Physiol.* 234(4):5086-5096.
62. He Q, Gao Y, Wang T, Zhou L, Zhou W, Yuan Z. Deficiency of Yes-Associated Protein Induces Cataract in Mice. *Aging Dis.* 10(2):293-306.
  63. Zheng BN, Ding CH, Chen SJ, Zhu K, Shao J, Feng J, Xu WP, Cai LY, Zhu CP, Duan W, Ding J, Zhang X, Luo C, Xie WF. Targeting PRMT5 Activity Inhibits the Malignancy of Hepatocellular Carcinoma by Promoting the Transcription of HNF4  $\alpha$ . *Theranostics.* 9(9):2606-2617.
  64. He Q, Xue S, Tan Y, Zhang L, Shao Q, Xing L, Li Y, Xiang T, Luo X, Ren G. Dual inhibition of Akt and ERK signaling induces cell senescence in triple-negative breast cancer. *Cancer Lett.* 448:94-104.
  65. Chen X, Wang L, Hou J, Li J, Chen L, Xia J, Wang Z, Xiao M, Wang Y. Study on the Dynamic Biological Characteristics of Human Bone Marrow Mesenchymal Stem Cell Senescence. *Stem Cells Int.* 2019:9271595.
  66. Zhu B, Zhang L, Liang C, Liu B, Pan X, Wang Y, Zhang Y, Zhang Y, Xie W, Yan B, Liu F, Yip HK, Yu XY, Li Y. Stem Cell-Derived Exosomes Prevent Aging-Induced Cardiac Dysfunction through a Novel Exosome/lncRNA MALAT1/NF- $\kappa$ B/TNF- $\alpha$  Signaling Pathway. *Oxid Med Cell Longev.* 2019:9739258.
  67. Guo Z, Tian Y, Guo Y, Li B, Liu X, Xie K, Song Y, Wang D. RAD6B Plays a Critical Role in Neuronal DNA Damage Response to Resist Neurodegeneration. *Front Cell Neurosci.* 13:392.
  68. Lin F, Zeng Z, Song Y, Li L, Wu Z, Zhang X, Li Z, Ke X, Hu X. YBX-1 mediated sorting of miR-133 into hypoxia/reoxygenation-induced EPC-derived exosomes to increase fibroblast angiogenesis and MEndoT. *Stem Cell Res Ther.* 10(1):263.
  69. Shen L, Chen Y, Cheng J, Yuan S, Zhou S, Yan W, Liu J, Luo A, Wang S. CCL5 secreted by senescent theca-interstitial cells inhibits preantral follicular development via granulosa cellular apoptosis. *J Cell Physiol.* 234(12):22554-22564.
  70. Kang L, Xiang Q, Zhan S, Song Y, Wang K, Zhao K, Li S, Shao Z, Yang C, Zhang Y. Restoration of Autophagic Flux Rescues Oxidative Damage and Mitochondrial Dysfunction to Protect against Intervertebral Disc Degeneration. *Oxid Med Cell Longev.* 2019:7810320.
  71. Zhang P, Wang Q, Nie L, Zhu R, Zhou X, Zhao P, Ji N, Liang X, Ding Y, Yuan Q, Wang Q. Hyperglycemia-induced inflamm-aging accelerates gingival senescence via NLRP4 phosphorylation. *J Biol Chem.* 294(49):18807-18819.
  72. Han Y, Wu P, Wang Z, Zhang Z, Sun S, Liu J, Gong S, Gao P, Iwakuma T, Molina-Vila MA, Chen BP, Zhang Y, Ji T, Mo Q, Chen P, Hu J, Wang S, Zhou J, Lu H, Gao Q. Ubiquinol-cytochrome C reductase core protein II promotes tumorigenesis by facilitating p53 degradation. *EBioMedicine.* 40:92-105.
  73. Guan L, Chen Y, Wang Y, Zhang H, Fan S, Gao Y, Jiao T, Fu K, Sun J, Yu A, Huang M, Bi H. Effects of carnitine palmitoyltransferases on cancer cellular senescence. *J Cell Physiol.* 234(2):1707-1719.
  74. He S, Lin J, Xu Y, Lin L, Feng J. A positive feedback loop between ZNF205-AS1 and EGR4 promotes non-small cell lung cancer growth. *J Cell Mol Med.* 23(2):1495-1508.
  75. Zheng G, Pan Z, Zhan Y, Tang Q, Zheng F, Zhou Y, Wu Y, Zhou Y, Chen D, Chen J, Wang X, Gao W, Xu H, Tian N, Zhang X. TFEB protects nucleus pulposus cells against apoptosis and senescence via restoring autophagic flux. *OSTEOARTHR CARTILAGE.* 27(2):347-357.
  76. Li Y, Xue Y, Xu X, Wang G, Liu Y, Wu H, Li W, Wang Y, Chen Z, Zhang W, Zhu Y, Ji W, Xu T, Liu L, Chen Q. A mitochondrial FUNDCl/HSC70 interaction organizes the proteostatic stress response at the risk of cell morbidity. *EMBO J.* 38(3). pii: e98786.
  77. Wang L, Xie Y, Yang W, Yang Z, Jiang S, Zhang C, Zhang G. Alfalfa polysaccharide prevents H2O2-induced oxidative damage in MEFs by activating MAPK/Nrf2 signaling pathways and suppressing NF- $\kappa$ B signaling pathways. *SCI REP-UK.* 9(1):1782.
  78. Guo L, Wang X, Yuan J, Zhu M, Fu X, Xu RH, Wu C, Wu Y. TSA restores hair follicle-inductive capacity of skin-derived precursors. *SCI REP-UK.* 9(1):2867.
  79. Xia Y, Wei K, Yang FM, Hu LQ, Pan CF, Pan XL, Wu WB, Wang J, Wen W, He ZC, Xu J, Xu XF, Zhu Q, Chen L. miR-1260b, mediated by YY1, activates KIT signaling by targeting SOCS6 to regulate cell proliferation and apoptosis in NSCLC. *Cell Death Dis.* 10(2):112.
  80. Fu GB, Huang WJ, Zeng M, Zhou X, Wu HP, Liu CC, Wu H, Weng J, Zhang HD, Cai YC, Ashton C, Ding M, Tang D, Zhang BH, Gao Y, Yu WF, Zhai B, He ZY, Wang HY, Yan HX. Expansion and differentiation of human hepatocyte-derived liver progenitor-like cells and their use for the study of hepatotropic pathogens. *Cell Res.* 29(1):8-22.
  81. Wei Z, Guo H, Qin J, Lu S, Liu Q, Zhang X, Zou Y, Gong Y, Shao C. Pan-senescence transcriptome analysis identified RRAD as a marker and negative regulator of cellular senescence. *FREE RADICAL BIO MED.* 130:267-277.
  82. Liu R, Liu J, Wang S, Wang Y, Zhang T, Liu Y, Geng X, Wang F. Combined treatment with emodin and a telomerase inhibitor induces significant telomere damage/dysfunction and cell death. *Cell Death Dis.* 10(7):527.
  83. Zhang L, Li F, Su X, Li Y, Wang Y, Fang R, Guo Y, Jin T, Shan H, Zhao X, Yang R, Shan H, Liang H. Melatonin prevents lung injury by regulating apelin 13 to improve mitochondrial dysfunction. *Exp Mol Med.* 51(7):73.
  84. Zuo R, Wang Y, Li J, Wu J, Wang W, Li B, Sun C, Wang Z, Shi C, Zhou Y, Liu M, Zhang C. Rapamycin Induced Autophagy Inhibits Inflammation-Mediated Endplate Degeneration by Enhancing Nrf2/Keap1 Signaling of Cartilage Endplate Stem Cells. *Stem Cells.* 37(6):828-840.
  85. He H, Yu B, Liu Z, Ye G, You W, Hong Y, Lian Q, Zhang Y, Li X. Vascular progenitor cell senescence in patients with Marfan syndrome. *J Cell Mol Med.* 23(6):4139-4152.
  86. Wang H, Hu Z, Wu J, Mei Y, Zhang Q, Zhang H, Miao D, Sun W. Sirt1 Promotes Osteogenic Differentiation and Increases Alveolar Bone Mass via Bmi1 Activation in Mice. *J Bone Miner Res.* 34(6):1169-1181.
  87. Guo X, Gu X, Hareshwaree S, Rong X, Li L, Chu M. Induced pluripotent stem cell-conditional medium inhibits H9C2 cardiomyocytes apoptosis via autophagy flux and Wnt/ $\beta$ -catenin pathway. *J Cell Mol Med.* 23(6):4358-4374.
  88. Gao S, Zhang X, Song Q, Liu J, Ji X, Wang P. POLD1 deficiency is involved in cognitive function impairment in AD patients and SAMP8 mice. *Biomed Pharmacother.* 114:108833.
  89. Yao S, Tan L, Chen H, Huang X, Zhao W, Wang Y. Potential Research Tool of Stem Cells from Human Exfoliated Deciduous Teeth: Lentiviral Bmi-1 Immortalization with EGFP Marker. *Stem Cells Int.* 2019:3526409.
  90. Hu M, Hong L, Liu C, Hong S, He S, Zhou M, Huang G, Chen Q. Electrical stimulation enhances neuronal cell activity mediated by Schwann cell derived exosomes. *SCI REP-UK.* 9(1):4206.
  91. Jin D, Guo J, Wu Y, Du J, Wang X, An J, Hu B, Kong L, Di W, Wang W. UBE2C, Directly Targeted by miR-548e-5p, Increases the Cellular Growth and Invasive Abilities of Cancer Cells Interacting with the EMT Marker Protein Zinc Finger E-box Binding Homeobox 1/2 in NSCLC. *Theranostics.* 9(7):2036-2055.
  92. Wang J, Zheng X, Qin Z, Wei L, Lu Y, Peng Q, Gao Y, Zhang X, Zhang X, Li Z, Fu Y, Liu P, Liu C, Yan Q, Xiong W, Li G, Lu J, Ma J. Epstein-Barr virus miR-BART3-3p promotes tumorigenesis by regulating the senescence pathway in gastric cancer. *J Biol Chem.* 294(13):4854-4866.
  93. Chen B, Sun Y, Zhang J, Zhu Q, Yang Y, Niu X, Deng Z, Li Q, Wang Y. Human embryonic stem cell-derived exosomes promote pressure ulcer healing in aged mice by rejuvenating senescent endothelial cells. *Stem Cell Res Ther.* 10(1):142.
  94. You W, Hong Y, He H, Huang X, Tao W, Liang X, Zhang Y, Li X. TGF- $\beta$  mediates aortic smooth muscle cell senescence in Marfan syndrome. *AGING-US.* 11(11):3574-3584.
  95. Sun J, Hu X, Gao Y, Tang Q, Zhao Z, Xi W, Yang F, Zhang W, Song Y, Song B, Wang T, Wang H. MYSM1-AR complex-mediated repression of Akt/c-Raf/GSK-3 $\beta$  signaling impedes castration-resistant prostate cancer growth. *AGING-US.* 11(22):10644-10663.
  96. Huang X, Zhang H, Liang X, Hong Y, Mao M, Han Q, He H, Tao W, Jiang G, Zhang Y, Li X. Adipose-Derived Mesenchymal Stem Cells Isolated from Patients with Abdominal Aortic Aneurysm Exhibit Senescence Phenomena. *Oxid Med Cell Longev.* 2019:1305049.
  97. Wang M, Yin X, Zhang S, Mao C, Cao N, Yang X, Bian J, Hao W, Fan Q, Liu H. Autoantibodies against AT1 Receptor Contribute to Vascular Aging and Endothelial Cell Senescence. *Aging Dis.* 10(5):1012-1025.

98. Liao N, Shi Y, Zhang C, Zheng Y, Wang Y, Zhao B, Zeng Y, Liu X, Liu J. Antioxidants inhibit cell senescence and preserve stemness of adipose tissue-derived stem cells by reducing ROS generation during long-term in vitro expansion. *Stem Cell Res Ther.* 10(1):306.
99. Yao J, Shi Z, Ma X, Xu D, Ming G. lncRNA GAS5/miR-223/NAMPT axis modulates the cell proliferation and senescence of endothelial progenitor cells through PI3K/AKT signaling. *J Cell Biochem.* 120(9):14518-14530.
100. Xu D, Zeng F, Han L, Wang J, Yin Z, Lv L, Guo L, Wang D, Xu Y, Zhou H. The synergistic action of phosphate and interleukin-6 enhances senescence-associated calcification in vascular smooth muscle cells depending on p53. *Mech Ageing Dev.* 182:111124.
101. Wang Y, Zuo R, Wang Z, Luo L, Wu J, Zhang C, Liu M, Shi C, Zhou Y. Kinsenoside ameliorates intervertebral disc degeneration through the activation of AKT-ERK1/2-Nrf2 signaling pathway. *AGING-US.* 11(18):7961-7977.
102. Wang R, Liu L, Liu H, Wu K, Liu Y, Bai L, Wang Q, Qi B, Qi B, Zhang L. Reduced NRF2 expression suppresses endothelial progenitor cell function and induces senescence during aging. *AGING-US.* 11(17):7021-7035.
103. Sun X, Thorne RF, Zhang XD, He M, Li J, Feng S, Liu X, Wu M. LncRNA GUARDIN suppresses cellular senescence through a LRP130-PGC1 $\alpha$ -FOXO4-p21-dependent signaling axis. *EMBO Rep.* 21(4):e48796.
104. Guo J, Jin D, Wu Y, Yang L, Du J, Gong K, Chen W, Dai J, Miao S, Xi S. The miR 495-UBE2C-ABCG2/ERCC1 axis reverses cisplatin resistance by downregulating drug resistance genes in cisplatin-resistant non-small cell lung cancer cells. *EBioMedicine.* 2018 Sep;35:204-221.
105. Zhiqiang Chen, Xueliang Zuo, Liyong Pu, Yao Zhang, Guoyong Han, Long Zhang, Jindao Wu, Xuehao Wang. circLARP4 Induces Cellular Senescence Through Regulating miR-761/RUNX3/p53/p21 Signaling in Hepatocellular Carcinoma. *Cancer Sci.* doi: 10.1111/cas.13901.
106. Zuo R, Liu M, Wang Y, Li J, Wang W, Wu J, Sun C, Li B, Wang Z, Lan W, Zhang C, Shi C, Zhou Y. BM-MSC-derived exosomes alleviate radiation-induced bone loss by restoring the function of recipient BM-MSCs and activating Wnt/ $\beta$ -catenin signaling. *Stem Cell Res Ther.* 10(1):30.
107. T Ning, J Shao, X Zhang, X Luo, X Huang, H Wu, S Xu, B Wu, D Ma. Ageing Affects the Proliferation and Mineralization of Rat Dental Pulp Stem Cells Under Inflammatory Conditions. *Int Endod J.* doi: 10.1111/iej.13205.
108. Mengdi Yang, Zhiyuan Jiang, Guangyu Yao, Zhiyu Wang, Jing Sun, Huanlong Qin, Hui Zhao. GALC Triggers Tumorigenicity of Colorectal Cancer via Senescent Fibroblasts. *Front Oncol.* doi: 10.3389/fonc.2020.00380.
109. Liangzhi Gong, Bi Chen, Juntao Zhang, Yongjin Sun, Ji Yuan, Xin Niu, Guowen Hu, Yu Chen, Zongping Xie, Zhifeng Deng, Qing Li, Yang Wang. Human ESC-sEVs alleviate age-related bone loss by rejuvenating senescent bone marrow-derived mesenchymal stem cells. *J Extracell Vesicles.* doi: 10.1080/20013078.2020.1800971.
110. Junyan Shen, Xiaoqi Zhu, Hailiang Liu. MiR-483 induces senescence of human adipose-derived mesenchymal stem cells through IGF1 inhibition. *AGING-US.* doi: 10.18632/aging.103818.
111. Xudan Lei, Ningning He, Lihong Zhu, Manqian Zhou, Kaiyue Zhang, Chen Wang, Haoyan Huang, Shang Chen, Yuhao Li, Qiang Liu, Zhibo Han, Zhikun Guo, Zhongchao Han, Zongjin Li. Mesenchymal Stem Cell-Derived Extracellular Vesicles Attenuate Radiation-Induced Lung Injury via miRNA-214-3p. *ANTIOXID REDOX SIGN.* doi: 10.1089/ars.2019.7965.
112. Li Zhang, Xiang Tong, Jizhen Huang, Man Wu, Shijie Zhang, Dongguang Wang, Sitong Liu, Hong Pan. Fisetin Alleviated Bleomycin-Induced Pulmonary Fibrosis Partly by Rescuing Alveolar Epithelial Cells From Senescence. *Front Pharmacol.* doi: 10.3389/fphar.2020.553690.
113. Chi Zhang, Yun Xie, Haicheng Chen, Linyan Lv, Jiahui Yao, Min Zhang, Kai Xia, Xin Feng, Yanqing Li, Xiaoyan Liang, Xiangzhou Sun, Chunhua Deng, Guihua Liu. FOXO4-DRI alleviates age-related testosterone secretion insufficiency by targeting senescent Leydig cells in aged mice. *AGING-US.* doi: 10.18632/aging.102682.
114. Jian Chen, Zi-Han Sun, Li-Ying Chen, Feng Xu, Yun-Pei Zhao, Gui-Qing Li, Ming Tang, You Li, Quan-You Zheng, Shu-Feng Wang, Xin-Hua Yang, Yu-Zhang Wu, Gui-Lian Xu. C5aR deficiency attenuates the breast cancer development via the p38/p21 axis. *AGING-US.* doi: 10.18632/aging.103468.
115. Yong Wang, Shouzhen Chen, Shuna Sun, Guangyi Liu, Lipeng Chen, Yangyang Xia, Jianfeng Cui, Wenfu Wang, Xuewen Jiang, Lei Zhang, Yaofeng Zhu, Yongxin Zou, Benkang Shi. Wogonin Induces Apoptosis and Reverses Sunitinib Resistance of Renal Cell Carcinoma Cells via Inhibiting CDK4-RB Pathway. *Front Pharmacol.* doi: 10.3389/fphar.2020.01152.
116. Kai Zheng, Jingxin Ma, Yifei Wang, Zhenan He, Kejun Deng. Sulforaphane Inhibits Autophagy and Induces Exosome-Mediated Paracrine Senescence via Regulating mTOR/TFE3. *Mol Nutr Food Res.* doi: 10.1002/mnfr.201901231.
117. Manying Zhang, Lin Liang, Junyu He, Zhengxi He, Chunxue Yue, Xi Jin, Mengxiang Gao, Songshu Xiao, Yanhong Zhou. Fra-1 Inhibits Cell Growth and the Warburg Effect in Cervical Cancer Cells via STAT1 Regulation of the p53 Signaling Pathway. *Front Cell Dev Biol.* doi: 10.3389/fcell.2020.579629.
118. Yangyang Xia, Shouzhen Chen, Jianfeng Cui, Yong Wang, Xiaochen Liu, Yangli Shen, Li Gong, Xuewen Jiang, Wenfu Wang, Yaofeng Zhu, Shuna Sun, Jiangxia Li, Yongxin Zou, Benkang Shi. Berberine suppresses bladder cancer cell proliferation by inhibiting JAK1-STAT3 signaling via upregulation of miR-17-5p. *Biochem Pharmacol.* doi: 10.1016/j.bcp.2021.114575.
119. Guangyan Kan, Ziyang Wang, Chunjie Sheng, Chen Yao, Yizhi Mao, Shuai Chen. Inhibition of DKC1 induces telomere-related senescence and apoptosis in lung adenocarcinoma. *J Transl Med.* doi: 10.1186/s12967-021-02827-0.
120. Haibiao Xie, Kaifang Ma, Kenan Zhang, Jingcheng Zhou, Lei Li, Wuping Yang, Yanqing Gong, Lin Cai, Kan Gong. Cell-cycle arrest and senescence in TP53-wild type renal carcinoma by enhancer RNA-P53-bound enhancer regions 2 (p53BER2) in a p53-dependent pathway. *Cell Death Dis.* doi: 10.1038/s41419-020-03229-8.
121. Xiaoying Luo, Yangqiu Bai, Shuli He, Suofeng Sun, Xiaoke Jiang, Zhiyu Yang, Di Lu, Peiru Wei, Yuan Liang, Cong Peng, Yaru Wang, Ruli Sheng, Shuangyin Han, Xiuling Li, Bingyong Zhang. Sirtuin 1 ameliorates defenestration in hepatic sinusoidal endothelial cells during liver fibrosis via inhibiting stress-induced premature senescence. *CELL PROLIFERAT.* doi: 10.1111/cpr.12991.
122. Rongrong Le, Yixin Huang, Yanping Zhang, Hu Wang, Jiaming Lin, Yu Dong, Ziyi Li, Mingyue Guo, Xiaochen Kou, Yanhong Zhao, Mo Chen, Qianshu Zhu, Anqi Zhao, Jiqing Yin, Jiatong Sun, Zhongqu Su, Kerong Shi, Yawei Gao, Jiayu Chen, Wenqiang Liu, Lan Kang, Yixuan Wang, Chong Li, Xiaoyu Liu, Rui Gao, Hong Wang, Zhenyu Ju, Shaorong Gao. Dcaf11 activates Zscan4-mediated alternative telomere lengthening in early embryos and embryonic stem cells. *Cell Stem Cell.* doi: 10.1016/j.stem.2020.11.018.
123. Li Tao, Wei Zhang, Yaxin Zhang, Mei Zhang, Yueying Zhang, Xing Niu, Qing Zhao, Zhenxing Liu, Yuyin Li, Aipo Diao. Caffeine promotes the expression of telomerase reverse transcriptase to regulate cellular senescence and aging. *Food Funct.* doi: 10.1039/d0fo03246h.
124. Meiling Zhou, Yumei Yuan, Zishu Lin, Bin Zhang, Wenluo Qin, Yongqiang Liu, Cong Yang, Qi Wang, Hong Zhang. Acupoint catgut embedding improves senescence in a rat model of ageing by regulating mitophagy via the PINK1 pathway. *J Cell Mol Med.* doi: 10.1111/jcmm.16295.
125. Juan Wang, Yuhong Wu, Yan Li, Yamei Wang, Fangrong Shen, Jinhua Zhou, Youguo Chen. Guanosine monophosphate synthase upregulation mediates cervical cancer progression by inhibiting the apoptosis of cervical cancer cells via the Stat3/P53 pathway. *Int J Oncol.* doi: 10.3892/ijo.2021.5183.
126. Linshan Xu, Yuyang Wang, Jianping Wang, Jianglong Zhai, Li Ren, Guoying Zhu. Radiation-Induced Osteocyte Senescence Alters Bone Marrow Mesenchymal Stem Cell Differentiation Potential via Paracrine Signaling. *Int J Mol Sci.* doi: 10.3390/ijms22179323.
127. Huigen Luo, Renjie Gu, Huiya Ouyang, Lihong Wang, Shanwei Shi, Yuna Ji, Baicheng Bao, Guiqing Liao, Baoshan Xu. Cadmium exposure induces osteoporosis through cellular senescence, associated with activation of NF- $\kappa$ B pathway and mitochondrial dysfunction. *Environ Pollut.* doi: 10.1016/j.envpol.2021.118043.
128. Xuan Zhang, Ying Dong, Wan-Chen Li, Bi-Xi Tang, Jia Li, Yi Zang. Roxithromycin attenuates bleomycin-induced pulmonary fibrosis by targeting senescent cells.

- Acta Pharmacol Sin. doi: 10.1038/s41401-021-00618-3.
129. Gen-Long Bai, Ping Wang, Xin Huang, Zi-Yue Wang, Di Cao, Chuan Liu, Yi-Yi Liu, Ruo-Lin Li, Ai-Jun Chen. Rapamycin Protects Skin Fibroblasts From UVA-Induced Photoaging by Inhibition of p53 and Phosphorylated HSP27. *Front Cell Dev Biol*. doi: 10.3389/fcell.2021.633331.
130. Lei Xu, Fenbo Ma, Jun Huang, Ka Li Frankie Leung, Chenghe Qin, William Weijia Lu, X Edward Guo, Bin Tang. Metformin Hydrochloride Encapsulation by Alginate Strontium Hydrogel for Cartilage Regeneration by Reliving Cellular Senescence. *Biomacromolecules*. doi: 10.1021/acs.biomac.0c01488.
131. Panpan Chen, Qianbin Zhang, Huizhen Zhang, Yue Gao, Yanying Zhou, Yixin Chen, Lihuan Guan, Tingying Jiao, Yingyuan Zhao, Min Huang, Huichang Bi. Carnitine palmitoyltransferase 1C reverses cellular senescence of MRC-5 fibroblasts via regulating lipid accumulation and mitochondrial function. *J Cell Physiol*. doi: 10.1002/jcp.29906.
132. Huan-Tian Zhang, Tao Gui, Ri-Xu Liu, Kui-Leung Tong, Chong-Jie Wu, Zhenyan Li, Xun Huang, Qiu-Tong Xu, Jie Yang, Wang Tang, Yuan Sang, Wanting Liu, Ning Liu, Ryan D Ross, Qing-Yu He, Zhen-Gang Zha. Sequential targeting of YAP1 and p21 enhances the elimination of senescent cells induced by the BET inhibitor JQ1. *Cell Death Dis*. doi: 10.1038/s41419-021-03416-1.
133. Xiao Yang, Zhiqian Chen, Chen Chen, Chen Han, Yifan Zhou, Xunlin Li, Haijun Tian, Xiaofei Cheng, Kai Zhang, An Qin, Tangjun Zhou, Jie Zhao. Bleomycin induces fibrotic transformation of bone marrow stromal cells to treat height loss of intervertebral disc through the TGF $\beta$ 1/Smad2/3 pathway. *Stem Cell Res Ther*. doi: 10.1186/s13287-020-02093-9.
134. Qian Zhou, Sun Hye Kim, Rolando Pérez-Lorenzo, Chang Liu, Man Huang, Gian Paolo Dotto, Bin Zheng, Xunwei Wu. Phenformin Promotes Keratinocyte Differentiation via the Calcineurin/NFAT Pathway. *J Invest Dermatol*. doi: 10.1016/j.jid.2020.05.114.
135. Xianmei Pan, Bo Wu, Xianglin Fan, Guanghui Xu, Caiwen Ou, Minsheng Chen. YAP accelerates vascular senescence via blocking autophagic flux and activating mTOR. *J Cell Mol Med*. doi: 10.1111/jcmm.15902.
136. Yi Qiao Hua, Yi Zeng, Jin Xu, Xiao Le Xu. Naringenin alleviates nonalcoholic steatohepatitis in middle-aged Apoe $^{-/-}$  mice: role of SIRT1. *Phytomedicine*. doi: 10.1016/j.phymed.2020.153412.
137. Jianyang Liu, Jialin He, Lite Ge, Han Xiao, Yan Huang, Liuwang Zeng, Zheng Jiang, Ming Lu, Zhiping Hu. Hypoxic preconditioning rejuvenates mesenchymal stem cells and enhances neuroprotection following intracerebral hemorrhage via the miR-326-mediated autophagy. *Stem Cell Res Ther*. doi: 10.1186/s13287-021-02480-w.
138. Zhi Chen, Liling Xiong, Huili Jin, Jiachao Yu, Xin Li, Huijia Fu, Li Wen, Hongbo Qi, Chao Tong, Richard Saffery, Mark D Kilby, Philip N Baker. Advanced maternal age causes premature placental senescence and malformation via dysregulated  $\alpha$ -Klotho expression in trophoblasts. *Aging Cell*. doi: 10.1111/accel.13417.
139. Mengzhu Lv, Simeng Zhang, Bo Jiang, Sunrun Cao, Yuqing Dong, Liu Cao, Shu Guo. Adipose-derived stem cells regulate metabolic homeostasis and delay aging by promoting mitophagy. *FASEB J*. doi: 10.1096/fj.202100332R.
140. Wuyuan Tao, Yimei Hong, Haiwei He, Qian Han, Mengmeng Mao, Bei Hu, Hao Zhang, Xiaoran Huang, Wei You, Xiaoting Liang, Yuelin Zhang, Xin Li. MicroRNA-199a-5p aggravates angiotensin II-induced vascular smooth muscle cell senescence by targeting Sirtuin-1 in abdominal aortic aneurysm. *J Cell Mol Med*. doi: 10.1111/jcmm.16485.
141. Yanzhuo Zhang, Zhe Lv, Yudong Liu, Huan Cao, Jianwang Yang, Baoshan Wang. PIN1 Protects Hair Cells and Auditory HEI-OC1 Cells against Senescence by Inhibiting the PI3K/Akt/mTOR Pathway. *Oxid Med Cell Longev*. doi: 10.1155/2021/9980444.
142. Jie Wen, Wenqian Chen, Baoxiang Zhao, Qiuping Xu, Chang Liu, Qun Zhang, Zhiwei Xie, Yonggan Yan, Jing Guo, Jun Huang, Junying Miao, Xunwei Wu. A carbazole compound, 9-ethyl-9H-carbazole-3-carbaldehyde, plays an antitumor function through reactivation of the p53 pathway in human melanoma cells. *Cell Death Dis*. doi: 10.1038/s41419-021-03867-6.
143. Yaojie Liu, Ying Liu, Yatu Guo, Lin Xu, Hao Wang. Phlorizin exerts potent effects against aging induced by D-galactose in mice and PC12 cells. *Food Funct*. doi: 10.1039/d0fo02707c.
144. Yu-Kun Yang, Ye Li, Yan-Ying Wang, Guang-Ping Ruan, Chuan Tian, Qiang Wang, Huan-Yu He, Gao-Hong Zhu, Dong Fang, Mao Wang, Xiang-Qing Zhu, Xing-Hua Pan. The effects of BMMSC treatment on lung tissue degeneration in elderly macaques. *Stem Cell Res Ther*. doi: 10.1186/s13287-021-02201-3.
145. Qiaoxia Dai, Xia Fan, Xue Meng, Shiyue Sun, Yue Su, Xiao Ling, Xiangjuan Chen, Kai Wang, Xiaozhen Dai, Chi Zhang, Sun Da, Guigui Zhang, Chunjie Gu, Hui Chen, Junhong He, Haiqi Hu, Lechu Yu, Xiaohong Pan, Yi Tan, Xiaoqing Yan. FGF21 promotes ischaemic angiogenesis and endothelial progenitor cells function under diabetic conditions in an AMPK/NAD $^{+}$ -dependent manner. *J Cell Mol Med*. doi: 10.1111/jcmm.16369.
146. Yuanyuan Guo, Xiong Jia, Yongzhi Cui, Yu Song, Siyuan Wang, Yongtao Geng, Rui Li, Weihang Gao, Dehao Fu. Sirt3-mediated mitophagy regulates AGEs-induced BMSCs senescence and senile osteoporosis. *Redox Biol*. doi: 10.1016/j.redox.2021.101915.
147. Bowei Ni, Wenbin Pei, Yunkun Qu, Rui Zhang, Xiangyu Chu, Yingguang Wang, Xiaojian Huang, Hongbo You. MCC950, the NLRP3 Inhibitor, Protects against Cartilage Degradation in a Mouse Model of Osteoarthritis. *Oxid Med Cell Longev*. doi: 10.1155/2021/4139048.
148. Dexiang Zhang, Pinxiang Lu, Kaihua Zhu, Haixia Wu, Yuedi Dai. TFPC2 Overcomes Senescence by Cooperating With SREBP2 to Activate Cholesterol Synthesis in Pancreatic Cancer. *Front Oncol*. doi: 10.3389/fonc.2021.724437.
149. Shen Li, Song Jiang, Qingyan Zhang, Bo Jin, Daoyuan Lv, Wenju Li, Min Zhao, Chunming Jiang, Chunsun Dai, Zhihong Liu. Integrin  $\beta$ 3 Induction Promotes Tubular Cell Senescence and Kidney Fibrosis. *Front Cell Dev Biol*. doi: 10.3389/fcell.2021.733831.
150. Zheng-Zhao Liu, Chun-Gu Hong, Wen-Bao Hu, Meng-Lu Chen, Ran Duan, Hong-Ming Li, Tao Yue, Jia Cao, Zhen-Xing Wang, Chun-Yuan Chen, Xiong-Ke Hu, Ben Wu, Hao-Ming Liu, Yi-Juan Tan, Jiang-Hua Liu, Zhong-Wei Luo, Yan Zhang, Shan-Shan Rao, Ming-Jie Luo, Hao Yin, Yi-Yi Wang, Kun Xia, Lang Xu, Si-Yuan Tang, Rong-Gui Hu, Hui Xie. Autophagy receptor OPTN (optineurin) regulates mesenchymal stem cell fate and bone-fat balance during aging by clearing FABP3. *Autophagy*. doi: 10.1080/15548627.2020.1839286.
151. Mingming Sun, Hao Sheng, Tingfeng Wu, Jiaqi Song, Huanran Sun, Yingzhi Wang, Jiyan Wang, Zhen Li, Huifang Zhao, Junzhen Tan, Yanping Li, Guo Chen, Qingrong Huang, Yuan Zhang, Bei Lan, Shuangping Liu, Changliang Shan, Shuai Zhang. PIKE-A promotes glioblastoma growth by driving PPP flux through increasing G6PD expression mediated by phosphorylation of STAT3. *Biochem Pharmacol*. doi: 10.1016/j.bcp.2021.114736.
152. Liling Xiong, Xuan Ye, Zhi Chen, Huijia Fu, Sisi Li, Ping Xu, Jiachao Yu, Li Wen, Rufeí Gao, Yong Fu, Hongbo Qi, Mark D Kilby, Richard Saffery, Philip N Baker, Chao Tong. Advanced Maternal Age-associated SIRT1 Deficiency Compromises Trophoblast Epithelial-Mesenchymal Transition through an Increase in Vimentin Acetylation. *Aging Cell*. doi: 10.1111/accel.13491.
153. Yizhen Huang, Zhenlei Zhang, Jianle Wang, Shuying Shen, Teng Yao, Yining Xu, Zizheng Chen, Bin Fang, Jianjun Ma. circSPG21 protects against intervertebral disc disease by targeting miR-1197/ATP1B3. *Exp Mol Med*. doi: 10.1038/s12276-021-00674-z.
154. Shuang Qiao, Wenhua Lu, Christophe Glorieux, Jiangjiang Li, Peiting Zeng, Ning Meng, Huiqin Zhang, Shijun Wen, Peng Huang. Wild-type IDH2 protects nuclear DNA from oxidative damage and is a potential therapeutic target in colorectal cancer. *Oncogene*. doi: 10.1038/s41388-021-01968-2.
155. Dong Yuan, Jing Luo, Yixuan Sun, Lijuan Hao, Jing Zheng, Zhu Yang. PCOS follicular fluid derived exosomal miR-424-5p induces granulosa cells senescence by targeting CDCA4 expression. *Cell Signal*. doi: 10.1016/j.cellsig.2021.110030.
156. Qian Li, Yuanjun Deng, Lele Liu, Chunjiang Zhang, Yang Cai, Tianjing Zhang, Min Han, Gang Xu. Sympathetic Denervation Ameliorates Renal Fibrosis via Inhibition of Cellular Senescence. *Front Immunol*. doi: 10.3389/fimmu.2021.823935.
157. Senyu Yao, Xiaoyue Wei, Wenrui Deng, Boyan Wang, Jianye Cai, Yinong Huang,

- Xiaofan Lai, Yuan Qiu, Yi Wang, Yuanjun Guan, Jiancheng Wang. Nestin-dependent mitochondria-ER contacts define stem Leydig cell differentiation to attenuate male reproductive ageing. *Nat Commun*. doi: 10.1038/s41467-022-31755-w.
158. Bo Jiang, Xuan Wu, Fang Meng, Limiao Si, Sunrun Cao, Yuqing Dong, Huayi Sun, Mengzhu Lv, Hongde Xu, Ning Bai, Qiqiang Guo, Xiaoyu Song, Yang Yu, Wendong Guo, Fei Yi, Tingting Zhou, Xiaoman Li, Yanling Feng, Zhuo Wang, Dan Zhang, Yi Guan, Mengtao Ma, Jingwei Liu, Xining Li, Weidong Zhao, Baohua Liu, Toren Finkel, Liu Cao. Progerin modulates the IGF-1R/Akt signaling involved in aging. *Sci Adv*. doi: 10.1126/sciadv.abo0322.
159. Li Lin, Wujiang Gao, Yumei Chen, Taoqiong Li, Chunli Sha, Lu Chen, Meiling Yang, Hong Wei, Yunpeng Chen, Xiaolan Zhu. Reactive oxygen species-induced SIAH1 promotes granulosa cells' senescence in premature ovarian failure. *J Cell Mol Med*. doi: 10.1111/jcmm.17264.
160. Yuqian An, Yi Li, Yujun Hou, Shichao Huang, Gang Pei. Alzheimer's Amyloid- $\beta$  Accelerates Cell Senescence and Suppresses the SIRT1/NRF2 Pathway in Human Microglial Cells. *Oxid Med Cell Longev*. doi: 10.1155/2022/3086010.
161. Zhilei Hu, Yu Wang, Xiaoxin Gao, Yuyao Zhang, Chenhao Liu, Yu Zhai, Xian Chang, Haiyin Li, Yueyang Li, Jinhui Lou, Changqing Li. Optineurin-mediated mitophagy as a potential therapeutic target for intervertebral disc degeneration. *Front Pharmacol*. doi: 10.3389/fphar.2022.893307.
162. Jiaqian Li, Di Huang, Bingxi Lei, Jingying Huang, Linbing Yang, Man Nie, Shicheng Su, Qiyi Zhao, Ying Wang. VLA-4 suppression by senescence signals regulates meningeal immunity and leptomeningeal metastasis. *Elife*. doi: 10.7554/eLife.83272.
163. Peng Zou, Xiaoping Zhang, Rui Zhang, Xin Chai, Yuanting Zhao, Erliang Li, Qian Zhang, Rongbao Yan, Junsong Yang, Bo Liao. Blockage of ERCC6 Alleviates Spinal Cord Injury Through Weakening Apoptosis, Inflammation, Senescence, and Oxidative Stress. *Front Mol Biosci*. doi: 10.3389/fmolb.2022.853654.
164. Zhicong Zhao, Zhenjun Zhao, Jianye Wang, He Zhang, Zhifeng Xi, Qiang Xia. ABCG6 Knockdown Fuels Cell Proliferation by Regulating PPAR $\alpha$  in Hepatocellular Carcinoma. *Front Oncol*. doi: 10.3389/fonc.2022.840287.
165. Jiani Huang, Wangshu Yu, Qin He, Xiaoying He, Ming Yang, Wei Chen, Wei Han. Autophagy facilitates age-related cell apoptosis—a new insight from senile cataract. *Cell Death Dis*. doi: 10.1038/s41419-021-04489-8.
166. Yulin Huang, Jiaming Yang, Xizhe Liu, Xiaoshuai Wang, Kai Zhu, Zemin Ling, Baozhu Zeng, Ningning Chen, Shaoyu Liu, Fuxin Wei. Cationic Polymer Brush-Modified Carbon Nanotube-Meditated eRNA LINC02569 Silencing Attenuates Nucleus Pulposus Degeneration by Blocking NF- $\kappa$ B Signaling Pathway and Alleviate Cell Senescence. *Front Cell Dev Biol*. doi: 10.3389/fcell.2021.837777.
167. Qi Gao, Fang Chen, Lijun Zhang, Ai Wei, Yongxiang Wang, Zhiwei Wu, Wangsen Cao. Inhibition of DNA methyltransferase aberrations reinstates antioxidant aging suppressors and ameliorates renal aging. *Aging Cell*. doi: 10.1111/acer.13526.
168. Yang Huang, Ying Liu, Siwei Yu, Wenzhe Li, Jinglan Li, Bo Zhao, Xin Hu, Haiying Jin. Biliverdin Reductase A Protects Lens Epithelial Cells against Oxidative Damage and Cellular Senescence in Age-Related Cataract. *Oxid Med Cell Longev*. doi: 10.1155/2022/5628946.
169. Changjian Lin, Lujie Ge, Luping Tang, Yuzhe He, Safwat Adel Abdo Moqbel, Kai Xu, Diana Ma, Xing Zhou, Jisheng Ran, Lidong Wu. Nitidine Chloride Alleviates Inflammation and Cellular Senescence in Murine Osteoarthritis Through Scavenging ROS. *Front Pharmacol*. doi: 10.3389/fphar.2022.919940.
170. Jiao Zhou, Honghan Chen, Jintao Du, Haoran Tai, Xiaojuan Han, Ning Huang, Xiaobo Wang, Hui Gong, Mingyao Yang, Hengyi Xiao. Glutamine Availability Regulates the Development of Aging Mediated by mTOR Signaling and Autophagy. *Front Pharmacol*. doi: 10.3389/fphar.2022.924081.
171. Wenxian Zhou, Yifeng Shi, Hui Wang, Linjie Chen, Caiyu Yu, Xufei Zhang, Lei Yang, Xiaolei Zhang, Aimin Wu. Exercise-induced FNDC5/irisin protects nucleus pulposus cells against senescence and apoptosis by activating autophagy. *Exp Mol Med*. doi: 10.1038/s12276-022-00811-2.
172. Yanqi Zhong, Yang Zhang, Weifang Liu, Yin Zhao, Li Zou, Xiaoxia Liu. TLR4 Modulates Senescence and Paracrine Action in Placental Mesenchymal Stem Cells via Inhibiting Hedgehog Signaling Pathway in Preeclampsia. *Oxid Med Cell Longev*. doi: 10.1155/2022/7202837.
173. Yi Li, Juan Lu, Yujun Hou, Shichao Huang, Gang Pei. Alzheimer's Amyloid- $\beta$  Accelerates Human Neuronal Cell Senescence Which Could Be Rescued by Sirtuin-1 and Aspirin. *Front Cell Neurosci*. doi: 10.3389/fncel.2022.906270.
174. Chao-Yue Sun, Milton Talukder, Di Cao, Cun-Wu Chen. Gilteritinib Enhances Anti-Tumor Efficacy of CDK4/6 Inhibitor, Abemaciclib in Lung Cancer Cells. *Front Pharmacol*. doi: 10.3389/fphar.2022.829759.
175. Yu Cao, Meng Zhang, Ye Li, Jingjing Lu, Wanhui Zhou, Xiaoshuang Li, Hao Shi, Bin Xu, Shize Li. O-GlcNAcylation of SIRT1 Protects against Cold Stress-Induced Skeletal Muscle Damage via Amelioration of Mitochondrial Homeostasis. *Int J Mol Sci*. doi: 10.3390/ijms232314520.
176. Nuan Lin, Jiazhe Lin, Torsten Plosch, Pingnan Sun, Xiaoling Zhou. An Oxidative Stress-Related Gene Signature in Granulosa Cells Is Associated with Ovarian Aging. *Oxid Med Cell Longev*. doi: 10.1155/2022/1070968.
177. Chaohua Kong, Dongchen Wang, Feng Wang, Yifei Lv, Wenying Zhou, Peng Ye, Yue Gu, Xiaomin Jiang, Linlin Zhu, Zhen Ge, Yuelin Chao, Shaoliang Chen. PCAF Accelerates Vascular Senescence via the Hippo Signaling Pathway. *Oxid Med Cell Longev*. doi: 10.1155/2022/1200602.
178. Hongyi Chen, Luo Cheng Lv, Ruoxu Liang, Weimin Guo, Zhaofu Liao, Yilin Chen, Kuikui Zhu, Ruijin Huang, Hui Zhao, Qin Pu, Ziqiang Yuan, Zhaohua Zeng, Xin Zheng, Shanshan Feng, Xufeng Qi, Dongqing Cai. miR-486 improves fibrotic activity in myocardial infarction by targeting SRSF3/p21-Mediated cardiac myofibroblast senescence. *J Cell Mol Med*. doi: 10.1111/jcmm.17539.
179. Xiaoman Zhu, Zhiyang Chen, Li Wang, Qingjian Ou, Zhong Feng, Honglei Xiao, Qi Shen, Yingao Li, Caixia Jin, Jing-Ying Xu, Furong Gao, Juan Wang, Jingfa Zhang, Jieping Zhang, Zhiguo Xu, Guo-Tong Xu, Lixia Lu, Haibin Tian. Direct conversion of human umbilical cord mesenchymal stem cells into retinal pigment epithelial cells for treatment of retinal degeneration. *Cell Death Dis*. doi: 10.1038/s41419-022-05199-5.
180. Meng-Yun Cai, Chung-Yin Yip, Kewu Pan, Yan Zhang, Renee Wan-Yi Chan, Wood Yee Chan, Wing-Hung Ko. Role of Carbon Monoxide in Oxidative Stress-Induced Senescence in Human Bronchial Epithelium. *Oxid Med Cell Longev*. doi: 10.1155/2022/5199572.
181. Jiyang Hou, Ruoxiang Ma, Shisheng Zhu, Yaping Wang. Revealing the Therapeutic Targets and Mechanism of Ginsenoside Rg1 for Liver Damage Related to Anti-Oxidative Stress Using Proteomic Analysis. *Int J Mol Sci*. doi: 10.3390/ijms231710045.

注：更多使用本产品的文献请参考产品网页。

Version 2024.03.12