




碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

DiO (细胞膜绿色荧光探针)

产品编号	产品名称	包装
C1038	DiO (细胞膜绿色荧光探针)	10mg

产品简介:

- DiO即DiOC18(3), 全称为3,3'-dioctadecyloxacarbocyanine perchlorate, 是最常用的细胞膜荧光探针之一, 呈现绿色荧光。DiO是一种亲脂性膜染料, 进入细胞膜后可以侧向扩散逐渐使整个细胞的细胞膜被染色。
- DiO在进入细胞膜之前荧光非常弱, 仅当进入到细胞膜后才可以被激发出很强的荧光。DiO被激发后可以发出绿色的荧光, DiO和磷脂双层膜结合后的激发光谱和发射光谱参考下图。其中, 最大激发波长为484nm, 最大发射波长为501nm。



- DiO的分子式为 $C_{53}H_{85}ClN_2O_6$, 分子量为881.72, CAS number为34215-57-1。
- DiO可以溶解于无水乙醇、DMSO和DMF, 其中在DMSO溶解度大于为10mg/ml。发现较难溶解时可以适当加热, 并用超声处理以促进溶解。
- DiO 被广泛用于正向或逆向的, 活的或固定的神经等细胞或组织的示踪剂或长期示踪剂(long-term tracer)。DiO通常不会明显影响细胞的生存力(viability)。DiO对于细胞膜染色的荧光强度通常要低于DiI, 有时对于某些经过固定的组织的染色效果欠佳。
- DiO除了最简单的细胞膜荧光标记外, 还可以用于检测细胞的融合和粘附, 检测发育或移植过程中细胞迁移, 通过FRAP(Fluorescence Recovery After Photobleaching)检测脂在细胞膜上的扩散, 检测细胞毒性和标记脂蛋白等。
- 用于细胞膜荧光标记时, DiO的常用浓度为1-30 μ M, 最常用的浓度为5-10 μ M。DiO可以直接染色活的细胞或组织, 染色时间通常为5-20分钟。对于固定的细胞或组织, 通常宜使用配制在PBS中的4%多聚甲醛进行固定, 使用其它不适当的固定液会导致荧光背景较高。

包装清单:

产品编号	产品名称	包装
C1038	DiO (细胞膜绿色荧光探针)	10mg
—	说明书	1份

保存条件:

4°C避光保存, 一年有效。配制的储存液-20°C避光保存, 半年有效。

注意事项:

- 荧光染料均存在淬灭问题, 请尽量注意避光, 以减缓荧光淬灭。
- 本产品仅限于专业人员的科学的研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用本产品的文献:

1. Ma Y, Zhuang Y, Xie X, Wang C, Wang F, Zhou D, Zeng J, Cai L. The role of surface charge density in cationic liposome-promoted dendritic cell maturation and vaccine-induced immune responses. *Nanoscale*. 2011 May;3(5):2307-14.
2. Jiang QY, Lai LH, Shen J, Wang QQ, Xu FJ, Tang GP. Gene delivery to

- tumor cells by cationic polymeric nanovectors coupled to folic acid and the cell-penetrating peptide octaarginine. *Biomaterials*. 2011 Oct;32(29):7253-62.
3. Shi YY, Wang YS, Zhang ZX, Cai Y, Zhou J, Hou HY, van Rooijen N. Monocyte/macrophages promote vasculogenesis in choroidal neovascularization in mice by stimulating SDF-1 expression in RPE cells. *Graefes Arch Clin Exp Ophthalmol*. 2011 Nov;249(11):1667-79.
 4. Zhao W, Zhang WP, Zhang ZL, He RL, Lin Y, Xie M, Wang HZ, Pang DW. Robust and highly sensitive fluorescence approach for point-of-care virus detection based on immunomagnetic separation. *Anal Chem*. 2012 Mar 6;84(5):2358-65.
 5. Xiong XY, Guo L, Gong YC, Li ZL, Li YP, Liu ZY, Zhou M. In vitro & in vivo targeting behaviors of biotinylated Pluronic F127/poly(lactic acid) nanoparticles through biotin-avidin interaction. *Eur J Pharm Sci*. 2012 Aug 15;46(5):537-44.
 6. Zhang Z, Xiong X, Wan J, Xiao L, Gan L, Feng Y, Xu H, Yang X. Cellular uptake and intracellular trafficking of PEG-b-PLA polymeric micelles. *Biomaterials*. 2012 Oct;33(29):7233-40.
 7. Zhang W, Dang S, Hong T, Tang J, Fan J, Bu D, Sun Y, Wang Z, Wisniewski T. A humanized single-chain antibody against beta 3 integrin inhibits pulmonary metastasis by preferentially fragmenting activated platelets in the tumor microenvironment. *Blood*. 2012 Oct 4;120(14):2889-98.
 8. Zhuang Y, Ma Y, Wang C, Hai L, Yan C, Zhang Y, Liu F, Cai L. PEGylated cationic liposomes robustly augment vaccine-induced immune responses: Role of lymphatic trafficking and biodistribution. *J Control Release*. 2012 Apr 10;159(1):135-42.
 9. Wu Y, Feng W, Zhang H, Li S, Wang D, Pan X, Hu S. Ca²⁺-regulatory proteins in cardiomyocytes from the right ventricle in children with congenital heart disease. *J Transl Med*. 2012 Apr 2;10:67.
 10. Wang ZG, Liu SL, Tian ZQ, Zhang ZL, Tang HW, Pang DW. Myosin-driven intercellular transportation of wheat germ agglutinin mediated by membrane nanotubes between human lung cancer cells. *ACS Nano*. 2012 Nov 27;6(11):10033-41.
 11. Song W, Wu K, Yan J, Zhang Y, Zhao L. MiR-148b laden titanium implant promoting osteogenic differentiation of rat bone marrow mesenchymal stem cells. *RSC Adv*. 2013 Apr;28:11292-300.
 12. Li P, Gao Y, Zhang J, Liu Z, Tan K, Hua X, Gong J. Renal interstitial permeability changes induced by microbubble-enhanced diagnostic ultrasound. *J Drug Target*. 2013 May;21(5):507-14.
 13. Luo J, Zhao X, Tan Z, Su Z, Meng F, Zhang M. Mesenchymal-like progenitors derived from human embryonic stem cells promote recovery from acute kidney injury via paracrine actions. *Cyotherapy*. 2013 Jun;15(6):649-62.
 14. Wu K, Xu J, Liu M, Song W, Yan J, Gao S, Zhao L, Zhang Y. Induction of osteogenic differentiation of stem cells via a lyophilized microRNA reverse transfection formulation on a tissue culture plate. *Int J Nanomedicine*. 2013;8:1595-607.
 15. Zhang B, Luo Z, Liu J, Ding X, Li J, Cai K. Cytochrome c end-capped mesoporous silica nanoparticles as redox-responsive drug delivery vehicles for liver tumor-targeted triple therapy in vitro and in vivo. *J Control Release*. 2014 Oct 28;192:192-201.
 16. He M, Wang Z, Cao Y, Zhao Y, Duan B, Chen Y, Xu M, Zhang L. Construction of Chitin/PVA Composite Hydrogels with Jellyfish Gel-Like Structure and Their Biocompatibility. *Biomacromolecules*. 2014 Sep 8;15(9):3358-65.
 17. Wang W, Yuan C, Wang S, Song X, Xu L, Yan R, Hasson IA, Li X. Transcriptional and proteomic analysis reveal recombinant galectins of Haemonchus contortus down-regulated functions of goat PBMC and modulation of several signaling cascades in vitro. *J Proteomics*. 2014 Feb 26;98:123-37.
 18. He M, Zhao Y, Duan J, Wang Z, Chen Y, Zhang L. Fast contact of solid-liquid interface created high strength multi-layered cellulose hydrogels with controllable size. *ACS Appl Mater Interfaces*. 2014 Feb 12;6(3):1872-8.
 19. Wu M, Zhang ZL, Chen G, Wen CY, Wu LL, Hu J, Xiong CC, Chen JJ, Pang DW. Rapid and Quantitative Detection of Avian Influenza A(H7N9) Virions in Complex Matrices Based on Combined Magnetic Capture and Quantum Dot Labeling. *Small*. 2015 Oct;11(39):5280-8.
 20. Wang X, Lan H, Li J, Su Y, Xu L. Muc1 promotes migration and lung metastasis of melanoma cells. *Am J Cancer Res*. 2015 Aug 15;5(9):2590-604.
 21. Zheng Y, Yin G, Le V, Zhang A, Chen S, Liang X, Liu J. Photodynamic-therapy Activates Immune Response by disrupting Immunity Homeostasis of Tumor Cells, which Generates Vaccine for Cancer Therapy. *Int J Biol Sci*. 2016 Jan 1;12(1):120-32.
 22. Zhou W, Zeng C, Liu R, Chen J, Li R, Wang X, Bai W, Liu X, Xiang T, Zhang L, Wan Y. Antiviral activity and specific modes of action of bacterial prodigiosin against *Bombyx mori* nucleopolyhedrovirus in vitro. *Appl Microbiol Biotechnol*. 2016 May;100(9):3979-88.
 23. Jiang Q, Hao S, Xiao X, Yao J, Ou B, Zhao Z, Liu F, Pan X, Luo B, Zhi H. Production and characterization of a novel long-acting Herceptin-targeted nanobubble contrast agentspecific for Her-2-positive breast cancers. *Breast Cancer*. 2016 May;23(3):445-55.
 24. Xu H, Li Z, Yu Y, Sizdahkhani S, Ho WS, Yin F, Wang L, Zhu G, Zhang M, Jiang L, Zhuang Z, Qin J . A dynamic in vivo-like organotypic blood-brain barrier model to probe metastatic brain tumors. *Sci Rep*. 2016 Nov 10;6:36670.
 25. Qiao H, Hu N, Bai J, Ren L, Liu Q, Fang L, Wang Z . Encapsulation of Nucleic Acids into Giant Unilamellar Vesicles by Freeze-Thaw: a Way Protocells MayForm. *Orig Life Evol Biosph*. 2016 Nov 2. [Epub ahead of print]
 26. Qi L, Yu Y, Chi X, Lu D, Song Y, Zhang Y, Zhang H . Depletion of Kindlin-2 induces cardiac dysfunction in mice. *Sci China Life Sci*. 2016 Nov;59(11):1123-1130.
 27. Zhong Q, Li W, Su X, Li G, Zhou Y, Kundu SC, Yao J, Cai Y. Degradation pattern of porous CaCO₃ and hydroxyapatite microspheres in vitro and in vivo for potential application in bone tissue engineering. *Colloids Surf B Biointerfaces*. 2016 Jul 1;143:56-63
 28. Kong X, Chen L, Ye P, Wang Z, Zhang J, Ye F, Chen S . The role of HYAL2 in LSS-induced glycocalyx impairment and the PKA-mediated decrease in eNOS-Ser-633 phosphorylation and nitric oxide production. *Mol Biol Cell*. 2016 Dec 15;27(25):3972-3979.
 29. Li J, Chen H, Xu Y, Hu J, Xie FQ, Yang C. Integrin endocytosis on elastic substrates mediates mechanosensing. *J Biomech*. 2016 Sep 6;49(13):2644-2654.
 30. Zhou W, Zeng C, Liu R, Chen J, Li R, Wang X, Bai W, Liu X, Xiang T, Zhang L, Wan Y . Antiviral activity and specific modes of action of bacterial prodigiosin against *Bombyx mori* nucleopolyhedrovirus in vitro. *Appl Microbiol Biotechnol*. 2016 May;100(9):3979-88.
 31. Li Z, Guo Y, Yu Y, Xu C, Xu H, Qin J . Assessment of metabolism-dependent drug efficacy and toxicity on a multilayer organs-on-a-chip. *Integr Biol (Camb)*. 2016 Oct 10;8(10):1022-1029.
 32. Chen L, Wu LL, Zhang ZL, Hu J, Tang M, Qi CB, Li N, Pang DW . Biofunctionalized magnetic nanospheres-based cell sorting strategy for efficient isolation, detection and subtype analyses of heterogeneous circulating hepatocellular carcinoma cells. *Biosens Bioelectron*. 2016 Nov 15;85:633-40.
 33. Yu J, Zhang J, Xing H, Sun Y, Yang Z, Yang T, Cai C, Zhao X, Yang L, Ding P. Novel guanidinylated bioresponsive poly(amidoamine)s

- designed for short hairpin RNA delivery. *Int J Nanomedicine*. 2016 Dec 8;11:6651-6666.
34. Yu J, Zhang J, Xing H, Yang Z, Cai C, Zhang C, Zhao X, Wei M, Yang L, Ding P. Guanidinylated bioresponsive poly(amido amine)s designed for intranuclear gene delivery. *Int J Nanomedicine*. 2016 Aug 17;11:4011-24.
 35. Zheng Y, Yin G, Le V, Zhang A, Chen S, Liang X, Liu J. Photodynamic-therapy Activates Immune Response by disrupting Immunity Homeostasis of Tumor Cells, which Generates Vaccine for Cancer Therapy. *Int J Biol Sci*. 2016 Jan 1;12(1):120-32.
 36. Fan H, Li Y, Yang J, Ye X. Effect of Hydrophobic Chain Length on the Stability and Guest Exchange Behavior of Shell-Sheddable Micelles Formed by Disulfide-Linked Diblock Copolymers. *J Phys Chem B*. 2017 Oct 19;121(41):9708-9717
 37. Qiao H, Hu N, Bai J, Ren L, Liu Q, Fang L, Wang Z. Encapsulation of Nucleic Acids into Giant Unilamellar Vesicles by Freeze-Thaw: a Way Protocells May Form. *ORIGINS LIFE EVOL B*. 2017 Dec;47(4):499-510
 38. Xia D, He Y, Li Q, Hu C, Huang W, Zhang Y, Wan F, Wang C, Gan Y. Transport mechanism of lipid covered saquinavir pure drug nanoparticles in intestinal epithelium. *J Control Release*. 2018 Jan 10;269:159-170
 39. Liu C, Zhou Z, Chen Y, Liu J, Wang Y, Liu H. Targeted delivery of garcinia glycosides by reconstituted high-density lipoprotein nano-complexes. *J Microencapsul*. 2018 Mar;35(2):115-120
 40. Wang H, Shen X, Tian G, Shi X, Huang W, Wu Y, Sun L, Peng C, Liu S, Huang Y, Chen X, Zhang F, Chen Y, Ding W, Lu Z. AMPKa2 deficiency exacerbates long-term PM2.5 exposure-induced lung injury and cardiac dysfunction. *FREE RADICAL BIO MED*. 2018 Jun;121:202-214
 41. Wang Y, Lei T, Yuan J, Wu Y, Shen X, Gao J, Feng W, Lu Z. GCN2 deficiency ameliorates doxorubicin-induced cardiotoxicity by decreasing cardiomyocyte apoptosis and myocardial oxidative stress. *Redox Biol*. 2018 Jul;17:25-34
 42. Jin Z, Wang P, Chen J, He L, Xiao L, Yong K, Deng S, Zhou L. A Tumor-Specific Tissue-Penetrating Peptide Enhances the Efficacy of Chemotherapy Drugs in Gastric Cancer. *Yonsei Med J*. 2018 Jul;59(5):595-601
 43. Tengda L, Shuping L, Mingli G, Jie G, Yun L, Weiwei Z, Anmei D. Serum exosomal microRNAs as potent circulating biomarkers for melanoma. *Melanoma Res*. 2018 Aug;28(4):295-303
 44. Zhang Z, Qian H, Huang J, Sha H, Zhang H, Yu L, Liu B, Hua D, Qian X. Anti-EGFR-iRGD recombinant protein modified biomimetic nanoparticles loaded with gambogic acid to enhance targeting and antitumor ability in colorectal cancer treatment. *INT J NANOMED*. 2018 Aug 31;13:4961-4975
 45. Chen H, Sha H, Zhang L, Qian H, Chen F, Ding N, Ji L, Zhu A, Xu Q, Meng F, Yu L, Zhou Y, Liu B. Lipid insertion enables targeted functionalization of paclitaxel-loaded erythrocyte membrane nanosystem by tumor-penetrating bispecific recombinant protein. *INT J NANOMED*. 2018 Sep 11;13:5347-5359
 46. Li M, Liu Y, Chen J, Liu T, Gu Z, Zhang J, Gu X, Teng G, Yang F, Gu N. Platelet bio-nanobubbles as microvascular recanalization nanoformulation for acute ischemic stroke lesion theranostics. *Theranostics*. 2018 Sep 9;8(18):4870-4883
 47. Liu J, Zhang F, Zhu L, Qu X, Chu D. Laser interference-based technique for dynamic measurement of single cell deformation manipulated by optical tweezers. *Electrophoresis*. 2018 Oct 10
 48. Wang GH, Ma KL, Zhang Y, Hu ZB, Liu L, Lu J, Chen PP, Lu CC, Ruan XZ, Liu BC. Platelet microparticles contribute to aortic vascular endothelial injury in diabetes via the mTORC1 pathway. *Acta Pharmacol Sin*. 2018 Nov 16
 49. Yang X, Shi G, Guo J, Wang C, He Y. Exosome-encapsulated antibiotic against intracellular infections of methicillin-resistant *Staphylococcus aureus*. *INT J NANOMED*. 2018 Nov 29;13:8095-8104
 50. Wang Y, Yin P, Chen YH, Yu YS, Ye WX, Huang HY, Ji ZC, Shen ZY. A functional variant of SMAD4 enhances macrophage recruitment and inflammatory response via TGF-β signal activation in Thoracic aortic aneurysm and dissection. *AGING-US*. 2018 Dec 7;10(12):3683-3701
 51. Li Y, Huang K, Liu L, Qu Y, Huang Y, Wu Y, Wei J. Effects of complement and serum IgG on rituximab-dependent natural killer cell-mediated cytotoxicity against Raji cells. *Oncol Lett*. 2019 Jan;17(1):339-347
 52. Zuo R, Liu M, Wang Y, Li J, Wang W, Wu J, Sun C, Li B, Wang Z, Lan W, Zhang C, Shi C, Zhou Y. BM-MSC-derived exosomes alleviate radiation-induced bone loss by restoring the function of recipient BM-MSCs and activating Wnt/β-catenin signaling. *Stem Cell Res Ther*. 2019 Jan 15;10(1):30
 53. Wu X, Dai H, Liu L, Xu C, Yin Y, Yi J, Bielec MD, Han Y, Li S. Citrate reduced oxidative damage in stem cells by regulating cellular redox signaling pathways and represent a potential treatment for oxidative stress-induced diseases. *Redox Biol*. 2019 Feb;21:101057
 54. Zheng Y, Tu C, Zhang J, Wang J. Inhibition of multiple myeloma-derived exosomes uptake suppresses the functional response in bone marrow stromal cell. *Int J Oncol*. 2019 Mar;54(3):1061-1070
 55. Wang Y, Zhang K, Qin X, Li T, Qiu J, Yin T, Huang J, McGinty S, Pontrelli G, Ren J, Wang Q, Wu W, Wang G. Biomimetic Nanotherapies: Red Blood Cell Based Core-Shell Structured Nanocomplexes for Atherosclerosis Management. *Adv Sci (Weinh)*. 2019 Apr 24;6(12):1900172.
 56. Zhuang M, Xin G, Wei Z, Li S, Xing Z, Ji C, Du J, Niu H, Huang W. Dihydrodiosgenin inhibits endothelial cell-derived factor VIII and platelet-mediated hepatocellular carcinoma metastasis. *Cancer Manag Res*. 2019 May 31;11:4871-4882
 57. Zhang J, Kong X, Wang Z, Gao X, Ge Z, Gu Y, Ye P, Chao Y, Zhu L, Li X, Chen S. AMP-activated protein kinase regulates glycocalyx impairment and macrophage recruitment in response to low shear stress. *FASEB J*. 2019 Jun;33(6):7202-7212
 58. Liu P, Chen N, Yan L, Gao F, Ji D, Zhang S, Zhang L, Li Y, Xiao Y. Preparation, characterisation and in vitro and in vivo evaluation of CD44-targeted chondroitin sulphate-conjugated doxorubicin PLGA nanoparticles. *CARBOHYD POLYM*. 2019 Jun 1;213:17-26
 59. Zhu X, Huang L, Zheng Y, Song Y, Xu Q, Wang J, Si K, Duan S, Gong W. Ultrafast optical clearing method for three-dimensional imaging with cellular resolution. *P NATL ACAD SCI USA*. 2019 Jun 4;116(23):11480-11489
 60. Xu C, Wang J, Yang J, Lei C, Hu J, Sun X. NSP2 forms viroplasms during *Dendrolimus punctatus* cytopivirus infection. *Virology*. 2019 Jul;533:68-76
 61. Li Y, Hao L, Liu F, Yin L, Yan S, Zhao H, Ding X, Guo Y, Cao Y, Li P, Wang Z, Ran H, Sun Y. Cell penetrating peptide-modified nanoparticles for tumor targeted imaging and synergistic effect of sonodynamic/HIFU therapy. *INT J NANOMED*. 2019 Jul 29;14:5875-5894
 62. Dong D, Zhu L, Hu J, Pang DW, Zhang ZL. Simple and rapid extracellular vesicles quantification via membrane biotinylation strategy coupled with fluorescent nanospheres-based lateral flow assay. *Talanta*. 2019 Aug 1;200:408-414
 63. Zhang Y, Deng H, Hu Y, Pan C, Wu G, Li Q, Tang Z. Adipose-derived mesenchymal stem cells stereotactic transplantation alleviate brain edema from intracerebral hemorrhage. *J Cell Biochem*. 2019 Sep;120(9):14372-14382
 64. Liu P, Shen H, Zhi Y, Si J, Shi J, Guo L, Shen SG. 3D bioprinting and in vitro study of bilayered membranous construct with human cells-laden alginate/gelatin composite hydrogels. *COLLOID SURFACE B*. 2019 Sep

1;181:1026-1034

65. Wang R,Yao X,Li T,Li X,Jin M,Ni Y,Yuan W,Xie X,Lu L,Li M.Reversible Thermoresponsive Hydrogel Fabricated from Natural Biopolymer for the Improvement of Critical Limb Ischemia by Controlling Release of Stem Cells.*Adv Healthc Mater.* 2019 Oct;8(20):e1900967
66. Pan Q,Zhang J,Li X,Han X,Zou Q,Zhang P,Luo Y,Jin Y.Preparation and pharmacokinetics of bifunctional epirubicin-loaded micelles.*Pharmazie.* 2019 Oct 1;74(10):577-582
67. Wang F,Zuo Z,Yang Z,Chen K,Fang J,Cui H,Shu G,Zhou Y,Geng Y,Ouyang P.Delayed Pulmonary Apoptosis of Diet-Induced Obesity Mice following Escherichia coli Infection through the Mitochondrial Apoptotic Pathway.*Oxid Med Cell Longev.* 2019 Oct 22;2019:1968539
68. Bai L,Li J,Li H,Song J,Zhou Y,Lu R,Liu B,Pang Y,Zhang P,Chen J,Liu X,Wu J,Liang C,Zhou J.Renoprotective effects of artemisinin and hydroxychloroquine combination therapy on IgA nephropathy via suppressing NF- κ B signaling and NLRP3 inflammasome activation by exosomes in rats.*Biochem Pharmacol.* 2019 Nov;169:113619
69. Yang W,Li Z,Qin R,Wang X,An H,Wang Y,Zhu Y,Liu Y,Cai S,Chen S,Sun T,Meng J,Yang C.YY1 Promotes Endothelial Cell-Dependent Tumor Angiogenesis in Hepatocellular Carcinoma by Transcriptionally Activating VEGFA.*Front Oncol.* 2019 Nov 14;9:1187
70. Zhang C,Xiong X,Li Y,Huang K,Liu L,Peng X,Weng W.Cytokine-induced killer cells/natural killer cells combined with anti-GD2 monoclonal antibody increase cell death rate in neuroblastoma SK-N-SH cells.*Oncol Lett.* 2019 Dec;18(6):6525-6535
71. Li L,Han B,Wang Y,Zhao J,Cao Y.Simple and universal signal labeling of cell surface for amplified detection of cancer cells via mild reduction.*Biosens Bioelectron.* 2019 Dec 1;145:111714
72. Cen J,Feng L,Ke H,Bao L,Li LZ,Tanaka Y,Weng J,Su L.Exosomal Thrombospondin-1 Disrupts the Integrity of Endothelial Intercellular Junctions to Facilitate Breast Cancer Cell Metastasis.*Cancers (Basel).* 2019 Dec 5;11(12). pii: E1946.
73. Fan JH,Fan GL,Yuan P,Deng FA,Liu LS,Zhou X,Yu XY,Cheng H,Li SY.A Theranostic Nanoprobe for Hypoxia Imaging and Photodynamic Tumor Therapy.*Front Chem.* 2019 Dec 20;7:868
74. Yang L,Zhai Y,Hao Y,Zhu Z,Cheng G.The Regulatory Functionality of Exosomes Derived from hUMSCs in 3D Culture for Alzheimer's Disease Therapy.*Small.* 2020 Jan;16(3):e1906273
75. Dayong Chen,Hai Qiao,Yiting Wang,Ling Zhou,Na Yin,Liaoqiong Fang,Zhibiao Wang.Adenomyosis-derived extracellular vesicles endow endometrial epithelial cells with an invasive phenotype through epithelial-mesenchymal transition.*Genes & Diseases.* 2020 Jan 23;7(4):636-648.;doi: 10.1016/j.gendis.2020.01.011.
76. Zhiping Yu,Yixuan Wang,Dan Xu,Lianhua Zhu,Ming Hu,Qiuli Liu,Weihua Lan,Jun Jiang,Luofu Wang.G250 Antigen-Targeting Drug-Loaded Nanobubbles Combined with Ultrasound Targeted Nanobubble Destruction: A Potential Novel Treatment for Renal Cell Carcinoma.*INT J NANOMED.* 2020 Jan 8;15:81-95.;doi: 10.2147/IJN.S230879
77. Xiaozhen Liu,Tao Jiang,Xuefei Li,Chao Zhao,JIayu Li,Fei Zhou,Limin Zhang,Sha Zhao,Yijun Jia,Jinpeng Shi,Guanghui Gao,Wei Li,Jing Zhao,Xiaoxia Chen,Chunxia Su,Shengxiang Ren,Caicun Zhou.Exosomes transmit T790M mutation-induced resistance in EGFR-mutant NSCLC by activating PI3K/AKT signalling pathway.*J Cell Mol Med.* 2020 Jan;24(2):1529-1540.;doi: 10.1111/jcmm.14838
78. Junling Gao,Juntao Yuan,Qiao'e Wang,Tong Lei,Xiyue Shen,Bingqing Cui,Fang Zhang,Wenjun Ding,Zhongbing Lu.Metformin protects against PM 2.5-induced lung injury and cardiac dysfunction independent of AMP-activated protein kinase α 2 Redox Biol. 2020 Jan;28:101345.;doi: 10.1016/j.redox.2019.101345
79. Cai X,Yang Q,Weng Q,Wang S.pH sensitive doxorubicin-loaded nanoparticle based on Radix pseudostellariae protein-polysaccharide conjugate and its improvement on HepG2 cellular uptake of doxorubicin.*Food Chem Toxicol.* 2020 Feb;136:111099
80. Lin M,Guo W,Zhang Z,Zhou Y,Chen J,Wang T,Zhong X,Lu Y,Yang Q,Wei Q,Han M,Xu D,Gao J.Reduced Toxicity of Liposomal Nitrogen Mustard Prodrug Formulation Activated by an Intracellular ROS Feedback Mechanism in Hematological Neoplasm Models.*MOL PHARMACOL.* 2020 Feb 3;17(2):499-506
81. Peitao Zhou,Qianqin Li,Shuwen Su,Wenhui Dong,Suyu Zong,Qiong Ma,Xi Yang,Daming Zuo,Shaoyi Zheng,Xianzhong Meng,Dingli Xu,Qingchun Zeng.Interleukin 37 Suppresses M1 Macrophage Polarization Through Inhibition of the Notch1 and Nuclear Factor Kappa B Pathways.*Front Cell Dev Biol.* 2020 Feb 14;8:56.;doi: 10.3389/fcell.2020.00056
82. Feng Liu,Junlin Sun,Wenqian Yu,Qunying Jiang,Min Pan,Zhen Xu,Fengye Mo,Xiaoqing Liu.Quantum dot-pulsed dendritic cell vaccines plus macrophage polarization for amplified cancer immunotherapy.*Biomaterials.* 2020 Feb 29;242:119928.;doi: 10.1016/j.biomaterials.2020.119928
83. Xixi Cai,Qian Yang,Qingxia Weng,Shaoyun Wang,pH sensitive doxorubicin-loaded nanoparticle based on Radix pseudostellariae protein-polysaccharide conjugate and its improvement on HepG2 cellular uptake of doxorubicin.*Food Chem Toxicol.* 2020 Feb;136:111099.;doi: 10.1016/j.fct.2019.111099
84. Di Wang,Changfu Hao,Lin Zhang,Jianhui Zhang,Suna Liu,Yiping Li,Yaqian Qu,Youliang Zhao,Ruojuan Huang,Jingjing Wei,Wu Yao.Exosomal miR-125a-5p derived from silica-exposed macrophages induces fibroblast transdifferentiation.*ECOTOX ENVIRON SAFE.* 2020 Apr 1;192:110253.;doi: 10.1016/j.ecoenv.2020.110253
85. Yi Liu,Jingtong Zhao,Jinlan Jiang,Fangfang Chen,Xuedong Fang.Doxorubicin Delivered Using Nanoparticles Camouflaged with Mesenchymal Stem Cell Membranes to Treat Colon Cancer.*INT J NANOMED.* 2020 Apr 23;15:2873-2884.;doi: 10.2147/IJN.S242787
86. Xin Chang,Lei Xing,Yi Wang,Tian-Jiao Zhou,Li-Jun Shen,Hu-Lin Jiang.Nanoengineered immunosuppressive therapeutics modulating M1/M2 macrophages into the balanced status for enhanced idiopathic pulmonary fibrosis therapy.*Nanoscale.* 2020 Apr 30;12(16):8664-8678.;doi: 10.1039/d0nr00750a
87. Gen-Que Fu,Sheng-Chun Zhang,Ge-Gu Chen,Xiang Hao,Jing Bian,Feng Peng.Xylan-based hydrogels for potential skin care application.*Int J Biol Macromol.* 2020 Apr 30;158:244-250.;doi: 10.1016/j.ijbiomac.2020.04.235
88. Xianmou Fan,Shaobing Wang,Yan Fang,Peiyuan Li,Weikang Zhou,Zhengchao Wang,Mingfeng Chen,Haiqing Liu.Tough polyacrylamide-tannic acid-kaolin adhesive hydrogels for quick hemostatic application.*MAT SCI ENG C-MATER.* 2020 Apr;109:110649.;doi: 10.1016/j.msec.2020.110649
89. Zhiping Yu,Ming Hu,Zhouquan Li,Dan Xu,Lianhua Zhu,Yanli Guo,Qiuli Liu,Weihua Lan,Jun Jiang,Luofu Wang.Anti-G250 nanobody-functionalized nanobubbles targeting renal cell carcinoma cells for ultrasound molecular imaging.*Nanotechnology.* 2020 May 15;31(20):205101.;doi: 10.1088/1361-6528/ab7040
90. Xiaohong Yang,Timon Chengyi Liu,Shaojie Liu,Weicong Zhu,Honglin Li,Peihong Liang,Suihui Ye,Shuliang Cui.Promoted Viability and Differentiated Phenotype of Cultured Chondrocytes With Low Level Laser Irradiation Potentiate Efficacious Cells for Therapeutics.*Front Bioeng Biotechnol.* 2020 May 29;8:468.;doi: 10.3389/fbioe.2020.00468.
91. Hai Qiao,Zixin Wei,Yiting Wang,Na Hu,Sineng Sun,Jin Bai,Liaoqiong Fang,Zhibiao Wang.Focused characteristics and effects of light

- reflected from spherical lipid membrane of giant unilamellar vesicles COLLOID SURFACE B. 2020 May;189:110828.;doi: 10.1016/j.colsurfb.2020.110828
92. Yuchen Dai,Yi Yang,Baohong Liu,Liang Qiao.Water-in-oil microcompartments for the study of biomimetic drug metabolism J COLLOID INTERF SCI. 2020 Jun 1;569:378-385.;doi: 10.1016/j.jcis.2020.02.096
93. Hao Shen,Emile Nibona,Gongyu Xu,Md Abdullah Al Hafiz,Xiaomei Ke,Xiaoting Liang,Qiting Yao,Xueping Zhong,Qingchun Zhou,Haobin Zhao.Identification, expression pattern, and immune response of Tim-1 and Tim-4 in embryos and adult medaka (*Oryzias latipes*) J Exp Zool B Mol Dev Evol. 2020 Jun;334(4):235-244.;doi: 10.1002/jez.b.22939
94. Honglian Li,Ruirui Lu,Yu Pang,Jicheng Li,Yiwen Cao,Hongxin Fu,Guoxing Fang,Qiuhe Chen,Bihao Liu,Junbiao Wu,Yuan Zhou,Jiuyao Zhou,Zhen-Wu-Tang Protects IgA Nephropathy in Rats by Regulating Exosomes to Inhibit NF- κ B/NLRP3 Pathway Front Pharmacol. 2020 Jul 16;11:1080.;doi: 10.3389/fphar.2020.01080
95. Jinheng Wang,Guangmeng Li,Chenggong Tu,Xiaoming Chen,Bin Yang,Yongliang Huo,Yi Li,Ai-Zheng Chen,Ping Lan,Yu Shrike Zhang,Maoxin Xie.High-throughput single-cell analysis of exosome mediated dual drug delivery, in vivo fate and synergistic tumor therapy Nanoscale. 2020 Jul 2;12(25):13742-13756.;doi: 10.1039/d0nr02344b
96. Feifei Li,Lei Zhan,Qian Dong,Qiwei Wang,Yuanhe Wang,Xiaoyan Li,Yong Zhang,Jingdong Zhang.Tumor-Derived Exosome-Educated Hepatic Stellate Cells Regulate Lactate Metabolism of Hypoxic Colorectal Tumor Cells via the IL-6/STAT3 Pathway to Confer Drug Resistance ONCOTARGETS THER. 2020 Aug 7;13:7851-7864.;doi: 10.2147/OTT.S253485
97. Yuan Liang,Hengli Yang,Qiaoying Li,Ping Zhao,Han Li,Yuxin Zhang,Wenbin Cai,Xiaojun Ma,Yunyou Duan.Novel biomimetic dual-mode nanodroplets as ultrasound contrast agents with potential ability of precise detection and photothermal ablation of tumors CANCER CHEMOTH PHARM. 2020 Sep;86(3):405-418.;doi: 10.1007/s00280-020-04124-x
98. Zhen Wang,Chenghui Yang,Lili Li,Zhigang Zhang,Jun Pan,Ke Su,Wuzhen Chen,Jinfan Li,Fuming Qiu,Jian Huang.CD62L dim Neutrophils Specifically Migrate to the Lung and Participate in the Formation of the Pre-Metastatic Niche of Breast Cancer Front Oncol. 2020 Oct 15;10:540484.;doi: 10.3389/fonc.2020.540484
99. Xingyu Zhang,Jun Chen,Qin Jiang,Xiaoquan Ding,Yunxia Li,Chen Chen,Wuli Yang,Shiyi Chen.Highly biosafe biomimetic stem cell membrane-disguised nanovehicles for cartilage regeneration J Mater Chem B. 2020 Oct 7;8(38):8884-8893.;doi: 10.1039/d0tb01686a
100. Ning Zhang,Yanan Song,Zheyong Huang,Jing Chen,Haipeng Tan,Hongbo Yang,Mengkang Fan,Qiyu Li,Qiaozi Wang,Jinfeng Gao,Zhiqing Pang,Juying Qian,Junbo Ge.Monocyte mimics improve mesenchymal stem cell-derived extracellular vesicle homing in a mouse MI/RI model Biomaterials. 2020 Oct;255:120168.;doi: 10.1016/j.biomaterials.2020.120168
101. Pengwei Luan,Xinyue Ding,Jiazheng Xu,Lixian Jiang,Yulan Xu,Yuying Zhu,Ruixiang Li,Jiange Zhang.Salvianolate reduces neuronal apoptosis by suppressing OGD-induced microglial activation Life Sci. 2020 Nov 1;260:118393.;doi: 10.1016/j.lfs.2020.118393
102. Xiaofen Xie,Xinying Li,Jinfeng Lei,Xi Zhao,Yongbo Lyu,Changdao Mu,Defu Li,Liming Ge,Yongbin Xu.Oxidized starch cross-linked porous collagen-based hydrogel for spontaneous agglomeration growth of adipose-derived stem cells MAT SCI ENG C-MATER. 2020 Nov;116:111165.;doi: 10.1016/j.msec.2020.111165
103. Chang-Chun Chen, Jing Chen, Wen-Liang Wang, Liang Xie, Chuan-Qiang Shao, Yan-Xiang Zhang. Inhibition of the P53/P21 Pathway Attenuates the Effects of Senescent Nucleus Pulpous Cell-Derived Exosomes on the Senescence of Nucleus Pulpous Cells. Orthop Surg. doi: 10.1111/os.12886.
104. Guangshu Wang,Zhenyu Wang,Ning Sun,Yurong Cai,Xiaogang Yang. Reactive oxygen species-responsive silk sericin microcapsules used for antioxidative stress damage. Microsc Res Tech. doi: 10.1002/jemt.23620.
105. Yubing Han, Zhimin Zhang, Wenjie Liu, Yuanfa Yao, Yingke Xu, Xu Liu, Cuifang Kuang, Xiang Hao. A Labeling Strategy for Living Specimens in Long-Term/Super-Resolution Fluorescence Imaging. Front Chem. doi: 10.3389/fchem.2020.601436.
106. Yongjin Sun, Xu Li, Xiaoxu Yang, Bi Chen, Wenzhi Zhang. Small Extracellular Vesicles Derived from Adipocytes Attenuate Intervertebral Disc Degeneration in Rats by Rejuvenating Senescent Nucleus Pulpous Cells and Endplate Cells by Delivering Exogenous NAMP. Oxid Med Cell Longev. doi: 10.1155/2021/9955448.
107. Fei Yan, Xuexi Zhang, Rongying Tan, Mingchen Li, Zengtuan Xiao, Hao Wang, Zhenfa Zhang, Zhenyi Ma, Zhe Liu. Autophagic flux in cancer cells at the invasive front in the tumor-stroma border. Aging (Albany NY). doi: 10.18632/aging.203406.
108. Li-Juan Zhao, Ying-Ying Li, Yu-Tong Zhang, Qi-Qi Fan, Hong-Mei Ren, Cheng Zhang, Adil Mardinoglu, Wen-Chao Chen, Jing-Ru Pang, Dan-Dan Shen, Jun-Wei Wang, Long-Fei Zhao, Jian-Ying Zhang, Zhen-Ya Wang, Yi-Chao Zheng, Hong-Min Liu. Lysine demethylase LSD1 delivered via small extracellular vesicles promotes gastric cancer cell stemness. EMBO Rep. doi: 10.15252/embr.202050922.
109. Shiji Lin, Dehui Wang, Lijuan Zhang, Yakang Jin, Zhigang Li, Elmar Bonaccorso, Zili You, Xu Deng, Longquan Chen. Macrodrop-Impact-Mediated Fluid Microdispensing. Adv Sci (Weinh). doi: 10.1002/advs.202101331.
110. Mingda Li, Jie Jia, Shanshan Li, Baocheng Cui, Jiao Huang, Zhaoming Guo, Kun Ma, Li Wang, Changhao Cui. Exosomes derived from tendon stem cells promote cell proliferation and migration through the TGF β signal pathway. Biochem Biophys Res Commun. doi: 10.1016/j.bbrc.2020.12.057.
111. Ning Zhang, Nianrong Sun, Chunhui Deng. Rapid isolation and proteome analysis of urinary exosome based on double interactions of Fe3O4@TiO2-DNA aptamer. Talanta. doi: 10.1016/j.talanta.2020.121571.
112. Jian Wang, Yuanting Liu, Ying Li, Xuan Zheng, Jianhui Gan, Zhaoyuan Wan, Jun Zhang, Yan Liu, Yaqi Wang, Wanning Hu, Yufeng Li, Yankun Liu. Exosomal - miR - 10a derived from colorectal cancer cells suppresses migration of human lung fibroblasts, and expression of IL-6, IL-8 and IL-1 β . Mol Med Rep. doi: 10.3892/mmr.2020.11723.
113. Xinglei Qin, Min Lu, Gang Li, Yajun Zhou, Zhaoyang Liu. Downregulation of tumor-derived exosomal miR-34c induces cancer-associated fibroblast activation to promote cholangiocarcinoma progress. Cancer Cell Int. doi: 10.1186/s12935-020-01726-6.
114. Xuechun Li, Yang Yu, Renyue Wei, Yimei Li, Jiawei Lv, Zhonghua Liu, Yu Zhang. In vitro and in vivo study on angiogenesis of porcine induced pluripotent stem cell-derived endothelial cells. Differentiation. doi: 10.1016/j.diff.2021.05.003.
115. Shuang-Shuang Dong, Dan-Dan Dong, Zhang-Fu Yang, Gui-Qi Zhu, Dong-Mei Gao, Jie Chen, Yan Zhao, Bin-Bin Liu. Exosomal miR-3682-3p Suppresses Angiogenesis by Targeting ANGPT1 via the RAS-MEK1/2-ERK1/2 Pathway in Hepatocellular Carcinoma. Front Cell Dev Biol. doi: 10.3389/fcell.2021.633358.
116. Zichu Yin, Han Liu, Minjie Lin, Wenjiao Xie, Xiaogang Yang, Yurong

- Cai. Controllable performance of a dopamine-modified silk fibroin-based bio-adhesive by doping metal ions. *Biomed Mater.* doi: 10.1088/1748-605X/ac0087.
117. Yan Zhang, Yifei Yin, Wei Zhang, Hongyan Li, Taixia Wang, Haohao Yin, Liping Sun, Chunxia Su, Kun Zhang, Huixiong Xu. Reactive oxygen species scavenging and inflammation mitigation enabled by biomimetic prussian blue analogues boycott atherosclerosis. *J Nanobiotechnology.* doi: 10.1186/s12951-021-00897-2.
118. Panmiao Liu, Zhongde Mu, Muhuo Ji, Xiaojiang Liu, Hanwen Gu, Yi Peng, Jianjun Yang, Zhuoying Xie, Fuyin Zheng. Robust Carbonated Structural Color Barcodes with Ultralow Ontology Fluorescence as Biomimic Culture Platform. *Research (Wash D C).* doi: 10.34133/2021/9851609.
119. Guangzheng Yang, Fei Jiang, Yuezhi Lu, Sihan Lin, Chang Liu, Anshuo Li, David L Kaplan, Shilei Zhang, Yue He, Cui Huang, Wenjie Zhang, Xinquan Jiang. Rapid construction and enhanced vascularization of microtissue using a magnetic control method. *Biofabrication.* doi: 10.1088/1758-5090/abe4c2.
120. Baixiong Zhao, Jun Zheng, Yang Qiao, Yongquan Wang, Yang Luo, Dinglin Zhang, Qiyan Cai, Yang Xu, Zhansong Zhou, Wenhao Shen. Prostatic fluid exosome-mediated microRNA-155 promotes the pathogenesis of type IIIA chronic prostatitis. *Transl Androl Urol.* doi: 10.21037/tau-21-139.
121. Xiuping Cao, Xinxin Shang, Yingshu Guo, Xiaofei Zheng, Wenxin Li, Di Wu, Li Sun, Shanliang Mu, Chuanen Guo. Lysosomal escaped protein nanocarriers for nuclear-targeted siRNA delivery. *Anal Bioanal Chem.* doi: 10.1007/s00216-021-03297-5.
122. Fei Xiao, Bin Zuo, Bo Tao, Chuandong Wang, Yang Li, Jianping Peng, Chao Shen, Yiming Cui, Junfeng Zhu, Xiaodong Chen. Exosomes derived from cyclic mechanical stretch-exposed bone marrow mesenchymal stem cells inhibit RANKL-induced osteoclastogenesis through the NF- κ B signaling pathway. *Ann Transl Med.* doi: 10.21037/atm-21-1838.
123. Yanxing Wang, Miao Feng, Bi Lin, Xiangrong Peng, Zhan Wang, Ruichan Lv. MET-targeted NIR II luminescence diagnosis and up-conversion guided photodynamic therapy for triple-negative breast cancer based on a lanthanide nanoprobe. *Nanoscale.* doi: 10.1039/d1nr05847a.
124. Lizhong Sun, Libang He, Wei Wu, Li Luo, Mingyue Han, Yifang Liu, Shijie Shi, Kaijing Zhong, Jiaoqiao Yang, Jiyao Li. Fibroblast membrane-camouflaged nanoparticles for inflammation treatment in the early stage. *Int J Oral Sci.* doi: 10.1038/s41368-021-00144-2.
125. Zhongjuan Xu, Junjun Cao, Zhe Zhao, Yong Qiao, Xingzhi Liu, Junjie Zhong, Bin Wang, Guangli Suo. A functional extracellular matrix biomaterial enriched with VEGFA and bFGF as vehicle of human umbilical cord mesenchymal stem cells in skin wound healing. *Biomed Mater.* doi: 10.1088/1748-605X/ac37b0.
126. Yinuo Cen, Yue Lou, Junjun Wang, Shicheng Wang, Peng Peng, Aili Zhang, Ping Liu. Supplementation with Serum-Derived Extracellular Vesicles Reinforces Antitumor Immunity Induced by Cryo-Thermal Therapy. *Int J Mol Sci.* doi: 10.3390/ijms222011021.
127. Yu Fu, Leilei Wang, Wei Liu, Lan Yang, Lin Li, Luyao Wang, Xun Sun, Zhi-Rong Zhang, Qing Lin, Ling Zhang. OX40L blockade cellular nanovesicles for autoimmune diseases therapy. *J Control Release.* doi: 10.1016/j.jconrel.2021.08.008.
128. Xiaopeng Hao, Huixia Fan, Jian Yang, Jinfu Tang, Junhui Zhou, Yuyang Zhao, Luqi Huang, Yong Xia. Network Pharmacology Research and Dual-omic Analyses Reveal the Molecular Mechanism of Natural Product Nodosin Inhibiting Muscle-Invasive Bladder Cancer in Vitro and in Vivo. *J Nat Prod.* doi: 10.1021/acs.jnatprod.2c00400.
129. Xiang Li, Erhui Jiang, Hui Zhao, Yang Chen, Yuming Xu, Chunyu Feng, Ji Li, Zhengjun Shang. Glycometabolic reprogramming-mediated proangiogenic phenotype enhancement of cancer-associated fibroblasts in oral squamous cell carcinoma: role of PGC-1 α /PFKFB3 axis. *Br J Cancer.* doi: 10.1038/s41416-022-01818-2.
130. Yawen Nie, Xiao-Qian Xie, Lingxi Zhou, Qijie Guan, Yilin Ren, Yong Mao, Jin-Song Shi, Zheng-Hong Xu, Yan Geng. Desulfovibrio fairfieldensis-Derived Outer Membrane Vesicles Damage Epithelial Barrier and Induce Inflammation and Pyroptosis in Macrophages. *Cells.* doi: 10.3390/cells12010089.
131. Yupeng Wang, Zhou Li, Wei Wu, Ying Liu, Yu Xiao, Dongdong Qi, Guangming Zhao, Meijuan Zhou, Hua Wang, Jing Liu, Zhiqi Song. TRPA1 promotes melanosome phagocytosis in keratinocytes via PAR-2/CYLD axis. *J Dermatol Sci.* doi: 10.1016/j.jdermsci.2022.05.005.
132. Yuzhou Li, Yiru Fu, He Zhang, Xu Wang, Tao Chen, Yanqiu Wu, Xinxin Xu, Sheng Yang, Ping Ji, Jinlin Song. Natural Plant Tissue with Bioinspired Nano Amyloid and Hydroxyapatite as Green Scaffolds for Bone Regeneration. *Adv Health Mater.* doi: 10.1002/adhm.202102807.
133. Hongxiao Cui, Yajun Sun, Hua Lin, Yan Zhao, Xin Zhao. The Outer Membrane Vesicles of *Salmonella enterica* Serovar Typhimurium Activate Chicken Immune Cells through Lipopolysaccharides and Membrane Proteins. *Pathogens.* doi: 10.3390/pathogens11030339.
134. Yao Liu, Yiping Lu, Bo Ning, Xiaomin Su, Binru Yang, Haiqing Dong, Bo Yin, Zhiqing Pang, Shun Shen. Intravenous Delivery of Living *Listeria monocytogenes* Elicits Gasdmermin-Dependent Tumor Pyroptosis and Motivates Anti-Tumor Immune Response. *ACS Nano.* doi: 10.1021/acsnano.1c09818.
135. Chenyang Cui, Tingting Guo, Shuai Zhang, Mingyan Yang, Jiaqi Cheng, Jiajia Wang, Jie Kang, Wenjie Ma, Yuanru Nian, Zhaowei Sun, Haibo Weng. Bacteria-derived outer membrane vesicles engineered with over-expressed pre-miRNA as delivery nanocarriers for cancer therapy. *Nanomedicine.* doi: 10.1016/j.nano.2022.102585.

Version 2024.03.12