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1. Ke Song, Xinyan Yang, Geng An, Xinyu Xia, Jiexiang Zhao, Xiaoheng Xu, Cong Wan, Tianyuan Liu, Yi Zheng, Shaofang Ren, Mei Wang, Gang Chang, Shane J F Cronin, Josef M Penninger, Tao Jing, Xianghong Ou, Shuan Rao, Zhaoting Liu, Xiao-Yang Zhao. Targeting APLN/APJ restores blood-testis barrier and improves spermatogenesis in murine and human diabetic models. Nat Commun. 2022 Nov 28;13(1):7335.
2. Bin Qiao, Xinye Song, Weiyi Zhang, Ming Xu, Bowen Zhuang, Wei Li, Huanling Guo, Wenxin Wu, Guangliang Huang, Minru Zhang, Xiaoyan Xie, Nan Zhang, Yong Luan, Chunyang Zhang. Intensity-adjustable pain management with prolonged duration based on phase-transitional nanoparticles-assisted ultrasound imaging-guided nerve blockade. J Nanobiotechnology. 2022 Nov 24;20(1):498.
3. Qing Zhang, Li Chen, Liyi Huang, Hongxin Cheng, Lu Wang, Lin Xu, Danrong Hu, Chengqi He, Chenying Fu, Quan Wei. CD44 promotes angiogenesis in myocardial infarction through regulating plasma exosome uptake and further enhancing FGFR2 signaling transduction. Mol Med. 2022 Dec 3;28(1):145.
4. Hui Zheng, Jialing Ji, Tangming Zhao, E Wang, Aiqing Zhang. Exosome-encapsulated miR-26a attenuates aldosterone-induced tubulointerstitial fibrosis by inhibiting the CTGF/SMAD3 signaling pathway. Int J Mol Med. 2023 Feb;51(2):11.

Version 2024.05.21