

## ER-Tracker Red (内质网红色荧光探针)

产品编号	产品名称	包装
C1041S	ER-Tracker Red (内质网红色荧光探针)	20 $\mu$ l
C1041M	ER-Tracker Red (内质网红色荧光探针)	100 $\mu$ l
C1041L	ER-Tracker Red (内质网红色荧光探针)	500 $\mu$ l

### 产品简介:

- ER-Tracker Red是一种具有细胞膜通透性的内质网(endoplasmic reticulum, ER)红色荧光探针, 对内质网有高度选择性, 可以用于活细胞内质网特异性荧光染色。
- ER-Tracker Red为采用Molecular Probes公司的BODIPY TR进行了荧光标记的glibenclamide。Glibenclamide即glyburide, 中文名为格列本脲, 是一种2型糖尿病治疗药物, 可以结合主要定位在内质网上的包含ATP敏感的钾离子通道( $K_{ATP}$  channel)的磺脲类(sulphonylurea)受体。因此荧光标记的glibenclamide就可以用作内质网特异性的荧光探针。ER-Tracker Red适用于活细胞内质网的荧光染色, 但不适合用于固定细胞内质网的荧光染色。
- ER-Tracker Red的分子式为 $C_{44}H_{42}BClF_2N_6O_7S_2$ , 分子量为915.23, ER-Tracker Red的化学结构式参考图1。




图1. ER-Tracker Red的化学结构式。

- ER-Tracker Red对于细胞的毒性极低。而传统的DiOC6(3)对ER染色的同时也对细胞有一定的毒性。
- ER-Tracker Red呈红色荧光, 检测时的最大激发波长为587nm, 最大发射波长为615nm。ER-Tracker Red的激发光谱和发射光谱参考图2。




图2. ER-Tracker Red的激发光谱和发射光谱。

- ER-Tracker Red按照后附的使用说明染色后, 用甲醛等固定后染色效果可以被部分保留。ER-Tracker Red在活细胞中的染色效果参考图2。




图3. ER-Tracker Red(内质网红色荧光探针)对于NRK-52E细胞(大鼠肾小管上皮细胞)的染色效果。Mito-Tracker Green染色的NRK-52E细胞其线粒体呈现绿色荧光(图B), ER-Tracker Red染色的NRK-52E细胞其内质网呈现红色荧光(图C), 绿色荧光、红色荧光及细胞核蓝色荧光的叠加(merge)效果见图D。其中细胞核使用Hoechst 33342 (C1027)染色。本图仅作参考, 实际检测效果会因实验条件、检测仪器等的不同而存在差异。

- 提供了ER-Tracker Red稀释液, 使ER-Tracker Red的使用更加便捷。
- 按照1:1000的比例稀释, 每10 $\mu$ l可以配制10ml ER-Tracker Red工作液; 按照1:3000的比例稀释, 每10 $\mu$ l可以配制30ml ER-Tracker Red工作液。

#### 包装清单:

产品编号	产品名称	包装
C1041S-1	ER-Tracker Red (1mM)	20 $\mu$ l
C1041S-2	ER-Tracker Red稀释液	60ml
—	说明书	1份

产品编号	产品名称	包装
C1041M-1	ER-Tracker Red (1mM)	100 $\mu$ l
C1041M-2	ER-Tracker Red稀释液	250ml
—	说明书	1份

产品编号	产品名称	包装
C1041L-1	ER-Tracker Red (1mM)	500 $\mu$ l
C1041L-2	ER-Tracker Red稀释液	500ml $\times$ 2
—	说明书	1份

#### 保存条件:

-20 $^{\circ}$ C保存, 半年有效。ER-Tracker Red需-20 $^{\circ}$ C避光保存。

#### 注意事项:

- ER-Tracker Red (1mM)在4 $^{\circ}$ C、冰浴等较低温度情况下会凝固而粘在离心管管底、管壁或管盖内, 可以20-25 $^{\circ}$ C水浴温育片刻至全部融解后使用。对于微量的液体, 每次使用前先离心数秒钟, 使液体充分沉降到管底。
- 荧光染料均存在淬灭问题, 请尽量注意避光, 以减缓荧光淬灭。
- 需自备盖玻片和载玻片(可以向碧云天订购)。
- 格列本脲的药理学活性可能会影响内质网的一些功能; 某些特殊细胞中磺脲类受体的可变表达可能会导致非内质网特异性染色。
- ER-Tracker Red适用于活细胞内质网的荧光染色, 但不适用于固定细胞内质网的荧光染色。如果经ER-Tracker Red染色后的细胞需要进行固定操作, 使用4%多聚甲醛在37 $^{\circ}$ C固定2分钟。
- ER-Tracker Red染色后的细胞不能用Triton X-100通透, Triton X-100通透处理会导致ER-Tracker Red的荧光染色消失。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

#### 使用说明:

- ER-Tracker Red工作液的配制:
  - 取少量ER-Tracker Red按照1:1000的比例加入到ER-Tracker Red稀释液中。例如取1 $\mu$ l ER-Tracker Red加入到1ml ER-Tracker Red稀释液中。混匀后即为ER-Tracker Red工作液。
  - ER-Tracker Red工作液使用前需37 $^{\circ}$ C预温育。  
注: 工作液中ER-Tracker Red的浓度可以根据实际情况进行适当调整, 推荐的稀释比例调整范围为1:1000-1:3000。为降低背景, 在染色效果可以接受的范围内, 建议尽量使用较低浓度的ER-Tracker Red。
- 内质网的荧光标记:
  - 去除细胞培养液, 用适量的溶液如HBSS with Ca $^{2+}$  & Mg $^{2+}$  (Hanks' Balanced Salt Solution with Ca $^{2+}$  & Mg $^{2+}$ )洗涤生长在盖玻片上的细胞。注: HBSS with Ca $^{2+}$  & Mg $^{2+}$  (C0219)可以向碧云天订购; 对于悬浮细胞的染色可以参考贴壁细胞的染色方法进行。
  - 去除洗涤液, 加入步骤1配制好的并37 $^{\circ}$ C预温育的ER-Tracker Red染色工作液, 与细胞37 $^{\circ}$ C共孵育15-30分钟。
  - 去除ER-Tracker Red染色工作液, 用细胞培养液洗涤细胞1-2次。
  - 随后通常用荧光显微镜或激光共聚焦显微镜进行观察。此时可观察到内质网呈明亮的强荧光染色。
  - 如果经ER-Tracker Red染色后的细胞需要进行固定, 可以使用4%甲醛37 $^{\circ}$ C固定2分钟。固定后用适当的洗涤液洗涤2-3次, 每次5分钟, 随后可以进行复染或滴加适当的抗荧光淬灭封片液, 最后封片观察。注意: ER-Tracker Red染色的细胞不能用Triton X-100通透, Triton X-100通透处理会导致ER-Tracker Red的荧光染色消失。

相关产品:

产品编号	产品名称	包装
C1002	DAPI	5mg/ml×0.2ml
C1005/C1006	DAPI 染色液	10ml/50ml
C1011	Hoechst 33258	10mg
C1017/C1018	Hoechst 33258 染色液	10ml/50ml
C1022	Hoechst 33342	10mg
C1025/C1026	Hoechst 33342 染色液	10ml/50ml
C1027/C1028/C1029	Hoechst 33342 活细胞染色液(100X)	0.1ml/0.5ml/3ml
C1033	Actin-Tracker Green (微丝绿色荧光探针)	0.2ml
C1036	DiI (细胞膜红色荧光探针)	10mg
C1038	DiO (细胞膜绿色荧光探针)	10mg
C1039-10mg	DiD (细胞膜远红外荧光探针)	10mg
C1041S/M/L	ER-Tracker Red (内质网红色荧光探针)	20μl/100μl/500μl
C1042S/M/L	ER-Tracker Green (内质网绿色荧光探针)	20μl/100μl/500μl
C1043	Golgi-Tracker Red (高尔基体红色荧光探针)	1mg
C1045S	Golgi-Tracker Green (高尔基体绿色荧光探针)	1mg
C1046	Lyso-Tracker Red (溶酶体红色荧光探针)	50μl
C1047S	Lyso-Tracker Green (溶酶体绿色荧光探针)	50μl
C1048	Mito-Tracker Green (线粒体绿色荧光探针)	50μg
C1049B-50μg	Mito-Tracker Red CMXRos (线粒体红色荧光探针)	50μg
C1049B-250μg	Mito-Tracker Red CMXRos (线粒体红色荧光探针)	50μg×5
C1050	Tubulin-Tracker Red (抗体法微管红色荧光探针)	40μl
C1051S	Tubulin-Tracker Green (抗体法微管绿色荧光探针)	40μl
C1991S	细胞膜红色荧光染色试剂盒(DiI)	200-1000次
C1993S	细胞膜绿色荧光染色试剂盒(DiO)	200-1000 次
C1995S	细胞膜远红外荧光染色试剂盒(DiD)	200-1000 次
C2005	JC-1	1mg
C2007	Rhodamine 123	5mg

使用本产品的文献:

- Zhang S, Wu W, Wu Y, Zheng J, Suo T, Tang H, Tang J. RNF152, a novel lysosome localized E3 ligase with pro-apoptotic activities. *Protein Cell* . 2010 Jul;1(7):656-63.
- Shao J, Dai Y, Zhao W, Xie J, Xue J, Ye J, Jia L. Intracellular distribution and mechanisms of actions of photosensitizer Zinc(II)-phthalocyanine solubilized in Cremophor EL against human hepatocellular carcinoma HepG2 cells. *Cancer Lett* .2013 Mar 1;330(1):49-56.
- Wu S, Wang F, Huang J, Fang Q, Shen Z, Ye G. Molecular and cellular analyses of a ryanodine receptor from hemocytes of *Pieris rapae*. *Dev Comp Immunol* . 2013 Sep;41(1):1-10.
- Wu L, Yang X, Duan X, Cui L, Li G. Exogenous expression of marine lectins DIFBL and SpRBL induces cancer cell apoptosis possibly through PRMT5-E2F-1 pathway. *SCI REP-UK* . 2014 Mar 28;4:4505.
- Pan C, Hu YF, Yi HS, Song J, Wang L, Pan MH, Lu C. Role of Bmbuffy in hydroxycamptothecine-induced apoptosis in BmN-SWU1 cells of the silkworm, *Bombyx mori*. *BIOCHEM BIOPH RES CO* . 2014 May 2;447(2):237-43.
- Yang X, Wu L, Duan X, Cui L, Luo J, Li G. Adenovirus carrying gene encoding Haliotis discus discus sialic acid binding lectin induces cancer cell apoptosis. *Mar Drugs* . 2014 Jun 30;12(7):3994-4004.
- Qin C, He B, Dai W, Lin Z, Zhang H, Wang X, Wang J, Zhang X, Wang G, Yin L, Zhang Q. The impact of a chlorotoxin-modified liposome system on receptor MMP-2 and the receptor-associated protein CIC-3. *Biomaterials* . 2014 Jul;35(22):5908-20.
- Zhang BC, Zhang J, Xiao ZZ, Sun L. Rock bream (*Oplegnathus fasciatus*) viperin is a virus-responsive protein that modulates innate immunity and promotes resistance against megalocytivirus infection. *Dev Comp Immunol* . 2014 Jul;45(1):35-42.
- Zhang Y, Zhou L, Zou J, Wang M, Qi J, Qi Z. Palmitoylation of STREX domain confers cerebroside sensitivity to the BKCa channel. *BBA-BIOMEMBRANES*. 2014 Oct;1838(10):2451-9.
- Ao D, Guo HC, Sun SQ, Sun DH, Fung TS, Wei YQ, Han SC, Yao XP, Cao SZ, Liu DX, Liu XT. Viroporin Activity of the Foot-and-Mouth Disease Virus Non-Structural 2B Protein. *PLoS One* . 2015 May 6;10(5):e0125828.
- Lu J, Cao Y, Cheng K, Xu B, Wang T, Yang Q, Yang Q, Feng X, Xia Q. Berberine regulates neurite outgrowth through AMPK-dependent pathways by lowering energy status. *Exp Cell Res*. 2015 Jun 10;334(2):194-206.
- Zhang B, Han H, Fu S, Yang P, Gu Z, Zhou Q, Cao Z. Dehydroeffusol inhibits gastric cancer cell growth and tumorigenicity by selectively inducing tumor-suppressive endoplasmic reticulum stress and a moderate apoptosis. *Biochem Pharmacol* . 2016 Mar 15;104:8-18.
- Ma R, Zhang J, Liu X, Li L, Liu H, Rui R, Gu L, Wang Q. Involvement of Rab6a in organelle rearrangement and cytoskeletal organization during mouse oocytematuration. *SCI REP-UK* . 2016 Mar 31;6:23560.
- Xu Y, Xu J, Shan W, Liu M, Cui Y, Li L, Liu C, Huang Y. The transport mechanism of integrin αvβ3 receptor targeting nanoparticles in Caco-2 cells. *INT J PHARMACOL*. 2016 Mar 16;500(1-2):42-53.
- Li G, Gao Y, Cui L, Wu L, Yang X, Chen J. *Anguilla japonica* lectin 1

- delivery through adenovirus vector induces apoptotic cancer cell death through interaction with PRMT5. *J Gene Med.* 2016 Apr;18(4-6):65-74.
16. Zhang ZZ, Yuan K, Yue HT, Yuan FH, Bi HT, Weng SP, He JG, Chen YH. Identification and functional characterization of an endoplasmic reticulum oxidoreductin 1- $\alpha$  gene in *Litopenaeus vannamei*. *Dev Comp Immunol.* 2016 Apr;57:10-9.
  17. Bi HT, Yuan FH, Yuan K, Weng SP, He JG, Chen YH. Identification and functional characterization of a glucose regulated protein 94 gene in *Litopenaeus vannamei* and its responsiveness in WSSV infection. *Mol Immunol.* 2016 May;73:29-36.
  18. Chen X, Zhang X, Wang HY, Chen Z, Wu FG. Subcellular Fate of a Fluorescent Cholesterol-Poly(ethylene glycol) Conjugate: An Excellent Plasma Membrane Imaging Reagent. *Langmuir.* 2016 Oct 4;32(39):10126-10135.
  19. Sun W, Li L, Yang QQ, Zhang ZR, Huang Y. Two birds, one stone: dual targeting of the cancer cell surface and subcellular mitochondria by the galectin-3-binding peptide G3-C12. *Acta Pharmacol Sin.* 2017 Jan 9.
  20. Wang H, Niu H, Zhai Y, Lu M. Characterization of BiP Genes from Pepper (*Capsicum annuum* L.) and the Role of CaBiP1 in Response to Endoplasmic Reticulum and Multiple Abiotic Stresses. *Front Plant Sci.* 2017 Jun 28;8:1122.
  21. Dong W, Lv H, Guo K, Wang T, Ouyang Y, Jin M, Zhang Y. Classical Swine Fever Virus Infection and Its NS4A Protein Expression Induce IL-8 Production through MAVS Signaling Pathway in Swine Umbilical Vein Endothelial Cells. *Front Microbiol.* 2018 Jan 12;8:2687.
  22. Zheng S, Huang C, Zhao X, Zhang Y, Liu S, Zhu Q. A hydrophobic organelle probe based on aggregation-induced emission: Nanosuspension preparation and direct use for endoplasmic reticulum imaging in living cells. *SPECTROCHIM ACTA A.* 2018 Jan 15;189:231-238.
  23. La X, Zhang L, Li H, Li Z, Song G, Yang P, Yang Y. Ajuba receptor mediates the internalization of tumor-secreted GRP78 into macrophages through different endocytosis pathways. *ONCOTARGET.* 2018 Jan 9;9(21):15464-15479.
  24. Wang B, Wen A, Feng C, Niu L, Xiao X, Luo L, Shen C, Zhu J, Lei J, Zhang X. The in vivo anti-fibrotic function of calcium sensitive receptor (CaSR) modulating poly(p-dioxanone-co-l-phenylalanine) prodrug. *Acta Biomater.* 2018 Jun;73:180-189.
  25. E S, Mao QX, Yuan XL, Kong XL, Chen XW, Wang JH. Targeted imaging of the lysosome and endoplasmic reticulum and their pH monitoring with surface regulated carbon dots. *Nanoscale.* 2018 Jul 9;10(26):12788-12796.
  26. Gong X, Sun R, Gao Z, Han W, Liu Y, Zhao L, Jing L, Yao X, Sun X. Tubeimoside 1 Acts as a Chemotherapeutic Synergist via Stimulating Macropinocytosis. *Front Pharmacol.* 2018 Sep 26;9:1044.
  27. Ling J, Ma Z, Liu L, Yin J, Su J, Shen F, Xie L, Hu S. Identification of a crucial tryptophan residue in ADAMTS13 required for its secretion and enzymatic activity. *CLIN EXP PHARMACOL P.* 2018 Nov;45(11):1181-1186.
  28. Guo C, Ma R, Liu X, Xia Y, Niu P, Ma J, Zhou X, Li Y, Sun Z. Silica nanoparticles induced endothelial apoptosis via endoplasmic reticulum stress-mitochondrial apoptotic signaling pathway. *Chemosphere.* 2018 Nov;210:183-192.
  29. Xu A, Tang Y, Lin W. Endoplasmic reticulum-targeted two-photon turn-on fluorescent probe for nitroreductase in tumor cells and tissues. *SPECTROCHIM ACTA A.* 2018 Nov 5;204:770-776.
  30. Lian YY, He HH, Zhang CZ, Li XC, Chen YH. Functional characterization of a matrix metalloproteinase 2 gene in *Litopenaeus vannamei*. *FISH SHELLFISH IMMUN.* 2019 Jan;84:404-413.
  31. Li X, Du X, Ni J. Zn<sup>2+</sup> Aggravates Tau Aggregation and Neurotoxicity. *Int J Mol Sci.* 2019 Jan 23;20(3).
  32. Yan JJ, Wang XY, Wang MZ, Pan DH, Yang RL, Xu YP, Wang LZ, Yang M. Self-Assembling Nonconjugated Poly(amide-imide) into Thermoresponsive Nanovesicles with Unexpected Red Fluorescence for Bioimaging. *Biomacromolecules.* 2019 Mar 11;20(3):1455-1463.
  33. Xu D, Wu L, Jiang X, Yang L, Cheng J, Chen H, Hua R, Geng G, Yang L, Li Q. SIRT2 Inhibition Results in Meiotic Arrest, Mitochondrial Dysfunction, and Disturbance of Redox Homeostasis during Bovine Oocyte Maturation. *Int J Mol Sci.* 2019 Mar 18;20(6).
  34. Jiang YJ, Lin M, Yang T, Li RS, Huang CZ, Wang J, Li YF. Nitrogen and phosphorus doped polymer carbon dots as a sensitive cellular mapping probe of nitrite. *J Mater Chem B.* 2019 Mar 28;7(12):2074-2080.
  35. Bao XZ, Dai F, Wang Q, Jin XL, Zhou B. Developing glutathione-activated catechol-type diphenylpolyenes as small molecule-based and mitochondria-targeted prooxidative anticancer theranostic prodrugs. *FREE RADICAL BIO MED.* 2019 Apr;134:406-418.
  36. Chen P, Zhao M, Chen Q, Fan L, Gao F, Zhao L. Absorption Characteristics of Chitobiose and Chitopentaose in the Human Intestinal Cell Line Caco-2 and Everted Gut Sacs. *J AGR FOOD CHEM.* 2019 Apr 24;67(16):4513-4523.
  37. Tan X, Zhang Y, Wang Q, Ren T, Gou J, Guo W, Yin T, He H, Zhang Y, Tang X. Cell-penetrating peptide together with PEG-modified mesostructured silica nanoparticles promotes mucous permeation and oral delivery of therapeutic proteins and peptides. *BIOMATER SCI-UK.* 2019 Jun 25;7(7):2934-2950.
  38. Lu S, Zhang J, Lian X, Sun L, Meng K, Chen Y, Sun Z, Yin X, Li Y, Zhao J, Wang T, Zhang G, He QY. A hidden human proteome encoded by 'non-coding' genes. *Nucleic Acids Res.* 2019 Sep 5;47(15):8111-8125.
  39. Liang XD, Zhang YN, Liu CC, Chen J, Chen XN, Sattar Baloch A, Zhou B. U18666A inhibits classical swine fever virus replication through interference with intracellular cholesterol trafficking. *Vet Microbiol.* 2019 Nov;238:108436.
  40. Dong C, Tu W, He M, Fu J, Kobayashi A, Konishi T, Shao C. Role of Endoplasmic Reticulum and Mitochondrion in Proton Microbeam Radiation-Induced Bystander Effect. *Radiat Res.* 2020 Jan;193(1):63-72.
  41. Xia Q, Wang X, Liu Y, Shen Z, Ge Z, Huang H, Li X, Wang Y. An endoplasmic reticulum-targeted two-photon fluorescent probe for bioimaging of HClO generated during sleep deprivation. *SPECTROCHIM ACTA A.* 2020 Mar 15;229:117992.
  42. Zhaoyang Wang, Jiaoyan Jia, Lu Wang, Feng Li, Yiliang Wang, Yuzhou Jiang, Xiaowei Song, Shurong Qin, Kai Zheng, Ju Ye, Zhe Ren, Yifei Wang, Shuhua Qi. Anti-HSV-1 activity of Aspergillipeptide D, a cyclic pentapeptide isolated from fungus *Aspergillus* sp. SCSIO 41501. *Virology J.* 2020 Mar 19;17(1):41.
  43. Ru-Jia Xie, Xiao-Xia Hu, Lu Zheng, Shuang Cai, Yu-Si Chen, Yi Yang, Ting Yang, Bing Han, Qin Yang. Calpain-2 activity promotes aberrant endoplasmic reticulum stress-related apoptosis in hepatocytes. *WORLD J GASTROENTERO.* 2020 Apr 7;26(13):1450-1462.
  44. Ming-Hong Sun, Xiao-Han Li, Yao Xu, Yi Xu, Zhen-Nan Pan, Shao-Chen Sun. Citrinin exposure disrupts organelle distribution and functions in mouse oocytes. *Environ Res.* 2020 Jun;185:109476.
  45. Zhi-Fei Liu, Jian-Fei Ji, Xiao-Feng Jiang, Tong Shao, Dong-Dong Fan, Xin-Hang Jiang, Ai-Fu Lin, Li-Xin Xiang, Jian-Zhong Shao. Characterization of cGAS homologs in innate and adaptive mucosal immunities in zebrafish gives evolutionary insights into cGAS-STING pathway. *FASEB J.* 2020 Jun;34(6):7786-7809.
  46. Zhen-Nan Pan, Meng-Hao Pan, Ming-Hong Sun, Xiao-Han Li, Yu Zhang, Shao-Chen Sun. RAB7 GTPase regulates actin dynamics for

- DRP1-mediated mitochondria function and spindle migration in mouse oocyte meiosis. *FASEB J.* 2020 Jul;34(7):9615-9627.
47. Zhiyong Xiong, Changfei Yuan, Jian Shi, Wei Xiong, Yu Huang, Wen Xiao, Hongmei Yang, Ke Chen, Xiaoping Zhang. Restoring the epigenetically silenced PCK2 suppresses renal cell carcinoma progression and increases sensitivity to sunitinib by promoting endoplasmic reticulum stress. *Theranostics.* 2020 Sep 15;10(25):11444-11461.
  48. Wei Xin, Min Zhang, Yang Yu, Songlin Li, Cui Ma, Junting Zhang, Yuan Jiang, Yiyang Li, Xiaodong Zheng, Lixin Zhang, Xijuan Zhao, Xuzhong Pei, Daling Zhu. BCAT1 binds the RNA-binding protein ZNF423 to activate autophagy via the IRE1-XBP-1-RIDD axis in hypoxic PSMCs. *Cell Death Dis.* 2020 Sep 16;11(9):764.
  49. Fan Xiao, Meixin Gao, Junru Yang, Lingling He, Hongshan Wei. FAM172A Deletion May Enhance Hepatic Steatosis by Promoting ER Stress. *DIGEST DIS SCI.* 2020 Sep 18.
  50. Xiuge Gao, Xiangchun Ruan, Hui Ji, Lin Peng, Yawei Qiu, Dan Yang, Xinhao Song, Chunlei Ji, Dawei Guo, Shanxiang Jiang. Maduramicin triggers methuosis-like cell death in primary chicken myocardial cells. *Toxicol Lett.* 2020 Oct 15;333:105-114.
  51. Yongnan Li, Shuodong Wu. Curcumin inhibits the proteolytic process of SREBP-2 by first inhibiting the expression of S1P rather than directly inhibiting SREBP-2 expression. *Food Sci Nutr.* 2020 Nov 8;9(1):209-216.
  52. Yajun Lu, Yuxiao Yang, Siqi Yang, Qianfeng Xia. Immunomodulatory action of excretory-secretory products of *Angiostrongylus cantonensis* in a mouse tumour model. *Parasitol Res.* 2020 Nov;119(11):3705-3718.
  53. Yaxian Zheng, Liyun Xing, Liqiang Chen, Rui Zhou, Jiawei Wu, Xi Zhu, Lian Li, Yucheng Xiang, Ruinan Wu, Ling Zhang, Yuan Huang. Tailored elasticity combined with biomimetic surface promotes nanoparticle transcytosis to overcome mucosal epithelial barrier. *Biomaterials.* 2020 Dec;262:120323.
  54. Yan Zhang, Ye Wang, Xie'an Feng, Shuo Zhang, Xueqiang Xu, Lingyu Li, Shudong Niu, Yingnan Bo, Chao Wang, Zhen Li, Guoliang Xia, Hua Zhang. Oocyte-derived microvilli control female fertility by optimizing ovarian follicle selection in mice. *Nat Commun.* 2021 May 5;12(1):2523.
  55. Yi Li, Yu-Fan Feng, Xiao-Tian Liu, Yu-Chen Li, Hui-Min Zhu, Meng-Ru Sun, Ping Li, Baolin Liu, Hua Yang. Songorine promotes cardiac mitochondrial biogenesis via Nrf2 induction during sepsis. *Redox Biol.* 2021 Jan;38:101771.
  56. Lijie Zhou, Zhengshuai Song, Junyi Hu, Lilong Liu, Yaxin Hou, Xiaoping Zhang, Xiong Yang, Ke Chen. ACSS3 represses prostate cancer progression through downregulating lipid droplet-associated protein PLIN3. *Theranostics.* 2021 Jan 1;11(2):841-860.
  57. Yiqing Tan, Ran Sun, Lei Liu, Dejuan Yang, Qin Xiang, Li Li, Jun Tang, Zhu Qiu, Weiyan Peng, Yuanyuan Wang, Lin Ye, Guosheng Ren, Tingxiu Xiang. Tumor suppressor DRD2 facilitates M1 macrophages and restricts NF- $\kappa$ B signaling to trigger pyroptosis in breast cancer. *Theranostics.* 2021 Mar 5;11(11):5214-5231.
  58. Runyu Zhang, Hailiang Deng, Yuxing Lin, Xing Wang, Bing He, Wenbing Dai, Hua Zhang, Ying Zheng, Qiang Zhang, Xueqing Wang. A common strategy to improve transmembrane transport in polarized epithelial cells based on sorting signals: Guiding nanocarriers to TGN rather than to the basolateral plasma membrane directly. *J Control Release.* 2021 Nov 10;339:430-444.
  59. Xiaoyu Wang, Yu Zhuang, Yukun Fang, Huabin Cao, Caiying Zhang, Chenghong Xing, Xiaoquan Guo, Guyue Li, Ping Liu, Guoliang Hu, Fan Yang. Endoplasmic reticulum stress aggravates copper-induced apoptosis via the PERK/ATF4/CHOP signaling pathway in duck renal tubular epithelial cells. *Environ Pollut.* 2021 Mar 1;272:115981.
  60. Guanzheng Liu, Jiefeng Yu, Runqiu Wu, Lin Shi, Xu Zhang, Wanhong Zhang, Xiaomin Zhong, Yifeng Wang, Huan Li, Yang Shen, Changyong Wu, Rutong Yu, Mingshan Niu, Xuejiao Liu. GRP78 determines glioblastoma sensitivity to UBA1 inhibition-induced UPR signaling and cell death. *Cell Death Dis.* 2021 Jul 23;12(8):733.
  61. Qi Yang, Rongjuan Pei, Yun Wang, Yuan Zhou, Min Yang, Xinwen Chen, Jizheng Chen. ADAM15 Participates in Tick-Borne Encephalitis Virus Replication. *J Virol.* 2021 Jan 28;95(4):e01926-20.
  62. Sijing Hu, Hao Zhang, Yunqiang Liu, Mohan Liu, Jingjing Li, Shun Yao Liao. Whole-exome sequencing of de novo genetic variants in a Chinese family with a sporadic case of congenital nonsyndromic hearing loss. *F1000Res.* 2021 Feb 2;10:61.
  63. Yan Hui Yang, Heng Yang, Rui Fang Li, Cui Xiang Li, Lei Zeng, Chao Jie Wang, Na Li, Zhuang Luo. A *Rehmannia glutinosa* cinnamate 4-hydroxylase promotes phenolic accumulation and enhances tolerance to oxidative stress. *Plant Cell Rep.* 2021 Feb;40(2):375-391.
  64. Xueling Su, Xi Yan, Xing Chen, Mei Guo, Yuxian Xia, Yueqing Cao. Calcofluor white hypersensitive proteins contribute to stress tolerance and pathogenicity in entomopathogenic fungus, *Metarhizium acridum*. *Pest Manag Sci.* 2021 Apr;77(4):1915-1924.
  65. Kangyong Chen, Jiayin Tian, Junya Wang, Zhao Jia, Qin Zhang, Wenji Huang, Xin Zhao, Zhipeng Gao, Qian Gao, Jun Zou. Lipopolysaccharide-induced TNF  $\alpha$  factor (LITAF) promotes inflammatory responses and activates apoptosis in zebrafish *Danio rerio*. *Gene.* 2021 May 15;780:145487.
  66. Qi Peng, Yan Liu, Xuehua Kong, Jie Xian, Lin Ye, Li Yang, Shuliang Guo, Yan Zhang, Lan Zhou, Tingxiu Xiang. The Novel Methylation Biomarker SCARA5 Sensitizes Cancer Cells to DNA Damage Chemotherapy Drugs in NSCLC. *Front Oncol.* 2021 Jun 4;11:666589.
  67. Ya-Yun Liu, Xiao-Dong Liang, Chun-Chun Liu, Yan Cheng, Huan Chen, Abdul Sattar Baloch, Jin Zhang, Yun Young Go, Bin Zhou. Fatty Acid Synthase Is Involved in Classical Swine Fever Virus Replication by Interaction with NS4B. *J Virol.* 2021 Aug 10;95(17):e0078121.
  68. Jinming Zhang, Wenshan Zhong, Yuanyuan Liu, Weimou Chen, Ye Lu, Zhaojin Zeng, Yujie Qiao, Haohua Huang, Xuan Wan, Wei Li, Xiaojing Meng, Fei Zou, Shaoxi Cai, Hangming Dong. Extracellular HSP90 $\alpha$  Interacts With ER Stress to Promote Fibroblasts Activation Through PI3K/AKT Pathway in Pulmonary Fibrosis. *Front Pharmacol.* 2021 Aug 23;12:708462.
  69. Rongjin Luo, Shuai Li, Gaocai Li, Saideng Lu, Weifeng Zhang, Hui Liu, Jie Lei, Liang Ma, Wencan Ke, Zhiwei Liao, Bingjin Wang, Yu Song, Kun Wang, Yukun Zhang, Cao Yang. FAM134B-Mediated ER-phagy Upregulation Attenuates AGEs-Induced Apoptosis and Senescence in Human Nucleus Pulposus Cells. *Oxid Med Cell Longev.* 2021 Aug 5;2021:3843145.
  70. Jinglun Zhang, Qian Wu, Chengcheng Yin, Xiaoshi Jia, Zifan Zhao, Xiaoxin Zhang, Guohua Yuan, Hao Hu, Qin Zhao. Sustained calcium ion release from bioceramics promotes CaSR-mediated M2 macrophage polarization for osteoinduction. *J Leukoc Biol.* 2021 Sep;110(3):485-496.
  71. Lin Meng, Hongmei Hu, Zhiqiang Liu, Luyao Zhang, Qingrui Zhuan, Xue Li, Xiangwei Fu, Shien Zhu, Yunpeng Hou. The Role of Ca<sup>2+</sup> in Maturation and Reprogramming of Bovine Oocytes: A System Study of Low-Calcium Model. *Front Cell Dev Biol.* 2021 Oct 26;9:746237.
  72. Zhihao Liu, Xiyun Bian, Wenbo Gao, Jing Su, Chuanrui Ma, Xiaolin Xiao, Tian Yu, Han Zhang, Xiaozhi Liu, Guanwei Fan. Rg3 promotes the SUMOylation of SERCA2a and corrects cardiac dysfunction in heart failure. *Pharmacol Res.* 2021 Oct;172:105843.
  73. Yi-Ming Li, Shao-Yang Zhao, Huan-Huan Zhao, Bao-Hua Wang,

- Sai-Mei Li. Procyanidin B2 Alleviates Palmitic Acid-Induced Injury in HepG2 Cells via Endoplasmic Reticulum Stress Pathway. *Evid Based Complement Alternat Med.* 2021 Dec 16;2021:8920757.
74. Lin Lei, Wenwen Gao, Juan J Loor, Ahmad Aboragah, Zhiyuan Fang, Xiliang Du, Min Zhang, Yuxiang Song, Guowen Liu, Xinwei Li. Reducing hepatic endoplasmic reticulum stress ameliorates the impairment in insulin signaling induced by high levels of  $\beta$ -hydroxybutyrate in bovine hepatocytes. *J Dairy Sci.* 2021 Dec;104(12):12845-12858.
  75. Shaojuan Song, Xin Xia, Jiajia Qi, Xiaopei Hu, Qian Chen, Jiang Liu, Ning Ji, Hang Zhao. Silmitasertib-induced macropinocytosis promoting DDP intracellular uptake to enhance cell apoptosis in oral squamous cell carcinoma. *Drug Deliv.* 2021 Dec;28(1):2480-2494.
  76. Lin-Lin Hu, Hong-Ge Li, Xiao-Mei Li, Yi Xu, Ya-Qin Pang, Bin Wang, Jun-Li Wang, Shao-Chen Sun. Nonylphenol exposure-induced oocyte quality deterioration could be reversed by melatonin supplementation in mice. *Environ Pollut.* 2022 Jul 15;305:119317.
  77. Chen Xie, Liu Yang, Guixian Jia, Kang Yan, Shizhong Zhang, Guodong Yang, Changai Wu, Yingping Gai, Chengchao Zheng, Jinguang Huang. Maize HEAT UP-REGULATED GENE 1 plays vital roles in heat stress tolerance. *J Exp Bot.* 2022 Oct 18;73(18):6417-6433.
  78. Hao Zhao, Yue Xu, Xinlei Song, Qingchen Zhang, Yajie Wang, Haiyan Yin, Xiaohui Bai, Jianfeng Li. Cisplatin induces damage of auditory cells: Possible relation with dynamic variation in calcium homeostasis and responding channels. *Eur J Pharmacol.* 2022 Jan 5;914:174662.
  79. Jianzhao Liao, Zhuoying Hu, Quanwei Li, Hongji Li, Weijin Chen, Haihua Huo, Qingyue Han, Hui Zhang, Jianying Guo, Lianmei Hu, Jiaqiang Pan, Ying Li, Zhaoxin Tang. Endoplasmic Reticulum Stress Contributes to Copper-Induced Pyroptosis via Regulating the IRE1 $\alpha$ -XBP1 Pathway in Pig Jejunal Epithelial Cells. *J Agric Food Chem.* 2022 Feb 2;70(4):1293-1303.
  80. Meng-Meng Shan, Yuan-Jing Zou, Zhen-Nan Pan, Hao-Lin Zhang, Yi Xu, Jia-Qian Ju, Shao-Chen Sun. Kinesin motor KIF1C1 is required for tubulin acetylation and actin-dependent spindle migration in mouse oocyte meiosis. *Development.* 2022 Mar 1;149(5):dev200231.
  81. Lin Zhu, Bin Li, Dongqin Chen, Ning Chen, Le Xu, Qinjin Li, Xiaodong Chen. sSTEAP4 regulates cellular homeostasis and improves high-fat-diet-caused oxidative stress in hepatocytes. *Life Sci.* 2022 May 1:296:120438.
  82. Fangjiao Lv, Yang Xu, Dean W Gabriel, Xue Wang, Ning Zhang, Wenxing Liang. Quantitative Proteomic Analysis Reveals Important Roles of the Acetylation of ER-Resident Molecular Chaperones for Conidiation in *Fusarium oxysporum*. *Mol Cell Proteomics.* 2022 May;21(5):100231.
  83. Yu-Shen Luo, Quan-Kuo He, Ming-Xin Sun, Feng-Xin Qiao, Yue-Cen Liu, Chang-Long Xu, Zhi-Ran Xu, Si-Cheng Zhao, Hai-Long Wang, Zhong-Quan Qi, Yu Liu. Acrylonitrile exposure triggers ovarian inflammation and decreases oocyte quality probably via mitochondrial dysfunction induced apoptosis in mice. *Chem Biol Interact.* 2022 Jun 1:360:109934.
  84. Shuanghang Liu, Tingxian Deng, Liping Hua, Xinzhe Zhao, Hanxiao Wu, Peihao Sun, Mingxiao Liu, Shujun Zhang, Liguang Yang, Aixin Liang. Novel functional mutation of the PDIA3 gene affects milk composition traits in Chinese Holstein cattle. *J Dairy Sci.* 2022 Jun;105(6):5153-5166.
  85. Yunyi Liu, Juan Li, Hailong Ou, Dan Qi, Bei Hu, Yuxi Xu, Jian Hu, Yi Xiong, Luling Xia, Jason H Huang, Xiaoxiao Hu, Erxi Wu. Identification of new aptamer BC-3 targeting RPS7 from rapid screening for bladder carcinoma. *Genes Dis.* 2022 Aug 3;10(5):2137-2150.
  86. Huidi Xie, Yang Shi, Ying Zhou, Hongfang Liu. TMBIM6 promotes diabetic tubular epithelial cell survival and albumin endocytosis by inhibiting the endoplasmic reticulum stress sensor, IRE1 $\alpha$ . *Mol Biol Rep.* 2022 Oct;49(10):9181-9194.
  87. Yue Li, Ran Yao, Miao Ren, Ke Yuan, Yuwei Du, Yuan He, Haiquan Kang, Shengnan Yuan, Wen Ju, Jianlin Qiao, Kailin Xu, Lingyu Zeng. Liposomes trigger bone marrow niche macrophage "foam" cell formation and affect hematopoiesis in mice. *J Lipid Res.* 2022 Oct;63(10):100273.
  88. Yan Li, Hong-Ying Li, Jun Shao, Lingpeng Zhu, Tian-Hua Xie, Jiping Cai, Wenjuan Wang, Meng-Xia Cai, Zi-Li Wang, Yong Yao, Ting-Ting Wei. GRP75 Modulates Endoplasmic Reticulum-Mitochondria Coupling and Accelerates Ca<sup>2+</sup>-Dependent Endothelial Cell Apoptosis in Diabetic Retinopathy. *Biomolecules.* 2022 Nov 29;12(12):1778.
  89. Han Zhang, Yunfeng Tian, Xiaoya Yuan, Fei Xie, Siqi Yu, Jiayou Cai, Bin Sun, Changliang Shan, Weicheng Zhang. Site-directed late-stage diversification of macrocyclic nannocystins facilitating anticancer SAR and mode of action studies. *RSC Med Chem.* 2022 Dec 22;14(2):299-312.
  90. Zhao-Ran Wang, Ting-Jian Zhang, Qiu-Yin Wang, En-Yu Xu, Xu Zhang, Zhen-Hao Zhang, Peng-Fei Lu, Hai-Yang Zhao, Lin Wang, Fan-Hao Meng. (E)-2-styrylanthracene-9,10-dione derivatives as novel fluorescent probes: synthesis, photophysical properties and application in mitochondria imaging. *Spectrochim Acta A Mol Biomol Spectrosc.* 2023 Feb 5:286:121988.

Version 2024.03.12