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1. Yingying Shi, Yichao Lu, Chunqi Zhu, Zhenyu Luo, Xiang Li, Yu Liu, Mengshi Jiang, Xu Liu, Lihua Luo, Yongzhong Du, Jian You. Targeted regulation of lymphocytic ER stress response with an overall immunosuppression to alleviate allograft rejection. *Biomaterials*. 2021 May;272:120757.
2. Liu Chen, Zheng Ni, Jionggang Hua, Weicheng Ye, Keshu Liu, Tao Yun, Yinchu Zhu, Cun Zhang. Simultaneous tracking of capsid VP26, envelope protein gC localization in living cells infected with double fluorescent duck enteritis virus. *Virus Res*. 2021 May;297:198393.
3. Qi Peng, Yan Liu, Xuehua Kong, Jie Xian, Lin Ye, Li Yang, Shuliang Guo, Yan Zhang, Lan Zhou, Tingxiu Xiang. The Novel Methylation Biomarker SCARA5 Sensitizes Cancer Cells to DNA Damage Chemotherapy Drugs in NSCLC. *Front Oncol*. 2021 Jun 4;11:666589.
4. Caixia Jin, Qingjian Ou, Jie Chen, Tao Wang, Jieping Zhang, Zhe Wang, Yuanyuan Wang, Haibin Tian, Jing-Ying Xu, Furong Gao, Juan Wang, Jiao Li, Lixia Lu, Guo-Tong Xu. Chaperone-mediated autophagy plays an important role in regulating retinal progenitor cell homeostasis. *Stem Cell Res Ther*. 2022 Apr 1;13(1):136.
5. Delong Jiao, Jing Wang, Wenting Yu, Ke Zhang, Ning Zhang, Lingyan Cao, Xinquan Jiang, Yuxing Bai. Biocompatible reduced graphene oxide stimulated BMSCs induce acceleration of bone remodeling and orthodontic tooth movement through promotion on osteoclastogenesis and angiogenesis. *Bioact Mater*. 2022 Feb 6;15:409-425.
6. Hong Wang, Guoxin Jing, Jintong Niu, Li Yang, Youyuan Li, Yi Gao, Huichao Wang, Xiaorong Xu, Yechang Qian, Shilong Wang. A mitochondria-anchored supramolecular photosensitizer as a pyroptosis inducer for potent photodynamic therapy and enhanced antitumor immunity. *J Nanobiotechnology*. 2022 Dec 3;20(1):513.

Version 2024.03.12