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1. Zhu K, Guo C, Lai H, Yang W, Wang C. Novel hyperbranched polyamidoamine nanoparticle based gene delivery: transfection, cytotoxicity and in vitro evaluation. *Int J Pharm.* 2012 Feb 28;423(2):378-83.
2. Ou J, Pan F, Geng P, Wei X, Xie G, Deng J, Pang X, Liang H. Silencing fibronectin extra domain A enhances radiosensitivity in nasopharyngeal carcinomas involving an FAK/Akt/JNK pathway. *Int J Radiat Oncol Biol Phys.* 2012 Mar 15;82(4):e685-91.
3. Shi S, Deng YZ, Zhao JS, Ji XD, Shi J, Feng YX, Li G, Li JJ, Zhu D, Koeffler HP, Zhao Y, Xie D. RACK1 promotes non-small-cell lung cancer tumorigenicity through activating sonic hedgehog signaling pathway. *J Biol Chem.* 2012 Mar 9;287(11):7845-58.
4. Zhu K, Lai H, Guo C, Xu D, Wang C. Novel vascular endothelial growth factor gene delivery system-manipulated mesenchymal stem cells repair infarcted myocardium. *Exp Biol Med (Maywood).* 2012 Jun 1;237(6):678-87.
5. Ou J, Luan W, Deng J, Sa R, Liang H. αV Integrin Induces Multicellular Radioresistance in Human Nasopharyngeal Carcinoma via Activating SAPK/JNK Pathway. *PLoS One.* 2012;7(6):e38737.
6. Zhu K, Guo C, Xia Y, Lai H, Yang W, Wang Y, Song D, Wang C. Transplantation of novel vascular endothelial growth factor gene delivery system manipulated skeletal myoblasts promote myocardial repair. *Int J Cardiol.* 2013 Oct 3;168(3):2622-31.
7. Li H, Qin J, Jin G, Zou L, Shi J, Han X, Cheng X, Zhang X. Overexpression of Lhx8 inhibits cell proliferation and induces cell cycle arrest in PC12 cell line. *In Vitro Cell Dev Biol Anim.* 2015 Apr;51(4):329-35.
8. Xu R, Bi C, Song J, Wang L, Ge C, Liu X, Zhang M. Upregulation of miR-142-5p in atherosclerotic plaques and regulation of oxidized low-density lipoprotein-induced apoptosis in macrophages. *Mol Med Rep.* 2015 May;11(5):3229-34.
9. Peilin Huang, Baozhen Zeng, Zhuoxian Mai, Juntao Deng, Yueping Fang, Wenhua Huang, Hongwu Zhang, Jinying Yuan, Yen Wei and Wuyi Zhou. Novel drug delivery nanosystems based on outside bifunctionalized mesoporous silica yolk-shell magnetic nanostars used as nanocarriers for curcumin. *J. Mater. Chem. B.*

- 2016;4(1):46-56.
10. Dong SM, Cui JH, Zhang W, Zhang XW, Kou TC, Cai QC, Xu S, You S1, Yu DS, Ding L, Lai JH, Li M, Luo KJ. Inhibition of translation initiation factor eIF4A is required for apoptosis mediated by *Microplitisbicoloratus* bracovirus. *ARCH INSECT BIOCHEM.* 2017 Nov;96(3).
 11. Mai Z, Chen J, Hu Y, Liu F, Fu B, Zhang H, Dong X, Huang W, Zhou W. Novel functional mesoporous silica nanoparticles loaded with Vitamin E acetate as smartplatforms for pH responsive delivery with high bioactivity. *J COLLOID INTERF SCI.* 2017 Dec 15;508:184-195.
 12. Hua J,Chen H,Chen Y,Zheng G,Li F,Qu J,Ma X,Hou L. MITF acts as an anti-oxidant transcription factor to regulate mitochondrial biogenesis and redox signaling in retinal pigment epithelial cells. *Exp Eye Res.* 2018 May;170:138-147.
 13. Zhang T,Ni S,Luo Z,Lang Y,Hu J,Lu H. The protective effect of microRNA-21 in neurons after spinal cord injury. *Spinal Cord.* 2019 Feb;57(2):141-149.
 14. Tian HY,Hu Y,Zhang P,Xing WX,Xu C,Yu D,Yang Y,Luo K,Li M. *Spodoptera litura* cyclophilin A is required for *Microplitis bicoloratus* bracovirus-induced apoptosis during insect cellular immune response. *ARCH INSECT BIOCHEM.* 2019 Mar;100(3).
 15. Shi Q,Dai J,Huang L. microRNA-29a functions as a tumor suppressor in nasopharyngeal carcinoma 5-8F cells through targeting VEGF. *Iran J Basic Med Sci.* 2019 May;22(5):541-546.
 16. Cui JH,Dong SM,Chen CX,Xiao W,Cai QC,Zhang LD,He HJ,Zhang W,Zhang XW,Liu T,Ding L,Yang Y,Lai JH,Zhu QS,Luo KJ. *Microplitis bicoloratus* bracovirus modulates innate immune suppression through the eIF4E-eIF4A axis in the insect *Spodoptera litura*. *Dev Comp Immunol.* 2019 Jun;95:101-107.
 17. Shen L,Bing T,Zhang N,Wang L,Wang J,Liu X,Shangguan D. A Nucleus-Targeting DNA Aptamer for Dead Cell Indication. *ACS Sens.* 2019 Jun 28;4(6):1612-1618.
 18. Wang K,Chu D,Wu J,Zhao M,Zhang M,Li B,Du W,Du J,Guo R. Cinobufagin induced cell apoptosis and protective autophagy through the ROS/MAPK signaling pathway. *Life Sci.* 2019 Jul 10:116642.
 19. Liu C,Xing H,Guo C,Yang Z,Wang Y,Wang Y. MiR-124 reversed the doxorubicin resistance of breast cancer stem cells through STAT3/HIF-1 signaling pathways. *Cell Cycle.* 2019 Sep;18(18):2215-2227.
 20. Qi Z,Chen M,Song Y,Wang X,Li B,Chen ZF,Tsang SY,Cai Z. Acute exposure to triphenyl phosphate inhibits the proliferation and cardiac differentiation of mouse embryonic stem cells and zebrafish embryos. *J Cell Physiol.* 2019 Nov;234(11):21235-21248.
 21. Zhao HY,Zhang ST,Cheng X,Li HM,Zhang L,He H,Qin JB,Zhang WY,Sun Y,Jin GH. Long non-coding RNA GAS5 promotes PC12 cells differentiation into Tuj1-positive neuron-like cells and induces cell cycle arrest. *Neural Regen Res.* 2019 Dec;14(12):2118-2125.
 22. Yang F,Zhu B,Liu J,Liu Y,Jiang C,Sheng Q,Qiu J,Nie Z. The effect of acetylation on the protein stability of BmApoLp-III in the silkworm, *Bombyx mori*. *Insect Mol Biol.* 2020 Feb;29(1):104-111.
 23. Wu C,Miao C,Tang Q,Zhou X,Xi P,Chang P,Hua L,Ni H. MiR-129-5p promotes docetaxel resistance in prostate cancer by down-regulating CAMK2N1 expression. *J Cell Mol Med.* 2020 Feb;24(3):2098-2108.
 24. Jingru Zhang,Ling Wang,Jing Jiang,Zhiwei Qiao. Elevation of microRNA-512-5p inhibits MUC1 to reduce radioresistance in cervical cancer. *Cell Cycle.* 2020 Mar;19(6):652-665.
 25. Yan Zhang,Zhen Li,Zhi Lan. Silencing UNC5B antisense lncRNA 1 represses growth and metastasis of human Colon cancer cells via raising miR-622. *ARTIF CELL NANOMED B.* 2020 Dec;48(1):60-67.
 26. Yan Cao,Rongyang Li,Wei Jian Li,Honglin Liu,Yafei Cai. Melatonin Attenuates Peroxynitrite-Induced Meiosis Dysfunction in Porcine Oocytes. *Reprod Sci.* 2021 May;28(5):1281-1289.
 27. Chen Jing,Li Yan,Zhang Wei,Zhang Shoumin,Yin Guangwen,Zhang Jiangan,Jia Xuesong,Chen Hongxiang,Duan Ziyu,Li Jianguo. Exogenous delivery of microRNA-134 (miR-134) using α -tocopherol-based PEGylated liposome for effective treatment in skin squamous cell carcinoma. *DRUG DELIV TRANSL RE.* 2021 Jun;11(3):1000-1008.
 28. Hua Zhang, Yanfang Wang, Xiaolei Bian, Hang Yin. MicroRNA-194 acts as a suppressor during abdominal aortic aneurysm via inhibition of KDM3A-mediated BNIP3. *Life Sci.* 2021 Jul 15;277:119309.
 29. Wei Liu, Shu Tan, Xiaoxu Bai, Shihong Ma, Xiuwei Chen. Long non-coding RNA LINC01215 promotes epithelial-mesenchymal transition and lymph node metastasis in epithelial ovarian cancer through RUNX3 promoter methylation. *Transl Oncol.* 2021 Aug;14(8):101135.
 30. Ruixuan Liu, Zhibin Peng, Yubo Zhang, Rui Li, Yansong Wang. Upregulation of miR-128 inhibits neuronal cell apoptosis following spinal cord injury via FasL downregulation by repressing ULK1. *Mol Med Rep.* 2021 Sep;24(3):667.
 31. Xian Long, Yongpan Huang, Jianbin He, Xiang Zhang, Yan Zhou, Yingmin Wei, Ying Tang, Lijing Liu. Upregulation of miR-335 exerts protective effects against sepsis-induced myocardial injury. *Mol Med Rep.* 2021 Nov;24(5):806.
 32. Yuting Zhang, Peifang Qin, Xingfeng Xu, Mao Li, Haitao Huang, Jianguo Yan, Yali Zhou. Mediator Complex Subunit 19 Promotes the Development of Hepatocellular Carcinoma by Regulating the AKT/mTOR Signaling Pathway. *Front Oncol.* 2022 Jan 3;11:792285.
 33. Jing Xia, Changbin Chu, Wanqing Li, Hong Chen, Wenhua Xie, Rui Cheng, Kai Hu, Xi Li. Mitochondrial Protein UCP1 Inhibits the Malignant Behaviors of Triple-negative Breast Cancer through Activation of Mitophagy and Pyroptosis. *Int J Biol Sci.* 2022 Apr 18;18(7):2949-2961.
 34. Ming Zhai, Shiyu Gong, Peipei Luan, Yefei Shi, Wenxin Kou, Yanxi Zeng, Jiayun Shi, Guanye Yu, Jiayun Hou, Qing Yu, Weixia Jian, Jianhui Zhuang, Mark W Feinberg, Wenhui Peng. Extracellular traps from activated vascular smooth muscle cells drive the progression of atherosclerosis. *Nat Commun.* 2022 Dec 6;13(1):7500.

Version 2024.03.12