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- Dane Huang, Chao Zhao, Ruyue Li, Bingyi Chen, Yuting Zhang, Zhejun Sun, Junkang Wei, Huihao Zhou, Qiong Gu, Jun Xu. Identification of a binding site on soluble RANKL that can be targeted to inhibit soluble RANK-RANKL interactions and treat osteoporosis. *Nat Commun.* 2022 Sep 12;13(1):5338.
- Xuan Wang, Yingqi Liu, Chencheng Xue, Yan Hu, Yuanyuan Zhao, Kaiyong Cai, Menghuan Li, Zhong Luo. A protein-based cGAS-STING nanoagonist enhances T cell-mediated anti-tumor immune responses. *Nat Commun.* 2022 Sep 28;13(1):5685.
- Jin-Fei Lin, Pei-Shan Hu, Yi-Yu Wang, Yue-Tao Tan, Kai Yu, Kun Liao, Qi-Nian Wu, Ting Li, Qi Meng, Jun-Zhong Lin, Ze-Xian Liu, Heng-Ying Pu, Huai-Qiang Ju, Rui-Hua Xu, Miao-Zhen Qiu. Phosphorylated NFS1 weakens oxaliplatin-based chemosensitivity of colorectal cancer by preventing PANoptosis. *Signal Transduct Target Ther.* 2022 Feb 28;7(1):54.
- Xiawei Dong, Jing Ye, Yihan Wang, Hongjie Xiong, Hui Jiang, Hongbing Lu, Xiaohui Liu, Xuemei Wang. Ultra-Small and Metabolizable Near-Infrared Au/Gd Nanoclusters for Targeted FL/MRI Imaging and Cancer Theranostics. *Biosensors (Basel).* 2022 Jul 24;12(8):558.
- Huifang Yin, Zhenglei Wang, Shoushen Yang, Xintian Zheng, Yinli Bao, Weiming Lin, Cuiqin Huang, Longxin Qiu. Taurine inhibits necroptosis helps to alleviate inflammatory and injury induced by Klebsiella infection. *Vet Immunol Immunopathol.* 2022 Aug;250:110444.
- Long-Long Zhang, Ya-Jun Liu, Yong-Hong Chen, Zhuang Wu, Bo-Ran Liu, Qian-Yi Cheng, Ke-Qin Zhang, Xue-Mei Niu. Modulating Activity Evaluation of Gut Microbiota with Versatile Toluquinol. *Int J Mol Sci.* 2022 Sep 14;23(18):10700.
- Xiaojin Li, Meng Qi, Kai He, Haiyan Liu, Wenlan Yan, Lizhuo Zhao, Yanyan Jia, Lei He, Chaochao Lv, Min Zhang, Zhiguo Wei, Wenchao Yan, Tianqi Wang, Fuchang Yu, Weifeng Qian. Neospora caninum inhibits tumor development by activating the immune response and destroying tumor cells in a B16F10 melanoma model. *Parasit Vectors.* 2022 Sep 23;15(1):332.
- Zhuqing Tian, Longbing Yang, Mingjiao Huang, Chaoqin Sun, Mingming Chen, Wenjing Zhao, Jian Peng, Guo Guo. Antitumor Activity and Mechanism of Action of the Antimicrobial Peptide AMP-17 on Human Leukemia K562 Cells. *Molecules.* 2022 Nov 21;27(22):8109.
- Weisin Chen, Miersalijiang Yasen, Hanquan Wang, Chenyang Zhuang, Zixiang Wang, Shunyi Lu, Libo Jiang, Hong Lin. Celecoxib activates autophagy

by inhibiting the mTOR signaling pathway and prevents apoptosis in nucleus pulposus cells. BMC Pharmacol Toxicol. 2022 Dec 1;23(1):90.

Version 2024.03.12