碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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1. Ge X,Hua H,Wang P,Liu J,Zhang Y,Ding G,Zhu C,Huang S,Jia Z,Zhang A. Inhibition of mitochondrial complex I by rotenone protects against acetaminophen-induced liver injury. Am J Transl Res. 2019 Jan 15;11(1):188-198. eCollection 2019.
2. Li J,Du H,Zhang M,Zhang Z,Teng F,Zhao Y,Zhang W,Yu Y,Feng L,Cui X,Zhang M,Lu T,Guan F,Chen L. Amorphous solid dispersion of Berberine mitigates apoptosis via iPLA2β/Cardiolipin/Opa1 pathway in db/db mice and in Palmitate-treated MIN6 β-cells. Int J Biol Sci. 2019 Jun 2;15(7):1533-1545.
3. He X,Zhu Y,Zhang Y,Geng Y,Gong J,Geng J,Zhang P,Zhang X,Liu N,Peng Y,Wang C,Wang Y,Liu X,Wan L,Gong F,Wei C,Zhong H. RNF34 functions in immunity and selective mitophagy by targeting MAVS for autophagic degradation. EMBO J. 2019 Jul 15;38(14):e100978.
4. Liu Q,Fang J,Chen P,Die Y,Wang J,Liu Z,Liu X. Chicoric acid improves neuron survival against inflammation by promoting mitochondrial function and energy metabolism. Food Funct. 2019 Sep 1;10(9):6157-6169.
5. Xu G,Shen H,Nibona E,Wu K,Ke X,Al Hafiz MA,Liang X,Zhong X,Zhou Q,Qi C,Zhao H. Fundc1 is necessary for proper body axis formation during embryogenesis in zebrafish. SCI REP-UK. 2019 Dec 11;9(1):18910.
6. Du F,Bai L,He M,Zhang WY,Gu YY,Yin H,Liu YJ. Design, synthesis and biological evaluation of iridium(III) complexes as potential antitumor agents. J Inorg Biochem. 2019 Dec;201:110822.
7. Tao Shu-Ya,Zhang Qiu-Yang,Li Jing-Jing,Yao Jin,Yan Biao. Suppression of pathological ocular neovascularization by a small molecule, SU1498. Biomed Pharmacother. 2020 Aug;128:110248.
8. Fang Li,Jie Zhang,Rui Liao,Yongchun Duan,Lili Tao,Yuwei Xu,Anbao Chen. Mesenchymal stem cell-derived extracellular vesicles prevent neural stem cell hypoxia injury via promoting miR-210-3p expression. Mol Med Rep. 2020 Nov;22(5):3813-3821.

9. Jinghui Zhai, Huan Gao, Shuo Wang, Sixi Zhang, Xiaoyu Qu, Yueming Zhang, Lina Tao, Jingmeng Sun, Yanqing Song, Li Fu. Ginsenoside Rg3 attenuates cisplatin-induced kidney injury through inhibition of apoptosis and autophagy-inhibited NLRP3. *J Biochem Mol Toxicol.* 2021 Nov;35(11):e22896.
10. Kai Meng, Chengbo Fang. Knockdown of Tripartite motif-containing 22 (TRIM22) relieved the apoptosis of lens epithelial cells by suppressing the expression of TNF receptor-associated factor 6 (TRAF6). *Bioengineered.* 2021 Dec;12(1):7213-7222.
11. Yuanyang Li, Yufan Zhang, Xiangzhong Zhou, Xianghong Lei, Xinhang Li, Liping Wei. Dynamic observation of 5-fluorouracil-induced myocardial injury and mitochondrial autophagy in aging rats. *Exp Ther Med.* 2021 Dec;22(6):1451.

Version 2024.03.12