碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-1683301或800-8283301  
 订货e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

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- Li P, Yin Y, Yu Q, Yang Q. Lactobacillus acidophilus S-layer protein-mediated inhibition of Salmonella-induced apoptosis in Caco-2 cells. *Biochem Biophys Res Commun.* 2011 May 27; 409(1):142-7.
- Zhang C, Li C, Jia S, Yao P, Yang Q, Zhang Y. High-mobility group box 1 inhibition alleviates lupus-like disease in BXSB mice. *Scand J Immunol.* 2014 May; 79(5):333-7.
- Luo B, Li B, Wang W, Liu X, Xia Y, Zhang C, Zhang M, Zhang Y, An F. NLRP3 gene silencing ameliorates diabetic cardiomyopathy in a type 2 diabetes rat model. *PLoS One.* 2014 Aug 19; 9(8):e104771.
- Juan Liu, Yanhong Xu, Zheng Zhang and Jianhe Wei. Hydrogen peroxide promotes programmed cell death and salicylic acid accumulation during the induced production of sesquiterpenes in cultured cell suspensions of *Aquilaria sinensis* Functional Plant Biology. 2014 Nov; 42(4):337-46.
- Wang L, Chen K, Wan X, Wang F, Guo Z, Mo Z. NLRP3 inflammasome activation in mesenchymal stem cells inhibits osteogenic differentiation and enhances adipogenic differentiation. *Biochem Biophys Res Commun.* 2017 Mar 18; 484(4):871-877.
- Qiu Z, Lei S, Zhao B, Wu Y, Su W, Liu M, Meng Q, Zhou B, Leng Y, Xia ZY. NLRP3 Inflammasome Activation-Mediated Pyroptosis Aggravates Myocardial Ischemia/Reperfusion Injury in Diabetic Rats. *Oxid Med Cell Longev.* 2017; 2017:9743280.
- Jin Q, Li R, Hu N, Xin T, Zhu P, Hu S, Ma S, Zhu H, Ren J, Zhou H. DUSP1 alleviates cardiac ischemia/reperfusion injury by suppressing the Mff-required mitochondrial fission and Bnip3-related mitophagy via the JNK pathways. *Redox Biol.* 2018 Apr; 14:576-587.
- Zhai Y, Meng X, Ye T, Xie W, Sun G, Sun X. Inhibiting the NLRP3 Inflammasome Activation with MCC950 Ameliorates Diabetic Encephalopathy in db/db Mice. *Molecules.* 2018 Feb 27; 23(3). pii: E522.
- Fan S, Yuan J, Deng S, Chen Y, Xie B, Wu K, Zhu M, Xu H, Huang Y, Yang J, Zhang Y, Chen J, Zhao M. Activation of Interleukin-1 $\beta$  Release by the Classical Swine Fever Virus Is Dependent on the NLRP3 Inflammasome, Which Affects Virus Growth in Monocytes. *Front Cell Infect Microbiol.* 2018 Jul 2; 8:225.
- Wu D, Han R, Deng S, Liu T, Zhang T, Xie H, Xu Y. Protective Effects of Flagellin A N/C Against Radiation-Induced NLR Pyrin Domain Containing 3 Inflammasome-Dependent Pyroptosis in Intestinal Cells. *Int J Radiat Oncol Biol Phys.* 2018 May 1; 101(1):107-117.
- Zhang X, Luan J, Chen W, Fan J, Nan Y, Wang Y, Liang Y, Meng G, Ju D. Mesoporous silica nanoparticles induced hepatotoxicity via NLRP3 inflammasome activation and caspase-1-dependent pyroptosis. *Nanoscale.* 2018 May 17; 10(19):9141-9152.
- Yuan J, Zhu M, Deng S, Fan S, Xu H, Liao J, Li P, Zheng J, Zhao M, Chen J. Classical swine fever virus induces pyroptosis in the

- peripheral lymphoid organs of infected pigs. *Virus Res.* 2018 May 2;250:37-42.
13. Li J, Zhou Y, Liu Y, Dai B, Zhang YH, Zhang PF, Shi XL. Sorafenib inhibits caspase-1 expression through suppressing TLR4/stat3/SUMO1 pathway in hepatocellular carcinoma. *Cancer Biol Ther.* 2018 Oct 2;1-8.
  14. Xu M, Li F, Wang M, Zhang H, Xu L, Adcock IM, Chung KF, Zhang Y. Protective effects of VGX-1027 in PM2.5-induced airway inflammation and bronchial hyperresponsiveness. *Eur J Pharmacol.* 2019 Jan 5;842:373-383.
  15. Liu LB, Shen HF, Cha W, Jin ZJ, Xia HJ, Liu JJ, Hu JF. SXBX pill suppresses homocysteine-induced vascular smooth muscle cells dedifferentiation by inhibiting NLRP3 inflammasomes activation via ERK/p38 MAPK pathways. *Am J Transl Res.* 2019 Feb 15;11(2):806-818. eCollection 2019.
  16. Yang G, Wang J, Luo T, Zhang X. White spot syndrome virus infection activates Caspase 1-mediated cell death in crustacean. *Virology.* 2019 Feb;528:37-47.
  17. Qiu Z, He Y, Ming H, Lei S, Leng Y, Xia ZY. Lipopolysaccharide (LPS) Aggravates High Glucose- and Hypoxia/Reoxygenation-Induced Injury through Activating ROS-Dependent NLRP3 Inflammasome-Mediated Pyroptosis in H9C2 Cardiomyocytes. *J Diabetes Res.* 2019 Feb 17;2019:8151836.
  18. Chen Y, Ma K. NLR4 inflammasome activation regulated by TNF- $\alpha$  promotes inflammatory responses in nonalcoholic fatty liver disease. *BIOCHEM BIOPH RES CO.* 2019 Apr 9;511(3):524-530.
  19. Ye P, Wang X, Ge S, Chen W, Wang W, Han X. Long-term cigarette smoking suppresses NLRP3 inflammasome activation in oral mucosal epithelium and attenuates host defense against *Candida albicans* in a rat model. *Biomed Pharmacother.* 2019 May;113:108597.
  20. Xu M, Zhang Y, Wang M, Zhang H, Chen Y, Adcock IM, Chung KF, Mo J, Zhang Y, Li F. TRPV1 and TRPA1 in Lung Inflammation and Airway Hyperresponsiveness Induced by Fine Particulate Matter (PM2.5). *Oxid Med Cell Longev.* 2019 Jun 2;2019:7450151.
  21. Zhang X, Dong P, Xu L, Tian Y, Sun H, Shi H, Duan Z, Chen L, Ren F. The different expression of caspase-1 in HBV-related liver disease and acts as a biomarker for acute-on-chronic liver failure. *BMC Gastroenterol.* 2019 Aug 20;19(1):148.
  22. Wu Z, Li Y, Liu Q, Liu Y, Chen L, Zhao H, Guo H, Zhu K, Zhou N, Chai TC, Shi B. Pyroptosis engagement and bladder urothelial cell-derived exosomes recruit mast cells and induce barrier dysfunction of bladder urothelium after uropathogenic *E. coli* infection. *AM J PHYSIOL-CELL PH.* 2019 Sep 1;317(3):C544-C555.
  23. Wang S, Fan J, Mei X, Luan J, Li Y, Zhang X, Chen W, Wang Y, Meng G, Ju D. Interleukin-22 Attenuated Renal Tubular Injury in Aristolochic Acid Nephropathy via Suppressing Activation of NLRP3 Inflammasome. *Front Immunol.* 2019 Sep 24;10:2277.
  24. Shen T, Li X, Looor JJ, Zhu Y, Du X, Wang X, Xing D, Shi Z, Fang Z, Li X, Liu G. Hepatic nuclear factor kappa B signaling pathway and NLR family pyrin domain containing 3 inflammasome is over-activated in ketotic dairy cows. *J Dairy Sci.* 2019 Nov;102(11):10554-10563.
  25. Wu G, Zhu Q, Zeng J, Gu X, Miao Y, Xu W, Lv T, Song Y. Extracellular mitochondrial DNA promote NLRP3 inflammasome activation and induce acute lung injury through TLR9 and NF- $\kappa$ B. *J Thorac Dis.* 2019 Nov;11(11):4816-4828.
  26. Zhong X, Liu M, Yao W, Du K, He M, Jin X, Jiao L, Ma G, Wei B, Wei M. Epigallocatechin-3-Gallate Attenuates Microglial Inflammation and Neurotoxicity by Suppressing the Activation of Canonical and Noncanonical Inflammasome via TLR4/NF- $\kappa$ B Pathway. *Mol Nutr Food Res.* 2019 Nov;63(21):e1801230.
  27. Yu LM, Zhang WH, Han XX, Li YY, Lu Y, Pan J, Mao JQ, Zhu LY, Deng JJ, Huang W, Liu YH. Hypoxia-Induced ROS Contribute to Myoblast Pyroptosis during Obstructive Sleep Apnea via the NF- $\kappa$ B/HIF-1 $\alpha$  Signaling Pathway. *Oxid Med Cell Longev.* 2019 Dec 11;2019:4596368.
  28. Liu Y, He H, Fan L, Yuan J, Huang H, Yang W, Wang L, Mo Z, Wang F. Compound C attenuates NLRP3 inflammasome despite AMPK knockdown in LPS plus palmitate-induced THP-1 cells. *N-S ARCH PHARMACOL.* 2020 Jan;393(1):67-76.
  29. Feng Yao, Zhen Jin, Zihan Zheng, Xiaohan Lv, Lingxuan Ren, Jianjun Yang, Danli Chen, Bo Wang, Wei Yang, Lifang Chen, Weirong Wang, Jianli Gu, Rong Lin. HDAC11 promotes both NLRP3/caspase-1/GSDMD and caspase-3/GSDME pathways causing pyroptosis via ERG in vascular endothelial cells. *Cell Death Discov.* 2022 Mar 12;8(1):112.
  30. Jingying Gao, Lixia Xia, Yuanyuan Wei. Oxymatrine inhibits the pyroptosis in rat insulinoma cells by affecting nuclear factor kappa B and nuclear factor (erythroid-derived 2)-like 2 protein/heme oxygenase-1 pathways. *Korean J Physiol Pharmacol.* 2022 May 1;26(3):165-174.
  31. Jinjian Li, Kun Yang, Xinghui Pan, Hui Peng, Chenting Hou, Jie Xiao, Qing Wang. Long Noncoding RNA MIAT Regulates Hyperosmotic Stress-Induced Corneal Epithelial Cell Injury via Inhibiting the Caspase-1-Dependent Pyroptosis and Apoptosis in Dry Eye Disease. *J Inflamm Res.* 2022 Jun 2;15:3269-3283.
  32. Ying Su, Shan Lu, Chenjian Hou, Kehan Ren, Meili Wang, Xiaoli Liu, Shanyu Zhao, Xiuping Liu. Mitigation of liver fibrosis via hepatic stellate cells mitochondrial apoptosis induced by metformin. *Int Immunopharmacol.* 2022 Jul;108:108683.
  33. Yang Li, Cheng Wen, Jialin Zhong, Junqi Ling, Qianzhou Jiang. *Enterococcus faecalis* OG1RF induces apoptosis in MG63 cells via caspase-3/-8/-9 without activation of caspase-1/GSDMD. *Oral Dis.* 2022 Oct;28(7):2026-2035.

Version 2024.03.12