



碧云天生物技术/Beyotime Biotechnology
订货热线: 400-168-3301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

细胞线粒体分离试剂盒

产品编号	产品名称	包装
C3601	细胞线粒体分离试剂盒	50-100次

产品简介:

- 细胞线粒体分离试剂盒(Cell Mitochondria Isolation Kit)是用于快速便捷地分离培养细胞线粒体的试剂盒。
- 本试剂盒在分离线粒体的同时可以获得去除线粒体的细胞浆蛋白, 可用于研究细胞色素c等线粒体蛋白向胞浆的释放。
- 使用本试剂盒分离获得的线粒体纯度较高, 并且绝大部分分离获得的线粒体都含有完整的内膜和外膜, 并具有线粒体的生理功能。因此本试剂盒分离得到的线粒体可以用于线粒体的生理功能等方面的研究。例如可以使用碧云天的C2006线粒体膜电位检测试剂盒(JC-1)测定分离得到的线粒体的膜电位。
- 本试剂盒分离得到的线粒体也可以被试剂盒中的线粒体裂解液或其它适当裂解液裂解后用于SDS-PAGE、Western、双向电泳等蛋白分析。
- 本试剂盒提供了线粒体制备过程中匀浆程度的重要判断指标, 即台盼蓝染色, 使分离得到的线粒体的质量更加有保证。台盼蓝染色液为选用试剂, 在实验条件成熟后可以不必使用。
- 本试剂盒提供了蛋白酶抑制剂PMSF, 使匀浆时蛋白酶的活性被适当抑制, 这样在获得线粒体的同时还可以获得有时有用的去除了大量线粒体的蛋白, 在进行蛋白或酶活分析时可以作为对照。
- 如果每个样品的细胞数量为2000-5000万, 本试剂盒可以处理50-100个样品。

包装清单:

产品编号	产品名称	包装
C3601-1	线粒体分离试剂	125ml
C3601-2	台盼蓝染色液	10ml
C3601-3	线粒体储存液	15ml
C3601-4	线粒体裂解液	15ml
C3601-5	PMSF(晶体)	For 1.5ml 100mM
C3601-6	PMSF(溶剂)	1.5ml
—	说明书	1份

保存条件:

-20°C保存, 一年有效。其中台盼蓝染色液也可以4°C保存, PMSF(晶体)和PMSF(溶剂)在配制成100mM PMSF溶液前可以室温保存。

注意事项:

- 试剂盒中的试剂对于不同的实验目的不必全部使用。
- 如果不是用于制备线粒体蛋白样品, 线粒体分离试剂和线粒体裂解液中不必加入PMSF。
- 如果用于制备线粒体蛋白样品, 线粒体分离试剂和线粒体裂解液中需添加PMSF。PMSF一定要在线粒体分离试剂或线粒体裂解液加入到样品中前2-3分钟内加入, 以免PMSF在水溶液中很快失效。
- 分离线粒体的所有步骤均需在冰上或4°C进行, 所用溶液需冰浴或4°C预冷。
- 通常在分离线粒体时前后两次离心速度选取600g和11,000g, 如果希望纯度更高, 但对线粒体的得率要求不高, 前后两次离心速度可以采用1000g和3500g。
- 使用本试剂盒中的线粒体储存液稀释或储存的线粒体样品应及时使用, 以免线粒体膜电位受影响。如果不能及时使用, 建议在-80°C保存。冻存后的线粒体样品不推荐用于膜电位的检测, 但可以用线粒体蛋白或核酸的相关检测。
- 如果出现本试剂盒中线粒体储存液不够用的情况, 可以单独订购线粒体储存液(C3609)。
- 台盼蓝与PMSF对人体有毒, 操作时请特别小心, 并注意有效防护以避免直接接触人体或吸入体内。
- Bradford法蛋白浓度测定试剂盒(P0006)、Bradford蛋白浓度测定试剂盒(去垢剂兼容型)(P0006C)和BCA蛋白浓度测定试剂盒(P0009/P0010/P0010S/P0011/P0012/P0012S)可以向碧云天订购。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用说明:

1. **准备溶液:** 室温融解试剂盒中的各种溶液, 融解后立即置于冰上并混匀。第一次使用时, 把1.5ml PMSF(溶剂)加入到

PMSF(晶体)中, 溶解并混匀, 即得到1.5ml 100mM PMSF。配制好的100mM PMSF溶液-20°C保存。如果最终目的是制备线粒体蛋白样品, 根据样品数量, 取适量线粒体分离试剂备用, 在线粒体分离试剂加入到细胞样品中前数分钟内加入PMSF, 使PMSF的最终浓度为1mM。同时, 根据样品数量取适量线粒体裂解液备用, 在线粒体裂解液加入到线粒体样品中前数分钟内加入PMSF, 使PMSF的最终浓度为1mM。

2. 收集细胞:

对于贴壁细胞: 用PBS洗一遍, 用胰酶细胞消化液(Trypsin-EDTA Solution)消化细胞, 100-200g, 室温离心5-10分钟收集细胞。胰酶细胞消化液(C0201)可以向碧云天订购。

对于悬浮细胞: 直接离心收集细胞。

3. **洗涤细胞:** 用冰浴预冷的PBS轻轻重悬细胞沉淀, 取少量细胞用于计数, 剩余细胞600g, 4°C离心5分钟沉淀细胞。弃上清。

4. **预处理:** 加入1-2.5ml线粒体分离试剂或临用前添加了PMSF的线粒体分离试剂至2000-5000万细胞中, 轻轻悬浮细胞, 冰浴放置10-15分钟。

5. **匀浆:** 把细胞悬液转移到一适当大小的玻璃匀浆器中, 匀浆10-30下左右。

6. **匀浆效果的鉴定:** 不同细胞不同匀浆器所需的匀浆次数有所不同, 需自行优化。通常可以在匀浆10次后取约2微升细胞匀浆, 加入30-50微升台盼蓝染色液, 混匀后显微镜下观察台盼蓝染色阳性(蓝色)细胞的比例。如果阳性细胞比例不足50%, 增加5次匀浆。随后再同前取样进行台盼蓝染色鉴定。当阳性比例超过50%时即可停止匀浆进入下一步。请勿过度匀浆, 过度匀浆会导致线粒体的机械损伤。同时记录对于该细胞的匀浆次数, 通常在后续实验时不必再摸索匀浆次数。

7. 把细胞匀浆在600g, 4°C离心10分钟。

注: 如需获得纯度更高的线粒体, 可以将此步骤的离心速度改为1000g, 其它离心条件不变; 获得更高纯度线粒体的缺点是相同数量细胞的线粒体抽提得率会下降。

8. 小心把上清转移到另一离心管中, 在11,000g, 4°C离心10分钟。

注: 如需获得纯度更高的线粒体, 可以将此步骤的离心速度改为3500g, 其它离心条件不变; 获得更高纯度线粒体的缺点是相同数量细胞的线粒体抽提得率会下降。

9. 小心去除上清。沉淀即为分离得到的细胞线粒体。

注: 如果希望获得去除线粒体的细胞浆蛋白, 在本步骤需收集上清, 并且在收集上清时注意勿触及沉淀。随后把收集的上清12,000g, 4°C离心10分钟。取上清即为去除线粒体的细胞浆蛋白。该细胞浆蛋白可以使用BCA法或Bradford法测定蛋白浓度。

10. 线粒体的使用:

a. 如果用于完整线粒体的功能或酶活性研究, 初始2000-5000万细胞分离得到的线粒体样品中可以加入150-200微升线粒体储存液, 重悬线粒体。

b. 如果用于线粒体的蛋白分析, 初始2000-5000万细胞分离得到的线粒体样品中可以加入150-200微升临用前添加了PMSF的线粒体裂解液裂解线粒体。裂解后的线粒体可以用于PAGE、Western、IP以及线粒体中的一些酶活性的测定等。裂解后的蛋白样品可以使用BCA法或去垢剂兼容型Bradford法(P0006C)测定蛋白浓度。

c. 如果用于双向电泳, 请使用适当的用于双向电泳的裂解液处理线粒体。

相关产品:

产品编号	产品名称	包装
C2006	线粒体膜电位检测试剂盒(JC-1)	>100次
C3601	细胞线粒体分离试剂盒	50-100次
C3606	组织线粒体分离试剂盒	50-100次
C3609	线粒体储存液	50ml

使用本产品的文献:

- Shen J, Huang C, Jiang L, Gao F, Wang Z, Zhang Y, Bai J, Zhou H, Chen Q. Enhancement of cisplatin induced apoptosis by suberoylanilide hydroxamic acid in human oral squamous cell carcinoma cell lines. *Biochem Pharmacol.* 2007 Jun 15;73(12):1901-9.
- Zhang WC, Peng YJ, He WQ, Lv N, Chen C, Zhi G, Chen HQ, Zhu MS. Identification and functional characterization of an aggregation domain in long myosin light chain kinase. *FEBS J.* 2008 May;275(10):2489-500.
- Yu K, He Y, Yeung LW, Lam PK, Wu RS, Zhou B. DE-71-induced apoptosis involving intracellular calcium and the Bax-mitochondria-caspase protease pathway in human neuroblastoma cells in vitro. *Toxicol Sci.* 2008;104(2):341-51.
- Chen K, Zhang Q, Wang J, Liu F, Mi M, Xu H, Chen F, Zeng K. Taurine protects transformed rat retinal ganglion cells from hypoxia-induced apoptosis by preventing mitochondrial dysfunction. *Brain Res.* 2009 Jul 7;1279:131-8.
- Xu WN, Liu WB, Liu ZP. Trichlorfon-induced apoptosis in hepatocyte primary cultures of *Carassius auratus gibelio*. *Chemosphere.* 2009 Nov;77(7):895-901.
- Zhang Z, Wang S, Qiu H, Duan C, Ding K, Wang Z. Waltonitone induces human hepatocellular carcinoma cells apoptosis in vitro and in vivo. *Cancer Lett.* 2009 Dec 28;286(2):223-31.
- Zhang Y, Zhou L, Bao YL, Wu Y, Yu CL, Huang YX, Sun Y, Zheng LH, Li YX. Butyrate induces cell apoptosis through activation of JNK MAP kinase pathway in human colon cancer RKO cells. *Chem Biol Interact.* 2010;185(3):174-81.
- Huang XJ, Zhang HH, Wang X, Huang LL, Zhang LY, Yan CX, Liu Y, Yuan ZH. ROS mediated cytotoxicity of porcine adrenocortical cells induced by QdNOs derivatives in vitro. *Chem Biol Interact.* 2010;185(3):227-34.
- Li Z, Wei H, Liu X, Hu S, Cong X, Chen X. PA rescues ER stress-associated apoptosis in hypoxia and serum deprivation-stimulated mesenchymal stem cells. *J Cell Biochem.* 2010 Nov 1;111(4):811-20.
- Xiong Y, Fang JH, Yun JP, Yang J, Zhang Y, Jia WH, Zhuang SM. Effects of microRNA-29 on apoptosis, tumorigenicity, and prognosis of hepatocellular carcinoma. *Hepatology.* 2010 Mar;51(3):836-45.

11. Lu B, Hu M, Liu K, Peng J. Cytotoxicity of berberine on human cervical carcinoma HeLa cells through mitochondria, deathreceptor and MAPK pathways, and in-silico drug-target prediction. *Toxicol In Vitro*. 2010 Sep;24(6):1482-90.
12. Liu H, Xiao Y, Xiong C, Wei A, Ruan J. Apoptosis induced by a new flavonoid in human hepatoma HepG2 cells involves reactive oxygenspecies-mediated mitochondrial dysfunction and MAPK activation. *Eur J Pharmacol*. 2011 Mar 11;654(3):209-16.
13. Zhang Y, Li RJ, Ying X, Tian W, Yao HJ, Men Y, Yu Y, Zhang L, Ju RJ, Wang XX, Zhou J, Chen JX, Li N, Lu WL. Targeting therapy with mitosomal daunorubicin plus amlodipine has the potential to circumventintrinsic resistant breast cancer. *Mol Pharm*. 2011 Feb 7;8(1):162-75.
14. Shu D, Qing Y, Tong Q, He Y, Xing Z, Zhao Y, Li Y, Wei Y, Huang W, Wu X. Deltonin isolated from *Dioscorea zingiberensis* inhibits cancer cell growth through inducingmitochondrial apoptosis and suppressing Akt and mitogen activated protein kinase signals. *Biol Pharm Bull*. 2011;34(8):1231-9.
15. Liu W, Fan Z, Han Y, Lu S, Zhang D, Bai X, Xu W, Li J, Wang H. Curcumin attenuates peroxynitrite-induced neurotoxicity in spiral ganglion neurons. *Neurotoxicology*. 2011 Jan;32(1):150-7.
16. Li L, Lu N, Dai Q, Wei L, Zhao Q, Li Z, He Q, Dai Y, Guo Q. GL-V9, a newly synthetic flavonoid derivative, induces mitochondrial-mediated apoptosis andG2/M cell cycle arrest in human hepatocellular carcinoma HepG2 cells. *Eur J Pharmacol*. 2011 Nov 16;670(1):13-21.
17. Wang Y, Deng L, Zhong H, Wang Y, Jiang X, Chen J. Natural plant extract tubeimoside I promotes apoptosis-mediated cell death in cultured humanhepatoma (HepG2) cells. *Biol Pharm Bull*. 2011;34(6):831-8.
18. Fu K, Ren H, Wang Y, Fei E, Wang H, Wang G. DJ-1 inhibits TRAIL-induced apoptosis by blocking pro-caspase-8 recruitment to FADD. *Oncogene*. 2012 Mar 8;31(10):1311-22.
19. Zhang L, Yao HJ, Yu Y, Zhang Y, Li RJ, Ju RJ, Wang XX, Sun MG, Shi JF, Lu WL. Mitochondrial targeting liposomes incorporating daunorubicin and quinacrine for treatment of relapsed breast cancer arising from cancer stem cells. *Biomaterials*. 2012 Jan;33(2):565-82.
20. Shi ZY, Li YQ, Kang YH, Hu GQ, Huang-fu CS, Deng JB, Liu B. Piperonal ciprofloxacin hydrazone induces growth arrest and apoptosis of humanhepatocarcinoma SMMC-7721 cells. *Acta Pharmacol Sin*. 2012 Feb;33(2):271-8.
21. Cao H, Hu Y, Wang P, Zhou J, Deng Z, Wen J. Down-regulation of Notch receptor signaling pathway induces caspase-dependent and caspase-independent apoptosis in lung squamous cell carcinoma cells. *APMIS*. 2012 Jun;120(6):441-50.
22. Gong K, Xie J, Yi H, Li W. CS055 (Chidamide/HBI-8000), a novel histone deacetylase inhibitor, induces G1 arrest, ROS-dependent apoptosis and differentiation in human leukaemia cells. *Biochem J*. 2012 May 1;443(3):735-46.
23. Wang W, Lu J, Zhu F, Wei J, Jia C, Zhang Y, Zhou L, Xie H, Zheng S. Pro-apoptotic and anti-proliferative effects of mitofusin-2 via Bax signaling in hepatocellularcarcinoma cells. *Med Oncol*. 2012 Mar;29(1):70-6.
24. Li P, Wang K, Zhang J, Zhao L, Liang H, Shao C, Sutherland LC. The 3p21.3 tumor suppressor RBM5 resensitizes cisplatin-resistant human non-small cell lungcancer cells to cisplatin. *Cancer Epidemiol*. 2012 Oct;36(5):481-9.
25. Sun WH, Liu F, Chen Y, Zhu YC. Hydrogen sulfide decreases the levels of ROS by inhibiting mitochondrial complex IV andincreasing SOD activities in cardiomyocytes under ischemia/reperfusion. *Biochem Biophys Res Commun*. 2012 May 4;421(2):164-9.
26. Ge XH, Zhu GJ, Geng DQ, Zhang ZJ, Liu CF. Erythropoietin attenuates 6-hydroxydopamine-induced apoptosis via glycogen synthase kinase3β-mediated mitochondrial translocation of Bax in PC12 cells. *Neurol Sci*. 2012 Dec;33(6):1249-56.
27. Yu Y, Wang ZH, Zhang L, Yao HJ, Zhang Y, Li RJ, Ju RJ, Wang XX, Zhou J, Li N, Lu WL. Mitochondrial targeting topotecan-loaded liposomes for treating drug-resistant breast cancer and inhibitinginvasive metastases of melanoma. *Biomaterials*. 2012 Feb;33(6):1808-20.
28. Du W, Hong J, Wang YC, Zhang YJ, Wang P, Su WY, Lin YW, Lu R, Zou WP, Xiong H, Fang JY. Inhibition of JAK2/STAT3 signalling induces colorectal cancer cell apoptosis via mitochondrial pathway. *J Cell Mol Med*. 2012 Aug;16(8):1878-88.
29. Qin Y, Zhang S, Gong W, Li J, Jia J, Quan Z. Adenovirus-mediated gene transferoftissuefactor pathway inhibitor-2 inhibits gallbladder carcinoma growth invitro and in vivo. *Cancer Sci*. 2012 Apr;103(4):723-30.
30. Xiong Y, Lu QJ, Zhao J, Wu GY. Metformin inhibits growth of hepatocellular carcinoma cells by inducing apoptosis via mitochondrion-mediatedpathway. *Asian Pac J Cancer Prev*. 2012;13(7):3275-9.
31. Liu W, Fan Z, Han Y, Zhang D, Li J, Wang H. Intranuclear localization of apoptosis-inducing factor and endonuclease G involves in peroxynitrite-inducedapoptosis of spiral ganglion neurons. *Neurol Res*. 2012 Dec;34(10):915-22.
32. Tu JB, Dong Q, Hu XY, Jiang F, Ma RZ, He LY, Yang ZQ. Proteomic analysis of mitochondria from infantile hemangioma endothelial cells treated with sodium morrhuateand its liposomal formulation. *J Biochem Mol Toxicol*. 2012 Sep;26(9):374-80.
33. Gong K, Chen C, Zhan Y, Chen Y, Huang Z, Li W. Autophagy-related gene 7 (ATG7) and reactive oxygen species/extracellular signal-regulated kinase regulate tetrandrine-induced autophagy in human hepatocellular carcinoma. *J Biol Chem*. 2012 Oct 12;287(42):35576-88.
34. Hu M, Xu L, Yin L, Qi Y, Li H, Xu Y, Han X, Peng J, Wan X. Cytotoxicity of dioscin in human gastric carcinoma cells through death receptor and mitochondrialpathways. *J Appl Toxicol*. 2013 Aug; 33(8):712-22.
35. Shen XL, Zhang Y, Xu W, Liang R, Zheng J, Luo Y, Wang Y, Huang K. An iTRAQ-based mitoproteomics approach for profiling the nephrotoxicity mechanisms of ochratoxin A in HEK293 cells. *J Proteomics*. 2013 Jan 14;78:398-415.
36. Xiao H, Wang J, Yuan L, Xiao C, Wang Y, Liu X. Chicoric acid induces apoptosis in 3T3-L1 preadipocytes through ROS-mediated PI3K/Akt and MAPK signalingpathways. *J Agric Food Chem*. 2013 Feb 20;61(7):1509-20.
37. He B, Lin P, Jia Z, Du W, Qu W, Yuan L, Dai W, Zhang H, Wang X, Wang J, Zhang X, Zhang Q. The transport mechanisms of polymer nanoparticles in Caco-2 epithelial cells. *Biomaterials*. 2013 Aug;34(25):6082-98.
38. Li N, Zhang CX, Wang XX, Zhang L, Ma X, Zhou J, Ju RJ, Li XY, Zhao WY, Lu WL. Development of targeting lonidamine liposomes that circumvent drug-resistant cancer by acting on mitochondrialsignaling pathways. *Biomaterials*. 2013 Apr;34(13):3366-80.
39. Xu S, Pi H, Chen Y, Zhang N, Guo P, Lu Y, He M, Xie J, Zhong M, Zhang Y, Yu Z, Zhou Z. Cadmium induced Drp1-dependent mitochondrial fragmentation by disturbing calcium homeostasis in itshepatotoxicity. *Cell Death Dis*. 2013 Mar 14;4:e540.
40. Sun JP, Shi ZY, Liu SM, Kang YH, Hu GQ, Huangfu CS, Deng JB, Liu B. Trimethoxy-benzaldehyde levofloxacin hydrazone inducing the growth arrest and apoptosis of humanhepatocarcinoma cells. *Cancer Cell Int*. 2013 Jul 2;13(1):67.
41. Li J, He M, Shen B, Yuan D, Shao C. Alpha particle-induced bystander effect is mediated by ROS via a p53-dependent SCO2 pathway in hepatomacells. *Int J Radiat Biol*. 2013 Dec;89(12):1028-34.
42. Wu J, Yang J, Liu Q, Wu S, Ma H, Cai Y. Lanthanum induced primary neuronal apoptosis through mitochondrial dysfunction modulated by Ca²⁺ and Bcl-2 family. *Biol Trace Elem Res*. 2013 Apr;152(1):125-34.
43. Wang Y, Gao D, Chen Z, Li S, Gao C, Cao D, Liu F, Liu H, Jiang Y. Acridone derivative 8a induces oxidative stress-mediated apoptosis in CCRF-CEM leukemia cells: application of metabolomics in mechanistic studies of antitumor agents. *PLoS One*. 2013 May 7;8(5):e63572.
44. Zhou J, Zhao WY, Ma X, Ju RJ, Li XY, Li N, Sun MG, Shi JF, Zhang CX, Lu WL. The anticancer efficacy of paclitaxel liposomes modified with mitochondrial targeting conjugate in resistant lungcancer. *Biomaterials*. 2013 May;34(14):3626-38.
45. Cui Z, Lin D, Cheng F, Luo L, Kong L, Xu J, Hu J, Lan F. The role of the WWOX gene in leukemia and its mechanisms of action. *Oncol Rep*. 2013 Jun;29(6):2154-62.
46. Zhou D, Wei A, Cao C, Ruan J. DICO, a novel nonaromatic B-ring flavonoid, induces G2/M cell cycle arrest and apoptosis in human hepatomacells. *Food Chem Toxicol*. 2013 Jul;57:322-9.
47. Chang G, Zhang D, Liu J, Zhang P, Ye L, Lu K, Duan Q, Zheng A, Qin S. Exenatide protects against hypoxia/reoxygenation-induced apoptosis by improving mitochondrial function in H9c2 cells. *Exp Biol Med (Maywood)*. 2014 Apr;239(4):414-22.
48. Liu D, Yi B, Liao Z, Tang L, Yin D, Zeng S, Yao J, He M. 14-3-3γ protein attenuates lipopolysaccharide-induced cardiomyocytes injury through the Bcl-2 family/mitochondria pathway. *Int Immunopharmacol*. 2014 Aug;21(2):509-15.
49. Yan X, Tian J, Wu H, Liu Y, Ren J, Zheng S, Zhang C, Yang C, Li Y,

- Wang S. Ginsenoside rb1 protects neonatal rat cardiomyocytes from hypoxia/ischemia induced apoptosis and inhibits activation of the mitochondrial apoptotic pathway. *Evid Based Complement Alternat Med*. 2014;2014:149195.
50. Qin Y, Chu B, Gong W, Wang J, Tang Z, Shen J, Quan Z. The inhibitory effects of deleted in liver cancer 1 gene on gallbladder cancer growth through induction of cell cycle arrest and apoptosis. *J Gastroenterol Hepatol*. 2014 May;29(5):964-72.
 51. Zhang L, Zhou L, Du J, Li M, Qian C, Cheng Y, Peng Y, Xie J, Wang D. Induction of apoptosis in human multiple myeloma cell lines by ebselen via enhancing the endogenous reactive oxygen species production. *Biomed Res Int*. 2014;2014:696107.
 52. Zhou Y, Zhao W, Xie G, Huang M, Hu M, Jiang X, Zeng D, Liu J, Zhou H, Chen H, Wang GH, Zhang XK. Induction of Nur77-dependent apoptotic pathway by a coumarin derivative through activation of JNK and p38 MAPK. *Carcinogenesis*. 2014 Dec;35(12):2660-9.
 53. Zhang Y, Liu Z, Ng TB, Chen Z, Qiao W, Liu F. Purification and characterization of a novel antitumor protein with antioxidant and deoxyribonuclease activity from edible mushroom *Pholiota nameko*. *Biochimie*. 2014 Apr;99:28-37.
 54. Gao M, Gao L, Tao Y, Hou J, Yang G, Wu X, Xu H, Tompkins VS, Han Y, Wu H, Zhan F, Shi J. Proteasome inhibitor carfilzomib interacts synergistically with histone deacetylase inhibitor vorinostat in Jurkat T-leukemia cells. *Acta Biochim Biophys Sin (Shanghai)*. 2014 Jun;46(6):484-91.
 55. Wang YJ, Guo XL, Li SA, Zhao YQ, Liu ZC, Lee WH, Xiang Y, Zhang Y. Prohibitin is involved in the activated internalization and degradation of protease-activated receptor 1. *Biochim Biophys Acta*. 2014 Jul;1843(7):1393-401.
 56. Li Y, Zhang B, Huang K, He X, Luo Y, Liang R, Luo H, Shen XL, Xu W. Mitochondrial proteomic analysis reveals the molecular mechanisms underlying reproductive toxicity of zearalenone in MLTC-1 cells. *Toxicology*. 2014 Oct 3;324:55-67.
 57. He Y, Wu YT, Huang C, Meng XM, Ma TT, Wu BM, Xu FY, Zhang L, Lv XW, Li J. Inhibitory effects of long noncoding RNA MEG3 on hepatic stellate cells activation and liver fibrogenesis. *Biochim Biophys Acta*. 2014 Nov;1842(11):2204-15.
 58. Weng H, Huang H, Dong B, Zhao P, Zhou H, Qu L. Inhibition of miR-17 and miR-20a by oridonin triggers apoptosis and reverses chemoresistance by derepressing BIM-S. *Cancer Res*. 2014 Aug 15;74(16):4409-19.
 59. Cheng L, Xia TS, Wang YF, Zhou W, Liang XQ, Xue JQ, Shi L, Wang Y, Ding Q. The apoptotic effect of D Rhamnose β -Hederin, a novel oleanane-type triterpenoid saponin on breast cancer cells. *PLoS One*. 2014 Mar 6;9(6):e90848.
 60. Hao Y, Ren J, Liu C, Li H, Liu J, Yang Z, Li R, Su Y. Zinc Protects Human Kidney Cells from Depleted Uranium-induced Apoptosis. *Basic Clin Pharmacol Toxicol*. 2014 Mar;114(3):271-80.
 61. Zhou J, Xi C, Wang W, Fu X, Jinqiang L, Qiu Y, Jin J, Xu J, Huang Z. Triptolide-induced oxidative stress involved with Nrf2 contribute to cardiomyocyte apoptosis through mitochondrial dependent pathways. *Toxicol Lett*. 2014 Nov 4;230(3):454-66.
 62. Wang YH, Xu XJ, Zhang LF, Li HL. Mimic of manganese superoxide dismutase induces apoptosis in human acute myeloid leukemia cells. *Leuk Lymphoma*. 2014 May;55(5):1166-75.
 63. Yan J, Liu XH, Han MZ, Wang YM, Sun XL, Yu N, Li T, Su B, Chen ZY. Blockage of GSK3 β -mediated Drp1 phosphorylation provides neuroprotection in neuronal and mouse models of Alzheimer's disease. *Neurobiol Aging*. 2015 Jan;36(1):211-27.
 64. Zhang Q, Cheng G, Qiu H, Zhu L, Ren Z, Zhao W, Zhang T, Liu L. The p53-inducible gene 3 involved in flavonoid-induced cytotoxicity through the reactive oxygen species-mediated mitochondrial apoptotic pathway in human hepatoma cells. *Food Funct*. 2015 May;6(5):1518-25.
 65. Guo Z, Li CS, Wang CM, Xie YJ, Wang AL. CSE/H2S system protects mesenchymal stem cells from hypoxia and serum deprivation induced apoptosis via mitochondrial injury, endoplasmic reticulum stress and PI3K/Akt activation pathways. *Mol Med Rep*. 2015 Aug;12(2):2128-34.
 66. Jiang L, Li L, He X, Yi Q, He B, Cao J, Pan W, Gu Z. Overcoming drug-resistant lung cancer by paclitaxel loaded dual-functional liposomes with mitochondria targeting and pH-response. *Biomaterials*. 2015 Jun;52:126-39.
 67. Wang J, Gao Q, Wang D, Wang Z, Hu C. Metformin inhibits growth of lung adenocarcinoma cells by inducing apoptosis via the mitochondria-mediated pathway. *Oncol Lett*. 2015 Sep;10(3):1343-1349.
 68. Cai XY, Xia Y, Yang SH, Liu XZ, Shao ZW, Liu YL, Yang W, Xiong LM. Ropivacaine- and bupivacaine-induced death of rabbit annulus fibrosus cells in vitro: involvement of the mitochondrial apoptotic pathway. *Osteoarthritis Cartilage*. 2015 Oct;23(10):1763-75.
 69. Lu M, Sun L, Zhou J, Zhao Y, Deng X. Dihydroartemisinin-Induced Apoptosis is Associated with Inhibition of Sarco/Endoplasmic Reticulum Calcium ATPase Activity in Colorectal Cancer. *Cell Biochem Biophys*. 2015 Sep;73(1):137-45.
 70. Wen X, Zhou J, Zhang D, Li J, Wang Q, Feng N, Zhu H, Song Y, Li H, Bai C. Denatonium inhibits growth and induces apoptosis of airway epithelial cells through mitochondrial signaling pathways. *Respir Res*. 2015 Feb 5;16:13.
 71. Huang G, Mao J, Ji Z, Ailati A. Stachyose-induced apoptosis of Caco-2 cells via the caspase-dependent mitochondrial pathway. *Food Funct*. 2015 Mar;6(3):765-71.
 72. Sun W, Zhang X, Ding X, Li H, Geng M, Xie Z, Wu H, Huang M. Lactate dehydrogenase B is associated with the response to neoadjuvant chemotherapy in oral squamous cell carcinoma. *PLoS One*. 2015 May 14;10(5):e0125976.
 73. Cao L, Chen J, Li M, Qin YY, Sun M, Sheng R, Han F, Wang G, Qin ZH. Endogenous level of TIGAR in brain is associated with vulnerability of neurons to ischemic injury. *Neurosci Bull*. 2015 Oct;31(5):527-40.
 74. Zhang BF, Xing L, Cui PF, Wang FZ, Xie RL, Zhang JL, Zhang M, He YJ, Lyu JY, Qiao JB, Chen BA, Jiang HL. Mitochondria apoptosis pathway synergistically activated by hierarchical targeted nanoparticles co-delivering siRNA and lisdamine. *Biomaterials*. 2015 Aug;61:178-89.
 75. Li X, Liu X, Xu Y, He Y, Liu J, Xie M. Expression profile of apoptotic and proliferative proteins in hypoxic HUVEC treated with statins. *Int J Oncol*. 2015 Feb;46(2):677-84.
 76. Gong Y, Ni ZH, Zhang X, Chen XH, Zou ZM. Valproic Acid Enhances the Anti-tumor Effect of (-)-gossypol to Burkitt Lymphoma Namalwa Cells. *Biomed Environ Sci*. 2015 Oct;28(10):773-7.
 77. Zheng G, Lyu J, Liu S, Huang J, Liu C, Xiang D, Xie M, Zeng Q. Silencing of uncoupling protein 2 by small interfering RNA aggravates mitochondrial dysfunction in cardiomyocytes under septic conditions. *Int J Mol Med*. 2015 Jun;35(6):1525-36.
 78. Deng S, Tang S, Zhang S, Zhang C, Wang C, Zhou Y, Dai C, Xiao X. Furazolidone induces apoptosis through activating reactive oxygen species-dependent mitochondrial signaling pathway and suppressing PI3K/Akt signaling pathway in HepG2 cells. *Food Chem Toxicol*. 2015 Jan;75:173-86.
 79. Lai W, Zheng Z, Zhang X, Wei Y, Chu K, Brown J, Hong G, Chen L. Salidroside-Mediated Neuroprotection is Associated with Induction of Early Growth Response Genes (Egrs) Across a Wide Therapeutic Window. *Neurotox Res*. 2015 Aug;28(2):108-21.
 80. Yang J, Yao S. JNK-Bcl-2/Bcl-xL-Bax/Bak Pathway Mediates the Crosstalk between Matrine-Induced Autophagy and Apoptosis via Interplay with Beclin 1. *Int J Mol Sci*. 2015 Oct 27;16(10):25744-58.
 81. Yan J, Liu XH, Han MZ, Wang YM, Sun XL, Yu N, Li T, Su B, Chen ZY. Blockage of GSK3 β -mediated Drp1 phosphorylation provides neuroprotection in neuronal and mouse models of Alzheimer's disease. *Neurobiol Aging*. 2015 Jan;36(1):211-27.
 82. Chu B, Liu F, Li L, Ding C, Chen K, Sun Q, Shen Z, Tan Y, Tan C, Jiang Y. A benzimidazole derivative exhibiting antitumor activity blocks EGFR and HER2 activity and upregulates DR5 in breast cancer cells. *Cell Death Dis*. 2015 Mar 12;6:e1686.
 83. Yang P, Fu S, Cao Z, Liao H, Huo Z, Pan Y, Zhang G, Gao A, Zhou Q. Oroxin B selectively induces tumor-suppressive ER stress and concurrently inhibits tumor-adaptive ER stress in B-lymphoma cells for effective anti-lymphoma therapy. *Toxicol Appl Pharmacol*. 2015 Oct 15;288(2):269-79.
 84. Wu ZJ, Yu J, Fang QJ, Wang RX, Lian JB, He RL, Jiao HX, Lin MJ. Sodium Ferulate Prevents Daunorubicin-Induced Apoptosis in H9c2 Cells via Inhibition of the ERKs Pathway. *Cell Physiol Biochem*. 2015;36(6):2121-36.
 85. Wang J, Cai J, Huang Y, Ke Q, Wu B, Wang S, Han X, Wang T, Wang Y, Li W, Lao C, Song W, Xiang AP. Nestin regulates proliferation and invasion of gastrointestinal stromal tumor cells by altering mitochondrial dynamics. *Oncogene*. 2016 Jun;34:3139-50.
 86. Li M, Song M, Ren LM, Xiu CY, Liu JY, Zhu YZ, Li YF. AIC13 induces lymphocyte apoptosis in rats through the mitochondria-caspase dependent pathway. *Environ Toxicol*. 2016 Apr;31(4):385-94.
 87. Zhao FL, Qiao PF, Yan N, Gao D, Liu MJ, Yan Y. Hydrogen Sulfide

- Selectively Inhibits γ -Secretase Activity and Decreases Mitochondrial A β Production in Neurons from APP/PS1 Transgenic Mice. *Neurochem Res*. 2016 May;41(5):1145-59.
88. Song E, Tang S, Xu J, Yin B, Bao E, Hartung J. Lenti-siRNA Hsp60 promote bax in mitochondria and induces apoptosis during heat stress. *Biochem Biophys Res Commun*. 2016 Dec 2;481(1-2):125-131.
 89. Dong WQ, Chao M, Lu QH, Chai WL, Zhang W, Chen XY, Liang ES, Wang LB, Tian HL, Chen YG, Zhang MX. Prohibitin overexpression improves myocardial function in diabetic cardiomyopathy. *Oncotarget*. 2016 Jan 5;7(1):66-80.
 90. Li Q, Tang Y, Qin J, Yi L, Yang Y, Wang J, He J, Su Q, Tan H. Subcellular localization of DJ-1 in human HL-60 leukemia cells in response to diallyl disulfide treatment. *Mol Med Rep*. 2016 Nov;14(5):4666-4672.
 91. Kong F, Wang H, Guo J, Peng M, Ji H, Yang H, Liu B, Wang J, Zhang X, Li S. Hsp70 suppresses apoptosis of BRL cells by regulating the expression of Bcl-, cytochrome C, and caspase 8/3. *In Vitro Cell Dev Biol Anim*. 2016 May;52(5):568-75.
 92. Tan G, Liu D, Pan F, Zhao J, Li T, Ma Y, Shen B, Lyu J. His-87 ligand in mitoNEET is crucial for the transfer of iron sulfur clusters from mitochondria to cytosolic aconitase. *Biochem Biophys Res Commun*. 2016 Jan 29;470(1):226-32.
 93. Jin F, Qiao C, Luan N, Li H. Lentivirus-mediated PHLDA2 overexpression inhibits trophoblast proliferation, migration and invasion, and induces apoptosis. *Int J Mol Med*. 2016 Apr;37(4):949-57.
 94. Feng B, Qiu L, Ye C, Chen L, Fu Y, Sun W. Exposure to a 50-Hz magnetic field induced mitochondrial permeability transition through the ROS/GSK-3 β signaling pathway. *Int J Radiat Biol*. 2016;92(3):148-55.
 95. Li Z, Liu X, Wang L, Wang Y, Du C, Xu S, Zhang Y, Wang C, Yang C. The role of PGC-1 α and MRP1 in lead-induced mitochondrial toxicity in testicular Sertoli cells. *Toxicology*. 2016 Apr 29;355-356:39-48.
 96. Li W, Wang Q, Du M, Ma X, Wu L, Guo F, Zhao S, Huang F, Wang H, Qin G. Effects of overexpressing FoxO1 on apoptosis in glomeruli of diabetic mice and in podocytes cultured in high glucose medium. *Biochem Biophys Res Commun*. 2016 Sep 16;478(2):612-7.
 97. Wang XD, Li CY, Jiang MM, Li D, Wen P, Song X, Chen JD, Guo LX, Hu XP, Li GQ, Zhang J, Wang CH, He ZD. Induction of apoptosis in human leukemia cells through an intrinsic pathway by cathachunine, a unique alkaloid isolated from *Catharanthus roseus*. *Phytomedicine*. 2016 Jun 1;23(6):641-53.
 98. Wei C, Li H, Wang Y, Peng X, Shao H, Li H, Bai S, Xu C. Exogenous spermine inhibits hypoxia/ischemia-induced myocardial apoptosis via regulation of mitochondrial permeability transition pore and associated pathways. *Exp Biol Med (Maywood)*. 2016 Aug;241(14):1505-15.
 99. Wu L, Bao K, Song R, Wang D, Zhang L, Wang W, Zhang W, Bin W. Development of Novel Erythromycin Derivatives with Inhibitory Activity against Proliferation of Tumor Cells. *PLoS One*. 2016 Jul 22; 11(7): e0159503.
 100. Yu H, Guan Q, Guo L, Zhang H, Pang X, Cheng Y, Zhang X, Sun Y. Gypenosides alleviate myocardial ischemia-reperfusion injury via attenuation of oxidative stress and preservation of mitochondrial function in rat heart. *Cell Stress Chaperones*. 2016 May;21(3):429-37.
 101. Teng H, Huang Q, Chen L. Inhibition of cell proliferation and triggering of apoptosis by agrimonolide through MAP kinase (ERK and p38) pathways in human gastric cancer AGS cells. *Food Funct*. 2016 Nov 9; 7(11): 4605-4613.
 102. Yang Y, Wang W, Xiong Z, Kong J, Qiu Y, Shen F, Huang Z. Activation of SIRT3 attenuates triptolide-induced toxicity through closing mitochondrial permeability transition pore in cardiomyocytes. *Toxicol In Vitro*. 2016 Aug;34:128-37.
 103. Zhang GM, Deng MT, Zhang YL, Fan YX, Wan YJ, Nie HT, Wang ZY, Wang F, Lei ZH. Effect of PGC-1 α overexpression or silencing on mitochondrial apoptosis of goat luteinized granulosa cells. *J Bioenerg Biomembr*. 2016 Oct;48(5):493-507.
 104. Zhao FL, Qiao PF, Yan N, Gao D, Liu MJ, Yan Y. Hydrogen Sulfide Selectively Inhibits γ -Secretase Activity and Decreases Mitochondrial A β Production in Neurons from APP/PS1 Transgenic Mice. *Neurochem Res*. 2016 May;41(5):1145-59.
 105. Huang JY, Yuan YH, Yan JQ, Wang YN, Chu SF, Zhu CG, Guo QL, Shi JG, Chen NH. 20C, a dibenzyl compound isolated from *Gastrodia elata*, protects PC12 cells against rotenone-induced apoptosis via activation of the Nrf2/ARE/HO-1 signaling pathway. *Acta Pharmacol Sin*. 2016 Jun;37(6):731-40.
 106. Wo L, Lu D, Gu X. Knockdown of miR-182 promotes apoptosis via regulating RIP1 ubiquitination in TNF- α -treated triple-negative breast cancer cells. *Tumour Biol*. 2016 Oct;37(10):13733-13742.
 107. Gong L, Di C, Xia X, Wang J, Chen G, Shi J, Chen P, Xu H, Zhang W. AKT/mTOR signaling pathway is involved in salvianolic acid B-induced autophagy and apoptosis in hepatocellular carcinoma cells. *Int J Oncol*. 2016 Dec;49(6):2538-2548.
 108. Gao X, Xu H, Shang J, Yuan L, Zhang Y, Wang L, Zhang W, Luan X, Hu G, Chu H, Zhu T, Jia G. Ozonized carbon black induces mitochondrial dysfunction and DNA damage. *Environ Toxicol*. 2017 Mar;32(3):944-955.
 109. Zhang Y, Wu JY, Weng LH, Li XX, Yu LJ, Xu Y. Valproic acid protects against MPP $^{+}$ -mediated neurotoxicity in SH-SY5Y Cells through autophagy. *Neurosci Lett*. 2017 Jan 18;638:60-68.
 110. Yang M, Wang B, Gao J, Zhang Y, Xu W, Tao L. Spinosad induces programmed cell death involves mitochondrial dysfunction and cytochrome C release in *Spodoptera frugiperda* Sf9 cells. *Chemosphere*. 2017 Feb; 169:155-161.
 111. Li L, Sun W, Li L, Liu Y, Wu L, Wang F, Zhou Z, Zhang Z, Huang Y. A pH-responsive sequential-disassembly nanohybrid for mitochondrial targeting. *Nanoscale*. 2017 Jan 7;9(1):314-325.
 112. Sun W, Li L, Yang QQ, Zhang ZR, Huang Y. Two birds, one stone: dual targeting of the cancer cell surface and subcellular mitochondria by the galectin-3-binding peptide G3-C12. *Acta Pharmacol Sin*. 2017 Jan 9. doi: 10.1038/aps.2016.137. [Epub ahead of print]

Version 2017.06.15