碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

Annealing Buffer for DNA Oligos(5X)

产品编号	产品名称	包装
D0251	Annealing Buffer for DNA Oligos(5X)	1ml

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1. Qu B, Han X, Tang Y, Shen N. A Novel Vector-Based Method for Exclusive Overexpression of Star-Form MicroRNAs. *PLoS One*. 2012;7(7):e41504.
2. Jiang W, Liang P, Wang B, Fang J, Lang J, Tian G, Jiang J, Zhu TF. Optimized DNA extraction and metagenomic sequencing of airborne microbial communities. *Nat Protoc*. 2015 May;10(5):768-79.
3. Huang G, Gao Q, Zhao Y, Dong Z, Li T, Guan X, Jiang J. A novel siRNA validation system for functional screening of effective RNAi targets in mammalian cells and development of a derivative lentivirus delivery system. *Gene*. 2015 Mar 10;558(2):278-86.
4. Yuanyou Li, Ziyang Yin, Yewei Dong, Shuqi Wang, Óscar Monroig, Douglas R Tocher, Cuihong You. Ppar γ Is Involved in the Transcriptional Regulation of Liver LC-PUFA Biosynthesis by Targeting the $\Delta 6\Delta 5$ Fatty Acyl Desaturase Gene in the Marine Teleost *Siganus Canaliculatus* *MAR BIOTECHNOL*. 2019 Feb;21(1):19-29.
5. Tian M, Qu J, Li P, Bao Y, Liu J, Ding C, Wang S, Li T, Qi J, Yu S. Identification of novel genes essential for *Brucella abortus* to establish infection by signature-tagged mutagenesis. *Vet Microbiol*. 2019 Mar;230:130-137.
6. Fang H, Dong Y, Yue X, Hu J, Jiang S, Xu H, Wang Y, Su M, Zhang J, Zhang Z, Wang N, Chen X. The B-box zinc finger protein MdbBX20 integrates anthocyanin accumulation in response to ultraviolet radiation and low temperature. *Plant Cell Environ*. 2019 Jul;42(7):2090-2104.
7. Chen H, Wang J, Cui J, Wang C, Liang S, Liu H, Wen J. Negative regulation of bleomycins biosynthesis by ArsR/SmtB family repressor BlmR in *Streptomyces verticillus*. *APPL MICROBIOL BIOT*. 2019 Aug;103(16):6629-6644.
8. Chang Y, Deng Y, Li T, Wang J, Wang T, Tan F, Li X, Tian K. Visual detection of porcine reproductive and respiratory syndrome virus using CRISPR-Cas13a. *Transbound Emerg Dis*. 2020 Mar;67(2):564-571.
9. Yaokang Wu, Yanfeng Liu, Xueqin Lv, Jianghua Li, Guocheng Du, Long Liu. CAMERS-B: CRISPR/Cpf1 assisted multiple-genes editing and regulation system for *Bacillus subtilis* *Biotechnol Bioeng*. 2020 Jun;117(6):1817-1825.
10. Tiyu Ding, Sumathi Tomes, Andrew P Gleave, Hengtao Zhang, Andrew P Dare, Blue Plunkett, Richard V Espley, Zhiwei Luo, Ruiping Zhang, Andrew C Allan, Zhe Zhou, Huan Wang, Mengmeng Wu, Haiqing Dong, Chonghuai Liu, Jihong Liu, Zhenli Yan, Jia-Long Yao. microRNA172 targets APETALA2 to regulate flavonoid biosynthesis in apple (*Malus domestica*). *Hortic Res*. 2022 Jan 18;9:uhab007.
11. Xueyang Chen, Xingming Wang, Yuxin Yang, Chun Fang, Jing Liu, Xiongyan Liang, Yuying Yang. Enhanced pathogenicity by up-regulation of A20 after avian leukemia subgroup a virus infection. *Front Vet Sci*. 2022 Nov 14;9:1031480.
12. Liu-Qing Yang, Han-Yin Hu, Yao Han, Ze-Yi Tang, Jie Gao, Qi-Yin Zhou, Yi-Xuan Liu, Hao-Sa Chen, Tu-Nan Xu, Lei Ao, Ying Xu, Xuan Che, Ya-Bo Jiang, Chun-Wei Xu, Xian-Chao Zhang, Yu-Xin Jiang, Michal Heger, Xiao-Min Wang, Shu-Qun Cheng, Wei-Wei Pan. CpG-binding protein CFP1 promotes ovarian cancer cell proliferation by regulating BST2 transcription. *Cancer Gene Ther*. 2022 Dec;29(12):1895-1907.
13. Minggui Lin, Xingyun Liu, Xinhui Zhang, Huimin Wang, Yu Fang, Xiaoting Wu, Anqi Yin, Wanqing Yang, Dong Zhang, Miaomiao Li, Ling Zhang, Songcheng Ying. Sp1 Controls the Basal Level of Interleukin-34 Transcription. *Immunol Invest*. 2023 Feb;52(2):224-240.
14. Jiaxuan Chen, Haitao Chen, Haoming Mai, Shuang Lou, Mengqi Luo, Haisheng Xie, Bin Zhou, Jinlin Hou, De-Ke Jiang. A functional variant of CD40 modulates clearance of hepatitis B virus in hepatocytes via regulation of the ANXA2/CD40/BST2 axis. *Hum Mol Genet*. 2023 Apr 6;32(8):1334-1347.

Version 2024.03.12