碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

DNase I

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- Zhu W, Xie K, Xu Y, Wang L, Chen K, Zhang L, Fang J. CRISPR/Cas9 produces anti-hepatitis B virus effect in hepatoma cells and transgenic mouse. *Virus Res.* 2016 Jun 2;217:125-32.
- Zhu G, Chen J, Tian J, Ge L, Xing A, Tang G. Expression of NLRC4 in children with septicemia and mechanisms of NLRC4 in vitro cytokine secretion. *Mol Med Rep.* 2016 Jul;14(1):509-14.
- Sun XY, Duan ZJ, Liu Z, Tang SX, Li Y, He SC, Wang QM, Chang QY. Inhibition of P-glycoprotein, multidrug resistance-associated protein 2 and cytochrome P450 3A4 improves the oral absorption of octreotide in rats with portal hypertension. *Exp Ther Med.* 2016 Dec;12(6):3716-3722.
- Wang H, Wu J. 17 β -estradiol suppresses hyperoxia-induced apoptosis of oligodendrocytes through paired-immunoglobulin-like receptor B. *Mol Med Rep.* 2016 Mar;13(3):2892-8.
- Hui H, Rao W, Zhang L, Xie Z, Peng C, Su N, Wang K, Wang L, Luo P, Hao YL, Zhang S, Fei Z. Inhibition of Na(+)-K(+)-2Cl(-) Cotransporter-1 attenuates traumatic brain injury-induced neuronal apoptosis via regulation of Erk signaling. *Neurochem Int.* 2016 Mar;94:23-31.
- Yang S, Li SS, Yang XM, Yin DH, Wang L. Embelin prevents LMP1-induced TRAIL resistance via inhibition of XIAP in nasopharyngeal carcinoma cells. *Oncol Lett.* 2016 Jun;11(6):4167-4176.
- Feng Q, Hu ZY, Liu XQ, Zhang X, Lan X, Geng YQ, Chen XM, He JL, Wang YX, Ding YB. Stomatin-like protein 2 is involved in endometrial stromal cell proliferation and differentiation during decidualization in mice and humans. *Reprod Biomed Online.* 2016 Nov 23. pii: S1472-6483(16)30614-9.
- Qiao C, Liu J, Yang J, Li Y, Weng J, Shao Y, Zhang X. Enhanced non-inflammasome mediated immune responses by mannoseylated zwitterionic-based cationic liposomes for HIV DNA vaccines. *Biomaterials.* 2016 Apr;85:1-17.
- Liu X, Tian F, Wang S, Wang F, Xiong L. Astrocyte Autophagy Flux Protects Neurons Against Oxygen-Glucose Deprivation and Ischemic/Reperfusion Injury. *Rejuven Res.* 2017 Dec 22. doi: 10.1089/rej.2017.1999. [Epub ahead of print]
- Li J, Mao Q, He J, She H, Zhang Z, Yin C. Human umbilical cord mesenchymal stem cells improve the reserve function of perimenopausal ovary via a paracrine mechanism. *Stem Cell Res Ther.* 2017 Mar 9;8(1):55.
- Zhang Q, Cheng G, Qiu H, Wang Y, Wang J, Xu H, Zhang T, Liu L, Tao Y, Ren Z. Expression of prostate stem cell antigen is downregulated during flavonoid-induced cytotoxicity in prostate cancer cells. *Exp Ther Med.* 2017 Aug;14(2):1795-1801.
- Peng C, Rao W, Zhang L, Gao F, Hui H, Wang K, Dai S, Yang Y, Luo P, Ma Y, Ma W, Yu X, Fei Z. Mitofusin 2 Exerts a Protective Role in Ischemia Reperfusion Injury Through Increasing Autophagy. *Cell Physiol Biochem.* 2018;46(6):2311-2324.
- Pan Z, Fan Z, Ma J, Liu H, Shen L, He B, Zhang M. Profiling and functional characterization of circulation lncRNAs that are associated with coronary atherosclerotic plaque stability. *Am J Transl Res.* 2019 Jun 15;11(6):3801-3815. eCollection 2019.
- Zhong J, Xu W. Characterization of DNA hydroxymethylation in the hypothalamus of elderly mice with post-operative cognitive dysfunction. *Exp Ther Med.* 2019 Nov;18(5):4002-4010.
- Zhao RZ, Jiang S, Ru NY, Jiao B, Yu ZB. Comparison of hypoxic effects induced by chemical and physical hypoxia on cardiomyocytes. *Can J Physiol Pharmacol.* 2019 Oct;97(10):980-988.
- Zhong J, Xu W. Characterization of DNA hydroxymethylation in the hypothalamus of elderly mice with post-operative cognitive dysfunction. *Exp Ther Med.* 2019 Nov;18(5):4002-4010.

17. Jiang Zhong, Wei Xu. Characterization of DNA Hydroxymethylation in the Hypothalamus of Elderly Mice With Post-Operative Cognitive Dysfunction *Exp Ther Med.* 2019 Nov; 18(5):4002-4010.; doi: 10.3892/etm.2019.8056
18. Jin Liu, Yong-Ming Zhu, Yi Guo, Liang Lin, Zhan-Xiang Wang, Feng Gu, Xin-Yi Dong, Ming Zhou, Yi-Fan Wang, Hui-Ling Zhang. Inhibition of GSK3 β and RIP1K Attenuates Glial Scar Formation Induced by Ischemic Stroke via Reduction of Inflammatory Cytokine Production *Front Pharmacol.* 2020 Jun 12; 11:812.; doi: 10.3389/fphar.2020.00812
19. Wen Li, Zhuo Luo, Chang-Yu Yan, Xiao-Hua Wang, Zheng-Jie He, Shu-Hua Ouyang, Chang Yan, Li-Fang Liu, Qing-Qing Zhou, Han-Lu Mu, Hai-Biao Gong, Wen-Jun Duan, Lei Liang, Hiroshi Kurihara, Du Feng, Yi-Fang Li, Rong-Rong He. Autophagic degradation of PML promotes susceptibility to HSV-1 by stress-induced corticosterone *Theranostics.* 2020 Jul 11; 10(20):9032-9049.; doi: 10.7150/thno.46921
20. Jian Tan, Zhiguo Wu, Jun Liu, Wenting Zhang, Wanqiu Yuan, Hong Peng. MicroRNA-203-mediated inhibition of doublecortin underpins cardioprotection conferred by sevoflurane in rats after myocardial ischaemia-reperfusion injury *J Cell Mol Med.* 2020 Sep; 24(17):9825-9838.; doi: 10.1111/jcmm.15566.

Version 2021.09.01