BeyoRT™ II cDNA第一链合成试剂盒(RNase H-)

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1. Bao W, Gu Y, Ta L, Wang K, Xu Z. Induction of autophagy by the MG 132 proteasome inhibitor is associated with endoplasmic reticulum stress in MCF 7 cells. *Mol Med Rep*. 2016 Jan;13(1):796-804.
2. Zhang CG, Yang SD, Zhu WJ, You BG, Liu Y, Yuan ZQ, Chen WL, Li JZ, Zhou XF, Liu C, Zhang XN. Distinctive polymer micelle designed for siRNA delivery and reversal of MDR1 gene-dependent multidrugresistance. *J Biomed Mater Res B Appl Biomater*. 2016 Jul 13. doi: 10.1002/jbm.b.33748. [Epub ahead of print]
3. Li D, Chen J, Ye J, Zhai X, Song J, Jiang C, Wang J, Zhang H, Jia X, Zhu F. Anti-inflammatory effect of the six compounds isolated from *Nauclea officinalis* Pierre ex Pitard, and molecular mechanism of strictosamide via suppressing the NF-κB and MAPK signaling pathway in LPS-induced RAW 264.7 macrophages. *J Ethnopharmacol*. 2017 Jan 20;196:66-74.
4. Shaoyuan Yang, Hairun Pei, Xiaoying Zhang, Qiang Wei, Jia Zhu, Jimin Zheng, Zongchao Jia. Characterization of DicB by partially masking its potent inhibitory activity of cell division *Open Biol*. 2016 Jul;6(7):160082.;doi: 10.1098/rsob.160082
5. Li D, Chen J, Ye J, Zhai X, Song J, Jiang C, Wang J, Zhang H, Jia X, Zhu F. Anti-inflammatory effect of the six compounds isolated from *Nauclea officinalis* Pierre ex Pitard, and molecular mechanism of strictosamide via suppressing the NF-κB and MAPK signaling pathway in LPS-induced RAW 264.7 macrophages. *J Ethnopharmacol*. 2017 Jan 20;196:66-74
6. Jiang L, Tang Z. Expression and regulation of the ERK1/2 and p38 MAPK signaling pathways in periodontal tissue remodeling of orthodontic tooth movement. *Mol Med Rep*. 2018 Jan;17(1):1499-1506
7. Cao L, Zhang X, Yuan S, Cheng K, Zhang X. Autophagy induced by enterovirus 71 regulates the production of IL-6 through the p38MAPK and ERK signaling pathways. *MICROB PATHOGENESIS*. 2019 Jun;131:120-127
8. Zhuang L, Qu H, Cong J, Dai H, Liu X. MiR-181c affects estrogen-dependent endometrial carcinoma cell growth by targeting PTEN. *Endocr J*. 2019 Jun 28;66(6):523-533
9. Su J, Wang J, Luo J, Li H. Ultrasound-mediated destruction of vascular endothelial growth factor (VEGF) targeted and paclitaxel loaded microbubbles for inhibition of human breast cancer cell MCF-7 proliferation. *MOL CELL PROBE*. 2019 Aug;46:101415
10. Zhang F, Fan B, Mao L. Radiosensitizing effects of Cyclocarya paliurus polysaccharide on hypoxic A549 and H520 human non-small cell lung carcinoma cells. *Int J Mol Med*. 2019 Oct;44(4):1233-1242
11. Haiming Cai, Ji Tian Wei, Haokun Shen, Jiayi Li, Qin Fan, Zengjue Zhao, Jinbo Deng, Feiping Ming, Min Zeng, Miaopeng Ma, Peijing Zhao, Qianyi Liang, Junhao Jia, Shuxia Zhang, Linghua Zhang. Molecular cloning, characterization and expression profiles of Annexin family (ANXA1~A6) in yellow catfish (*Pelteobagrus fulvidraco*) and ANX regulation by CpG ODN responding to bacterial infection *FISH SHELLFISH IMMUN*. 2020 Apr;99:609-630.;doi: 10.1016/j.fsi.2020.02.032
12. Jiaoqing Li, Qin Fan, Haiming Cai, Jinbo Deng, Feiping Ming, Jiayi Li, Min Zeng, Miaopeng Ma, Peijing Zhao, Qianyi Liang, Junhao Jia, Shuxia Zhang, Linghua Zhang. Identification of RBP4 from bighead carp (*Hypophthalmichthys nobilis*) / silver carp (*Hypophthalmichthys molitrix*) and effects of CpG ODN on RBP4 expression under *A. hydrophila* challenge *FISH SHELLFISH IMMUN*. 2020 May;100:476-488.;doi: 10.1016/j.fsi.2020.03.036
13. Lei Liu, Hexi Zhang, Yijun Shi, Lijian Pan. Prostaglandin E1 Improves Cerebral Microcirculation Through Activation of Endothelial NOS and GRP78 *J Mol Neurosci*. 2020 Dec;70(12):2041-2048.;doi: 10.1007/s12031-020-01610-y
14. Xiangrong Cui, Hongping Liang, Chonghua Hao, Xuan Jing. Liraglutide preconditioning attenuates myocardial ischemia/ reperfusion injury via homer1 activation *AGING-US*. 2021 Feb 1;13(5):6625-6633.;doi: 10.18632/aging.202429

Version 2021.09.01