碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

BeyoRT™ II cDNA合成试剂盒(with gDNA Eraser)

产品编号	产品名称	包装
D7170S	BeyoRT™ II cDNA合成试剂盒(with gDNA Eraser)	20次
D7170M	BeyoRT™ II cDNA合成试剂盒(with gDNA Eraser)	100次
D7170L	BeyoRT™ II cDNA合成试剂盒(with gDNA Eraser)	500次

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1. Zhang ZZ,Zhang HZ,Zhang ZY.3D printed poly(ϵ -caprolactone) scaffolds function with simvastatin-loaded poly(lactic-co-glycolic acid) microspheres to repair load-bearing segmental bone defects. *Exp Ther Med*. 2019 Jan;17(1):79-90
2. Peng Y,Xing SN,Tang HY,Wang CD,Yi FP,Liu GL,Wu XM.Influence of glucose transporter 1 activity inhibition on neuroblastoma in vitro. *Gene*. 2019 Mar 20;689:11-17
3. Fenbo M,Xingyu X,Bin T.Strontium chondroitin sulfate/silk fibroin blend membrane containing microporous structure modulates macrophage responses for guided bone regeneration. *CARBOHYD POLYM*. 2019 Jun 1;213:266-275
4. Gu Y,Zhang J,Zhang X,Liang G,Xu T,Niu W.Three-dimensional Printed Mg-Doped β -TCP Bone Tissue Engineering Scaffolds: Effects of Magnesium Ion Concentration on Osteogenesis and Angiogenesis In Vitro. *Tissue Eng Regen Med*. 2019 Jun 17;16(4):415-429
5. Deng Y,Ma F,Ruiz-Ortega LI,Peng Y,Tian Y,He W,Tang B.Fabrication of strontium Eucommia ulmoides polysaccharides and in vitro evaluation of their osteoimmunomodulatory property. *Int J Biol Macromol*. 2019 Nov 1;140:727-735.
6. Zhi-Yun Liu, Chun-Guang Zhao. miR-532-3p inhibits the progression of tongue squamous cell carcinoma by targeting podoplanin. *Chin Med J (Engl)*. 2021 Dec 8;134(24):2999-3008.
7. Jingjing Lu, Bingqi Zhu, Fangmei Zhou, Xinghong Ding, Chaodong Qian, Zhishan Ding, Xiaoqing Ye. Polysaccharides From the Aerial Parts of *Tetrastigma Hemsleyanum* Diels et Gilg Induce Bidirectional Immunity and Ameliorate LPS-Induced Acute Respiratory Distress Syndrome in Mice. *Front Pharmacol*. 2022 Mar 11;13:838873.
8. Panfeng Chen, Xiaoping Li, Xi Yu, Min Yang. Ginsenoside Rg1 Suppresses Non-Small-Cell Lung Cancer via MicroRNA-126-PI3K-AKT-mTOR Pathway. *Evid Based Complement Alternat Med*. 2022 Jul 1;2022:1244836.
9. Wendan He, Xianlong Xie, Chenxi Li, Huang Ding, Jishi Ye. Adenosine A2A Receptor Antagonist Improves Cognitive Impairment by Inhibiting Neuroinflammation and Excitatory Neurotoxicity in Chronic Periodontitis Mice. *Molecules*. 2022 Sep 23;27(19):6267.
10. Meng Yu, Weiyan Wang, Jingye Dang, Binghua Liu, Junling Xu, Jingjing Li, Yang Liu, Libo He, Yuling Ying, Jiaxiu Cai, Guohua Cheng, Ke Liu. Hydrogen sulfide protects retinal pigment epithelium cells against ferroptosis through the AMPK- and p62-dependent non-canonical NRF2-KEAP1 pathway. *Exp Cell Res*. 2023 Jan 1;422(1):113436.

Version 2024.03.12