碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

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1. Yu K, Zhang T. Metagenomic and metatranscriptomic analysis of microbial community structure and gene expression of activated sludge. *PLoS One*. 2012;7(5):e38183.
2. Ren H, Zhu F, Zheng C, Sun X, Wang W, Shu H. Transcriptome analysis reveals genes related to floral development in chrysanthemum responsive to photoperiods. *Biochem Genet*. 2013 Feb;51(1-2):20-32.
3. Cai L, Yu K, Yang Y, Chen BW, Li XD, Zhang T. Metagenomic exploration reveals high levels of microbial arsenic metabolism genes in activated sludge and coastal sediments. *Appl Microbiol Biotechnol*. 2013 Nov;97(21):9579-88.
4. Su M, Chao G, Liang M, Song J, Wu K. Anticytoproliferative effect of Vitamin C on rat hepatic stellate cell. *Am J Transl Res*. 2016 Jun 15;8(6):2820-5.
5. Kambura AK, Mwirichia RK, Kasili RW, Karanja EN, Makonde HM, Boga HI. Bacteria and Archaea diversity within the hot springs of Lake Magadi and Little Magadi in Kenya. *BMC Microbiol*. 2016 Jul 7;16(1):136.
6. He Y, Ma T, Zhang X. The Mechanism of Synchronous Precise Regulation of Two Shrimp White Spot Syndrome Virus Targets by a Viral MicroRNA. *Front Immunol*. 2017 Nov 27;8:1546
7. Zhang HZ, Li YY, An T, Huang FX, Wang MQ, Liu CX, Mao JJ, Zhang LS. Comparative Transcriptome and iTRAQ Proteome Analyses Reveal the Mechanisms of Diapause in *Aphidius gifuensis* Ashmead (Hymenoptera: Aphidiidae). *Front Physiol*. 2018 Nov 30;9:1697
8. Gong Y, Ju C, Zhang X. Shrimp miR-1000 Functions in Antiviral Immunity by Simultaneously Triggering the Degradation of Two Viral mRNAs. *Front Immunol*. 2018 Dec 18;9:2999
9. Zheng JM, Kong YY, Li YY, Zhang W. MagT1 regulated the odontogenic differentiation of BMMSCs induced by TGC-CM via ERK signaling pathway. *Stem Cell Res Ther*. 2019 Jan 31;10(1):48
10. Li C, Liu Q, Kong F, Guo D, Zhai J, Su M, Sun D. Circulation and genetic diversity of Feline coronavirus type I and II from clinically healthy and FIP-suspected cats in China. *Transbound Emerg Dis*. 2019 Mar;66(2):763-775
11. Zhang S, Zhang X. Shrimp miRNA suppresses the stemness of human cancer stem cells via the PIN1 pathway. *FASEB J*. 2019 Oct;33(10):10767-10779
12. Kong Y, Hu X, Zhong Y, Xu K, Wu B, Zheng J. Magnesium-enriched microenvironment promotes odontogenic differentiation in human dental pulp stem cells by activating ERK/BMP2/Smads signaling. *Stem Cell Res Ther*. 2019 Dec 10;10(1):378

Version 2021.09.01