碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

BeyoFast™ SYBR Green qPCR Mix (2X)

产品编号	产品名称	包装
D7260-1ml	BeyoFast™ SYBR Green qPCR Mix (2X)	1ml
D7260-5ml	BeyoFast™ SYBR Green qPCR Mix (2X)	5ml
D7260-25ml	BeyoFast™ SYBR Green qPCR Mix (2X)	25ml

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- Liu Y, Long YH, Wang SQ, Zhang YY, Li YF, Mi JS, Yu CH, Li DY, Zhang JH, Zhang XJ.JMJD6 regulates histone H2A.X phosphorylation and promotes autophagy in triple-negative breast cancer cells via a novel tyrosine kinase activity.Oncogene. 2019 Feb 38(7):980-997.
- Xiu Liu, Haiyi Fei, Cuiyu Yang, Jianmin Wang, Xiaohong Zhu, Anran Yang, Zhan Shi, Xiaoying Jin, Fei Yang, Dan Wu, Lingling Jiang, Songying Zhang.Trophoblast-Derived Extracellular Vesicles Promote Preeclampsia by Regulating Macrophage PolarizationHypertension. 2022 Oct;79(10):2274-2287.
- Ma Y, Yang H, Huang J.Icariin ameliorates dexamethasone-induced bone deterioration in an experimental mouse model via activation of microRNA-186 inhibition of cathepsin K.Mol Med Rep. 2018 Jan 17(1):1633-1641.
- Zhong L, Sun S, Yao S, Han X, Gu M, Shi J..Histone deacetylase 5 promotes the proliferation and invasion of lungcancer cells.. 2018 Oct
- Wei Y, Shao J, Shen H, Wang Y, Cao G, Chen W, Zhang J, Yin L.The inhibitory role of recombinant P-selectin glycoprotein ligand immunoglobulin G on portal vein thrombosis based on a novel rat model.Exp Ther Med. 2019 May 17(5):3589-3597.
- Gu Y, Zhang J, Zhang X, Liang G, Xu T, Niu W.Three-dimensional Printed Mg-Doped β -TCP Bone Tissue Engineering Scaffolds: Effects of Magnesium Ion Concentration on Osteogenesis and Angiogenesis In Vitro.Tissue Eng Regen Med. 2019 Jun 17 16(4):415-429.
- Chen Z, Du Y, Liu X, Chen H, Weng X, Guo J, Wang M, Wang X, Wang L.EZH2 inhibition suppresses bladder cancer cell growth and metastasis via the JAK2/STAT3 signaling pathway.Oncol Lett. 2019 Jul 18(1):907-915.
- Lou K, Huang P, Ma H, Wang X, Xu H, Wang W.Orlistat increases arsenite tolerance in THP-1 derived macrophages through the up-regulation of ABCA1.Drug Chem Toxicol. 2019 Oct 31:1-9 .
- Lingchen Kong, Peng Wu, Jianzhong Li.miR-331 inhibits CLDN2 expression and may alleviate the vascular endothelial injury induced by sepsisExp Ther Med. 2020 Aug;20(2):1343-1352.
- Chengjiang Qiu, Kairui Liu, Xuguang Li, Weirun Chen, Sheng Zhang, Youxing Huang .Restoration of Intestinal Mucosa in Euphorbia kansui-treated Severe Acute Pancreatitis Rats based on HMGB1/MFG-E8 ExpressionCurr Pharm Biotechnol. 2021;22(12):1672-1682.
- Wendan He, Xianlong Xie, Chenxi Li, Huang Ding, Jishi Ye.Adenosine A2A Receptor Antagonist Improves Cognitive Impairment by Inhibiting Neuroinflammation and Excitatory Neurotoxicity in Chronic Periodontitis MiceMolecules. 2022 Sep 23;27(19):6267.
- Binghua Yang, Yuxia Fan, Meilan Chen, Lingyun Tang, Xiaodi Tang, Hongyang Li, Aiping Gu, Renlong Liang, Yi Wu.Identification and validation of a CCL18-related signature for prediction of overall survival in patients with uveal melanomaExp Eye Res. 2023 May;230:109448.
- Xinyu Yao, Xiaoyi Xu, Kunhua Hu, Zhaoshou Yang, Shaodong Deng.BANF1 promotes glutamate-induced apoptosis of HT-22 hippocampal neuronsMol Biol Rep. 2023 Nov;50(11):9441-9452.

Version 2024.03.12