碧云天生物技术/Beyotime Biotechnology
订货热线：400-168-3301或800-8283301
订货e-mail：order@beyotime.com
技术咨询：info@beyotime.com
网址：http://www.beyotime.com

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1. Mao XM, Zhou Z, Hou XP, Guan WJ, Li YQ. Reciprocal regulation between SigK and differentiation programs in *Streptomyces coelicolor*. *J Bacteriol.* 2009 Nov;191(21):6473-81.
2. Xiang D, Zhang J, Chen Y, Guo Y, Schalow A, Zhang Z, Hu X, Yu H, Zhao M, Zhu S, Lu H, Wu M, Yu Y, Moldenhauer A, Han W. Expressions and purification of a mature form of recombinant human Chemerin in *Escherichia coli*. *Protein Expr Purif.* 2010 Feb;69(2):153-8.
3. Liu WX, Hu S, Qiao ZJ, Chen WY, Liu LT, Wang FK, Hua RH, Bu ZG, Li XR. Expression, purification, and improved antigenic specificity of a truncated recombinant bp26protein of *Brucella melitensis* M5-90: a potential antigen for differential serodiagnosis of brucellosis in sheep and goats. *Biotechnol Appl Biochem.* 2011 Jan-Feb;58(1):32-8.
4. Li G, Hu Z, Yin H, Zhang Y, Huang X, Wang S, Li W. A novel dendritic nanocarrier of polyamidoamine-polyethylene glycol-cyclic RGD for "smart" small interfering RNA delivery and in vitro antitumor effects by human ether-à-go-go-related gene silencing in anaplastic thyroid carcinoma cells. *Int J Nanomedicine.* 2013;8:1293-306.
5. Wang F, Wang L, Xu Z, Liang G. Identification and analysis of multi-protein complexes in placenta. *PLoS One.* 2013 Apr 29;8(4):e62988.
6. Lin K, Gao Z, Shang B, Sui S, Fu Q. Osthole suppresses the proliferation and accelerates the apoptosis of human glioma cells via the upregulation of microRNA-16 and downregulation of MMP-9. *Mol Med Rep.* 2015 Sep;12(3):4592-7.
7. Ao N, Liu Y, Bian X, Feng H, Liu Y. Ubiquitin-specific peptidase 22 inhibits colon cancer cell invasion by suppressing the signal transducer and activator of transcription 3/matrix metalloproteinase 9 pathway. *Mol Med Rep.* 2015 Aug;12(2):2107-13.
8. Deng ZY, Hu MM, Xin YF, Gang C. Resveratrol alleviates vascular inflammatory injury by inhibiting inflammasome activation in rats with hypercholesterolemia and vitamin D2 treatment. *Inflamm Res.* 2015 May;64(5):321-32.
9. Zhang Y, Xiang X, Lu Q, Zhang L, Ma F, Wang L. Adhesions of extracellular surface-layer associated proteins in *Lactobacillus M5-L* and *Q8-L*. *J Dairy Sci.* 2016 Feb;99(2):1011-8.
10. Deng ZY, Shan WG, Wang SF, Hu MM, Chen Y. Effects of astaxanthin on blood coagulation, fibrinolysis and platelet aggregation in hyperlipidemic rats. *Pharm Biol.* 2017 Dec;55(1):663-672.

Version 2017.03.22