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1. Song Y,Zhu JS,Hua R,Du L,Huang ST,Stackman RW Jr,Zhang G,Zhang YM. mall-Conductance Ca²⁺-Activated K⁺ Channel 2 in the Dorsal Horn of Spinal Cord Participates in Visceral Hypersensitivity in Rats. *Front Pharmacol.* 2018 Aug 3;9:840.
2. Cai JY,Xu TT,Wang Y,Chang JJ,Li J,Chen XY,Chen X,Yin YF,Ni XJ. Histone deacetylase HDAC4 promotes the proliferation and invasion of glioma cells. *Int J Oncol.* 2018 Dec;53(6):2758-2768.
3. Lu Y,Liao S,Tu W,Yang B,Liu S,Pei X,Tao D,Lu Y,Ma Y,Yang Y,Liu Y. DNA demethylation facilitates the specific transcription of the mouse X-linked Tsga8 gene in round spermatids. *Biol Reprod.* 2018 Dec 12.
4. Wang Y,Fan L,Cui C,Wang Y,Liang T. EZH2 inhibition promotes methyl jasmonate-induced apoptosis of human colorectal cancer through the Wnt/ β -catenin pathway. *Oncol Lett.* 2018 Jul;16(1):1231-1236.
5. Wang D,Lou X,Jiang XM,Yang C,Liu XL,Zhang N. Quercetin protects against inflammation, MMP-2 activation and apoptosis induction in rat model of cardiopulmonary resuscitation through modulating Bmi-1 expression. *Mol Med Rep.* 2018 Jul;18(1):610-616.
6. Ci Y,Zhang Y,Liu Y,Lu S,Cao J,Li H,Zhang J,Huang Z,Zhu X,Gao J,Han M. Myricetin suppresses breast cancer metastasis through down-regulating the activity of matrix metalloproteinase (MMP)-2/9. *Phytother Res.* 2018 Jul;32(7):1373-1381.
7. Zhang Q,Xin H,Fen T. Function of microRNA-141 in human breast cancer through cytotoxic CD4⁺ T cells regulated by MAP4K4 expression. *Mol Med Rep.* 2018 Jun;17(6):7893-7901.
8. Wang X,Zou P,He Y,Meng K,Quan F,Zhang Y. Effect of luteinizing hormone on goat theca cell apoptosis and steroidogenesis through activation of the PI3K/AKT pathway. *Anim Reprod Sci.* 2018 Mar;190:108-118.
9. Yang YZ,Zhang YF,Yang L,Xu J,Mo XM,Peng W. miR-760 mediates hypoxia-induced proliferation and apoptosis of human pulmonary artery smooth muscle cells via targeting TLR4. *Int J Mol Med.* 2018 Nov;42(5):2437-2446.
10. Yan L,Guo N,Cao Y,Zeng S,Wang J,Lv F,Wang Y,Cao X. miRNA-145 inhibits myocardial infarction-induced apoptosis through autophagy via Akt3/mTOR signaling pathway in vitro and in vivo. *Int J Mol Med.* 2018 Sep;42(3):1537-1547.
11. Xin Y,Li C,Guo Y,Xiao R,Zhang H,Zhou G. RNA-Seq analysis reveals a negative role of MSMO1 with a synergized NSDHL expression during adipogenesis of 3T3-L1. *BIOSCI BIOTECH BIOCH.* 2019 Apr;83(4):641-652.
12. Zhao GS,Zhang Q,Cao Y,Wang Y,Lv YF,Zhang ZS,Zhang Y,Tan QL,Chang Y,Quan ZX,Jiang DM,Guo QN. High expression of ID1 facilitates metastasis in human osteosarcoma by regulating the sensitivity of anoikis via PI3K/AKT depended suppression of the intrinsic apoptotic signaling pathway. *Am J Transl Res.* 2019 Apr 15;11(4):2117-2139.
13. Zhang J,Li H,Yi D,Lai C,Wang H,Zou W,Cao B. Knockdown of vascular cell adhesion molecule 1 impedes transforming growth factor beta 1-mediated proliferation, migration, and invasion of endometriotic cyst stromal cells. *Reprod Biol Endocrinol.* 2019 Aug 23;17(1):69.
14. Gong L,Chang H,Xu H. LncRNA MALAT1 knockdown alleviates oxygen-glucose deprivation and reperfusion induced cardiomyocyte apoptotic death by regulating miR-122. *Exp Mol Pathol.* 2019 Dec;111:104325.
15. Sun H,Sun Y. Lidocaine inhibits proliferation and metastasis of lung cancer cell via regulation of miR-539/EGFR axis. *Artif Cells Nanomed Biotechnol.* 2019 Dec;47(1):2866-2874.
16. Wang W,Hao Y,Li F. Notoginsenoside R1 alleviates high glucose-evoked damage in RSC96 cells through down-regulation of miR-503. *Artif Cells Nanomed Biotechnol.* 2019 Dec;47(1):3947-3954.
17. Luo J,Pan J,Jin Y,Li M,Chen M. MiR-195-5p Inhibits Proliferation and Induces Apoptosis of Non-Small Cell Lung Cancer Cells by Targeting CEP55. *Onco Targets Ther.* 2019 Dec 24;12:11465-11474.
18. Gao Z,Zhao GS,Lv Y,Peng D,Tang X,Song H,Guo QN. Anoikis-resistant human osteosarcoma cells display significant angiogenesis by activating the Src kinase-mediated MAPK pathway. *Oncol Rep.* 2019 Jan;41(1):235-245.
19. Sang Y,Zhang R,Sun L,Chen KK,Li SW,Xiong L,Peng Y,Zeng L,Huang G. MORF4L1 suppresses cell proliferation, migration and invasion by increasing p21 and E-cadherin expression in nasopharyngeal carcinoma. *Oncol Lett.* 2019 Jan;17(1):294-302.
20. Wang X,Meng K,He Y,Wang H,Zhang Y,Quan F. Melatonin Stimulates STAR Expression and Progesterone Production via Activation of the PI3K/AKT Pathway in Bovine Theca Cells. *Int J Biol Sci.* 2019 Jan 1;15(2):404-415.
21. Xiao R,Li C,Wang C,Cao Y,Zhang L,Guo Y,Xin Y,Zhang H,Zhou G. Adipogenesis associated Mth938 domain containing (AAMDC) protein expression is regulated by alternative polyadenylation and microRNAs. *FEBS Lett.* 2019 Jul;593(14):1724-1734.
22. Liu X,Chu Y,Wang D,Weng Y,Jia Z. MAPK-mediated upregulation of fibrinogen-like protein 2 promotes proliferation, migration, and invasion of colorectal cancer cells. *Cell Biol Int.* 2019 Jul 9;
23. Jiang J,Li Y,Jiang Z. Effects of LDOC1 on colorectal cancer cells via downregulation of the Wnt/ β -catenin signaling pathway. *Oncol Rep.* 2019 Jun;41(6):3281-3291.
24. Li Y,Niu S,Xi D,Zhao S,Sun J,Jiang Y,Liu J. Differences in Lipopolysaccharides-Induced Inflammatory Response Between Mouse Embryonic Fibroblasts and Bone Marrow-Derived Macrophages. *J Interferon Cytokine Res.* 2019 Jun;39(6):375-382.
25. Ouyang X,Zhu Z,Yang C,Wang L,Ding G,Jiang F. Epinephrine increases malignancy of breast cancer through p38 MAPK signaling pathway in depressive disorders. *Int J Clin Exp Pathol.* 2019 Jun 1;12(6):1932-1946.
26. Zhang K,Zhang Q,Deng J,Li J,Li J,Wen L,Ma J,Li C. ALK5 signaling pathway mediates neurogenesis and functional recovery after cerebral ischemia/reperfusion in rats via Gadd45b. *Cell Death Dis.* 2019 May 1;10(5):360.
27. Shan W,Cheng C,Huang W,Ding Z,Luo S,Cui G,Lu W,Liu F,Xu J,He W,Yin Z. Angiopoietin-like 2 upregulation promotes human chondrocyte injury via NF- κ B and p38/MAPK signaling pathway. *J Bone Miner Metab.* 2019 Nov;37(6):976-986.
28. Xi J,Luo X,Wang Y,Li J,Guo L,Wu G,Li Q. Tetrahydrocurcumin protects against spinal cord injury and inhibits the oxidative stress response by regulating FOXO4 in model rats. *Exp Ther Med.* 2019 Nov;18(5):3681-7.
29. Li S,Zhang S,Chen J. c-Myc induced upregulation of long non-coding RNA SNHG16 enhances progression and carcinogenesis in oral squamous cell carcinoma. *Cancer Gene Ther.* 2019 Nov;26(11-12):400-10.
30. Li Y,Jin S,Zhao X,Luo H,Li R,Li D,Xiao T. Sequence and expression analysis of the cytoplasmic pattern recognition receptor melanoma differentiation-associated gene 5 from the barbel chub *Squaliobarbus curriculus*. *Fish Shellfish Immunol.* 2019 Nov;94:485-496.
31. Du Y,Qian B,Gao L,Tan P,Chen H,Wang A,Zheng T,Pu S,Xia X,Fu W. Aloin Preconditioning Attenuates Hepatic Ischemia/Reperfusion Injury via Inhibiting TLR4/MyD88/NF- κ B Signal Pathway In Vivo and In Vitro. *Oxid Med Cell Longev.* 2019 Nov 20;2019:3765898.
32. Yang Y,Li S,Cao J,Li Y,Hu H,Wu Z. RRM2 Regulated By LINC00667/miR-143-3p Signal Is Responsible For Non-Small Cell Lung Cancer Cell Progression. *Onco Targets Ther.* 2019 Nov 20;12:9927-9939.
33. Liu Z,Yu Z,Chang H,Wang Y,Xiang H,Zhang X,Yu B. Strontium-containing α -calcium sulfate hemihydrate promotes bone repair via the TGF- β /Smad signaling pathway. *Mol Med Rep.* 2019 Oct;20(4):3555-3564.
34. Lou K,Huang P,Ma H,Wang X,Xu H,Wang W. Orlistat increases arsenite tolerance in THP-1 derived macrophages through the up-regulation of ABCA1. *Drug Chem Toxicol.* 2019 Oct 31:1-9;
35. Zhou Z,Zhou X,Dong Y,Li M,Xu Y. Formononetin ameliorates high glucose-induced endothelial dysfunction by inhibiting the JAK/STAT signaling pathway. *Mol Med Rep.* 2019 Sep;20(3):2893-2901.
36. Ji X,Wei X,Qian J,Mo X,Kai G,An F,Lu Y. 2',4'-Dihydroxy-6'-methoxy-3',5'-dimethylchalcone induced apoptosis and G1 cell cycle arrest through PI3K/AKT pathway in BEL-7402/5-FU cells. *Food Chem Toxicol.* 2019 Sep;131:110533.
37. Yang X,Zhang C,Tie H,Luo J,Wang Y,Wu Q. miR-760 exerts an antioncogenic effect in esophageal squamous cell carcinoma by negatively driving fat metabolism via targeting c-Myc. *J Cell Biochem.* 2020 Apr;121(4):2950-2961.

Version 2019.08.11