碧云天生物技术/Beyotime Biotechnology  
订货热线: 400-168-3301或800-8283301  
订货e-mail: order@beyotime.com  
技术咨询: info@beyotime.com  
网址: http://www.beyotime.com

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1. Wang MZ, Tian WR, Tian ZJ, Qu PP, Gao SS, Zhang ZH, Zhang QY, Hou CX and Li HH. Effect of RNA Interference Hsp72 Gene Expression on Development of Mouse Preimplantation Embryos. Agricultural Sciences in China. 2009;8(2):223-30.
2. Liu J, Zou L, Wang J, Schuler C, Zhao Z, Li X, Zhang J, Liu Y. Hydrostatic pressure promotes Wnt10b and Wnt4 expression dependent and independent on ERK signaling in early-osteoinduced MSCs. Biochem Biophys Res Commun. 2009 Feb 6;379(2):505-9.
3. Zhang L, Liang Y, Liu Y, Xiong CL. The role of brain-derived neurotrophic factor in mouse oocyte maturation in vitro involves activation of protein kinase B. Theriogenology. 2010;73(8):1096-103.
4. Chu YX, Zhang Y, Zhang YQ, Zhao ZQ. Involvement of microglial P2X7 receptors and downstream signaling pathways in long-term potentiation of spinal nociceptive responses. Brain Behav Immun.

2010;24(7):1176-89.

5. X. Zhang, J. Wang, Q. Zhou, Y. Xu, S. Pu, J. Wu, Y. Xue, Y. Tian, J. Lu, W. Jiang, D Du. BDNF-activated astrocytes produce mechanical allodynia in neuropathic pain. *Neuroscience*.2011 Dec;199(29):452-60.
6. Zhang Y, Xu Y, Zhu Q, Zhao F, Luo J, Zhang X, Wang X. Upregulation of dysbindin in temporal lobe epileptic foci of human and experimental animals. *Synapse*. 2012 Jul;66(7):622-9.
7. Zheng S, Wang G, Lin S. Heat shock effects and population survival in the polar dinoflagellate *Polarella glacialis*. *Journal of Experimental Marine Biology and Ecology*. 2012 Dec;138:100-8.
8. Cheng B, Wu X, Sun S, Wu Q, Mei C, Xu Q, Wu J, He P. MAPK-PPAR $\alpha$ / $\gamma$  signal transduction pathways are involved in Chlamydia pneumoniae-induced macrophage-derived foam cell formation. *Microb Pathog*. 2014 Apr-May;69-70:1-8.
9. Zhu W, Mao Z, Zhu C, Li M, Cao C, Guan Y, Yuan J, Xie G, Guan X. Adolescent exposure to cocaine increases anxiety-like behavior and induces morphologic and neurochemical changes in the hippocampus of adult rats. *Neuroscience*. 2016 Jan 28;313:174-83.
10. Tang Y, Lv P, Sun Z, Han L, Zhou W. 14-3-3 $\beta$  Promotes Migration and Invasion of Human Hepatocellular Carcinoma Cells by Modulating Expression of MMP2 and MMP9 through PI3K/Akt/NF- $\kappa$ B Pathway. *PLoS One*. 2016 Jan 5;11(1):e0146070.
11. Lv L, Liu HG, Dong SY, Yang F, Wang QX, Guo GL, Pan YF, Zhang XH. Upregulation of CD44v6 contributes to acquired chemoresistance via the modulation of autophagy in coloncancer SW480 cells. *Tumour Biol*. 2016 Jul;37(7):8811-24.
12. Shi B, Li Y, Wang X, Yang Y, Li D, Liu X, Yang X. Silencing of hypoxia inducible factor-1 $\alpha$  by RNA interference inhibits growth of SK-NEP-1 Wilms tumour cells in vitro, and suppresses tumourigenesis and angiogenesis in vivo. *Clin Exp Pharmacol Physiol*. 2016 Jun;43(6):626-33.
13. Liu C, Huang J, Li H, Yang Z, Zeng Y, Liu J, Hao Y, Li R. Ghrelin accelerates wound healing through GHS-R1a-mediated MAPK-NF- $\kappa$ B/GR signaling pathways in combined radiation and burn injury in rats. *Sci Rep*. 2016 Jun 7;6:27499.
14. Fang K, Dong H, Jiang S, Li F, Wang D, Yang D, Gong J, Huang W, Lu F. Diosgenin and 5-Methoxypsoralen Ameliorate Insulin Resistance through ER- $\alpha$ /PI3K/Akt-Signaling Pathways in HepG2 Cells. *Evid Based Complement Alternat Med*. 2016;2016:7493694.
15. Cai M, Wang H, Li JJ, Zhang YL, Xin L, Li F, Lou SJ. The signaling mechanisms of hippocampal endoplasmic reticulum stress affecting neuronal plasticity-related protein levels in high fat diet-induced obese rats and the regulation of aerobic exercise. *Brain Behav Immun*. 2016 Oct;57:347-59.
16. Xu Y, Liu J, He M, Liu R, Belegu V, Dai P, Liu W, Wang W, Xia QJ, Shang FF, Luo CZ, Zhou X, Liu S, McDonald J, Liu J, Zuo YX, Liu F, Wang TH. Mechanisms of PDGF siRNA-mediated inhibition of bone cancer pain in the spinal cord. *Sci Rep*. 2016 Jun 10;6:27512.
17. Pang M, Teng Y, Huang J, Yuan Y, Lin F, Xiong C. Substrate stiffness promotes latent TGF- $\beta$ 1 activation in hepatocellular carcinoma. *Biochem Biophys Res Commun*. 2017 Jan 29;483(1):553-558.
18. Zhao Y, Yu Y, Zhang Y, He L, Qiu L, Zhao J, Liu M, Zhang J. Letrozole regulates actin cytoskeleton polymerization dynamics in a SRC-1 dependent manner in the hippocampus of mice. *J Steroid Biochem Mol Biol*. 2017 Mar;167:86-97.
19. Tang Y, Liu L, Xu D, Zhang W, Zhang Y, Zhou J, Huang W. Interaction between astrocytic colony stimulating factor and its receptor on microglia mediates central sensitization and behavioral hypersensitivity in chronic post ischemic pain model. *Brain Behav Immun*. 2018 Feb;68:248-260.
20. Tao L, Nie Y, Wang G, Ding Y, Ding J, Xiong F, Tang S, Wang Y, Zhou B, Zhu H. All increases endothelial nitric oxide synthase phosphorylation in rabbits with atherosclerosis. *Mol Med Rep*. 2018 Feb;17(2):2619-2625.
21. Chen C, Wang YY, Wang YX, Cheng MQ, Yin JB, Zhang X, Hong ZP. Gentiopicroside ameliorates bleomycin-induced pulmonary fibrosis in mice via inhibiting inflammatory and fibrotic process. *Biochem Biophys Res Commun*. 2018 Jan 22;495(4):2396-2403.
22. Yi L, Wu Q, Chen N, Song G, Wang C, Zou Q, Zhang Z. Valproate Plays a Protective Role against Migraine by Inhibiting Protein Kinase C Signalling in Nitroglycerin-treated Mice. *Basic Clin Pharmacol Toxicol*. 2018 Mar;122(3):310-316.
23. Wang Y, Li L, Wang Y, Zhu X, Jiang M, Song E, Song Y. New application of the commercial sweetener rebaudioside A as a hepatoprotective candidate: Induction of the Nrf2 signaling pathway. *Eur J Pharmacol*. 2018 Mar 5;822:128-137.
24. Zhu W, Ge X, Gao P, Li M, Guan Y, Guan X. Adolescent cocaine exposure induces prolonged synaptic modifications in medial prefrontal cortex of adult rats. *Brain Struct Funct*. 2018 May;223(4):1829-1838.
25. Peng H, Zhang K, Liu Z, Xu Q, You B, Li C, Cao J, Zhou H, Li X, Chen J, Cheng G, Shi R, Zhang G. VPO1 Modulates Vascular Smooth Muscle Cell Phenotypic Switch by Activating Extracellular Signal-regulated Kinase 1/2 (ERK 1/2) in Abdominal Aortic Aneurysms. *J Am Heart Assoc*. 2018 Sep 4;7(17):e010069.

-trans ret

Version 2019.07.17