碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-168-3301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

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- Lin X, Zhao J, Qian J, Mao Y, Pan J, Li L, Peng H, Luo Y, Yan J. Identification of immunodominant B- and T-cell combined epitopes in outer membrane lipoproteins LipL32 and LipL21 of *Leptospira interrogans*. *Clin Vaccine Immunol*. 2010 May;17(5):778-83.
- Lin X, Sun A, Ruan P, Zhang Z, Yan J. Characterization of conserved combined T and B cell epitopes in *Leptospira interrogans* major outer membrane proteins OmpL1 and LipL41. *BMC Microbiol*. 2011 Jan 26;11(1):21.
- Yang Y, Xia T, Zhi W, Wei L, Weng J, Zhang C, Li X. Promotion of skin regeneration in diabetic rats by electrospun core-sheath fibers loaded with basic fibroblast growth factor. *Biomaterials*. 2011 Jun;32(18):4243-54.
- Yang Y, Xia T, Chen F, Wei W, Liu C, He S, Li X. Electrospun fibers with plasmid bFGF polyplex loadings promote skin wound healing in diabetic rats. *Mol Pharm*. 2012 Jan 1;9(1):48-58.
- Sun Q, Xiong J, Lu J, Xu S, Li Y, Zhong XP, Gao GK, Liu HQ. Secretory TAT-peptide-mediated protein transduction of LIF receptor α -chain distal cytoplasmic motifs into human myeloid HL-60 cells. *Braz J Med Biol Res*. 2012 Jun 21;45(10):913-920.
- He S, Xia T, Wang H, Wei L, Luo X, Li X. Multiple release of polyplexes of plasmids VEGF and bFGF from electrospun fibrous scaffolds towards regeneration of mature blood vessels. *Acta Biomater*. 2012 Jul;8(7):2659-69.
- Wang H, Zhang Y, Xia T, Wei W, Chen F, Guo X, Li X. Synergistic Promotion of Blood Vessel Regeneration by Astragaloside IV and Ferulic Acid from Electrospun Fibrous Mats. *Mol Pharm*. 2013 Jun 3;10(6):2394-403.
- Liu Q, Liu L, Zhou J, Shin HD, Chen RR, Madzak C, Li J, Du G, Chen J. Biosynthesis of homoeriodictyol from eriodictyol by flavone 3'-O-methyltransferase from recombinant *Yarrowia lipolytica*: Heterologous expression, biochemical characterization, and optimal transformation. *J Biotechnol*. 2013 Sep 20;167(4):472-8.
- Zhao L, Liu Q, Yan S, Chen Z, Chen J, Li X. Multimeric immobilization of alcohol oxidase on electrospun fibers for valid tests of alcoholic saliva. *J Biotechnol*. 2013 Oct 10;168(1):46-54.
- Chen F, Wan H, Xia T, Guo X, Wang H, Liu Y, Li X. Promoted regeneration of mature blood vessels by electrospun fibers with loaded multiple pDNA-calcium phosphate nanoparticles. *Eur J Pharm Biopharm*. 2013 Nov;85(3 Pt A):699-710.
- Luo X, Jia G, Song H, Liu C, Wu G, Li X. Promoting antitumor activities of hydroxycamptothecin by encapsulation into acid-labile nanoparticles using electrospraying. *Pharm Res*. 2014 Jan;31(1):46-59.
- Cai C, Zhao D, Ma C, Zhang Y, Wu X, Wei G, He D. Connexin 43 expression in Sprague-Dawley rat seminiferous epithelium after in utero exposure to flutamide. *Syst Biol Reprod Med*. 2014 Oct;60(5):257-62.
- Liu Y, Lu J, Li H, Wei J, Li X. Engineering blood vessels through micropatterned coculture of vascular endothelial and smooth muscle cells on bilayered electrospun fibrous mats with pDNA inoculations. *Acta Biomater*. 2015 Jan;11:114-25.
- Peng Y, Fu ZZ, Guo CS, Zhang YX, Di Y, Jiang B, Li QW. Effects and Mechanism of Baicalin on Apoptosis of Cervical Cancer HeLa Cells In vitro. *Iran J Pharm Res*. 2015 Winter;14(1):251-61.

Version 2017.12.19