碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-168-3301或800-8283301
 订货e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

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1. Liu X, Zhang Z, Li P, Zhu L, Wang Y, Wang C Polypeptide from *Chlamys farreri* modulates UVB-induced activation of NF-kappaB signaling pathway and protection HaCaT cells from apoptosis. *Regul Pept.* 2009 Feb 25;153(1-3):49-55.

2. Liu HT, Li WM, Li XY, Xu QS, Liu QS, Bai XF, Yu C, Du YG. Chitosan oligosaccharides inhibit the expression of interleukin-6 in lipopolysaccharide-induced human umbilical vein endothelial cells through p38 and ERK1/2 protein kinases. *Basic Clin Pharmacol Toxicol.* 2010 May; 106(5):362-71.

3. Wang XH, Fu ZG, Zhao Y, Shen W, Wu XY, Wang CY Survivin: an overexpression protein with notable cellular localization and multiple roles in colon cancer *THE CHINESE-GERMAN JOURNAL OF CLINICAL ONCOLOGY.* 2010 Sep;9(9):519-23.

4. Wang X, Fu Z, Zhao Y, Wu X, Shen W. Profile of protein expression of the colon cancer cell line SW480 with survivin/shRNA. *Eur J Cancer Prev.* 2011

May;20(3):190-8.

5. Wang X, Zhao Y, Fu Z, He Y, Xiang D, Zhang L. Prelining Autogenic Endothelial Cells in Allogeneic Vessels Inhibits Thrombosis and Intimal Hyperplasia: An Efficacy Study in Dogs. *J Surg Res.* 2011 Jul;169(1):148-55.
6. Wang F, Liang YJ, Wu XP, Su XD, Fu LW. ABCG2-overexpressing S1-M1-80 cell xenografts in nude mice keep original biochemistry and cellbiological properties. *Chin J Cancer.* 2012 Mar;31(3):150-8.
7. Zhou L, Xue H, Wang Z, Ni J, Yao T, Huang Y, Yu C, Lu L. Angiotensin-(1-7) attenuates high glucose-induced proximal tubular epithelial-to-mesenchymaltransition via inhibiting ERK1/2 and p38 phosphorylation. *Life Sci.* 2012 Mar 10;90(11-12):454-62.
8. Nie J, Cao X, Zhou M, Zhang X, Zhang R, Niu L, Xia Y, Hou L, Wang C. *Musca domestica* pupae lectin induces apoptosis in HepG2 cells through a NF- κ B/p65-mediated caspasepathway. *Bull Cancer.* 2012 Apr 1;99(4):E49-54.
9. Meng X, Li F, Chen S, Tang C, Zhang W, Wang Z, Zhao S. Cloning and expression of neuron-specific enolase in the corpus luteum of dairy goats. *Gene.* 2012 Jul 25;503(2):222-8.
10. Nie C, Zhang MN, Zhao HW, Olsen TD, Jackman K, Hu LN, Ma WP, Chen XF, Wang J, Zhang Y, Gao TS, Uehara H, Ambati BK, Luo L. Correlation of in vivo and in vitro methods in measuring choroidal vascularization volumes using a subretinalinjection induced choroidal neovascularization model. *Chin Med J (Engl).* 2015 Jun 5;128(11):1516-22.

Version 2016.12.06