碧云天生物技术/Beyotime Biotechnology
 订货热线: 400-1683301 或 800-8283301
 订货 e-mail: order@beyotime.com
 技术咨询: info@beyotime.com
 网址: http://www.beyotime.com

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- Fengxiang Wang, Emily T Graham, Nathchar Naowarajna, Zhennan Shi, Yuqi Wang, Guanglei Xie, Lili Zhou, Wendy Salmon, Jie-Min Jia, Xi Wang, Yuwei Huang, Stuart L Schreiber, Yilong Zou . PALP: A rapid imaging technique for stratifying ferroptosis sensitivity in normal and tumor tissues in situ Cell Chem Biol. 2022 Jan 20;29(1):157-170.e6.
- Deng Z, Wang Z, Jin J, Wang Y, Bao N, Gao Q, Zhao J . SIRT1 protects osteoblasts against particle-induced inflammatory responses and apoptosis in aseptic prosthesis loosening. Acta Biomater. 2017 Feb;49:541-554.
- Wang F, Liang YJ, Wu XP, Su XD, Fu LW. . ABCG2-overexpressing S1-M1-80 cell xenografts in nude mice keep original biochemistry and cell biological properties. Chin J Cancer. 2012 Mar;31(3):150-8.
- Deng Z, Jin J, Wang Z, Wang Y, Gao Q, Zhao J . The metal nanoparticle-induced inflammatory response is regulated by SIRT1 through NF-κB deacetylation in aseptic loosening. INT J NANOMED. 2017 May 10;12:3617-3636.

5. Liu X, Zhang Z, Li P, Zhu L, Wang Y, Wang C. . Polypeptide from *Chlamys farreri* modulates UVB-induced activation of NF-kappaB signaling pathway and protection HaCaT cells from apoptosis. *REGUL PEPTIDES*. 2009 Feb 25;153(1-3):49-55.
6. Liu HT, Li WM, Li XY, Xu QS, Liu QS, Bai XF, Yu C, Du YG. . Chitosan oligosaccharides inhibit the expression of interleukin-6 in lipopolysaccharide-induced human umbilical vein endothelial cells through p38 and ERK1/2 protein kinases. *BASIC CLIN PHARMACOL*. 2010 May;106(5):362-71.
7. Wang XH, Fu ZG, Zhao Y, Shen W, Wu XY, Wang CY. . Survivin: an overexpression protein with notable cellular localization and multiple roles in colon cancer. 《中德临床肿瘤学杂志》. 2010 Sep;9(9):519-23.
8. Wang X, Fu Z, Zhao Y, Wu X, Shen W. . Profile of protein expression of the colon cancer cell line SW480 with survivin/shRNA. *Eur J Cancer Prev*. 2011 May;20(3):190-8.
9. Wang X, Zhao Y, Fu Z, He Y, Xiang D, Zhang L. . Prelining autogenic endothelial cells in allogeneic vessels inhibits thrombosis and intimal hyperplasia: an efficacy study in dogs. *J Surg Res*. 2011 Jul;169(1):148-55.
10. Zhou L, Xue H, Wang Z, Ni J, Yao T, Huang Y, Yu C, Lu L. . Angiotensin-(1-7) attenuates high glucose-induced proximal tubular epithelial-to-mesenchymal transition via inhibiting ERK1/2 and p38 phosphorylation. *Life Sci*. 2012 Mar 10;90(11-12):454-62.
11. Nie J, Cao X, Zhou M, Zhang X, Zhang R, Niu L, Xia Y, Hou L, Wang C. . *Musca domestica* pupae lectin induces apoptosis in HepG2 cells through a NF-kB/p65-mediated caspase pathway. *B CANCER*. 2012 Apr 1;99(4):E49-54.
12. Meng X, Li F, Chen S, Tang C, Zhang W, Wang Z, Zhao S. . Cloning and expression of neuron-specific enolase in the corpus luteum of dairy goats. *Gene*. 2012 Jul 25;503(2):222-8.
13. Liu Q, Wu C, Huang S, Wu Q, Zhou T, Liu X, Liu X, Hu X, Li L. . Decreased hyperpolarization-activated cyclic nucleotide-gated channels are involved in bladder dysfunction associated with spinal cord injury. *Int J Mol Med*. 2018 May;41(5):2609-2618
14. Liu Q, Wu C, Huang S, Wu Q, Zhou T, Liu X, Liu X, Hu X, Li L. . Decreased hyperpolarization-activated cyclic nucleotide-gated channels are involved in bladder dysfunction associated with spinal cord injury. *Int J Mol Med*. 2018 May;41(5):2609-2618.

Version 2024.05.21