碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-1683301 或 800-8283301  
 订货 e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

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5. Wang XH, Fu ZG, Zhao Y, Shen W, Wu XY, Wang CY. . Survivin: an overexpression protein with notable cellular localization and multiple roles in colon cancer. 《中德临床肿瘤学杂志》. 2010 Sep;9(9):519-23.
6. Wang X, Fu Z, Zhao Y, Wu X, Shen W. . Profile of protein expression of the colon cancer cell line SW480 with survivin/shRNA. Eur J Cancer Prev. 2011 May;20(3):190-8.
7. Wang X, Zhao Y, Fu Z, He Y, Xiang D, Zhang L. . Prelining Autogenic Endothelial Cells in Allogeneic Vessels Inhibits Thrombosis and Intimal Hyperplasia: An Efficacy Study in Dogs. J Surg Res. 2011 Jul;169(1):148-55.
8. Zhou L, Xue H, Wang Z, Ni J, Yao T, Huang Y, Yu C, Lu L. . Angiotensin-(1-7) attenuates high glucose-induced proximal tubular epithelial-to-mesenchymaltransition via inhibiting ERK1/2 and p38 phosphorylation. Life Sci. 2012 Mar 10;90(11-12):454-62.
9. Nie J, Cao X, Zhou M, Zhang X, Zhang R, Niu L, Xia Y, Hou L, Wang C. . Musca domestica pupae lectin induces apoptosis in HepG2 cells through a NF- $\kappa$ B/p65-mediated caspasepathway. B CANCER. 2012 Apr 1;99(4):E49-54.
10. Meng X, Li F, Chen S, Tang C, Zhang W, Wang Z, Zhao S. . Cloning and expression of neuron-specific enolase in the corpus luteum of dairy goats. Gene. 2012 Jul 25;503(2):222-8.
11. Zhao G, Song X, Kong X, Zhang N, Qu S, Zhu W, Yang Y, Wang Q. . Immunization with Toxoplasma gondii aspartic protease 3 increases survival time of infectedmice. Acta Trop. 2017 Jul;171:17-23.
12. Liu Q, Wu L, Huang S, Wu Q, Zhou T, Liu X, Liu X, Hu X, Li L. . Decreased hyperpolarization-activated cyclic nucleotide-gated channels are involved in bladderdysfunction associated with spinal cord injury. Int J Mol Med. 2018 May;41(5):2609-2618.
13. Liu Q, Wu C, Huang S, Wu Q, Zhou T, Liu X, Liu X, Hu X, Li L. . Decreased hyperpolarization-activated cyclic nucleotide-gated channels are involved in bladder dysfunction associated with spinal cord injury. Int J Mol Med. 2018 May
14. Jie Liu, Zequn Niu, Rui Zhang, Zhuo Peng, Liming Wang, Zhong Liu, Yanxia Gao, Honghong Pei, Longfei Pan . MALAT1 shuttled by extracellular vesicles promotes M1 polarization of macrophages to induce acute pancreatitis via miR-181a-5p/HMGB1 axis J Cell Mol Med. 2021 Oct;25(19):9241-9254.

Version 2024.05.21