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1. Guan KL, Yu W, Lin Y, Xiong Y, Zhao S. Generation of acetyllysine antibodies and affinity enrichment of acetylated peptides. *Nat Protoc.* 2010;5(9):1583-95.
2. Hu Y, Yang G, Liang B, Fang L, Ma G, Zhu Q, Chen S, Ye X. The fabrication of superlow protein absorption zwitterionic coating by electrochemically mediated atom transfer radical polymerization and its application. *Acta Biomater.* 2015 Feb;13:142-9.
3. Han S, Wang H, Sun Z, Zhao H, Zhang P. Surface modification of PS microtiter plate with chitosan oligosaccharides by ^{60}Co irradiation. *CARBOHYD POLYM.* 2017 Nov 15;176:135-139.
4. Cai Q, Zhang WJ, Zhu QQ, Chen Q. Influence of heat treatment on the structure and core IgE-binding epitopes of rAra h 2.02. *Food Chem.* 2016 Jul 1;202:404-8.
5. Shijie Dai, Haofei Hong, Kun Zhou, Kai Zhao, Yuntian Xie, Chen Li, Jie Shi, Zhifang Zhou, Lei Nie, Zhimeng Wu. Exendin 4-Hapten Conjugate Capable of Binding with Endogenous Antibodies for Peptide Half-life Extension and Exerting Long-Acting Hypoglycemic Activity. *J Med Chem.* 2021 Apr 22;64(8):4947-4959.
6. Cheng C, Deng L, Li R. The Immunogenicity and Anti-Tumor Efficacy of a Rationally Designed EGFR Vaccine. *CELL PHYSIOL BIOCHEM.* 2018;46(1):46-56.
7. Chen L, Tian Y, Zhan K, Chen A, Weng Z, Huang J, Li Y, Sun Y, Zheng H, Li Y. A humanized TCR retaining authentic specificity and affinity conferred potent anti-tumour cytotoxicity. *Immunology.* 2018 Sep;155(1):123-136.
8. Du G, Ye L, Zhang G, Dong Q, Liu K, Tian J. Human IL18-IL2 fusion protein as a potential antitumor reagent by enhancing NK cell cytotoxicity and IFN- γ production. *J CANCER RES CLIN.* 2012 Oct;138(10):1727-36.
9. Xiao G, Tang L, Yuan F, Zhu W, Zhang S, Liu Z, Geng Y, Qiu X, Zhang Y, Su L. Eicosapentaenoic acid enhances heat stress-impaired intestinal epithelial barrier function in Caco-2 cells. *PLoS One.* 2013 Sep 16;8(9):e73571.
10. Wang C, Wang Z. Studying the relationship between cell cycle and Alzheimer's disease by gold nanoparticle probes. *Anal Biochem.* 2015 Nov 15;489:32-7.
11. Rao XJ, Wu P, Shahzad T, Liu S, Chen L, Yang YF, Shi Q, Yu XQ. Characterization of a dual-CRD galectin in the silkworm *Bombyx mori*. *Dev Comp Immunol.* 2016 Jul;60:149-59.
12. Xu R, Zhou J, Zhou XD, Li Q, Perelman JM, Kolosov VP. Munc13-4 mediates human neutrophil elastase-induced airway mucin5AC hypersecretion by interacting with syntaxin2. *Mol Med Rep.* 2018 Jul;18(1):1015-1024.
13. Batool K, Alam I, Zhao G, Wang J, Xu J, Yu X, Huang E, Guan X, Zhang L. C-Type Lectin-20 Interacts with ALP1 Receptor to Reduce Cry Toxicity in *Aedes aegypti*. *TOXINS.* 2018 Sep 25;10(10). pii: E390.
14. Wang Z, Guo Z, Song T, Zhang X, He N, Liu P, Wang P, Zhang Z. Proteome-Wide Identification of On- and Off-Targets of Bcl-2 Inhibitors in Native Biological Systems by Using Affinity-Based Probes (AfBPs). *Chembiochem.* 2018 Nov 2;19(21):2312-2320.
15. Leiqiong Gao, Jing Zhou, Sen Yang, Lisha Wang, Xiangyu Chen, Yang Yang, Ren Li, Zhiwei Pan, Jing Zhao, Zhirong Li, Qizhao Huang, Jianfang Tang, Li Hu, Pinghuang Liu, Guozhong Zhang, Yaokai Chen, Lilin Ye. The dichotomous and incomplete adaptive immunity in COVID-19 patients with different disease severity. *Signal Transduct Target Ther.* 2021 Mar 8;6(1):113.
16. ZhenSong Jiang, Jianru Zhang. Mesenchymal stem cell-derived exosomes containing miR-145-5p reduce inflammation in spinal cord injury by regulating the TLR4/NF- κ B signaling pathway. *Cell Cycle.* 2021 May;20(10):993-1009.
17. Fanwu Gong, Hua-Xing Wei, Ji Qi, Huan Ma, Lianxin Liu, Jianping Weng, Xucai Zheng, Qiangsheng Li, Dan Zhao, Haopeng Fang, Liu Liu, Hongliang He, Cuichen Ma, Jinglong Han, Anyuan Sun, Baolong Wang, Tengchuan Jin, Bowei Li, Bofeng Li. Pulling-Force Spinning Top for Serum Separation Combined with Paper-Based Microfluidic Devices in COVID-19 ELISA Diagnosis. *ACS Sens.* 2021 Jul 23;6(7):2709-2719.
18. Yujie Zhao, Yuan Chen, Jiayi Wang, Xue Zhu, Ke Wang, Yue Li, Fanfan Zhou. Ginkgolide J protects human synovial cells SW982 via suppression of p38-dependent production of pro-inflammatory mediators. *Mol Med Rep.* 2021 Aug;24(2):555.

19. Ning Su, Chong-Bin Hu, Tong Shao, Chun-Yu Jin, Hao Li, Jian-Fei Ji, Lu-Lu Qin, Dong-Dong Fan, Ai-Fu Lin, Li-Xin Xiang, Jian-Zhong Shao. Functional role of CD40 and CD154 costimulatory signals in IgZ-mediated immunity against bacterial infection. *Fish Shellfish Immunol Rep.* 2021 Nov 14;2:100038.
20. Jun Zhao, Rubo Zhang, Ling Zhu, Huidan Deng, Fengqing Li, Lei Xu, Jianbo Huan, Xiangang Sun, Zhiwen Xu. Establishment of a peptide-based enzyme-linked immunosorbent assay for detecting antibodies against PRRSV M protein. *BMC Vet Res.* 2021 Nov 19;17(1):355.
21. Cheng Chen, Jiao Tian, Zhicheng He, Wenyong Xiong, Yingying He, Shubai Liu. Identified Three Interferon Induced Proteins as Novel Biomarkers of Human Ischemic Cardiomyopathy. *Int J Mol Sci.* 2021 Dec 4;22(23):13116.
22. Jie Gao, Lijun Wang, Chenyi Zhao, Yongkang Wu, Zhiyuan Lu, Yining Gu, Zongtao Ba, Xingyu Wang, Jian Wang, Ying Xu. Peony seed oil ameliorates neuroinflammation-mediated cognitive deficits by suppressing microglial activation through inhibition of NF- κ B pathway in presenilin 1/2 conditional double knockout mice. *J Leukoc Biol.* 2021 Dec;110(6):1005-1022.

Version 2022.07.03