碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-1683301 或 800-8283301  
 订货 e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

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1. Chang Y, Yao S, Chen Y, Huang J, Wu A, Zhang M, Xu F, Li F, Huang Y . Genetically-engineered protein prodrug-like nanoconjugates for tumor-targeting biomimetic delivery via a SHEATH strategy. *Nanoscale*. 2019 Jan 3 11(2):611-621.
2. Zhang Z, Xiao H, Zhang X, Zhou P . Conformation, allergenicity and human cell allergy sensitization of tropomyosin from *Exopalaemon modestus*: Effects of deglycosylation and Maillard reaction. *Food Chem*. 2019 Mar 15 276:520-527.
3. Du G, Ye L, Zhang G, Dong Q, Liu K, Tian J. . Human IL18-IL2 fusion protein as a potential antitumor reagent by enhancing NK cell cytotoxicity and IFN- $\gamma$  production. *J CANCER RES CLIN*. 2012 Oct;138(10):1727-36.
4. Li M, Zhang YX, Zhang Z, Zhou XY, Zuo XL, Cong Y, Li YQ. . Endomicroscopy Will Track Injected Mesenchymal Stem Cells in Rat Colitis Models. *Inflamm Bowel Dis*. 2015 Sep;21(9):2068-77.
5. He G, Hu W, Li CM. . Spontaneous interfacial reaction between metallic copper and PBS to form cupric phosphate nanoflower and its enzyme hybrid with enhanced activity. *COLLOID SURFACE B*. 2015 Nov 1;135:613-8.
6. He S, Shen J, Hu N, Xu X, Li J. . DKK4 enhances resistance to chemotherapeutics 5-Fu and YN968D1 in colorectal cancer cells. *Oncol Lett*. 2017 Feb;13(2):587-592.
7. Ci Y, Yang Y, Xu C, Shi L . Vesicular stomatitis virus G protein transmembrane region is crucial for the hemi-fusion to full fusion transition. *SCI REP-UK*. 2018 Jul 13;8(1):10669.
8. Zhao D, Han K, Huang X, Zhang L, Wang H, Liu N, Tian Y, Liu Q, Yang J, Liu Y, Li Y . Screening and identification of B-cell epitopes within envelope protein of tembusu virus. *Virology*. 2018 Sep 17;15(1):142.
9. Yu W, Li Y, Xie B, Ma M, Chen C, Li C, Yu X, Wang Z, Wen K, Tang BZ, Shen J . An Aggregation-Induced Emission-Based Indirect Competitive Immunoassay for Fluorescence "Turn-On" Detection of Drug Residues in Foodstuffs. *Front Chem*. 2019 Apr 24 7:228.
10. Zhu H, Zhao L, Li Z, Wen B, Qiu C, Liu M, Xu Z, Hu S, Li H . Preparation and characterization of humanized nanobodies targeting the dimer interface of epidermal growth factor receptor (EGFR). *PROTEIN EXPRES PURIF*. 2019 May 157:57-62.
11. Zhang Z, Xiao H, Zhang X, Zhou P . Insight into the effects of deglycosylation and glycation of shrimp tropomyosin on in vivo allergenicity and mast cell function. *Food Funct*. 2019 Jul 17 10(7):3934-3941.
12. Yu Q, Liu M, Xiao H, Wu S, Qin X, Ke K, Li S, Mi H, Shi D, Li P . Development of novel aptamer-based enzyme-linked apta-sorbent assay (ELASA) for rapid detection of mariculture pathogen *Vibrio alginolyticus*. *J Fish Dis*. 2019 Nov 42(11):1523-1529.
13. Zhu D, Fu P, Huang R, Xiong L, Wang Y, He L, Liao L, Li Y, Zhu Z, Wang Y . Molecular characterization, tissue distribution and functional analysis of galectin 1-like 2 in grass carp (*Ctenopharyngodon idella*). *FISH SHELLFISH IMMUN*. 2019 Nov 94:455-463.
14. Zhang Z, Xiao H, Zhou P . Glycation by saccharides of different molecular sizes affected the allergenicity of shrimp tropomyosin via epitope loss and the generation of advanced glycation end products. *Food Funct*. 2019 Nov 1 10(11):7042-7051.
15. Jing Li, Guishan Tan, Yabo Cai, Ruihuan Liu, Xiaolin Xiong, Baohua Gu, Wei He, Bing Liu, Qingyun Ren, Jianping Wu, Bo Chi, Hang Zhang, Yanzhong Zhao, Yangrui Xu, Zhenxing Zou, Fenghua Kang, Kangping Xu . A novel Apigenin derivative suppresses renal cell carcinoma via directly inhibiting wild-type and mutant MET *Biochem Pharmacol*. 2021 Aug;190:114620.
16. Mingxu Fu, Yongyan Pei, Fang Lu, Huici Jiang, Yingying Bi, Jiajing Cheng, Jinlong Qin . Identification of Potential Hub Genes and miRNA-mRNA Pairs Related to the Progression and Prognosis of Cervical Cancer Through Integrated Bioinformatics Analysis *Front Genet*. 2021 Dec 22;12:775006.
17. Pan Liu, Zhengdong Zhang, Jinwu Wang, Xiao Zhang, Xiaoping Yu, Yao Li . Empagliflozin protects diabetic pancreatic tissue from damage by inhibiting the activation of the NLRP3/caspase-1/GSDMD pathway in pancreatic  $\beta$  cells: in vitro and in vivo studies *Bioengineered*. 2021 Dec;12(2):9356-9366.
18. Myagmarsend Enkhbat, Boya Zhong, Ray Chang, Jin Geng, Long-Sheng Lu, Yin-Ju Chen, Peng-Yuan Wang . Harnessing Focal Adhesions to Accelerate p53 Accumulation and Anoikis of A549 Cells Using Colloidal Self-Assembled Patterns (cSAPs) *ACS Appl Bio Mater*. 2022 Jan 17;5(1):322-333.
19. Xinlei Sun, Fujie Shi, Weiran Wang, Yunfei Wu, Shuang Qu, Jing Li, Hongwei Liang, Ke Zen . Myeloid-Specific Pyruvate-Kinase-Type-M2-Deficient Mice Are Resistant to Acute Lung Injury *Biomedicines*. 2022 May 21;10(5):1193.

Version 2024.03.12