BeyoMag™ Anti-MBP Magnetic Beads (Anti-MBP磁珠)

产品编号	产品名称	包装
P2123-0.5ml	BeyoMag™ Anti-MBP Magnetic Beads (Anti-MBP磁珠)	0.5ml
P2123-2ml	BeyoMag™ Anti-MBP Magnetic Beads (Anti-MBP磁珠)	2ml

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Characteristics	Description
Product content	10mg/ml magnetic beads in specific protective buffer
Beads size	~200nm
Magnetization	Superparamagnetic
Coupled antibody	Anti-MBP Mouse monoclonal antibody
Isotype	IgG1
M.W. of antibody	Approximately 150kDa
Antibody concentration	≥ 0.6mg MBP antibody per ml beads
Binding capacity	≥ 0.6mg MBP-tagged fusion protein per ml beads
Specificity	Met-MBP-Protein, MBP-Protein, Protein-MBP
Elution method	SDS-PAGE loading buffer or Acid elution. Note: If elute with SDS-PAGE loading buffer, the light (~25kDa) and heavy (~50kDa) chain of MBP antibody will be denatured and release from the beads.
Application	IP, Co-IP, Protein purification

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产品编号	产品名称	包装
P2123-0.5ml	BeyoMag™ Anti-MBP Magnetic Beads	0.5ml
P2123-2ml	BeyoMag™ Anti-MBP Magnetic Beads	2ml
—	说明书	1份

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Problem	Possible Causes	Solution
Very few or no MBP - tagged protein exists in the eluate.	Protein is not completely eluted.	Change elution methods.
	No target protein expressed.	Make sure the protein of interest contains the MBP-tag by Western blot or dot blot analyses.
	Very low protein expression level.	1. Use larger volume of cell lysate. 2. Optimize expression conditions to raise the protein expression level.
	Washes are too stringent.	Reduce the time and number of washes.
	Incubation times are inadequate.	Increase the incubation time.
	Interfering substance is present in sample.	Lysates containing high concentration of DTT, 2-mercaptoethanol, or other reducing agents may destroy antibody function, and must be avoided.
	Detection system is inadequate.	If Western blot detection is used: 1. Check primary and secondary antibodies using proper controls to confirm binding and reactivity. 2. Verify that the transfer was adequate by using prestained protein marker or staining the membrane with Ponceau S. 3. Use fresh detection substrate or try a different detection system.
Background is too high.	Proteins bind nonspecifically to the Anti-MBP monoclonal antibody, insufficient washing on magnetic beads, or the microcentrifuge tubes.	1. Pre-clear lysate with Mouse IgG Magnetic Beads (P2171) to remove nonspecific binding proteins. 2. After suspending beads for the final wash, transfer entire sample to a clean microcentrifuge tube before centrifugation.
	Washes are insufficient.	1. Increase the number of washes. 2. Prolong duration of the washes, incubating each wash for at least 15 minutes. 3. Increase the salt and/or detergent concentrations in the wash solutions. 4. Centrifuge at lower speed to avoid nonspecific trapping of denatured proteins.
Multiple protein bands found in the eluate.	The protein is not stable at room temperature.	Purify the target protein at lower temperature, such as 4°C.
	Protein degradation due to proteases activity during purification process.	Add protease inhibitors to cell lysate.
	Non - specific binding.	1. Prepare cell lysate again. 2. Add additional wash steps.

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1. Duplay P, Bedouelle H, Fowler A, Zabin I, Saurin W, et al. J Biol Chem. 1984. 259(16):10606-13.
2. Reuten R, Nikodemus D, Oliveira MB, Patel TR, Brachvogel B, et al. PLoS One. 2016. 11(3):e0152386.

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产品编号	产品名称	包装
P2102-1ml	BeyoMag™ Protein A磁珠	1ml
P2102-5ml	BeyoMag™ Protein A磁珠	5ml
P2105-1ml	BeyoMag™ Protein G磁珠	1ml
P2105-5ml	BeyoMag™ Protein G磁珠	5ml
P2108-1ml	BeyoMag™ Protein A+G磁珠	1ml

P2108-5ml	BeyoMag™ Protein A+G磁珠	5ml
P2115-0.5ml	BeyoMag™ Anti-Flag Magnetic Beads (Anti-Flag磁珠)	0.5ml
P2115-2ml	BeyoMag™ Anti-Flag Magnetic Beads (Anti-Flag磁珠)	2ml
P2118-0.5ml	BeyoMag™ Anti-Myc Magnetic Beads (Anti-Myc磁珠)	0.5ml
P2118-2ml	BeyoMag™ Anti-Myc Magnetic Beads (Anti-Myc磁珠)	2ml
P2121-0.5ml	BeyoMag™ Anti-HA Magnetic Beads (Anti-HA磁珠)	0.5ml
P2121-2ml	BeyoMag™ Anti-HA Magnetic Beads (Anti-HA磁珠)	2ml
P2123-0.5ml	BeyoMag™ Anti-MBP Magnetic Beads (Anti-MBP磁珠)	0.5ml
P2123-2ml	BeyoMag™ Anti-MBP Magnetic Beads (Anti-MBP磁珠)	2ml
P2135-0.5ml	BeyoMag™ Anti-His Magnetic Beads (Anti-His磁珠)	0.5ml
P2135-2ml	BeyoMag™ Anti-His Magnetic Beads (Anti-His磁珠)	2ml
P2138-0.5ml	BeyoMag™ Anti-GST Magnetic Beads (Anti-GST磁珠)	0.5ml
P2138-2ml	BeyoMag™ Anti-GST Magnetic Beads (Anti-GST磁珠)	2ml
P2151-200µl	BeyoMag™ Streptavidin Magnetic Beads (链霉亲和磁珠)	200µl
P2151-1ml	BeyoMag™ Streptavidin Magnetic Beads (链霉亲和磁珠)	1ml
P2151-5ml	BeyoMag™ Streptavidin Magnetic Beads (链霉亲和磁珠)	5ml
P2171-1ml	BeyoMag™ Mouse IgG Magnetic Beads (小鼠IgG磁珠)	1ml
P2171-5ml	BeyoMag™ Mouse IgG Magnetic Beads (小鼠IgG磁珠)	5ml
P2173-1ml	BeyoMag™ Rabbit IgG Magnetic Beads (兔IgG磁珠)	1ml
P2173-5ml	BeyoMag™ Rabbit IgG Magnetic Beads (兔IgG磁珠)	5ml

Version 2023.08.08