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Characteristics	Description
Product content	50% settled gel in specific protective buffer
Agarose structure	6% cross-linked agarose
Average particle size	45-165 μ m
Coupled protein	Streptavidin
M.W. of protein	~55kDa (Streptavidin)
Protein concentration	\geq 2mg streptavidin per ml settled gel
Binding capacity	Per ml settled gel: \geq 1mg biotinylated antibody or dsDNA, \geq 200nmol free biotin, \geq 8nmol biotinylated oligonucleotides or peptides
Specificity	Biotinylated ligands
Elution method	Elution with acid or SDS-PAGE loading buffer for single-use applications
Application	Purification of biotin-labeled proteins and nucleic acids, IP, Co-IP, DNA-protein pull-downs

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Buffer	Components	Application
TBS	Tris Buffered Saline (e.g. ST661/ST665)	Prewash
Binding & Washing Buffer I (2X)	10mM Tris-HCl (pH7.5), 1mM EDTA, 2M NaCl, 0.01%-0.1% Tween-20	for nucleic acid
Washing Buffer II (1X)	PBS (pH7.4), 0.05% Tween-20, with or without 0.01%-0.1% BSA	for antibody/protein
DNA/RNA Elution Buffer	95% formamide, 10mM EDTA, pH8.2	for nucleic acid
Acid Elution Buffer	0.1M Glycine (pH2.8) or 8M Guanidine·HCl (pH1.5)	for antibody/protein
Neutralization solution	1M Tris-HCl (pH8.5-9.5)	for antibody/protein
SDS-PAGE Loading Buffer	SDS-PAGE Loading Buffer (e.g. P0015)	for antibody/protein

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1. Lim KH, Huang H, Pralle A, Park S. Biotechnol Bioeng. 2013. 110(1):57-67.

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产品编号	产品名称	包装
P2015-2ml	Protein A Agarose (Fast Flow, 抗体纯化用)	2ml
P2015-10ml	Protein A Agarose (Fast Flow, 抗体纯化用)	10ml
P2015-50ml	Protein A Agarose (Fast Flow, 抗体纯化用)	50ml
P2015-200ml	Protein A Agarose (Fast Flow, 抗体纯化用)	200ml
P2017-2ml	Protein G Agarose (Fast Flow, 抗体纯化用)	2ml
P2017-10ml	Protein G Agarose (Fast Flow, 抗体纯化用)	10ml
P2017-50ml	Protein G Agarose (Fast Flow, 抗体纯化用)	50ml
P2017-200ml	Protein G Agarose (Fast Flow, 抗体纯化用)	200ml
P2019-2ml	Protein A+G Agarose (Fast Flow, 抗体纯化用)	2ml
P2019-10ml	Protein A+G Agarose (Fast Flow, 抗体纯化用)	10ml
P2019-50ml	Protein A+G Agarose (Fast Flow, 抗体纯化用)	50ml
P2019-200ml	Protein A+G Agarose (Fast Flow, 抗体纯化用)	200ml
P2159-1ml	Streptavidin Agarose(链霉亲和素琼脂糖凝胶)	1ml
P2159-5ml	Streptavidin Agarose(链霉亲和素琼脂糖凝胶)	5ml
P2159-20ml	Streptavidin Agarose(链霉亲和素琼脂糖凝胶)	20ml
P2151-200μl	BeyoMag™ Streptavidin Magnetic Beads (链霉亲和素磁珠)	200μl
P2151-1ml	BeyoMag™ Streptavidin Magnetic Beads (链霉亲和素磁珠)	1ml
P2151-5ml	BeyoMag™ Streptavidin Magnetic Beads (链霉亲和素磁珠)	5ml

Version 2023.05.22