碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-1683301或800-8283301  
 订货e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

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++++, Strong Binding  
 ++~+++, Medium Binding  
 +, Weak Binding  
 +/-, Weak or No Binding  
 -, No Binding

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Characteristics	Description
Product content	10mg/ml magnetic beads in specific protective buffer
Beads size	~200nm
Magnetization	Superparamagnetic
Coupled protein	Recombinant Protein A and Protein G
M.W. of protein	~25kDa (Protein A/Protein G)
Antibody concentration	≥0.6mg Protein A+G per ml beads
Binding capacity	≥ 0.7mg human IgG per ml beads
Specificity	Antibodies from many different species, including mouse, human, rabbit, cow, goat and sheep
Elution method	Elution with acid, competing peptide or SDS-PAGE loading buffer
Application	IP, Co-IP, Protein purification

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产品编号	产品名称	包装
P2179S-1	Lysis Buffer	50ml
P2179S-2	TBS (10X)	10ml
P2179S-3	Protease Inhibitor Cocktail (100X)	0.5ml
P2179S-4	BeyoMag™ Protein A+G Magnetic Beads	0.4ml
P2179S-5	Normal Mouse IgG (1mg/ml)	10μl
P2179S-6	Normal Rabbit IgG (1mg/ml)	10μl
P2179S-7	Acid Elution Buffer	2ml
P2179S-8	Neutralization Buffer	0.2ml
P2179S-9	SDS-PAGE Sample Loading Buffer (5X)	0.4ml
—	说明书	1份

产品编号	产品名称	包装
P2179M-1	Lysis Buffer	250ml
P2179M-2	TBS (10X)	30ml
P2179M-3	Protease Inhibitor Cocktail (100X)	2.5ml
P2179M-4	BeyoMag™ Protein A+G Magnetic Beads	2ml
P2179M-5	Normal Mouse IgG (1mg/ml)	50μl
P2179M-6	Normal Rabbit IgG (1mg/ml)	50μl
P2179M-7	Acid Elution Buffer	10ml
P2179M-8	Neutralization Buffer	1ml
P2179M-9	SDS-PAGE Sample Loading Buffer (5X)	2ml
—	说明书	1份

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Steps	Solution required	Volume per assay	Volume per assay
Cell lysis and sample preparation	Lysis Buffer with Protease Inhibitor Cocktail	100μl	500μl
Preparation of magnetic beads	TBS	~0.5ml	~1.5ml
Immunoprecipitation	Magnetic Beads	4μl	20μl
Wash for beads-Ab complex (3 times)	TBS	100μl each time	500μl each time
Wash for beads-Ab-Ag complex (3 times)	Lysis Buffer with Protease Inhibitor Cocktail	100μl each time	500μl each time
Acid elution and neutralization (optional)	Acid Elution Buffer	20μl	100μl
	Neutralization Buffer	2μl	10μl
SDS-PAGE sample loading buffer elution (optional)	SDS-PAGE Sample Loading Buffer (1X)	20μl	100μl

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Problem	Possible Causes	Solution
Very few or no target protein exists in the eluate.	Protein is not completely eluted.	Change elution methods.
	No target protein expressed.	Make sure the protein of interest contains the target protein by Western blot or dot blot analyses.
	Very low protein expression level.	1. Use larger volume of cell lysate. 2. Optimize expression conditions to raise the protein expression level.
	Washes are too stringent.	Reduce the time and number of washes.
	Incubation times are inadequate.	Increase the incubation time.
	Interfering substance is present in sample.	Lysates containing high concentration of DTT, 2-mercaptoethanol, or other reducing agents may destroy antibody function, and must

		be avoided.
	Detection system is inadequate.	If Western blot detection is used: 1. Check primary and secondary antibodies using proper controls to confirm binding and reactivity. 2. Verify that the transfer was adequate by using prestained protein marker or staining the membrane with Ponceau S. 3. Use fresh detection substrate or try a different detection system.
Background is too high.	Proteins bind nonspecifically to the antibody, insufficient washing on magnetic beads, or the microcentrifuge tubes.	1. Pre-clear lysate with Normal IgG to remove nonspecific binding proteins. 2. After suspending beads for the final wash, transfer entire sample to a clean microcentrifuge tube before separation.
	Washes are insufficient.	1. Increase the number of washes. 2. Prolong duration of the washes, incubating each wash for at least 15 minutes. 3. Increase the salt and/or detergent concentrations in the wash solutions. 4. Centrifuge at lower speed to avoid nonspecific trapping of denatured proteins.
Multiple protein bands found in the eluate.	The protein is not stable at room temperature.	Purify the target protein at lower temperature, such as 4°C.
	Protein degradation due to proteases activity during purification process.	Add protease inhibitors to cell lysate.
	Non-specific binding.	1. Prepare cell lysate again. 2. Add additional wash steps.

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