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Characteristics	Description
Product content	10mg/ml magnetic beads in specific protective buffer
Beads size	~200nm
Magnetization	Superparamagnetic
Coupled antibody	Anti-V5 mouse monoclonal antibody
Isotype	IgG2b
M.W. of antibody	Approximately 150kDa
Antibody concentration	≥0.6mg V5 antibody per ml beads
Binding capacity	≥ 0.6mg V5-tagged fusion protein per ml beads
Specificity	V5-Protein, Protein-V5
Elution method	Acid, peptide competitive or SDS-PAGE loading buffer elution. Note: If elute with SDS-PAGE loading buffer, the light (~25kDa) and heavy (~50kDa) chain of V5 antibody will be denatured and release from the beads.
Application	IP, Co-IP, Protein purification

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产品编号	产品名称	包装
P2187M-1	Lysis Buffer	250ml
P2187M-2	TBS (10X)	15ml
P2187M-3	Protease Inhibitor Cocktail (100X)	2.5ml
P2187M-4	BeyoMag™ Anti-V5 Magnetic Beads	2ml
P2187M-5	BeyoMag™ Mouse IgG Magnetic Beads	0.4ml
P2187M-6	V5 Peptide (25X)	0.4ml
P2187M-7	Acid Elution Buffer	10ml
P2187M-8	Neutralization Buffer	1ml
P2187M-9	SDS-PAGE Sample Loading Buffer (5X)	2ml
—	说明书	1份

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Steps	Solution required	Volume per assay	Volume per assay
Cell lysis and sample preparation	Lysis Buffer with Protease Inhibitor Cocktail	100μl	500μl
Preparation of V5 Peptide elution buffer and magnetic beads	TBS	~0.52ml	~1.6ml
Immunoprecipitation	Magnetic Beads	4μl	20μl
Wash (3 times)	Lysis Buffer with Protease Inhibitor Cocktail	100μl each time	500μl each time
Peptide competitive elution (optional)	Peptide solution (1X)	20μl	100μl
Acid elution and neutralization (optional)	Acid Elution Buffer	20μl	100μl
	Neutralization Buffer	2μl	10μl

SDS-PAGE sample loading buffer elution (optional)	SDS-PAGE Sample Loading Buffer (1X)	20 μ l	100 μ l
---	-------------------------------------	------------	-------------

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Problem	Possible Causes	Solution
Very few or no tagged protein exists in the eluate.	Protein is not completely eluted.	Change elution methods.
	No target protein expressed.	Make sure the protein of interest contains the tagged protein by Western blot or dot blot analyses.
	Very low protein expression level.	1. Use larger volume of cell lysate. 2. Optimize expression conditions to raise the protein expression level.
	Washes are too stringent.	Reduce the time and number of washes.
	Incubation times are inadequate.	Increase the incubation time.
	Interfering substance is present in sample.	Lysates containing high concentration of DTT, 2-mercaptoethanol, or other reducing agents may destroy antibody function, and must be avoided.
	Detection system is inadequate.	If Western blot detection is used: 1. Check primary and secondary antibodies using proper controls to confirm binding and reactivity. 2. Verify that the transfer was adequate by using prestained protein marker or staining the membrane with Ponceau S. 3. Use fresh detection substrate or try a different detection system.
Background is too high.	1. Pre-clear lysate with Mouse IgG Magnetic Beads (P2171) to remove nonspecific binding proteins. 2. After suspending beads for the final wash, transfer entire sample to a clean microcentrifuge tube before centrifugation.	

	tubes.	
	Washes are insufficient.	<ol style="list-style-type: none"> 1. Increase the number of washes. 2. Prolong duration of the washes, incubating each wash for at least 15 minutes. 3. Choose other wash buffers. Increase the salt and/or detergent concentrations in the wash solutions. 4. Centrifuge at lower speed to avoid nonspecific trapping of denatured proteins.
Multiple protein bands found in the eluate.	The protein is not stable at room temperature.	Purify the target protein at lower temperature, such as 4°C.
	Protein degradation due to proteases activity during purification process.	Add protease inhibitors to cell lysate.
	Non-specific binding.	<ol style="list-style-type: none"> 1. Prepare cell lysate again. 2. Add additional wash steps.

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1. Lee C. Methods Mol Biol. 2007. 362:401-406.
2. Gevaert K, Vandekerckhove J. Electrophoresis. 2000. 21:1145-1154.
3. Kolodziej K E, Pourfarzad F, de Boer E, Krpic S, Grosveld F, Strouboulis J. BMC Mol Biol. 2009. 10:6.

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