碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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- of RAV5 in activating nitrate reductases and repressing catalase activity to improve disease resistance in cassava. *Plant Biotechnol J*. 2021 Apr;19(4):785-800.
6. Danrui Cui, Xiufang Xiong, Jianfeng Shu, Xiaoqing Dai, Yi Sun, Yongchao Zhao. FBXW7 Confers Radiation Survival by Targeting p53 for Degradation. *Cell Rep*. 2020 Jan 14;30(2):497-509.e4.;doi: 10.1016/j.celrep.2019.12.032.
 7. Xiaoli Cui, Ruijin Zhang, Ye Yang, Erzhong Wu, Yiheng Tang, Zhihua Zhao, Chao Li, Lei Yang, Xueyi Teng, Yanzhen Ye, Ya Cui, Feng Xu, Zewen Su, Dongpeng Wang, Dongdong Zhang, Yan Yang, Jianyuan Sun, Jianjun Luo, Shuli Zhang, Runsheng Chen, Jianzhong Jeff Xi. Identification and characterization of long non-coding RNA Carip in modulating spatial learning and memory. *Cell Rep*. 2022 Feb 22;38(8):110398.
 8. Qidong Li, Hanhui Yao, Yichun Wang, Yang Wu, Rick F Thorne, Youming Zhu, Mian Wu, Lianxin Liu. circPRKAA1 activates a Ku80/Ku70/SREBP-1 axis driving de novo fatty acid synthesis in cancer cells. *Cell Rep*. 2022 Nov 22;41(8):111707
 9. Yao-Shun Yang, Xi-Zheng Jia, Qian-Yun Lu, Sun-Li Cai, Xue-Ting Huang, Shu-Hua Yang, Chris Wood, Yue-Hong Wang, Jiao-Jiao Zhou, Yi-Ding Chen, Jin-Shu Yang, Wei-Jun Yang. Exosomal DEK removes chemoradiotherapy resistance by triggering quiescence exit of breast cancer stem cells. *Oncogene*. 2022 Apr;41(18):2624-2637.
 10. Ling C, Qiao GQ, Shuai BW, Olavarria K, Yin J, Xiang RJ, Song KN, Shen YH, Guo Y, Chen GQ. Engineering NADH/NAD⁺ ratio in *Halomonas bluephagenesis* for enhanced production of polyhydroxyalkanoates (PHA). *Metab Eng*. 2018 Sep;49:275-286.
 11. Xunxun Wu, Lian Liu, Qiuling Zheng, Haiping Hao, Hui Ye, Ping Li, Hua Yang. Protocatechuic aldehyde protects cardiomyocytes against ischemic injury via regulation of nuclear pyruvate kinase M2. *Acta Pharm Sin B*. 2021 Nov;11(11):3553-3566.
 12. Liu C, Jiang W, Tian X, Yang P, Xiao L, Li J, Qiu L, Tu H, Tan W. Identification of Vigilin as a Potential Ischemia Biomarker by Brain Slice-Based Systematic Evolution of Ligands by Exponential Enrichment. *Anal Chem*. 2019 May 21;91(10):6675-6681.
 13. Hairong Shi, Hongshun Wang, Cheng Zhang, Yajie Lu, Jun Yao, Zhibin Chen, Guangqian Xing, Qinjun Wei, Xin Cao. Mutations in OSBPL2 cause hearing loss associated with primary cilia defects via sonic hedgehog signaling. *JCI Insight*. 2022 Feb 22;7(4):e149626.
 14. Yu Yan, Peng Wang, Yi Lu, Yujing Bai, Yunxie Wei, Guoyin Liu, Haitao Shi. MeRAV5 promotes drought stress resistance in cassava by modulating hydrogen peroxide and lignin accumulation. *Plant J*. 2021 Aug;107(3):847-860.
 15. Xingchen Huo, Zhensheng Wang, Xun Xiao, Chunrong Yang, Jianguo Su. Oral Administration of Nanopeptide CMCS-20H Conspicuously Boosts Immunity and Precautionary Effect Against Bacterial Infection in Fish. *Front Immunol*. 2022 Jan 11;12:811616.
 16. Pu Z, Che Y, Zhang W, Sun H, Meng T, Xie H, Cao L, Hao H. Dual roles of IL-18 in colitis through regulation of the function and quantity of goblet cells. *Int J Mol Med*. 2019 Jun;43(6):2291-2302.
 17. Liu G, Zeng H, Li X, Wei Y, Shi H. Functional Analysis of MaWRKY24 in Transcriptional Activation of Autophagy-Related Gene 8f/g and Plant Disease Susceptibility to Soil-Borne *Fusarium oxysporum* f. sp. *cubense*. *Pathogens*. 2019 Nov 25;8(4). pii: E264.
 18. Chen Q, Zhou W, Qian C, Shen K, Zhu X, Zhou D, Sun Z, Lu W, Liu H, Li K, Xu T, Bao Q, Lu J. OXA-830, a Novel Chromosomally Encoded Extended-Spectrum Class D β -Lactamase in *Aeromonas simiae*. *Front Microbiol*. 2019 Nov 26;10:2732.
 19. Yanli Chang, Yujing Bai, Yunxie Wei, Haitao Shi. CAMTA3 negatively regulates disease resistance through modulating immune response and extensive transcriptional reprogramming in cassava. *Tree Physiol*. 2020 Oct 29;40(11):1520-1533.
 20. Guobei Yan, Na Liu, Junhua Wang, Jun Tian, Hongbin Liu, Suping Li, Wencao Liu, Xiaoliang Li, Kai Li, Hong Wang. Deubiquitylation and stabilization of Acf7 by ubiquitin carboxylterminal hydrolase 14 (USP14) is critical for NSCLC migration. *J Biosci*. 2021;46:19.
 21. Ning Zhang, Lingran Zhang, Linjie Li, Junyou Geng, Lei Zhao, Yan Ren, Zhongdong Dong, Feng Chen. Global Profiling of 2-hydroxyisobutyrilome in Common Wheat. *Genomics Proteomics Bioinformatics*. 2021 Feb 10;S1672-0229(21)00007-3.
 22. Guobei Yan, Na Liu, Jun Tian, Yanli Fu, Wei Wei, Jingjing Zou, Suping Li, Qing Wang, Kai Li, Junhua Wang. Deubiquitylation and stabilization of ARMC5 by ubiquitin-specific processing protease 7 (USP7) are critical for RCC proliferation. *J Cell Mol Med*. 2021 Mar;25(6):3149-3159.
 23. Rongcan Luo, Yu Fan, Jing Yang, Maosen Ye, Deng-Feng Zhang, Kun Guo, Xiao Li, Rui Bi, Min Xu, Lu-Xiu Yang, Yu Li, Xiaoqian Ran, Hong-Yan Jiang, Chen Zhang, Liwen Tan, Nengyin Sheng, Yong-Gang Yao. A novel missense variant in ACAA1 contributes to early-onset Alzheimer's disease, impairs lysosomal function, and facilitates amyloid- β pathology and cognitive decline. *Signal Transduct Target Ther*. 2021 Aug 31;6(1):325.
 24. Qian Ma, Jia Wang, Haopeng Yu, Junyi Sun, Jun Sun, Jiahao Chen, Yangyinhuai Yu, Peihang Fang, Qi Tian, Junjun Ding. Protocol to alter a protein's phase separation capacity to control cell fate transitions. *STAR Protoc*. 2021 Oct 14;2(4):100887.
 25. Yan Hu, Ya-Ping Liang, Hang-Yu Tian, Cui-Xian Xu, Dan Yu, Pan Zhang, Hui Ye, Ming Li. *Microplitis bicoloratus* bracovirus regulates cyclophilin A-apoptosis-inducing factor interaction to induce cell apoptosis in the insect immunosuppressive process. *Arch Insect Biochem Physiol*. 2022 May;110(1):e21877.
 26. Ming-Li Li, Sheng Wang, Penghui Xu, Hang-Yu Tian, Mixue Bai, Ya-Ping Zhang, Yong Shao, Zi-Jun Xiong, Xiao-Guang Qi, David N Cooper, Guojie Zhang, He Helen Zhu, Dong-Dong Wu. Functional genomics analysis reveals the evolutionary adaptation and demographic history of pygmy lorises. *Proc Natl Acad Sci U S A*. 2022 Oct 4;119(40):e2123030119.
 27. Yang Wu, Yichun Wang, Hanhui Yao, Heng Li, Fanzheng Meng, Qidong Li, Xiansheng Lin, Lianxin Liu. MNX1-AS1, a c-Myc induced lncRNA, promotes the Warburg effect by regulating PKM2 nuclear translocation. *J Exp Clin Cancer Res*. 2022 Dec 7;41(1):337.

Version 2023.12.06