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Characteristics	Description
Product content	25% slurry in specific protective buffer
Beads size	30~150 $\mu\text{m}$
Magnetization	Superparamagnetic
Coupled Ligand	GSH
Ligand density	30-50 $\mu\text{mol/ml}$ beads
Binding capacity	10-20mg GST-tagged fusion protein per ml beads (precipitation)
Specificity	Met-GST-tag-Protein, GST-tag-Protein, Protein-GST-tag
Application	Protein purification, IP, Co-IP

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Buffer	Components
Binding/Wash Buffer	1X PBS, pH7.4
Elution Buffer	50mM Tris-HCl, 10mM GSH, pH8.0

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Problem	Possible Causes	Solution
Very few or no GST-tagged protein exists in the eluate.	Protein is not completely eluted.	1. Increase the GSH concentration to 15mM or higher in the elution buffer. 2. Increase the elution or incubation time. 3. Add Triton X-100 (0.1%) or n-Octylglucoside (2%) or NaCl (0.1-0.2M) to the Elution Buffer.
	No target protein expressed.	Make sure the protein of interest contains the GST-tag by Western blot or dot blot analyses.
	Very low protein expression level.	1. Use larger volume of cell lysate. 2. Optimize expression conditions to raise the protein expression level.
	Washes are too stringent.	Reduce the time and number of washes.
	Incubation times are inadequate.	Increase the incubation time.
	Detection system is inadequate.	If Western blot detection is used: 1. Check primary and secondary antibodies using proper controls to confirm binding and reactivity. 2. Verify that the transfer was adequate by using pre-stained protein marker or staining the membrane with Ponceau S. 3. Use fresh detection substrate or try a different detection system.
The purity of the target protein is low.	The target protein has degraded.	Add appropriate protease inhibitors such as PMSF to the cell lysate and Binding/Wash Buffer.
	Host proteins, such as chaperonins, may interact with the fusion protein	1. Add DTT (5mM) in the Binding/Wash Buffer. 2. Add Chaperonin Buffer (2mM ATP, 10mM MgSO <sub>4</sub> , 50mM Tris-HCl) to the cell lysate and incubate at 37°C for 10 minutes prior to the purification.
	Washes are insufficient.	1. Increase the number of washes. 2. Prolong duration of the washes, incubating each wash for at least 15 minutes. 3. Add more stringent wash conditions. Detergent such as 1% Triton X-100, 1% Tween-20, 0.03% SDS, or 0.1% NP-40 may be used. 4. Centrifuge at lower speed to avoid nonspecific trapping of denatured proteins.
The binding capacity of the magnetic beads has declined.	Protein or other substances have aggregated on the beads.	Wash the beads with NaOH.
	The beads have been reused too many times.	Use new beads.

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- Einarson MB, Pugacheva EN, Orlinick JR. CSH Protoc. 2007. 2007:pdb.prot4757.
- Schäfer F, Seip N, Maertens B, Block H, Kubicek J. Methods Enzymol. 2015. 559:127-39.

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产品编号	产品名称	包装
P2138-0.5ml	BeyoMag™ Anti-GST Magnetic Beads (Anti-GST磁珠)	0.5ml
P2138-2ml	BeyoMag™ Anti-GST Magnetic Beads (Anti-GST磁珠)	2ml
P2250	BeyoGold™ GST-tag Purification Resin	1ml
P2251	BeyoGold™ GST-tag Purification Resin	10ml
P2253	BeyoGold™ GST-tag Purification Resin	100ml

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