碧云天生物技术/Beyotime Biotechnology
订货热线: 400-168-3301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

Beyozol (总RNA抽提试剂)

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1. Li H, Zhang YY, Tan HW, Jia YF, Li D. Therapeutic effect of tripteryne on adjuvant arthritis in rats. *J Ethnopharmacol.* 2008 Aug 13;118(3):479-84.
2. ZHANG ZY, LIU XJ, LIU T, YAN L, WANG YJ, WANG CB. Polypeptide from *Chlamys farreri* inhibits UVB-induced apoptosis of HaCaT cells via iNOS/NO and HSP90. *Chinese Journal of Oceanology and Limnology.* 2009;27(3):594-9.
3. Longxi P, Buwu F, Yuan W, Sinan G. Expression of p53 in the effects of artesunate on induction of apoptosis and inhibition of proliferation in rat primary hepatic stellate cells. *PLoS One.* 2011;6(10):e26500.
4. Zhao YQ, Feng HW, Jia T, Chen XM, Zhang H, Xu AT, Zhang HL, Fan XL. Antiproliferative effects of celecoxib in Hep-2 cells through telomerase inhibition and induction of apoptosis. *Asian Pac J Cancer Prev.* 2014;15(12):4919-23.
5. Wang X, Jiang Q, Wang W, Su L, Han Y, Wang C. Molecular mechanism of polypeptides from *Chlamys farreri* (PCF) s anti-apoptotic effect in UVA-exposed HaCaT cells involves HSF1/HSP70, JNK, XO, iNOS and NO/ROS. *J Photochem Photobiol B.* 2014 Jan 5;130:47-56.
6. Jin Q, Jiang S, Wu YL, Bai T, Yang Y, Jin X, Lian LH, Nan JX. Hepatoprotective effect of cryptotanshinone from *Salvia miltiorrhiza* in D-galactosamine/lipopolysaccharide-induced fulminant hepatic failure. *Phytomedicine.* 2014 Jan 15;21(2):141-7.
7. Ye J, Han Y, Chen X, Xie J, Liu X, Qiao S, Wang C. L-Carnitine attenuates H2O2-induced neuron apoptosis via inhibition of endoplasmic reticulum stress. *Neurochem Int.* 2014 Sep 9;78C:86-95.
8. Zhao S, Wang Z, Chen J, Chen J. Preparation of heparan sulfate-like polysaccharide and application in stem cell chondrogenic differentiation. *Carbohydr Res.* 2015 Jan 12;401:32-8.
9. Zhao S, Deng C, Wang Z, Teng L, Chen J. Heparan sulfate 6-O-sulfotransferase 3 is involved in bone marrow mesenchymal stromal cell osteogenic differentiation. *Biochemistry (Mosc).* 2015 Mar;80(3):379-89.
10. Yu B, Dai CQ, Chen J, Deng L, Wu XL, Wu S, Zhao CL, Jiang ZY, Chen XY. Dysbiosis of gut microbiota induced the disorder of helper T cells in influenza virus-infected mice. *Hum Vaccin Immunother.* 2015;11(5):1140-6.

Version 2016.12.07