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- Zhou F, Zhang L, Gong K, Lu G, Sheng B, Wang A, Zhao N, Zhang X, Gong Y. LEF-1 activates the transcription of E2F1. *Biochem Biophys Res Commun.* 2008 Jan 4;365(1):149-53.
- Rong JJ, Hu R, Qi Q, Gu HY, Zhao Q, Wang J, Mu R, You QD, Guo QL. Gambogic acid down-regulates MDM2 oncogene and induces p21(Waf1/CIP1) expression independent of p53. *Cancer Lett.* 2009 Oct 18; 284(1):102-12.
- Yuan D, Pan Y, Zhang J, Shao C. Role of nuclear factor-kappaB and P53 in radioadaptive response in Chang live cells. *Mutat Res.* 2010 Jun 1; 688(1-2):66-71.
- Zhang J, Chen Y, Xin XL, Li QN, Li M, Lin LP, Geng MY, Ding J. Oligomannuric acid sulfate blocks tumor growth by inhibiting NF-kappaB activation. *Acta Pharmacol Sin.* 2010 Mar;31(3):375-81.
- Zhou ZL, Luo ZG, Yu B, Jiang Y, Chen Y, Feng JM, Dai M, Tong LJ, Li Z, Li YC, Ding J, Miao ZH. Increased accumulation of hypoxia-inducible factor-1 α with reduced transcriptional activity mediates the antitumor effect of triptolide. *Mol Cancer.* 2010 Oct 11; 9:268.
- Jiang ZJ, Jin TT, Gao F, Liu JW, Zhong JJ, Zhao H. Effects of Ganoderic acid Me on inhibiting multidrug resistance and inducing apoptosis in multidrug resistant colon cancer cells. *Process Biochemistry.* 2011 Jun;46(6):1307-14.
- Xu Y, Jiang Z, Yin P, Li Q, Liu J. Role for Class I histone deacetylases in multidrug resistance. *Exp Cell Res.* 2012 Feb 1;318(3):177-86.
- Feng N, Xu B, Tao J, Li P, Cheng G, Min Z, Mi Y, Wang M, Tong N, Tang J, Zhang Z, Wu H, Zhang W, Wang Z, Hua L. A miR-125b binding site polymorphism in bone morphogenetic protein membrane receptor type IBgene and prostate cancer risk in China. *Mol Biol Rep.* 2012 Jan;39(1):369-73.
- Zhang J, Xin X, Chen Q, Xie Z, Gui M, Chen Y, Lin L, Feng J, Li Q, Ding J, Geng M. Oligomannuric acid sulfate sensitizes cancer cells to doxorubicin by inhibiting atypical activation of NF-kB via targeting of Mre11. *Int J Cancer.* 2012 Jan 15;130(2):467-77.
- Guo L, Wu WJ, Liu LD, Wang LC, Zhang Y, Wu LQ, Guan Y, Li QH. Herpes simplex virus 1 ICP22 inhibits the transcription of viral gene promoters by binding to and blocking therecruitment of P-TEFb. *PLoS One.* 2012;7(9):e45749.
- Zhou H, Chen S, Wang W, Wang Z, Wu X, Zhang Z. Nanog inhibits lipopolysaccharide-induced expression of pro-inflammatory cytokines by blocking NF-kBtranscriptional activity in rat primary microglial cells. *Mol Med Rep.* 2012 Mar;5(3):842-6.
- Zhou Z, Wang M, Zhao J, Wang L, Gao Y, Zhang H, Liu R, Song L. The increased transcriptional response and translocation of a Rel/NF-kB homologue in scallop *Chlamys farreri*during the immune stimulation. *Fish Shellfish Immunol.* 2013 May;34(5):1209-15.
- Deng X, Rui W, Zhang F, Ding W. PM2.5 induces Nrf2-mediated defense mechanisms against oxidative stress by activating PI3K/AKT signaling pathway in human lung alveolar epithelial A549 cells. *Cell Biol Toxicol.* 2013 Jun;29(3):143-57.
- Song S, Zhou F, Chen WR, Xing D. PDT-induced HSP70 externalization up-regulates NO production via TLR2 signal pathway in macrophages. *FEBS Lett.* 2013 Jan 16;587(2):128-35.
- Zeng W, Chang H, Ma M, Li Y. CCL20/CCR6 promotes the invasion and migration of thyroid cancer cells via NF-kappa B signaling-induced MMP-3 production. *Exp Mol Pathol.* 2014 Aug;97(1):184-90.
- Gao F, Wang W. MicroRNA-96 promotes the proliferation of colorectal cancer cells and targets tumor protein p53 inducible nuclear protein 1, forkhead box protein O1 (FOXO1) and FOXO3a. *Mol Med Rep.* 2015 Feb;11(2):1200-6.
- Gong K, Qu B, Liao D, Liu D, Wang C, Zhou J, Pan X. MiR-132 regulates osteogenic differentiation via downregulating Sirtuin1 in a peroxisome proliferator-activated receptor β / δ -dependent manner. *Biochem Biophys Res Commun.* 2016 Sep 9;478(1):260-7.
- Pan Y, Zhang J, Fu H, Shen L. miR-144 functions as a tumor suppressor in breast cancer through inhibiting ZEB1/2-mediated epithelialmesenchymal transition process.*Onco Targets Ther.* 2016 Oct 11;9:6247-6255.
- Yuan G, Hua B, Yang Y, Xu L, Cai T, Sun N, Yan Z, Lu C, Qian R. The Circadian Gene Clock Regulates Bone Formation Via PDIA3.*J Bone Miner Res.* 2016 Nov 24. doi: 10.1002/jbmr.3046. [Epub ahead of print]
- Xu X, Fan S, Zhou J, Zhang Y, Che Y, Cai H, Wang L, Guo L, Liu L, Li Q. The mutated tegument protein UL7 attenuates the virulence of herpes simplex virus 1 by reducing the modulation of α -4 gene transcription. *Virology.* 2016 Sep 13;13:152.
- Li T, Kong AN, Ma Z, Liu H, Liu P, Xiao Y, Jiang X, Wang L. Protein arginine methyltransferase 1 may be involved in pregnane x receptor-activated overexpression ofmultidrug resistance 1 gene during acquired multidrug resistant. *Oncotarget.* 2016 Apr 12;7(15):20236-48.
- Qu H, Zou Z, Pan Z, Zhang T, Deng N, Chen G, Wang Z. IL-7/IL-7 receptor axis stimulates prostate cancer cell invasion and migration via AKT/NF-kB pathway. *Int Immunopharmacol.* 2016 Nov;40:203-210.
- Tong S, Xia T, Fan K, Jiang K, Zhai W, Li JS, Wang SH, Wang JJ.14-3-3 ζ promotes lung cancer cell invasion by increasing the Snail protein expression through atypicalprotein kinase C (aPKC)/NF-kB signaling. *Exp Cell Res.* 2016 Oct 15;348(1):1-9.
- Qu B, Ma Y, Yan M, Gong K, Liang F, Deng S, Jiang K, Ma Z, Pan X. Sirtuin1 promotes osteogenic differentiation through downregulation of peroxisome proliferator-activatedreceptor γ in MC3T3-E1 cells. *Biochem Biophys Res Commun.* 2016 Sep 9;478(1):439-45.
- Liu J, Li Y, Luo M, Yuan Z, Liu J. MicroRNA-214 inhibits the osteogenic differentiation of human osteoblasts through the direct regulation ofbaculoviral IAP repeat-containing 7. *Exp Cell Res.* 2017 Jan 18. pii: S0014-4827(17)30016-2.
- Liu Z, Liao L, Cao Z, Chen X, Du J. Synergistic effect of APRIL knockdown and Jiedu Xiaozheng Yin, a Chinese medicinal recipe, on theinhibition of hepatocellular carcinoma cell proliferation. *Oncol Rep.* 2017 Feb;37(2):754-760.

Version 2017.03.09