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1. Yanan Gao, Zihao Liang, Bingyong Mao, Xudong Zheng, Jinjun Shan, Cuiyuan Jin, Shijia Liu, Narasaiah Kolliputi, Yugen Chen, Feng Xu, Liyun Shi . Gut microbial GABAergic signaling improves stress-associated innate immunity to respiratory viral infection J Adv Res. 2023 Jun 22;S2090-1232(23)00171-6.
2. Changwei Chen, Xin Chai, Xueping Hu, Shengying Lou, Dan Li, Tingjun Hou, Sunliang Cui . Discovery of 2-(1-(3-Chloro-4-cyanophenyl)-1 H-pyrazol-3-yl)acetamides as Potent, Selective, and Orally Available Antagonists Targeting the Androgen Receptor J Med Chem. 2022 Oct 13;65(19):13074-13093.
3. Yajie Hu, Yang Liu, Yuying Yang, Hongming Lv, Shuai Lian, Bin Xu, Shize Li . OGT upregulates myogenic IL-6 by mediating O-GlcNAcylation of p65 in mouse skeletal muscle under cold exposure J Cell Physiol. 2022 Feb;237(2):1341-1352.
4. Lingling Tang, Li Liu, Xianhong Sun, Po Hu, Hui Zhang, Bohan Wang, Xiaona Zhang, Jinjin Jiang, Xia Zhao, Xiaolu Shi . BMAL1/FOXA2-induced rhythmic fluctuations in IL-6 contribute to nocturnal asthma attacks Front Immunol. 2022 Nov 25;13:947067.
5. Bin Zhang, Zhi-Yi Liu, Rui Wu, Cheng-Ming Zhang, Kuan Cao, Wen-Gang Shan, Zhen Liu, Ming Ji, Zi-Lu Tian, Gautam Sethi, Heng-Liang Shi, Ren-Hao Wang . Transcriptional regulator CTR9 promotes hepatocellular carcinoma progression and metastasis via increasing PEG10 transcriptional activity Acta Pharmacol Sin. 2022 Aug;43(8):2109-2118.
6. Tiezhu Guo, Yueting Liu, Xinliang Ren, Wei Wang, Hanrui Liu . Promoting Role of Long Non-Coding RNA Small Nucleolar RNA Host Gene 15 (SNHG15) in Neuronal Injury Following Ischemic Stroke via the MicroRNA-18a/CXC Chemokine Ligand 13 (CXCL13)/ERK/MEK Axis MED SCI MONITOR. 2020 Aug 30;26:e923610.
7. Ligang Zhang, Yabin Dang, Ying Wang, Xin Fan . Nucleolar and spindle-associated protein 1 accelerates cellular proliferation and invasion in nasopharyngeal carcinoma by potentiating Wnt/β-catenin signaling via modulation of GSK-3β J Bioenerg Biomembr. 2020 Dec;52(6):441-451.
8. Xiaobin Zhou, Yao Lu, Pengfei Guo, Chenglin Zhou . Upregulation of microRNA-140-3p mediates dachshund family transcription factor 1 expression in immunoglobulin A nephropathy through cell cycle-dependent mechanisms Mol Med Rep. 2021 Feb;23(2):134.
9. Ruimin Zhao, Shiyang Wang, Junsong Liu, Chongwen Xu, Shaoqiang Zhang, Yuan Shao, Xiaoyi Duan . KLK11 acts as a tumor-inhibitor in laryngeal squamous cell carcinoma through the inactivation of Akt/Wnt/β-catenin signaling J Bioenerg Biomembr. 2021 Feb;53(1):85-96.

10. Jie Feng, Liyi Xie, Xiaoyang Yu, Chao Liu, Hongjuan Dong, Wanhong Lu, Ranran Kong . Perilipin 5 ameliorates high-glucose-induced podocyte injury via Akt/GSK-3 β /Nrf2-mediated suppression of apoptosis, oxidative stress, and inflammation *BIOCHEM BIOPH RES CO*. 2021 Mar 12;544:22-30.
11. Chengliang Li, Hong Xin, Yingpeng Shi, Jianjun Mu . Glutaredoxin 2 protects cardiomyocytes from hypoxia/reoxygenation-induced injury by suppressing apoptosis, oxidative stress, and inflammation via enhancing Nrf2 signaling *Int Immunopharmacol*. 2021 May;94:107428.
12. Chenrui Wu, Mengran Du, Renqiang Yu, Yuqi Cheng, Biying Wu, Jiayuan Yuan Fu, Weilin Tan, Qiang Zhou, Ehab Balawi, Z B Liao . A novel mechanism linking ferroptosis and endoplasmic reticulum stress via the circPtpn14/miR-351-5p/5-LOX signaling in melatonin-mediated treatment of traumatic brain injury *Free Radic Biol Med*. 2022 Jan;178:271-294.
13. Xianting Sun, Jianchen Fang, Fen Ye, Shuxian Zhang, Honghui Huang, Jian Hou, Ting Wang . Diffuse Large B-Cell Lymphoma Promotes Endothelial-to-Mesenchymal Transition via WNT10A/Beta-Catenin/Snail Signaling *Front Oncol*. 2022 Apr 12;12:871788.
14. Yu-Qi Cheng, Chen-Rui Wu, Meng-Ran Du, Qiang Zhou, Bi-Ying Wu, Jia-Yuan-Yuan Fu, Ehab Balawi, Wei-Lin Tan, Zheng-Bu Liao . CircLphn3 protects the blood-brain barrier in traumatic brain injury *Neural Regen Res*. 2022 Apr;17(4):812-818.
15. Ruosi Yao, Menghui Zhang, Jian Zhou, Linlin Liu, Yan Zhang, Jian Gao, Kailin Xu . Novel dual-targeting c-Myc inhibitor D347-2761 represses myeloma growth via blocking c-Myc/Max heterodimerization and disturbing its stability *Cell Commun Signal*. 2022 May 26;20(1):73.
16. Weixu Ma, Mingzhu Cao, Shilei Bi, Lili Du, Jingsi Chen, Haibin Wang, Yufei Jiang, Yixuan Wu, Yixin Liao, Shuangbo Kong, Jianqiao Liu . MAX deficiency impairs human endometrial decidualization through down-regulating OSR2 in women with recurrent spontaneous abortion *Cell Tissue Res*. 2022 May;388(2):453-469.
17. Jingwei Yan, Weijie Chen, Hao Zeng, Hao Cheng, Jinwei Suo, Chenliang Yu, Baoru Yang, Heqiang Lou, Lili Song, Jiasheng Wu . Unraveling the malate biosynthesis during development of *Torreya grandis* nuts *Curr Res Food Sci*. 2022 Nov 19;5:2309-2315.
18. Mengran Du, Chenrui Wu, Renqiang Yu, Yuqi Cheng, Zhaohua Tang, Biying Wu, Jiayuan Yuan Fu, Weilin Tan, Qiang Zhou, Ziyu Zhu, Ehab Balawi, Xuekang Huang, Jun Ma, Z B Liao . A novel circular RNA, circIgf2bp2, links neural plasticity and anxiety through targeting mitochondrial dysfunction and oxidative stress-induced synapse dysfunction after traumatic brain injury *Mol Psychiatry*. 2022 Nov;27(11):4575-4589.
19. Yang Liu, Hongzhao Shi, Yajie Hu, Ruizhi Yao, Peng Liu, Yuying Yang, Shize Li . RNA binding motif protein 3 (RBM3) promotes protein kinase B (AKT) activation to enhance glucose metabolism and reduce apoptosis in skeletal muscle of mice under acute cold exposure *Cell Stress Chaperones*. 2022 Nov;27(6):603-618.

Version 2024.03.12