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1. Li MX, Wang YR, Cao JM, Tian GZ, Zhao LH. To Explore the Double - labeling Method of Monitoring the GHRP Regulatory Function on [Ca²⁺]_i and NO on Real Time in Cardiomyocytes Under LSCM. *J Med Res.* 2007 Aug;36(8):22-5.
2. Cao X, You QD, Li ZY, Guo QL, Shang J, Yan M, Chern JW, Chen ML. Design and synthesis of 7-alkoxy-4-heteroaryl-amino-3-quinolone carbonitriles as dual inhibitors of c-Src kinase and nitric oxide synthase. *Bioorg Med Chem.* 2008 Jun 1;16(11):5890-8.
3. Cao X, You QD, Li ZY, Liu XR, Xu D, Guo QL, Shang J, Chern JW, Chen ML. The design, synthesis and biological evaluation of 7-alkoxy-4-heteroaryl-amino-3-cyanoquinolines as dual inhibitors of c-Src and iNOS. *Bioorg Med Chem Lett.* 2008 Dec 1;18(23):6206-9.
4. Zhao F, Wang L, Liu K. In vitro anti-inflammatory effects of arctigenin, a lignan from *Arctium lappa* L., through inhibition on iNOS pathway. *J Ethnopharmacol.* 2009 Apr 21;122(3):457-62.
5. Xiong J, Lu H, Lu K, Duan Y, An L, Zhu C. Cadmium decreases crown root number by decreasing endogenous nitric oxide, which is indispensable for crown root primordia initiation in rice seedlings. *Planta.* 2009 Sep;230(4):599-610.
6. Xue L, Xu F, Meng L, Wei S, Wang J, Hao P, Bian Y, Zhang Y, Chen Y. Acetylation-dependent regulation of mitochondrial ALDH2 activation by SIRT3 mediates acute ethanol-induced eNOS activation. *FEBS Lett.* 2012 Jan 20;586(2):137-42.
7. Zhuang XY, Zhang W, Pang XF, Wu WB, Qu WJ. Combined effect of total flavonoids from seed residues of *Hippophae rhamnoides* L. and zinc on advanced glycation end products-induced endothelial cell dysfunction. *Food Chemistry.* 2012 Aug;133(3):906-11.
8. Song S, Zhou F, Chen WR, Xing D. PDT-induced HSP70 externalization up-regulates NO production via TLR2 signal pathway in macrophages. *FEBS Lett.* 2013 Jan 16;587(2):128-35.
9. Kang H, Liu M, Fan Y, Deng X. A potential gravity-sensing role of vascular smooth muscle cell glycocalyx in altered gravitational stimulation. *Astrobiology.* 2013 Jul;13(7):626-36.
10. Wang H, Xiao W, Niu Y, Jin C, Chai R, Tang C, Zhang Y. Nitric oxide enhances development of lateral roots in tomato (*Solanum lycopersicum* L.) under elevated carbon dioxide. *Planta.* 2013 Jan;237(1):137-44.
11. Chen J, Zhang HQ, Hu LB, Shi ZQ. Microcystin-LR-induced phytotoxicity in rice crown root is associated with the cross-talk between auxin and nitric oxide. *Chemosphere.* 2013 Sep;93(2):283-93.
12. Yao J, Pan D, Zhao Y, Zhao L, Sun J, Wang Y, You QD, Xi T, Guo QL, Lu N. Wogonin prevents lipopolysaccharide-induced acute lung injury and inflammation in mice via peroxisome proliferator-activated receptor gamma-mediated attenuation of the nuclear factor-kappaB pathway. *Immunology.* 2014 Oct;143(2):241-57.
13. Zhou Y, Chen S, Ding R, Yao W, Gao X. Inflammatory modulation effect of glycopeptide from *Ganoderma capense* (Lloyd) Teng. *Mediators Inflamm.* 2014;2014:691285.
14. Tong L, Chen T, Chen Z, Zhang P, Pi H, Ruan H, Wu J. Anti-inflammatory activity of omphalocarpin isolated from *Radix Toddaliae Asiaticae*. *J Ethnopharmacol.* 2014 Sep 29;155(3):1553-60.
15. Zhu W, Yang B, Fu H, Ma L, Liu T, Chai R, Zheng Z, Zhang Q, Li G. Flavone inhibits nitric oxide synthase (NOS) activity, nitric oxide production and protein S-nitrosylation in breast cancer cells. *Biochem Biophys Res Commun.* 2015 Mar 13;458(3):590-5.
16. Jiang ZS, Pu ZC, Hao ZH. Carvacrol protects against spinal cord injury in rats via suppressing oxidative stress and the endothelial nitric oxide synthase pathway. *Mol Med Rep.* 2015 Oct;12(4):5349-54.
17. Zhan R, Han Q, Zhang C, Tian Z, Zhang J. Toll-Like receptor 2 (TLR2) and TLR9 play opposing roles in host innate immunity against *Salmonella enterica* serovar Typhimurium infection. 2015 Apr; 83(4):1641-9.
18. Li P, Zhang L, Zhang M, Zhou C, Lin N. Uric acid enhances PKC-dependent eNOS phosphorylation and mediates cellular ER stress: A mechanism for uric acid-induced endothelial dysfunction. *Int J Mol Med.* 2016 Apr;37(4):989-97.
19. Kuang DB, Zhou JP, Yu LY, Zeng WJ, Xiao J, Zhu GZ, Zhang ZL, Chen XP. DDAH1-V3 transcript might act as miR-21 sponge to maintain balance of DDAH1-V1 in cultured HUVECs. *Nitric Oxide.* 2016 Nov 30;60:59-68.
20. Liu X, Li J, Peng X, Lv B, Wang P, Zhao X, Yu B. Geraniin Inhibits LPS-Induced THP-1 Macrophages Switching to M1 Phenotype via SOCS1/NF-kB Pathway. *Inflammation.* 2016 Aug;39(4):1421-33.
21. Shi SH, Yang WT, Huang KY, Jiang YL, Yang GL, Wang CF, Li Y. β -glucans from *Coriolus versicolor* protect mice against *S. typhimurium* challenge by activation of macrophages. *Int J Biol Macromol.* 2016 May;86:352-61.
22. Liu X, Cai Y, Zhang X, Zhang H, Zheng X, Zhang Z. Carbamoyl Phosphate Synthetase Subunit MoCpa2 Affects Development and Pathogenicity by Modulating Arginine Biosynthesis in *Magnaporthe oryzae*. *Front Microbiol.* 2016 Dec 19;7:2023.
23. Zeng Y, Liu J. Role of glypican-1 in endothelial NOS activation under various steady shear stress magnitudes. *Exp Cell Res.* 2016 Nov 1;348(2):184-189.
24. Chen Q, Tong Q, Ge H, Li W, Liu J, Wang Y, Guo Z, Quan F, Zhang Y. Identification of SP110 in horse (*Equus caballus*): Isolation of novel splice variants and evidence of activation effects on macrophages. *Tuberculosis (Edinb).* 2016 Dec;101:85-94.

Version 2016.12.08