碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301或800-8283301
订货e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

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- Jin WS, Kong ZL, Shen ZF, Jin YZ, Zhang WK, Chen GF. Regulation of hypoxia inducible factor-1α expression by the alteration of redox status in HepG2cells. J Exp Clin Cancer Res. 2011 May;30:61.
- Zhang H, Chen T, Jiang J, Wong YS, Yang F, Zheng W. Selenium-containing allophycocyanin purified from selenium-enriched Spirulina platensisattenuates AAPH-induced oxidative stress in human erythrocytes through inhibition of ROS generation. J Agric Food Chem. 2011 Aug;59(16):8683-90.
- Min L, He S, Chen Q, Peng F, Peng H, Xie M. Comparative proteomic analysis of cellular response of human airway epithelial cells (A549) tobenzo(a)pyrene. Toxicol Mech Methods. 2011 Jun;21(5):374-82.
- Wang L, Hao J, Hu J, Pu J, Lü Z, Zhao L, Wang Q, Yu Q, Wang Y, Li G. Protective Effects of Ginsenosides against Bisphenol A-Induced Cytotoxicity in 15P-1 Sertoli Cells via Extracellular Signal-Regulated Kinase 1/2 Signalling and Antioxidant Mechanisms. Basic Clin Pharmacol Toxicol. 2012 Jul;111(1):42-9.
- Huang HC, Xu, K, Jiang ZF. Curcumin-Mediated Neuroprotection Against Amyloid-β-Induced Mitochondrial Dysfunction Involves the Inhibition of GSK-3β. Journal of Alzheimer's Disease. 2012;32(4):981-96.
- Lu H, Wan J, Jiang R, Xie J, Peng X, Zhang L. Sodium butyrate potentiates carbon tetrachloride-induced acute liver injury in mice. Toxicol Mech Methods. 2012 Oct;22(8):648-55.
- Shi H, Wang Y, Cheng Z, Ye T, Chan Z. Analysis of natural variation in bermudagrass (Cynodon dactylon) reveals physiological responses underlying drought tolerance. PLoS One. 2012;7(12):e53422.
- Wang L, Hao J, Hu J, Pu J, Lü Z, Zhao L, Wang Q, Yu Q, Wang Y, Li G. Protective effects of ginsenosides against Bisphenol A-induced cytotoxicity in 15P-1 Sertoli cells via extracellular signal-regulated kinase 1/2 signalling and antioxidant mechanisms. Basic Clin Pharmacol Toxicol. 2012 Jul;111(1):42-9.
- Yan F, Fang F, Jin H and Jin Y Influence of Tartary Buckwheat Extracts Supplementation on Oxidative Stress Induced by Acute Exhaustive Exercise in Rats International Journal of Pharmacology. 2013 Feb;56:173-8.
- Li CW, Wu XL, Zhao XN, Su ZQ, Chen HM, Wang XF, Zhang XJ, Zeng HF, Chen JN, Li YC, Su ZR. Anti-inflammatory property of the ethanol extract of the root and rhizome of Pogostemon cablin (Blanco) Benth. ScientificWorldJournal. 2013 Dec 9;2013:434151.
- Chen H, Zha J, Liang X, Li J, Wang Z. Effects of the human antiepileptic drug carbamazepine on the behavior, biomarkers, and heat shock proteins in the Asian clam Corbicula fluminea. Aquat Toxicol. 2014 Oct;155:1-8.
- Jiang D, Chen K, Lu X, Gao HJ, Qin ZH, Lin F. Exercise ameliorates the detrimental effect of chloroquine on skeletal muscles in mice via restoring autophagy flux. Acta Pharmacol Sin. 2014 Jan;35(1):135-42.
- Luo C, Li Y, Yang L, Wang X, Long J, Liu J. Superparamagnetic iron oxide nanoparticles exacerbate the risks of reactive oxygen species-mediated external stresses. Arch Toxicol. 2015 Mar;89(3):357-69.
- Zhang G, Zhu J, Zhou Y, Wei Y, Xi L, Qin H, Rao Z, Han M, Ma Y, Wu X. Hesperidin Alleviates Oxidative Stress and Upregulates the Multidrug Resistance Protein 2 in Isoniazid and Rifampicin-Induced Liver Injury in Rats. J Biochem Mol Toxicol. 2016 Jul;30(7):342-9.
- Zhu Z, Chen Y, Shi G, Zhang X. Selenium delays tomato fruit ripening by inhibiting ethylene biosynthesis and enhancing the antioxidant defense system. Food Chem. 2017 Mar 15;219:179-184.
- Zhu Z, Chen Y, Shi G, Zhang X. Selenium delays tomato fruit ripening by inhibiting ethylene biosynthesis and enhancing the antioxidant defense system. Food Chem. 2017 Mar 15;219:179-184.
- He JH, Yan M, Zuo H, Niu S, Yuan J, Weng SP, He J, Xu X. High reduced/oxidized glutathione ratio in infectious spleen and kidney necrosis virus-infected cells contributes to degradation of VP08R multimers. Vet Microbiol. 2017 Aug;207:19-24.
- Zhang S, Shao Q, Geng H, Su S. The effect of royal jelly on the growth of breast cancer in mice. Oncol Lett. 2017 Dec;14(6):7615-7621.

19. Ren A, Liu R, Miao ZG, Zhang X, Cao PF, Chen TX, Li CY, Shi L, Jiang AL, Zhao MW. Hydrogen-rich water regulates effects of ROS balance on morphology, growth and secondary metabolism via glutathione peroxidase in Ganoderma lucidum. *Environ Microbiol*. 2017 Feb;19(2):566-583.
20. Duan G, Shi M, Xie L, Xu M, Wang Y, Yan H, Zhuge Y, Zou X. Increased Glutamine Consumption in Cisplatin-Resistant Cells Has a Negative Impact on Cell Growth. *Sci Rep-Uk*. 2018 Mar 6;8(1):4067.
21. Dai G, Zheng D, Guo W, Yang J, Cheng AY. Cinobufagin Induces Apoptosis in Osteosarcoma Cells Via the Mitochondria-Mediated Apoptotic Pathway. *Cell Physiol Biochem*. 2018;46(3):1134-1147.
22. Xu Y, Sun D, Song C, Wang R, Dong X. MnTMPyP inhibits paraquat-induced pulmonary epithelial-like cell injury by inhibiting oxidative stress. *J Toxicol Sci*. 2018;43(9):545-555.
23. Zhang X, Högy P, Wu X, Schmid I, Wang X, Schulze WX, Jiang D, Fangmeier A. Physiological and Proteomic Evidence for the Interactive Effects of Post-Anthesis Heat Stress and Elevated CO₂ on Wheat. *Proteomics*. 2018 Dec;18(23):e1800262.
24. Xu S, Zhao X, Liu F, Cao Y, Wang B, Wang X, Yin M, Wang Q, Feng X. Crucial role of oxidative stress in bactericidal effect of parthenolide against *Xanthomonas oryzae* pv. *oryzae*. *Pest Manag Sci*. 2018 Dec;74(12):2716-2723.
25. Dai J, Gu L, Su Y, Wang Q, Zhao Y, Chen X, Deng H, Li W, Wang G, Li K. Inhibition of curcumin on influenza A virus infection and influenza pneumonia via oxidative stress, TLR2/4, p38/JNK MAPK and NF-κB pathways. *Int Immunopharmacol*. 2018 Jan;54:177-187.
26. Duan G, Shi M, Xie L, Xu M, Wang Y, Yan H, Zhuge Y, Zou X. Increased Glutamine Consumption in Cisplatin-Resistant Cells Has a Negative Impact on Cell Growth. *Sci Rep-Uk*. 2018 Mar 6;8(1):4067.
27. Wang J, Zhang Y. Neuroprotective effect of berberine agonist against impairment of learning and memory skills in severe traumatic brain injury via Sirt1/p38 MAPK expression. *Mol Med Rep*. 2018 May;17(5):6881-6886.
28. Wu D, Han R, Deng S, Liu T, Zhang T, Xie H, Xu Y. Protective Effects of Flagellin A N/C Against Radiation-Induced NLR Pyrin Domain Containing 3 Inflammasome-Dependent Pyroptosis in Intestinal Cells. *Int J Radiat Oncol*. 2018 May 1;101(1):107-117.
29. Wang Y, An W, Zhang F, Niu M, Liu Y, Shi R. Nebivolol ameliorates kidney damage in Zucker diabetic fatty rats by regulation of oxidative stress/NO pathway: Comparison with captopril. *Clin Exp Pharmacol P*. 2018 Nov;45(11):1135-1148.
30. Yao S, Yan W. Overexpression of Mst1 reduces gastric cancer cell viability by repressing the AMPK-Sir3 pathway and activating mitochondrial fission. *Oncotargets Ther*. 2018 Nov 29;11:8465-8479.
31. Lu C, Chen X, Wang Q, Xu X, Xu B. TNFα promotes glioblastoma A172 cell mitochondrial apoptosis via augmenting mitochondrial fission and repression of MAPK-ERK-YAP signaling pathways. *Oncotargets Ther*. 2018 Oct 18;11:7213-7227.
32. Xie Z, Zhou Y, Duan X, Yang L. Inhibitory effect of Tanshinone IIA on inverted formin-2 protects HaCaT cells against oxidative injury via regulating mitochondrial stress. *J Recept Signal Transduct Res*. 2019 Apr;39(2):134-145.
33. Zhang L, Li S, Wang R, Chen C, Ma W, Cai H. Anti-tumor effect of LATS2 on liver cancer death: Role of DRP1-mediated mitochondrial division and the Wnt/β-catenin pathway. *Biomed Pharmacother*. 2019 Jun;114:108825.
34. Zhang J, Sun L, Li W, Wang Y, Li X, Liu Y. Overexpression of macrophage stimulating 1 enhances the anti-tumor effects of IL-24 in esophageal cancer via inhibiting ERK-Mfn2 signaling-dependent mitophagy. *Biomed Pharmacother*. 2019 Jun;114:108844.
35. Mao M, Zhang T, Wang Z, Wang H, Xu J, Yin F, Wang G, Sun M, Wang Z, Hua Y, Cai Z. Glaucomalyxin A-induced oxidative stress inhibits the activation of STAT3 signaling pathway and suppresses osteosarcoma progression in vitro and in vivo. *Biochim Biophys Acta Mol Basis Dis*. 2019 Jun 1;1865(6):1214-1225.
36. Xu X, Huang L, Zhang Z, Tong J, Mi J, Wu Y, Zhang C, Yan H. Targeting non-oncogene ROS pathway by alantolactone in B cell acute lymphoblastic leukemia cells. *Life Sci*. 2019 Jun 15;227:153-165.
37. Zhang X, Li F, Cui Y, Liu S, Sun H. Mst1 overexpression combined with Yap knockdown augments thyroid carcinoma apoptosis via promoting MIEF1-related mitochondrial fission and activating the JNK pathway. *Cancer Cell Int*. 2019 May 22;19:143.

Version 2020.05.20