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
1. Wang M, Zhao X, Xiao Z, Yin X, Xing T, Xia G. A wheat superoxide dismutase gene TaSOD2 enhances salt resistance through modulating redoxhomeostasis by promoting NADPH oxidase activity. *Plant Mol Biol* 2016 May;91(1-2):115-30.
2. Luo T, Liu G, Long M, Yang J, Song R, Wang Y, Yuan Y, Bian J, Liu X, Gu J, Zou H, Liu Z. Treatment of cadmium-induced renal oxidative damage in rats by administration of alpha-lipoic acid. *ENVIRON SCI POLLUT R* 2017 Jan;24(2):1832-1844.
3. Guo X, Zhang G, Chen L, Khan AA, Gu B, Li B. Newborn Neurons Are Damaged In Vitro by a Low Concentration of Silver Nanoparticles Through the Inflammatory Oxidative Stress Pathway. *DNA Cell Biol* 2017 Dec;36(12):1062-1070.
4. Zhao H, Liu YJ, Liu ZR, Tang DD, Chen XW, Chen YH, Zhou RN, Chen SQ, Niu HX. Role of mitochondrial dysfunction in renal fibrosis promoted by hypochlorite-modified albumin in a remnant kidney model and protective effects of antioxidant peptide SS-31. *Eur J Pharmacol* 2017 Jun 5;804:57-67.
5. Wang W, He Q, Yan W, Sun J, Chen Z, Liu Z, Lu Z, Hou J, Shao Y, Zhou X, Wang A. High glucose enhances the metastatic potential of tongue squamous cell carcinoma via the PKM2 pathway. *ONCOTARGET* 2017 Dec 4;8(67):111770-111779.
6. Wang W, He Q, Sun J, Liu Z, Zhao L, Lu Z, Zhou X, Wang A. Pyruvate kinase M2 deregulation enhances the metastatic potential of tongue squamous cell carcinoma. *ONCOTARGET* 2017 Jul 17;8(40):68252-68262.
7. Lan C, Chen X, Zhang Y, Wang W, Wang WE, Liu Y, Cai Y, Ren H, Zheng S, Zhou L, Zeng C. Curcumin prevents strokes in stroke-prone spontaneously hypertensive rats by improving vascular endothelial function. *BMC CARDIOVASC DISOR* 2018 Mar 1;18(1):43.
8. Zeng X, Yang J, Hu O, Huang J, Ran L, Chen M, Zhang Y, Zhou X, Zhu J, Zhang Q, Yi L, Mi M. Dihydromyricetin Ameliorates Nonalcoholic Fatty Liver Disease by Improving Mitochondrial Respiratory Capacity and Redox Homeostasis Through Modulation of SIRT3 Signaling. *ANTIOXID REDOX SIGN* 2018 Feb 21.
9. Zeng X, Yang J, Hu O, Huang J, Ran L, Chen M, Zhang Y, Zhou X, Zhu J, Zhang Q, Yi L, Mi M. Dihydromyricetin Ameliorates Nonalcoholic Fatty Liver Disease by Improving Mitochondrial Respiratory Capacity and Redox Homeostasis Through Modulation of SIRT3 Signaling. *ANTIOXID REDOX SIGN* 2018 Feb 21.
10. Huang Y, Zhang Y, Lin Z, Han M, Cheng H. Altered serum copper homeostasis suggests higher oxidative stress and lower antioxidant capability in patients with chronic hepatitis B. *MEDICINE* 2018 Jun;97(24):e11137.
11. Lan C, Chen X, Zhang Y, Wang W, Wang WE, Liu Y, Cai Y, Ren H, Zheng S, Zhou L, Zeng C. Curcumin prevents strokes in stroke-prone spontaneously hypertensive rats by improving vascular endothelial function. *BMC CARDIOVASC DISOR* 2018 Mar 1;18(1):43.
12. Liu ZR, Chen SQ, Zou YW, Wu XY, Li HY, Wang XQ, Shi Y, Niu HX. Hypochlorite modified albumins promote cell death in the tubule interstitium in rats via mitochondrial damage in obstructive nephropathy and the protective effects of antioxidant peptides. *FREE RADICAL RES* 2018 May;52(5):616-628.
13. Zhang Q, Liu X, Li N, Zhang J, Yang J, Bu P. Sirtuin 3 deficiency aggravates contrast-induced acute kidney injury. *J Transl Med* 2018 Nov 16;16(1):313.
14. Wu J, Chu Z, Ruan Z, Wang X, Dai T, Hu X. Changes of Intracellular Porphyrin, Reactive Oxygen Species, and Fatty Acids Profiles During Inactivation of Methicillin-Resistant *Staphylococcus aureus* by Antimicrobial Blue Light. *Front Physiol* 2018 Nov 28;9:1658.
15. Li Y, Zhang J, Liu H, Yuan J, Yin Y, Wang T, Cheng B, Sun S, Guo Z. Curcumin ameliorates glyoxylate-induced calcium oxalate deposition and renal injuries in mice. *Phytomedicine* 61:152861. 2019 Aug
16. Zhang Y, Liu Y, Bi X, Hu C, Ding F, Ding W. Therapeutic Approaches in Mitochondrial Dysfunction, Inflammation, and Autophagy in Uremic Cachexia: Role of Aerobic Exercise. *MEDIAT INFLAMM* 2019;2789014. 2019 Aug 18
17. Wang Q, Ye S, Chen X, Xu P, Li K, Zeng S, Huang M, Gao W, Chen J, Zhang Q, Zhong Z, Liu Q. Mitochondrial NOS1 suppresses apoptosis in colon cancer cells through increasing SIRT3 activity. *BIOCHEM BIOPH RES CO* 515(4):517-523. 2019 Aug 6
18. Zhou Q, Wang X, Shao X, Wang H, Liu X, Ke X, Xiong C, Wei L, Zou H. tert-Butylhydroquinone Treatment Alleviates Contrast-Induced Nephropathy in Rats by Activating the Nrf2/Sirt3/SOD2 Signaling Pathway. *Oxid Med Cell Longev* 2019;4657651. 2019 Dec 18
19. Li DP, Chen YL, Jiang HY, Chen Y, Zeng XQ, Xu LL, Ye Y, Ke CQ, Lin G, Wang JY, Gao H. Phosphocreatine attenuates Gynura segetum-induced hepatocyte apoptosis via a SIRT3-SOD2-mitochondrial reactive oxygen species pathway. *DRUG DES DEV THER* 13:2081-2096. 2019 Jun 27
20. Fang J, Zhao X, Li S, Xing X, Wang H, Lazarovici P, Zheng W. Protective mechanism of artemisinin on rat bone marrow-derived mesenchymal stem cells against apoptosis induced by hydrogen peroxide via activation of c-Raf-Erk1/2-p90rsk-CREB pathway. *Stem Cell Res Ther* 10(1):312. 2019 Oct 26
21. Chen S, Fan B. Myricetin protects cardiomyocytes from LPS-induced injury. *Herz* 43(3):265-274. 2018 May
22. Jin-Peng Chen, Hui-Yong Xu, Lin Liao, Zheng Zhang. Resolvin D2 prevents inflammation and oxidative stress in the retina of streptozocin-induced diabetic mice. *Int J Clin Exp Pathol* 2020 Aug 1;13(8):1986-1994.
23. Weijie Yan, Wei Sun, Jiahui Fan, Haiqing Wang, Song Han, Junfa Li, Yanling Yin. Sirt1-ROS-TRAF6 Signaling-Induced Pyroptosis Contributes to Early Injury in Ischemic Mice. *Neurosci Bull* 2020 Aug;36(8):845-859.
24. Siqian Wang, Jingyuan Yang, Tingting Lin, Shengbing Huang, Jianfeng Ma, Xin Xu. Excessive production of mitochondrion-derived reactive oxygen species induced by titanium ions leads to autophagic cell death of osteoblasts via the SIRT3/SOD2 pathway. *Mol Med Rep* 2020 Jul;22(1):257-264.
25. Xuan Xu, Chu Chen, Wen-Jiang Lu, Yi-Ling Su, Jia-Yu Shi, Yu-Chen Liu, Li Wang, Chen-Xi Xiao, Xiang Wu, Qi Lu. Pyrroloquinoline quinone can prevent chronic heart failure by regulating mitochondrial function. *Cardiovasc Diagn Ther* 2020 Jun;10(3):453-469.
26. Shu Li, Qisheng Lin, Xinghua Shao, Xuying Zhu, Jingkui Wu, Bei Wu, Minfang Zhang, Wenyan Zhou, Yijun Zhou, Haijiao Jin, Zhen Zhang, Chaojun Qi, Jianxiao Shen, Shan Mou, Leyi Gu, Zhaohui Ni. Drp1-regulated PARK2-dependent mitophagy protects against renal fibrosis in unilateral ureteral obstruction. *Free Radic Biol Med* 2020 May 20;152:632-649.
27. Jun Pei, Shuyu Cai, Shang Song, Yuangao Xu, Mei Feng, Guangheng Luo, Yuanlin Wang, Fa Sun, Hua Shi, Shuxiong Xu. Normobaric

- hyperoxia plays a protective role against renal ischemia-reperfusion injury by activating the Nrf2/HO-1 signaling pathway *Biochem Biophys Res Commun* 2020 Oct 29;532(1):151-158.
28. Min Lu, Xinglei Qin, Jungong Yao, Yuanyuan Yang, Minghu Zhao, Lin Sun MiR-134-5p targeting XIAP modulates oxidative stress and apoptosis in cardiomyocytes under hypoxia/reperfusion-induced injury *IUBMB Life* 2020 Oct;72(10):2154-2166.
 29. Zilong Lu, Haojie Feng, Xiaokang Shen, Ruyuan He, Heng Meng, Weichen Lin, Qing Geng MiR-122-5p protects against acute lung injury via regulation of DUSP4/ERK signaling in pulmonary microvascular endothelial cells *Life Sci* 2020 Sep 1;256:117851.
 30. Wen-Wen Wang, Ruiyu Han, Hai-Jun He, Jia Li, Si-Yan Chen, Yingying Gu, Chenglong Xie Administration of quercetin improves mitochondria quality control and protects the neurons in 6-OHDA-lesioned Parkinson's disease models *Aging (Albany NY)* 2021 Apr 20;13(8):11738-11751.
 31. Yu Wei, Jiayao Fu, Wenjing Wu, Pengfei Ma, Le Ren, Zimei Yi, Junhua Wu Quercetin Prevents Oxidative Stress-Induced Injury of Periodontal Ligament Cells and Alveolar Bone Loss in Periodontitis *Drug Des Devel Ther* 2021 Aug 12;15:3509-3522.
 32. Dan Liu, Xing Jin, Guanzhen Yu, Mingsong Wang, Lei Liu, Wenjuan Zhang, Jia Wu, Fengying Wang, Jing Yang, Qin Luo, Lili Cai, Xi Yang, Xisong Ke, Yi Qu, Zhenye Xu, Lijun Jia, Wen-Lian Chen Oleonic acid blocks the purine salvage pathway for cancer therapy by inactivating SOD1 and stimulating lysosomal proteolysis *Mol Ther Oncolytics* 2021 Aug 28;23:107-123.
 33. Gang Shen, Yanmei Li, Fuyan Hong, Jing Zhang, Zhenzhen Fang, Wei Xiang, Weiwei Qi, Xia Yang, Guoquan Gao, Ti Zhou A role for Snail-MnSOD axis in regulating epithelial-to-mesenchymal transition markers expression in RPE cells *Biochem Biophys Res Commun* 2021 Dec 31;585:146-154.
 34. Min Xie, Menglin Cheng, Bojun Wang, Ming Jiao, Liangzhu Yu, Haili Zhu 2-Bromopalmitate attenuates inflammatory pain by maintaining mitochondrial fission/fusion balance and function *Acta Biochim Biophys Sin (Shanghai)* 2021 Jan 12;53(1):72-84.
 35. Shujuan Zhou, Lan Sun, Shanhu Qian, Yongyong Ma, Ruyue Ma, Yuqing Dong, Yifen Shi, Songfu Jiang, Haige Ye, Zhijian Shen, Shenghui Zhang, Jianping Shen, Kang Yu, Siqian Wang Iron overload adversely effects bone marrow haematogenesis via SIRT-SOD2-mROS in a process ameliorated by curcumin *Cell Mol Biol Lett* 2021 Jan 13;26(1):2.
 36. Bo Wang, Yang Li, Chao You miR-129-3p Targeting of MCU Protects Against Glucose Fluctuation-Mediated Neuronal Damage via a Mitochondrial-Dependent Intrinsic Apoptotic Pathway *Diabetes Metab Syndr Obes* 2021 Jan 15;14:153-163.
 37. Lingyan Jin, Feng Gao, Tiannan Jiang, Binghua Liu, Caiyao Li, Xinghua Qin, Qiangsun Zheng Hyper-O-GlcNAcylation impairs insulin response against reperfusion-induced myocardial injury and arrhythmias in obesity *Biochem Biophys Res Commun* 2021 Jun 18;558:126-133.
 38. Li Zhang, Xiufen Chen, Mingxiu Chang, Boning Jiao MiR-30c-5p/ATG5 Axis Regulates the Progression of Parkinson's Disease *Front Cell Neurosci* 2021 May 25;15:644507.
 39. Bin Ye, Liqun Lai Yu Shi An Chang Fang Ameliorates TNBS-Induced Colitis in Mice by Reducing Inflammatory Response and Protecting the Intestinal Mucosal Barrier *Evid Based Complement Alternat Med* 2021 May 4;2021:8870901.
 40. Nan Zhang, Hai-Han Liao, Hong Feng, Shan-Qi Mou, Wen-Jing Li, Xiahenazi Aiyasiding, Zheng Lin, Wen Ding, Zi-Ying Zhou, Han Yan, Si Chen, Qi-Zhu Tang Knockout of AMPK α 2 Blocked the Protection of Sestrin2 Overexpression Against Cardiac Hypertrophy Induced by Pressure Overload *Front Pharmacol* 2021 Nov 17;12:716884.
 41. Wenqi Liang, Chunli Zhao, Zhongrui Chen, Zijing Yang, Ke Liu, Shusheng Gong Sirtuin-3 Protects Cochlear Hair Cells Against Noise-Induced Damage via the Superoxide Dismutase 2/Reactive Oxygen Species Signaling Pathway *Front Cell Dev Biol* 2021 Nov 18;9:766512.
 42. Shanshan Dai, Bozhi Ye, Lingfeng Zhong, Yanghao Chen, Guangliang Hong, Guangju Zhao, Lan Su, Zhongqiu Lu GSDMD Mediates LPS-Induced Septic Myocardial Dysfunction by Regulating ROS-dependent NLRP3 Inflammasome Activation *Front Cell Dev Biol* 2021 Nov 8;9:779432.
 43. Shouzhu Xu, Jie Xu, Ting Hao, Yu Yan, Shihao Zhang, Aihong Li, Chuandao Shi, Qiling Liu, Jing Zhao Paeonol alleviates lipopolysaccharide-induced hepatocytes injury through alteration of mitochondrial function and NF- κ B translocation *Mol Med Rep* 2021 Nov;24(5):779.
 44. Meng Wang, Lingchen Wang, Yuan Zhou, Xiaoxuan Feng, Chaoyang Ye, Chen Wang Icariin attenuates renal fibrosis in chronic kidney disease by inhibiting interleukin-1 β /transforming growth factor- β -mediated activation of renal fibroblasts *Phytother Res* 2021 Nov;35(11):6204-6215.
 45. Baoyue Zhang, Wenwen Lian, Jun Zhao, Zhe Wang, Ailin Liu, Guanhua Du DL0410 Alleviates Memory Impairment in D-Galactose-Induced Aging Rats by Suppressing Neuroinflammation via the TLR4/MyD88/NF- κ B Pathway *Oxid Med Cell Longev* 2021 Oct 4;2021:6521146.
 46. Jun Pei, Moudong Wu, Shuyu Cai, Jinpu Peng, Xiong Zhan, Dan Wang, Wei Wang, Nini An The Protective Effect of Ursolic Acid on Unilateral Ureteral Obstruction in Rats by Activating the Nrf2/HO-1 Antioxidant Signaling Pathway *Comput Intell Neurosci* 2022 Aug 25;2022:3690524.
 47. Yingying Zha, Yan Jin, Xinxing Wang, Lin Chen, Xulai Zhang, Ming Wang Long-term maintenance of synaptic plasticity by Fullereneol Ameliorates lead-induced-impaired learning and memory in vivo *J Nanobiotechnology* 2022 Aug 1;20(1):348.
 48. Yingchao Wu, Dajin Pi, Yiliu Chen, Qian Zuo, Lizhu Lin, Mingzi Ouyang Yifei sanjie Pills Alleviate Chemotherapy-Related Fatigue by Reducing Skeletal Muscle Injury and Inhibiting Tumor Growth in Lung Cancer Mice *Evid Based Complement Alternat Med* 2022 Aug 22;2022:2357616.
 49. Jingkui Wu, Xinghua Shao, Jianxiao Shen, Qisheng Lin, Xuying Zhu, Shu Li, Jialin Li, Wenyan Zhou, Chaojun Qi, Zhaohui Ni Downregulation of PPAR α mediates FABP1 expression, contributing to IgA nephropathy by stimulating ferroptosis in human mesangial cells *Int J Biol Sci* 2022 Aug 29;18(14):5438-5458.
 50. Yang Yao, Rong Li, Dan Liu, Lihui Long, Na He Rosmarinic acid alleviates acetaminophen-induced hepatotoxicity by targeting Nrf2 and NEK7-NLRP3 signaling pathway *Ecotoxicol Environ Saf* 2022 Aug;241:113773.
 51. Lu Chen, Yanping Chen, Jianbao Huang, Jiyan Zhang LncRNA LINC00707 serves as a sponge of miR-382-5p to alleviate lipopolysaccharide (LPS)-induced WI-38 cell injury through upregulating NKAP in infantile pneumonia *Autoimmunity* 2022 Aug;55(5):328-338.
 52. Fei Song, Jiale Lin, Houjian Zhang, Yuli Guo, Yijie Mao, Zuguo Liu, Gang Li, Yan Wang Long-Term Sleep Deprivation-Induced Myocardial Remodeling and Mitochondrial Dysfunction in Mice Were Attenuated by Lipoic Acid and N-Acetylcysteine *Pharmaceuticals (Basel)* 2022 Dec 29;16(1):51.

53. Ying Wang, Yuerong Xu, Wangang Guo, Yexian Fang, Lang Hu, Runze Wang, Ran Zhao, Dong Guo, Bingchao Qi, Gaotong Ren, Jun Ren, Yan Li, Mingming Zhang Ablation of Shank3 alleviates cardiac dysfunction in aging mice by promoting CaMKII activation and Parkin-mediated mitophagy *Redox Biol* 2022 Dec;58:102537.
54. Yuqiong Li, Huiyu Wen, Yuwei Yang, Zhengwei Zhao, Haihui Gao, Hongbing Li, Meizhou Huang Potential prognostic markers of retained placenta in dairy cows identified by plasma metabolomics coupled with clinical laboratory indicators *Vet Q* 2022 Dec;42(1):199-212.
55. Binbing Liu, Yang Tian, Yuchen Li, Pei Wu, Yongzhi Zhang, Jiaolin Zheng, Huaizhang Shi ACEA Attenuates Oxidative Stress by Promoting Mitophagy via CB1R/Nrf1/PINK1 Pathway after Subarachnoid Hemorrhage in Rats *Oxid Med Cell Longev* 2022 Feb 24;2022:1024279.
56. Wenxi He, Haifeng Zhou, Xiaoyu He Aloperine protects beta-cells against streptozocin-induced injury to attenuate diabetes by targeting NOS1 *Eur J Pharmacol* 2022 Feb 5;916:174721.
57. Jinzhu Yang, Tiantian Wang, Gang Lin, Mingzhu Li, Yanjiao Zhang, Kangsen Mai The Assessment of Dietary Organic Zinc on Zinc Homeostasis, Antioxidant Capacity, Immune Response, Glycolysis and Intestinal Microbiota in White Shrimp (*Litopenaeus vannamei* Boone, 1931) *Antioxidants (Basel)* 2022 Jul 29;11(8):1492.
58. Qian Liu, Shujuan Guo, Yanli Huang, Xiuqun Wei, Li Liu, Fangjun Huo, Ping Huang, Yafei Wu, Weidong Tian Inhibition of TRPA1 Ameliorates Periodontitis by Reducing Periodontal Ligament Cell Oxidative Stress and Apoptosis via PERK/eIF2 α /ATF-4/CHOP Signal Pathway *Oxid Med Cell Longev* 2022 Jun 10;2022:4107915.
59. Haiyan Li, Dong-Hua Yang, Yanmei Zhang, Fuchun Zheng, Fenfei Gao, Jiajia Sun, Ganggang Shi Geniposide suppresses NLRP3 inflammasome-mediated pyroptosis via the AMPK signaling pathway to mitigate myocardial ischemia/reperfusion injury *Chin Med* 2022 Jun 17;17(1):73.
60. Mou Zhang, Haotian Shi, Ningning Li, Nana Wei, Yan Tian, Jinfeng Peng, Xiaochen Chen, Liyuan Zhang, Meixiang Zhang, Hansong Dong Aquaporin OsPIP2;2 links the H2O2 signal and a membrane-anchored transcription factor to promote plant defense *Plant Physiol* 2022 Mar 28;188(4):2325-2341.
61. Xiaoxiao Tao, Yaping Jiang, Xian Zheng, Xiaoxiao Ji, Feifei Peng Dihydromyricetin ameliorates oxygen-glucose deprivation and re-oxygenation-induced injury in HT22 cells by activating the Wnt/ β -catenin signaling pathway *Mol Med Rep* 2022 Mar;25(3):103.
62. Fei Song, Yi-Jie Mao, Yu Hu, Shan-Shan Zhao, Ruiying Wang, Wei-Yin Wu, Gui-Rong Li, Yan Wang, Gang Li Acacetin attenuates diabetes-induced cardiomyopathy by inhibiting oxidative stress and energy metabolism via PPAR- α /AMPK pathway *Eur J Pharmacol* 2022 May 5;922:174916.
63. Ji-Hua Shi, Dong-Jing Yang, Qiang Jin, Nuo Cheng, Yuan-Bin Shi, Yang Bai, Dong-Sheng Yu, Wen-Zhi Guo, Guang-Bo Ge, Shui-Jun Zhang Cytochrome P450 2E1 predicts liver functional recovery from donation after circulatory death using air-ventilated normothermic machine perfusion *Sci Rep* 2022 May 6;12(1):7446.
64. Fan Li, Xiaojing Wu, Hongli Liu, Mengqi Liu, Zhengkai Yue, Zhenyu Wu, Lei Liu, Fuchang Li Copper Depletion Strongly Enhances Ferroptosis via Mitochondrial Perturbation and Reduction in Antioxidative Mechanisms *Antioxidants (Basel)* 2022 Oct 22;11(11):2084.
65. Shi-Meng Liu, Ya-Rong Zhang, Yao Chen, Deng-Ren Ji, Jie Zhao, Su Fu, Mo-Zhi Jia, Yan-Rong Yu, Chao-Shu Tang, Wei Huang, Ye-Bo Zhou, Yong-Fen Qi Intermedin Alleviates Vascular Calcification in CKD through Sirtuin 3-Mediated Inhibition of Mitochondrial Oxidative Stress *Pharmaceuticals (Basel)* 2022 Oct 2;15(10):1224.
66. Zhi-Bin Dong, Yu-Jia Wang, Meng-Lin Cheng, Bo-Jun Wang, Hong Lu, Hai-Li Zhu, Ling Liu, Min Xie 2-Bromopalmitate decreases spinal inflammation and attenuates oxaliplatin-induced neuropathic pain via reducing Drp1-mediated mitochondrial dysfunction *PLoS One* 2022 Oct 31;17(10):e0275428.
67. Tianyi Wang, Tingting Zhou, Mingming Xu, Shuo Wang, Anqi Wu, Mingyang Zhang, You Lang Zhou, Jiahai Shi Platelet membrane-camouflaged nanoparticles carry microRNA inhibitor against myocardial ischaemia-reperfusion injury *J Nanobiotechnology* 2022 Oct 4;20(1):434.
68. Xun Yu, Cui Meng, Xiaocheng Tan, Yuwei Su, Zhiping Cao, Tzann-Shun Hwang, Lin Li RsmA3 modulates RpoS through the RetS-Gac-Rsm signalling pathway in response to H2 O2 stress in the phytopathogen *Pseudomonas syringae* *Environ Microbiol* 2022 Oct;24(10):4755-4770.
69. Jiwei Sun, Tiantian Chen, Baoying Zhao, Wenjie Fan, Yufeng Shen, Haojie Wei, Man Zhang, Wenhao Zheng, Jinfeng Peng, Jinyu Wang, Yifan Wang, Lihong Fan, Yingying Chu, Lili Chen, Cheng Yang Acceleration of Oral Wound Healing under Diabetes Mellitus Conditions Using Bioadhesive Hydrogel *ACS Appl Mater Interfaces* 2023 Jan 11;15(1):416-431.
70. Hao-Xi Zhao, Zheng Zhang, Fang Hu, Qi-Feng Wei, Yong-Sheng Yu, Hao-Dong Zhao Swimming exercise activates peroxisome proliferator-activated receptor-alpha and mitigates age-related renal fibrosis in rats *Mol Cell Biochem* 2023 May;478(5):1109-1116.

Version 2024.03.12