碧云天生物技术/Beyotime Biotechnology  
 订货热线: 400-1683301 或 800-8283301  
 订货 e-mail: order@beyotime.com  
 技术咨询: info@beyotime.com  
 网址: http://www.beyotime.com

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- Citrus medica L. Food Chem. 2015 Apr 15;173:54-60.
5. Xia Q, Wang L, Xu C, Mei J, Li Y. . Effects of germination and high hydrostatic pressure processing on mineral elements, amino acids and antioxidants in vitro bioaccessibility, as well as starch digestibility Food Chem. 2017 Jan 1;214:533-42.
  6. Xia Q, Wang L, Xu C, Mei J, Li Y. . Effects of germination and high hydrostatic pressure processing on mineral elements, amino acids and antioxidants in vitro bioaccessibility, as well as starch digestibility in brown rice (*Oryza sativa* L.). Food Chem. 2017 Jan 1;214:533-542.
  7. Kong Y, Wei ZF, Fu YJ, Gu CB, Zhao CJ, Yao XH, Efferth Thomas. . Negative-pressure cavitation extraction of cajanin stilbene acid and pinostrobin from pigeon pea [*Cajanus cajan* (L.) Millsp.] leaves and evaluation of antioxidant activity. Food Chemistry. 2011 Oct;128(3):596-605.
  8. Luo JG, Li L, Kong LY. . Preparative separation of phenylpropanoid glycerides from the bulbs of *Lilium lancifolium* by high-speed counter-current chromatography and evaluation of their antioxidant activities. Food Chemistry. 2012 Apr;131(1):1056-62.
  9. He TB, Huang YP, Huang Y, Wang XJ, Hu JM, Sheng J. . Structural elucidation and antioxidant activity of an arabinogalactan from the leaves of *Moringa oleifera*. Int J Biol Macromol. 2018 Jun;112:126-133.
  10. Yao C, Li G, Qian Y, Cai M, Yin H, Xiao L, Tang W, Guo F, Shi B. . Protection of Pentoxifylline against Testis Injury Induced by Intermittent Hypobaric Hypoxia. Oxid Med Cell Longev. 2016;2016:3406802.
  11. Jiayin Lu, Jiaqiang Huang, Shisu Zhao, Wenjiao Xu, Yaoxing Chen, Yuanyuan Li, Zixu Wang, Yanjun Dong, Renrong You, Jing Cao, Yulan Dong. . FOXO1 Is a Critical Switch Molecule for Autophagy and Apoptosis of Sow Endometrial Epithelial Cells Caused by Oxidative Stress Oxid Med Cell Longev. 2021 Dec 21;2021:1172273.
  12. Zhang JX, Yang JR, Chen GX, Tang LJ, Li WX, Yang H, Kong X. . Sesamin ameliorates arterial dysfunction in spontaneously hypertensive rats via downregulation of NADPH oxidase subunits and upregulation of eNOS expression. Acta Pharmacol Sin. 2013 Jul 5;34(7):912-20.
  13. Li W, Wu Y, Ren C, Lu Y, Gao Y, Zheng X, Zhang C. . The activity of recombinant human neuroglobin as an antioxidant and free radical scavenger. Proteins. 2011 Jan;79(1):115-25.
  14. Liu F, Li XL, Lin T, He DW, Wei GH, Liu JH, Li LS. . The cyclophosphamide metabolite, acrolein, induces cytoskeletal changes and oxidative stress in Sertoli cells. Mol Biol Rep. 2012 Jan;39(1):493-500.
  15. Feng Liu, Xu-Liang Li, Tao Lin, Da-Wei He, Guang-Hui Wei, Jun-Hong Liu, Lu-Sheng Li. . The cyclophosphamide metabolite, acrolein, induces cytoskeletal changes and oxidative stress in Sertoli cells Mol Biol Rep. 2012 Jan;39(1):493-500.
  16. Yang T, Luo F, Shen Y, An J, Li X, Liu X, Ying B, Liao Z, Dong J, Guo L, Wang T, Xu D, Chen L, Wen F. . Quercetin attenuates airway inflammation and mucus production induced by cigarette smoke in rats. Int Immunopharmacol. 2012 May;13(1):73-81.
  17. Wang L, Hao J, Hu J, Pu J, Lv Z, Zhao L, Wang Q, Yu Q, Wang Y, Li G. . Protective Effects of Ginsenosides against Bisphenol A-Induced Cytotoxicity in 15P-1 Sertoli Cells via Extracellular Signal-Regulated Kinase 1/2 Signalling and Antioxidant Mechanisms. BASIC CLIN PHARMACOL. 2012 Jul;111(1):42-9.
  18. Zhou Y, Yuan J, Wu J, Han X. . The toxic effects of microcystin-LR on rat spermatogonia in vitro. Toxicol Lett. 2012 Jul 7;212(1):48-56.
  19. Yang L, Qu M, Wang Y, Duan H, Chen P, Wang Y, Shi W, Danielson P, Zhou Q. . Trichostatin A inhibits transforming growth factor- $\beta$ -induced reactive oxygen species accumulation and myofibroblast differentiation via enhanced NF-E2-related factor 2-antioxidant response element signaling. Mol Pharmacol. 2013 Mar;83(3):671-80.
  20. Qingfeng Sheng, Zhibao Lv, Wei Cai, Huanlei Song, Linxi Qian, Xueli Wang. . Protective effects of hydrogen-rich saline on necrotizing enterocolitis in neonatal rats J Pediatr Surg. 2013 Aug;48(8):1697-706.
  21. Zou Y, Wang R, Guo H, Dong M. . Phytoestrogen  $\beta$ -Ecdysterone Protects PC12 Cells Against MPP+-Induced Neurotoxicity In Vitro: Involvement of PI3K-Nrf2-Regulated Pathway. Toxicol Sci. 2015 Sep;147(1):28-38.
  22. Zuoqing Yuan, Jianyong Zhang, Changchao Tu, Zhijing Wang, Wenpeng Xin. . The protective effect of blueberry anthocyanins against perfluorooctanoic acid-induced disturbance in planarian (*Dugesia japonica*) Ecotoxicol Environ Saf. 2016 May;127:170-4.
  23. Sun CL, Chen L, Xu J, Qu W, Guan L, Liu WY, Akihisa T, Feng F, Zhang J. . Melanogenesis-Inhibitory and Antioxidant Activities of Phenolics from *Periploca forrestii*. Chem Biodivers. 2017 Aug;14(8).
  24. Hu XF, Zhang Q, Zhang PP, Sun LJ, Liang JC, Morris-Natschke SL, Chen Y, Lee KH. . Evaluation of in vitro/in vivo anti-diabetic effects and identification of compounds from *Physalis alkekengi*. Fitoterapia. 2018 Feb 12. pii: S0367-326X(18)30015-7.
  25. Tang Y, Han L, Chen X, Xie M, Kong W, Wu Z. . Dietary Supplementation of Probiotic *Bacillus subtilis* Affects Antioxidant Defenses and Immune Response in Grass Carp Under *Aeromonas hydrophila* Challenge. Probiotics Antimicrob Proteins. 2018 Apr 13.
  26. Hua Q, Chen C, Tel Zur N, Wang H, Wu J, Chen J, Zhang Z, Zhao J, Hu G, Qin Y. . Metabolomic characterization of pitaya fruit from three red-skinned cultivars with different pulp colors. PLANT PHYSIOL BIOCH. 2018 May;126:117-125.
  27. Huang ZQ, Chen P, Su WW, Wang YG, Wu H, Peng W, Li PB. . Antioxidant Activity and Hepatoprotective Potential of Quercetin 7-Rhamnoside In Vitro and In Vivo. Molecules. 2018 May 16;23(5). pii: E1188.
  28. Zhao C, Xia T, Du P, Duan W, Zhang B, Zhang J, Zhu S, Zheng Y, Wang M, Yu Y. . Chemical Composition and Antioxidant Characteristic of Traditional and Industrial Zhenjiang Aromatic Vinegars during the Aging Process. Molecules. 2018 Nov 12;23(11). pii: E2949.
  29. Zhang Y, Geng J, Hong Y, Jiao L, Li S, Sun R, Xie Y, Yan C, Aa J, Wang G. . Orally Administered Crocin Protects Against Cerebral Ischemia/Reperfusion Injury Through the Metabolic Transformation of Crocetin by Gut Microbiota. Front Pharmacol. 2019 Apr 30 10:440.
  30. Huo L, Su Y, Xu G, Zhai L, Zhao J. . Sulforaphane Protects the Male Reproductive System of Mice from Obesity-Induced Damage: Involvement of Oxidative Stress and Autophagy. INT J ENV RES PUB HE. 2019 Oct 7 16(19). pii: E3759.
  31. Wenjin Guo, Juxiong Liu, Jingxuan Sun, Qian Gong, He Ma, Xingchi Kan, Yu Cao, Jianfa Wang, Shoupeng Fu. . Butyrate alleviates oxidative stress by regulating NRF2 nuclear accumulation and H3K9/14 acetylation via GPR109A in bovine mammary epithelial cells and mammary glands Free Radic Biol Med. 2020 May 20;152:728-742.
  32. Linyu Cao, Xulin Li, Ruiming Xu, Kaiyun Yao, Wanqi Yang, Hongliang Zhu, Guibin Wang, Jianjun Zhang. . DUOX2, a common modulator in preventive

- effects of monoamine-based antidepressants on water immersion restraint stress- and indomethacin- induced gastric mucosal damage *Eur J Pharmacol.* 2020 Jun 5;876:173058.
33. Lin Qiaoxia, Zhou Yujie, Yin Meng, Cheng Yizhu, Wei Yan, Hu Yinchun, Lian Xiaojie, Chen Weiyi, Huang Di . Hydroxyapatite/tannic acid composite coating formation based on Ti modified by TiO<sub>2</sub> nanotubes *Colloids Surf B Biointerfaces.* 2020 Dec;196:111304.
34. Fan Zhang, Huabo Chen, Jie Lan, Kaichao Song, Xianggen Wu . Preparation and in vitro/in vivo evaluations of novel ocular micelle formulations of hesperetin with glycyrrhizin as a nanocarrier *Exp Eye Res.* 2021 Jan;202:108313.
35. Jiyu Chen, Yanli Huang, Xiaochuan Hu, Xiaohong Bian, Sihui Nian . Gastrodin prevents homocysteine-induced human umbilical vein endothelial cells injury via PI3K/Akt/eNOS and Nrf2/ARE pathway *J Cell Mol Med.* 2021 Jan;25(1):345-357.
36. Yong Hou, Guangjie Zhang, Min Li, Bowen Li, Li Chen, Ying Tian, Shijun Liu, Bin Li, Junxing Dong . Antioxidant and anti-inflammatory constituents from *Flos populi* *Nat Prod Res.* 2021 Feb;35(4):570-578.
37. Yuzhen Hou, Meng Xin, Qiqi Li, Xianggen Wu . Glycyrrhizin micelle as a genistein nanocarrier: Synergistically promoting corneal epithelial wound healing through blockage of the HMGB1 signaling pathway in diabetic mice *Exp Eye Res.* 2021 Mar;204:108454.
38. Leizhou Guo, Mingming Zhao, Yin Tang, Jiahui Han, Yuan Gui, Jiaming Ge, Shijie Jiang, Qilin Dai, Wei Zhang, Min Lin, Zhengfu Zhou, Jin Wang . Modular Assembly of Ordered Hydrophilic Proteins Improve Salinity Tolerance in *Escherichia coli* *Int J Mol Sci.* 2021 Apr 25;22(9):4482.
39. Qianwen Ding, Zhen Zhang, Yu Li, Hongliang Liu, Qiang Hao, Yalin Yang, Einar Ringø, Rolf Erik Olsen, Jihong Liu Clarke, Chao Ran, Zhigang Zhou . Propionate induces intestinal oxidative stress via Sod2 propionylation in zebrafish *iScience.* 2021 May 5;24(6):102515.
40. Xintong Wang, Zixu Wang, Jing Cao, Yulan Dong, Yaoxing Chen . Melatonin ameliorates anxiety-like behaviors induced by sleep deprivation in mice: Role of oxidative stress, neuroinflammation, autophagy and apoptosis *Brain Res Bull.* 2021 Sep;174:161-172.
41. Yang Bai, Lin Wang, Lisheng Zhao, E Lingling, Shuo Yang, Shunyi Jia, Ning Wen . Antibacterial and Antioxidant Effects of Magnesium Alloy on Titanium Dental Implants *Comput Math Methods Med.* 2022 Jan 6:2022:6537676.
42. Han Wu, Wenzheng Shi, Haiyuan Huang, Luzhu Chen, Nian Li, Yufan Du, Nuo Pei, Jie Sun . The impact of Flavourzyme on flavor, antioxidant activity, and sensory attributes of salted grass carp by brine injection and brining *J Sci Food Agric.* 2022 Mar 15;102(4):1665-1673.
43. Xiaojun Diao, Qi Cui, Ning Tian, Zixian Zhou, Wenjing Xiang, Yanlin Jiang, Jungang Deng, Hongzhan Liao, Xiaohui Lin, Qinghua Li, Rujia Liao . Hemorrhage-Induced Sphingosine Kinase 1 Contributes to Ferroptosis-Mediated Secondary Brain Injury in Intracerebral Hemorrhage *Mol Neurobiol.* 2022 Mar;59(3):1381-1397.
44. Ziru Liu, Jiameng Qu, Fan Ke, Haotian Zhang, Yiwen Zhang, Qian Zhang, Qing Li, Kaishun Bi, Huarong Xu . Material Basis Elucidation and Quantification of Dandelion through Spectrum-Effect Relationship Study between UHPLC Fingerprint and Antioxidant Activity via Multivariate Statistical Analysis *Molecules.* 2022 Apr 20;27(9):2632.
45. Jingming Ren, Rong Huang, Yanjie Li, Ruiyang Chen, Hongqi Tian, Chenlu Liu . Radioprotective effects and mechanism of HL-003 on radiation-induced salivary gland damage in mice *Sci Rep.* 2022 May 19;12(1):8419.
46. Qi Liu, Zixu Wang, Jing Cao, Yulan Dong, Yaoxing Chen . Dim Blue Light at Night Induces Spatial Memory Impairment in Mice by Hippocampal Neuroinflammation and Oxidative Stress Antioxidants (Basel). 2022 Jun 22;11(7):1218.
47. Liu-Hong Shen, Lei Fan, Yue Zhang, Yu Shen, Zhe-Tong Su, Guang-Neng Peng, Jun-Liang Deng, Zhi-Jun Zhong, Xiao-Feng Wu, Shu-Min Yu, Sui-Zhong Cao, Xiao-Lan Zong . Antioxidant Capacity and Protective Effect of Cow Placenta Extract on D-Galactose-Induced Skin Aging in Mice *Nutrients.* 2022 Nov 3;14(21):4659.

Version 2024.03.12