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- Li CG, Zeng QZ, Chen MY, Xu LH, Zhang CC, Mai FY, Zeng CY, He XH, Ouyang DY. Evodiamine Augments NLRP3 Inflammasome Activation and Anti-bacterial Responses Through Inducing α-Tubulin Acetylation. *Front Pharmacol.* 2019 Mar 26;10:290.
- Yi Zhang, Qingyang Liang, Yanan Zhang, Lei Hong, Da Lei, Li Zhang. Olmesartan alleviates bleomycin-mediated vascular smooth muscle cell senescence via the miR-665/SDC1 axis. *Am J Transl Res.* 2020 Sep 15;12(9):5205-5220.
- Chenyang Xue, Wei Chen, Aiwu Yuan, Cheng Chen, Shuaihu Li, Kai Chen, Yang Zhao, Tian Xiao, Genze Shao, Yongdong Zou, Duo Zheng. Dezocine, An Opioid Analgesic, Exerts Antitumor Effects in Triple-Negative Breast Cancer by Targeting Nicotinamide Phosphoribosyltransferase. *Front Pharmacol.* 2021 Apr 12;12:600296.
- Yu Zhang, Zhengzhe Zhang, Yinguang Chen. Biochar Mitigates N₂O Emission of Microbial Denitrification through Modulating Carbon Metabolism and Allocation of Reducing Power. *Environ Sci Technol.* 2021 Jun 15;55(12):8068-8078.
- Rui Jia, Jinliang Du, Liping Cao, Wenrong Feng, Qin He, Pao Xu, Guojun Yin. Application of transcriptome analysis to understand the adverse effects of hydrogen peroxide exposure on brain function in common carp (*Cyprinus carpio*). *Environ Pollut.* 2021 Oct 1;286:117240.
- Wen-Jing Cui, Biliang Zhang, Ran Zhao, Li-Xue Liu, Jian Jiao, Ziding Zhang, Chang-Fu Tian. Lineage-Specific Rewiring of Core Pathways Predating Innovation of Legume Nodules Shapes Symbiotic Efficiency. *mSystems.* doi: 10.1128/mSystems.01299-20.
- Jun Li, Linlin Xie, Shulan Qian, Yuhang Tang, Mingjie Shen, Shanshan Li, Jie Wang, Lie Xiong, Jie Lu, Weihong Zhong. A Type VI Secretion System Facilitates Fitness, Homeostasis, and Competitive Advantages for Environmental Adaptability and Efficient Nicotine Biodegradation. *Appl Environ Microbiol.* doi: 10.1128/AEM.03113-20.
- Feng Peng, Jing Chen, Xiuxia Liu, Ye Li, Chunli Liu, Yankun Yang, Zhonghu Bai. The PhoPR two-component system responds to oxygen deficiency and regulates the pathways for energy supply in *Corynebacterium glutamicum*. *World J Microbiol Biotechnol.* doi:

- 10.1007/s11274-021-03131-1.
9. Shan He, Hanxiang Zhang, Yang Lu, Zhaosi Zhang, Xiang Zhang, Nian Zhou, Zhenming Hu. Nampt promotes osteogenic differentiation and lipopolysaccharide-induced interleukin-6 secretion in osteoblastic MC3T3-E1 cells. *Aging* (Albany NY). doi: 10.18632/aging.202434.
 10. Yuan-Yuan Cheng, Wen-Jing Wang, Shi-Ting Ding, Ming-Xing Zhang, Ai-Guo Tang, Ling Zhang, Dao-Bo Li, Bing-Bing Li, Guo-Zhi Deng, Chao Wu. Pyruvate accelerates palladium reduction by regulating catabolism and electron transfer pathway in *Shewanella oneidensis*. *Appl Environ Microbiol*. doi: 10.1128/AEM.02716-20.
 11. Tianle Xu, Xubin Lu, Abdelaziz Adam Idriss Arbab, Xinyue Wu, Yongjiang Mao, Juan J Looor, Zhangping Yang. Metformin acts to suppress β -hydroxybutyric acid-mediated inflammatory responses through activation of AMPK signaling in bovine hepatocytes. *J Anim Sci*. doi: 10.1093/jas/skab153.
 12. Qiang-Zhong Pi, Xiao-Wen Wang, Zhao-Lei Jian, Dan Chen, Cheng Zhang, Qing-Chen Wu. Melatonin Alleviates Cardiac Dysfunction Via Increasing Sirt1-Mediated Beclin-1 Deacetylation and Autophagy During Sepsis. *Inflammation*. doi: 10.1007/s10753-021-01413-2.
 13. Haoling Chen, Jun Kang, Fuping Zhang, Tong Yan, Wenguo Fan, Hongwen He, Fang Huang. SIRT4 regulates rat dental papilla cell differentiation by promoting mitochondrial functions. *Int J Biochem Cell Biol*. doi: 10.1016/j.biocel.2021.105962.
 14. Yuliang Wen, Jiang Hu, Jiqing Wang, Xiu Liu, Shaobin Li, Yuzhu Luo. Effect of glycolysis and heat shock proteins on hypoxia adaptation of Tibetan sheep at different altitude. *Gene*. doi: 10.1016/j.gene.2021.145893.
 15. Yu Qiu, Ying Yu, Xiao-Mei Qin, Tao Jiang, Yan-Fang Tan, Wen-Xian Ouyang, Zheng-Hui Xiao, Shuang-Jie Li. CircTLK1 modulates sepsis-induced cardiomyocyte apoptosis via enhancing PARP1/HMGB1 axis-mediated mitochondrial DNA damage by sponging miR-17-5p. *J Cell Mol Med*. doi: 10.1111/jcmm.16738.
 16. Han-Fang Zeng, Jie Xu, Xin-Ling Wang, Shu-Jie Li, Zhao-Yu Han. Nicotinamide mononucleotide alleviates heat stress-induced oxidative stress and apoptosis in BMECs through reducing mitochondrial damage and endoplasmic reticulum stress. *Ecotoxicol Environ Saf*. doi: 10.1016/j.ecoenv.2022.113441.
 17. Fei Guo, Taniya Fernando, Xiaoyong Zhu, Yingli Shi. The overexpression of neurokinin B-neurokinin 3 receptor system exerts direct effects on the ovary under PCOS-like conditions to interfere with mitochondrial function. *Am J Reprod Immunol*. doi: 10.1111/aji.13663.
 18. Dan Wang, Yanling Kuang, Zhicheng Wan, Pei Li, Jiangchao Zhao, Huiling Zhu, Yulan Liu. Aspartate Alleviates Colonic Epithelial Damage by Regulating Intestinal Stem Cell Proliferation and Differentiation via Mitochondrial Dynamics. *Mol Nutr Food Res*. doi: 10.1002/mnfr.202200168.
 19. Zhuan Feng, Yifei Qin, Fei Huo, Zhe Jian, Xiaomin Li, Jiejie Geng, Yong Li, Jiao Wu. NMN recruits GSH to enhance GPX4-mediated ferroptosis defense in UV irradiation induced skin injury. *Biochim Biophys Acta Mol Basis Dis*. doi: 10.1016/j.bbadis.2021.166287.
 20. Bao Guo, Xiaoyan Xu, Miaomiao Shao, Xu Yang, Gaofei He, Kangwei Qi, Jianxin Gu, Lan Wang. UDP-glucose 6-dehydrogenase lessens sorafenib sensitivity via modulating unfolded protein response. *Biochem Biophys Res Commun*. doi: 10.1016/j.bbrc.2022.05.048.
 21. Xihui Chen, Fangfang Liu, Bowen Li, Yufeng Wang, Lijuan Yuan, Anan Yin, Qi Chen, Weihong Hu, Yan Yao, Mengjie Zhang, YuanMing Wu, Kun Chen. Neuropathy-associated Fars2 deficiency affects neuronal development and potentiates neuronal apoptosis by impairing mitochondrial function. *Cell Biosci*. doi: 10.1186/s13578-022-00838-y.
 22. Changzheng Li, Binghuo Wu, Yishan Li, Jie Chen, Zhitao Ye, Xiaobin Tian, Jin Wang, Xi Xu, Shuai Pan, Yucan Zheng, Xiongwei Cai, Linjia Jiang, Meng Zhao. Amino acid catabolism regulates hematopoietic stem cell proteostasis via a GCN2-eIF2 α axis. *Cell Stem Cell*. doi: 10.1016/j.stem.2022.06.004.
 23. Qingyong Li, Jun Peng, Yuelian Luo, Jiaxin Zhou, Tailin Li, Lin Cao, Shuling Peng, Zhiyi Zuo, Zhi Wang. Far infrared light irradiation enhances A β clearance via increased exocytotic microglial ATP and ameliorates cognitive deficit in Alzheimer's disease-like mice. *J Neuroinflammation*. doi: 10.1186/s12974-022-02521-y.
 24. Yan-Yan Sun, Yi-Jun Wu. Tri-ortho-cresyl phosphate induces axonal degeneration in chicken DRG neurons by the NAD⁺ pathway. *Toxicol Lett*. doi: 10.1016/j.toxlet.2022.05.007.
 25. Yuan Liu, Dan Fang, Kangni Yang, Tianqi Xu, Chengrui Su, Ruichao Li, Xia Xiao, Zhiqiang Wang. Sodium dehydroacetate confers broad antibiotic tolerance by remodeling bacterial metabolism. *J Hazard Mater*. doi: 10.1016/j.jhazmat.2022.128645.
 26. Zheng-Hao Zhang, Zi-Guan Zhang, Min-Wei Chen, Ying Yang, Run-Jing Li, Jia-Jia Xu, Cui Yang, Yu-Ying Li, Hong-Wei Chen, Shi-Xiao Liu, Yan-Ling Li, Ping Luo, Yi-Jiang Liu, Wen-Bo Chen, Zhong-Gui Shan, Zheng-Rong Huang. Inhibition of GSDMD Activates Poly(ADP-ribosylation) and Promotes Myocardial Ischemia-Reperfusion Injury. *Oxid Med Cell Longev*. doi: 10.1155/2022/1115749.
 27. Yi Wang, Jie Li, Meng-Yue Wang, Zhi-Yong Pan, Zhi-Qiang Li, Ze-Fen Wang. Chronic microglial inflammation promotes neuronal lactate supply but impairs its utilization in primary rat astrocyte-neuron co-cultures. *Biochem Biophys Res Commun*. doi: 10.1016/j.bbrc.2022.03.122.
 28. Fan Yang, Yanbin Zhang, Sheng Liu, Jiheng Xiao, Yuxin He, Zengwu Shao, Yuhui Zhang, Xianyi Cai, Liming Xiong. Tunneling Nanotube-Mediated Mitochondrial Transfer Rescues Nucleus Pulposus Cells from Mitochondrial Dysfunction and Apoptosis. *Oxid Med Cell Longev*. doi: 10.1155/2022/3613319.
 29. Chao Zhao, Hu Zhang, Jingjing Zhou, Qiang Lu, Ying Zhang, Xiaojin Yu, Shizhi Wang, Ran Liu, Yuepu Pu, Lihong Yin. Metabolomics-based molecular signatures reveal the toxic effect of co-exposure to nitrosamines in drinking water. *Environ Res*. doi: 10.1016/j.envres.2021.111997.
 30. Xiaojun Zhang, Wensi Wu, Yuelian Luo, Zhi Wang. Transcranial Photobiomodulation Therapy Ameliorates Perioperative Neurocognitive Disorder Through Modulation of Mitochondrial Function in Aged Mice. *Neuroscience*. doi: 10.1016/j.neuroscience.2021.12.033.
 31. Kaili Ma, Lina Sun, Mingjing Shen, Xin Zhang, Zhen Xiao, Jiajia Wang, Xiaowei Liu, Kanqiu Jiang, F Xiao-Feng Qin, Feng Guo, Baojun Zhang, Lianjun Zhang. Functional assessment of the cell-autonomous role of NADase CD38 in regulating CD8⁺ T cell exhaustion. *iScience*. doi: 10.1016/j.isci.2022.104347.
 32. Rui Huang, Yuyuan Gao, Qingrui Duan, Qingxi Zhang, Peikun He, Jianing Chen, Guixian Ma, Limin Wang, Yuhu Zhang, Kun Nie, Lijuan Wang. Endothelial LRP1-ICD Accelerates Cognition-Associated Alpha-Synuclein Pathology and Neurodegeneration through PARP1 Activation in a Mouse Model of Parkinson's Disease. *Mol Neurobiol*. doi: 10.1007/s12035-022-03119-4.
 33. Li-Wei Liu, Yu Xie, Guan-Qun Li, Tao Zhang, Yu-Hang Sui, Zhong-Jie Zhao, Yang-Yang Zhang, Wen-Bo Yang, Xing-Long Geng, Dong-Bo Xue, Hua Chen, Yong-Wei Wang, Tian-Qi Lu, Li-Ren Shang, Zhi-Bo Li, Le Li, Bei Sun. Gut microbiota-derived nicotinamide mononucleotide alleviates acute pancreatitis by activating pancreatic SIRT3 signalling.

- Br J Pharmacol. doi: 10.1111/bph.15980.
34. Yu Cao, Meng Zhang, Ye Li, Jingjing Lu, Wanhui Zhou, Xiaoshuang Li, Hao Shi, Bin Xu, Shize Li. O-GlcNAcylation of SIRT1 Protects against Cold Stress-Induced Skeletal Muscle Damage via Amelioration of Mitochondrial Homeostasis. *Int J Mol Sci*. doi: 10.3390/ijms232314520.
 35. Yi Wu, Minhua Zong, Zhenhui Zhang, Yongtai Wu, Lin Li, Xia Zhang, Hong Wu, Bing Li. Selective transportation and energy homeostasis regulation of dietary advanced glycation end-products in human intestinal Caco-2 cells. *Food Chem*. doi: 10.1016/j.foodchem.2022.133284.
 36. Wenchang Zhang, Jing Guo, Qiang Chen. Role of PARP-1 in Human Cytomegalovirus Infection and Functional Partners Encoded by This Virus. *Viruses*. doi: 10.3390/v14092049.
 37. Mengyuan Li, Yimei Lai, Binfeng Chen, Chaohuan Guo, Mianjing Zhou, Siyuan Zhao, Shuyi Wang, Jin Li, Niansheng Yang, Hui Zhang. NAMPT is a metabolic checkpoint of IFN γ -producing CD4 $^+$ T cells in lupus nephritis. *Mol Ther*. doi: 10.1016/j.ymthe.2022.09.013.
 38. Guanghui Han, Weizhe Zhen, Yuan Dai, Hongni Yu, Dongyue Li, Tao Ma. Dihuang-Yinzi Alleviates Cognition Deficits via Targeting Energy-Related Metabolism in an Alzheimer Mouse Model as Demonstrated by Integration of Metabolomics and Network Pharmacology. *Front Aging Neurosci*. doi: 10.3389/fnagi.2022.873929.
 39. Meng Wang, Lingchen Wang, Yuan Zhou, Xiaoxuan Feng, Chaoyang Ye, Chen Wang. Shen Shuai II Recipe attenuates renal fibrosis in chronic kidney disease by improving hypoxia-induced the imbalance of mitochondrial dynamics via PGC-1 α activation. *Phytomedicine*. doi: 10.1016/j.phymed.2022.153947.
 40. Guotao Huang, Yong He, Li Hong, Min Zhou, Xiaohu Zuo, Zhihan Zhao. Restoration of NAD $^+$ homeostasis protects C2C12 myoblasts and mouse levator ani muscle from mechanical stress-induced damage. *Anim Cells Syst (Seoul)*. doi: 10.1080/19768354.2022.2106303.
 41. Wenwen Xue, Xin Li, Wuhao Li, Yixuan Wang, Chengfei Jiang, Lin Zhou, Jian Gao, Ying Yu, Yan Shen, Qiang Xu. Intracellular CYTL1, a novel tumor suppressor, stabilizes NDUFB1 to inhibit metabolic reprogramming in breast cancer. *Signal Transduct Target Ther*. doi: 10.1038/s41392-021-00856-1.
 42. Jiliang Xia, Jingyu Zhang, Xuan Wu, Wanqing Du, Yinghong Zhu, Xing Liu, Zhenhao Liu, Bin Meng, Jiaojiao Guo, Qin Yang, Yihui Wang, Qinglin Wang, Xiangling Feng, Guoxiang Xie, Yi Shen, Yanjuan He, Juanjuan Xiang, Minghua Wu, Gang An, Luguai Qiu, Wei Jia, Wen Zhou. Blocking glycine utilization inhibits multiple myeloma progression by disrupting glutathione balance. *Nat Commun*. doi: 10.1038/s41467-022-31248-w.
 43. Wenrui Ma, Jingjing Zhang, Shaowen Liu, Shiqiang Yan, Kehua Xu, Yu Shrike Zhang, Mieradilijiang Abudupataer, Yang Ming, Shichao Zhu, Bitao Xiang, Xiaonan Zhou, Shaman Luo, Hui Huang, Yuyi Tang, Shan Zhang, Zhuxin Xie, Nan Chen, Xiaoning Sun, Jun Li, Hao Lai, Chunsheng Wang, Kai Zhu, Weijia Zhang. Patient-derived microphysiological model identifies the therapeutic potential of metformin for thoracic aortic aneurysm. *EBioMedicine*. doi: 10.1016/j.ebiom.2022.104080.
 44. Qifan Wang, Jin Hu, Guoqiang Han, Peipei Wang, Sha Li, Jiwei Chang, Kexin Gao, Rong Yin, Yashu Li, Tong Zhang, Jihua Chai, Zhuying Gao, Tiantian Zhang, Ying Cheng, Chengli Guo, Jing Wang, Weidong Liu, Manman Cui, Yu Xu, Jinxuan Hou, Quan-Fei Zhu, Yu-Qi Feng, Haojian Zhang. PTIP governs NAD $^+$ metabolism by regulating CD38 expression to drive macrophage inflammation. *Cell Rep*. doi: 10.1016/j.celrep.2022.110603.
 45. Xixi Chen, Gaozan Tong, Junfu Fan, Yingjie Shen, Nan Wang, Wenjie Gong, Zijing Hu, Kunxuan Zhu, Xiaokun Li, Litai Jin, Weitao Cong, Jian Xiao, Zhongxin Zhu. FGF21 promotes migration and differentiation of epidermal cells during wound healing via SIRT1-dependent autophagy. *Br J Pharmacol*. doi: 10.1111/bph.15701.
 46. Xiuxiu Jin, Li Li, Qinlu Peng, Chunmei Gan, Li Gao, Siyu He, Shuangyan Tan, Wenchen Pu, Yu Liu, Yanqiu Gong, Yuqin Yao, Gang Wang, Xiaohui Liu, Meng Gong, Peng Lei, Huiyuan Zhang, Shiqian Qi, Heng Xu, Hongbo Hu, Biao Dong, Yong Peng, Dan Su, Lunzhi Dai. Glycyrrhetic acid restricts mitochondrial energy metabolism by targeting SHMT2. *iScience*. doi: 10.1016/j.isci.2022.104349.
 47. Weifan Lin, Xiangwan Lu, Hang Yang, Linxuan Huang, Wuheng Huang, Yuluan Tang, Situn Liu, Hua Wang, Yan Zhang. Metabolic heterogeneity protects metastatic mucosal melanomas cells from ferroptosis. *Int J Mol Med*. doi: 10.3892/ijmm.2022.5180.
 48. Zhen Zhang, Zichen Yu, Jinduo Wang, Yifa Yu, Lanxiao Li, Pengjie Sun, Xiaoguang Fan, Qingyang Xu. Metabolic engineering of *Escherichia coli* for efficient production of L-5-hydroxytryptophan from glucose. *Microb Cell Fact*. doi: 10.1186/s12934-022-01920-3.
 49. Jian Zhang, Yan Li, Jiaqi Liu, Fu Han, Jihong Shi, Gaofeng Wu, Kejia Wang, Kuo Shen, Ming Zhao, Xiaowen Gao, Chenyang Tian, Yunchuan Wang, Ke Tao, Dahai Hu. Adiponectin ameliorates hypertrophic scar by inhibiting Yes-associated protein transcription through SIRT1-mediated deacetylation of C/EBP β and histone H3. *iScience*. doi: 10.1016/j.isci.2022.105236.
 50. Shuai Zhao, Yun Hong, Yue-Yue Liang, Xiao-Lu Li, Jiang-Chuan Shen, Cong-Cong Sun, Ling-Luo Chu, Jie Hu, Hua Wang, De-Xiang Xu, Shi-Chen Zhang, Dou-Dou Xu, Tao Xu, Ling-Li Zhao. Compartmentalized regulation of NAD $^+$ by Di (2-ethyl-hexyl) phthalate induces DNA damage in placental trophoblast. *Redox Biol*. doi: 10.1016/j.redox.2022.102414.

Version 2024.04.18