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碧云天生物技术/Beyotime Biotechnology

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订货e-mail: order@beyotime.com

技术咨询: info@beyotime.com

网址: http://www.beyotime.com

## Proteinase K (20mg/ml)

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1. Zhao L, Wang H, Huo K, Cui L, Zhang W, Ni H, Zhang Y, Wu Z, Chu PK. Antibacterial nano-structured titania coating incorporated with silver nanoparticles. *Biomaterials*. 2011 Aug;32(24):5706-16.
2. Lu C, Jiang Q, Hu M, Tan C, Ji Y, Yu H, Hua Z. Preliminary Biological Evaluation of Novel <sup>99m</sup>Tc-Cys-Annexin A5 as a Apoptosis Imaging Agent. *Molecules*. 2013 Jun 10;18(6):6908-18.
3. Luo Z, Liu H, Sun X, Guo R, Cui R, Ma X, Yan M. RNA interference against discoidin domain receptor 2 ameliorates alcoholic liver disease in rats. *PLoS One*. 2013;8(2):e55860.
4. Huang Y, Li P, Xia S, Zhuo Y, Wu L. Proapoptotic effect and the mechanism of action of pingyangmycin on cavernous hemangiomas. *Exp Ther Med*. 2014 Feb;7(2):473-477.
5. Jiang L, He D, Yang D, Chen Z, Pan Q, Mao A, Cai Y, Li X, Xing H, Shi M, Chen Y, Bruce IC, Wang T, Jin L, Qi X, Hua D, Jin J, Ma X. MiR-489 regulates chemoresistance in breast cancer via epithelial mesenchymal transition pathway. *FEBS Lett*. 2014 May 29;588(11):2009-15.
6. Lv W, Cheng L, Li B. Development and evaluation of a novel TPGS-mediated paclitaxel-loaded PLGA-mPEG nanoparticle for the treatment of ovarian cancer. *Chem Pharm Bull (Tokyo)*. 2015;63(2):68-74.
7. Lu C, Jiang Q, Hu M, Tan C, Yu H, Hua Z. Preliminary biological evaluation of <sup>18</sup>F-FBEM-Cys-Annexin V a novel apoptosis imaging agent. *Molecules*. 2015 Mar 17;20(3):4902-14.
8. Chen Y, Wang L, Zhu Y, Chen Z, Qi X, Jin L, Jin J, Hua D, Ma X. Breast cancer resistance protein (BCRP)-containing circulating microvesicles contribute to chemoresistance in breast cancer. *Oncol Lett*. 2015 Dec;10(6):3742-8.
9. Li Z, Wang S, Zhao W, Sun Z, Yan H, Zhu J. Oxidized low-density lipoprotein upregulates microRNA-146a via JNK and NF-κB signaling. *Mol Med Rep*. 2016 Feb;13(2):1709-16.
10. Meng YQ, Zhang LF, Liu DY, Liu LW, Zhang Y, Zhao MJ. Synthesis and antitumor activity evaluation of novel ursolic acid derivatives. *J Asian Nat Prod Res*. 2016;18(3):280-8.
11. Du Z, Li S, Liu L, Yang Q, Zhang H, Gao C. NADPH oxidase 3-associated oxidative stress and caspase 3-dependent apoptosis in the cochlea of D-galactose-induced aged rats. *Mol Med Rep*. 2016 Jan;13(1):1056.