



碧云天生物技术/Beyotime Biotechnology
订货热线: 400-1683301 或 800-8283301
订货 e-mail: order@beyotime.com
技术咨询: info@beyotime.com
网址: http://www.beyotime.com

GSH (抗氧化剂)

产品编号	产品名称	包装
S0073	GSH (抗氧化剂)	1g

产品简介:

- GSH即glutathione, reduced, 也称-Glu-Cys-Gly或还原型谷胱甘肽，是一种细胞内天然存在的抗氧化剂。
- GSH为白色粉末，分子量307.3，分子式为C₁₀H₁₇N₃O₆S，纯度大于98%。
- 溶解于水、乙醇和DMF。

包装清单:

产品编号	产品名称	包装
S0073	GSH (抗氧化剂)	1g
—	说明书	1份

保存条件:

4°C保存，一年有效。

注意事项:

- 如果配制成溶液，分装后-20°C保存，半年有效。
- 本产品仅限于专业人员的科学的研究用，不得用于临床诊断或治疗，不得用于食品或药品，不得存放于普通住宅内。
- 为了您的安全和健康，请穿实验服并戴一次性手套操作。

使用本产品的文献:

1. Xu Y, Zhang L, Liu Y, Jin Z, Zhao Q, Yang F, Xiao D. . Sensitive and selective determination of GSH based on the ECL quenching of Ru(II) 1,10-phenanthroline-5,6-dione complex. Biosens Bioelectron. 2016 Mar;15;77:182-7.
2. QiuHong Shen, Yin Shi, Jiaqi Liu, Hua Su, Jingtao Huang, Yi Zhang, Chao Peng, Tianhua Zhou, Qiming Sun, Wei Wan, Wei Liu . Acetylation of STX17 (syntaxin 17) controls autophagosome maturation Autophagy. 2021 May;17(5):1157-1169.
3. Haitao Wang, Jianhua Zhang, Haoqiu Liu, Man Wang, Yan Dong, Yijun Zhou, Sek-Man Wong, Kai Xu, Qiufang Xu . A plant virus hijacks phosphatidylinositol-3,5-bisphosphate to escape autophagic degradation in its insect vector Autophagy. 2023 Apr;19(4):1128-1143.
4. Liu W, Wang L, Zheng C, Liu L, Wang J, Li D, Tan Y, Zhao X, He L, Shu W. . Microcystin-LR increases genotoxicity induced by aflatoxin B1 through oxidative stress and DNAbase excision repair genes in human hepatic cell lines. Environ Pollut. 2018 Feb;233:455-463.
5. Liu X, Zhang X, Ding Y, Zhou W, Tao L, Lu P, Wang Y, Hu R. . Nuclear Factor E2-Related Factor-2 Negatively Regulates NLRP3 Inflammasome Activity by Inhibiting Reactive Oxygen Species-Induced NLRP3 Priming. ANTIOXID REDOX SIGN. 2017 Jan 1;26(1):28-43.
6. Wang Y, Fang W, Huang Y, Hu F, Ying Q, Yang W, Xiong B. . Reduction of selenium-binding protein 1 sensitizes cancer cells to selenite via elevating extracellular glutathione: a novel mechanism of cancer-specific cytotoxicity of selenite. FREE RADICAL BIO MED. 2015 Feb;79:186-96.
7. Wang L, Hao J, Hu J, Pu J, Lü Z, Zhao L, Wang Q, Yu Q, Wang Y, Li G. . Protective Effects of Ginsenosides against Bisphenol A-Induced Cytotoxicity in 15P-1 Sertoli Cells via Extracellular Signal-Regulated Kinase 1/2 Signalling and Antioxidant Mechanisms. BASIC CLIN PHARMACOL. 2012 Jul;111(1):42-9.
8. Zhou F, Chen S, Xiong J, Li Y, Qu L. . Luteolin Reduces Zinc-Induced Tau Phosphorylation at Ser262/356 in an ROS-Dependent Manner in SH-SY5Y Cells. Biol Trace Elem Res. 2012 Nov;149(2):273-9 .
9. Shi K, Wang D, Cao X, Ge Y. . Endoplasmic reticulum stress signaling is involved in mitomycin C (MMC)-induced apoptosis in human fibroblasts via PERK pathway. PLoS One. 2013;8(3):e59330.
10. Zhu L, Yuan C, Zhang D, Ma Y, Ding X, Zhu G. . BHV-1 induced oxidative stress contributes to mitochondrial dysfunction in MDBK cells. Vet Res. 2016 Mar 22;47:47.
11. Li SW, Liu CM, Guo J, Marcondes AM, Deeg J, Li X, Guan F. . Iron overload induced by ferric ammonium citrate triggers reactive oxygen species-mediated apoptosis via both extrinsic and intrinsic pathways in human hepatic cells. Hum Exp Toxicol. 2016 Jun;35(6):598-607.
12. Liu X, Li J, Peng X, Lv B, Wang P, Zhao X, Yu B. . Geraniin Inhibits LPS-Induced THP-1 Macrophages Switching to M1 Phenotype via SOCS1/NF-κB Pathway. Inflammation. 2016 Aug;39(4):1421-33.
13. Zhang Y, Cheng Z, Wang C, Ma H, Meng W, Zhao Q. . Neuroprotective Effects of Kukoamine a against Radiation-induced Rat Brain Injury through Inhibition of Oxidative Stress and Neuronal Apoptosis. Neurochem Res. 2016 Oct;41(10):2549-2558.
14. Xie Q, Wei M, Zhang B, Kang X, Liu D, Zheng W, Pan X, Quan Y, Liao D, Shen J. . MicroRNA-33 regulates the NLRP3 inflammasome signaling pathway in macrophages. Mol Med Rep. 2018 Feb;17(2):3318-3327.

15. Ding X, Sun W, Chen J . IL-2 augments the sorafenib-induced apoptosis in liver cancer by promoting mitochondrial fission and activating the JNK/TAZ pathway. *Cancer Cell Int.* 2018 Nov 9;18:176.
16. Zhang L, Li S, Wang R, Chen C, Ma W, Cai H . Cytokine augments the sorafenib-induced apoptosis in Huh7 liver cancer cell by inducing mitochondrial fragmentation and activating MAPK-JNK signalling pathway. *Biomed Pharmacother.* 2019 Feb 110:213-223.
17. Liang X, Wang Z, Gao M, Wu S, Zhang J, Liu Q, Yu Y, Wang J, Liu W . Cyclic stretch induced oxidative stress by mitochondrial and NADPH oxidase in retinal pigment epithelial cells. *BMC Ophthalmol.* 2019 Mar 18 19(1):79.
18. Huo L, Su Y, Xu G, Zhai L, Zhao J . Sulforaphane Protects the Male Reproductive System of Mice from Obesity-Induced Damage: Involvement of Oxidative Stress and Autophagy. *INT J ENV RES PUB HE.* 2019 Oct 7 16(19). pii: E3759.
19. Xiaojin Luan, Yidan Yan, Qianwen Zheng, Min Wang, Wanyin Chen, Jun Yu, Jie Fang . Excessive reactive oxygen species induce apoptosis via the APPL1-Nrf2/HO-1 antioxidant signalling pathway in trophoblasts with missed abortion *Life Sci.* 2020 Aug 1;254:117781.
20. Lei Chen, Miao Hao, Jingmin Yan, Lin Sun, Guihua Tai, Hairong Cheng, Yifa Zhou . Citrus-derived DHCP inhibits mitochondrial complex II to enhance TRAIL sensitivity via ROS-induced DR5 upregulation *J Biol Chem.* 2021 Jan-Jun;296:100515.
21. Peifu Feng, Ling Li, Jing Dai, Lingli Zhou, Jing Liu, Jinfeng Zhao, Xiaodong Li, Neng Ling, Siyuan Qiu, Lin Zhang, Tiantian Xie, Yinglei Chen, Michael J Donovan, Tianhuan Peng, Jianhui Song, Mao Ye . The regulation of NONO by USP11 via deubiquitination is linked to the proliferation of melanoma cells *J Cell Mol Med.* 2021 Feb;25(3):1507-1517.
22. Yuchen Xia, Guihua Wang, Manli Jiang, Xuetong Liu, Yan Zhao, Yinghui Song, Binyuan Jiang, Demao Zhu, Ling Hu, Zhao Zhang, Ting Cao, Ji Ming Wang, Jinyue Hu . A Novel Biological Activity of the STAT3 Inhibitor Stattic in Inhibiting Glutathione Reductase and Suppressing the Tumorigenicity of Human Cervical Cancer Cells via a ROS-Dependent Pathway *Onco Targets Ther.* 2021 Jul 5;14:4047-4060.
23. Fei Han, Shijie Li, Yankun Yang, Zhonghu Bai . Interleukin-6 promotes ferroptosis in bronchial epithelial cells by inducing reactive oxygen species-dependent lipid peroxidation and disrupting iron homeostasis *Bioengineered.* 2021 Dec;12(1):5279-5288.
24. Qing Su, Xiao-Jing Yu, Xiao-Min Wang, Bo Peng, Juan Bai, Hong-Bao Li, Ying Li, Wen-Jie Xia, Li-Yan Fu, Kai-Li Liu, Jin-Jun Liu, Yu-Ming Kang . Na+/K+-ATPase Alpha 2 Isoform Elicits Rac1-Dependent Oxidative Stress and TLR4-Induced Inflammation in the Hypothalamic Paraventricular Nucleus in High Salt-Induced Hypertension Antioxidants (Basel). 2022 Jan 31;11(2):288.
25. Hai-Fan Yu, Zhan-Qing Yang, Ming-Yue Xu, Ji-Cheng Huang, Zhan-Peng Yue, Bin Guo . Yap is essential for uterine decidualization through Rrm2/GSH/ROS pathway in response to Bmp2 *Int J Biol Sci.* 2022 Mar 6;18(6):2261-2276.
26. Kun Zeng, Wenqun Xi, Yuanjiao Qiao, Xiaosheng Huang, Xinhua Liu . Paeoniflorin inhibits epithelial mesenchymal transformation and oxidative damage of lens epithelial cells in diabetic cataract via sirtuin 1 upregulation *Bioengineered.* 2022 Mar;13(3):5903-5914.
27. Min Ling, Qing Liu, Yufei Wang, Xuetong Liu, Manli Jiang, Jinyue Hu . LCS-1 inhibition of superoxide dismutase 1 induces ROS-dependent death of glioma cells and degrades PARP and BRCA1 *Front Oncol.* 2022 Aug 1:12:937444.

Version 2024.03.12