



碧云天网站



微信公众号

碧云天生物技术/Beyotime Biotechnology

订货热线: 400-168-3301或800-8283301

订货e-mail: order@beyotime.com

技术咨询: info@beyotime.com

网址: http://www.beyotime.com

PD 98,059 (MAPKK抑制剂)

产品编号	产品名称	包装
S1805	PD 98,059 (MAPKK抑制剂)	2mg

产品简介:

- PD 98,059也写作PD98059, PD 98059, PD098059或PD 098059, 是一种常用的MAPKK抑制剂, 即MAPK kinase抑制剂或MEK抑制剂。PD 98,059可以通透细胞, 选择性抑制MEK1(一种MAPK kinase), 从而抑制MAP kinase (Erk1/2即p44/42 MAPK)的磷酸化和激活。
- PD 98,059分子量为267.28, 分子式为 $C_{16}H_{13}NO_3$, CAS Number:167869-21-8。本产品纯度大于98%。
- 本PD 98,059为进口分装, 用DMSO配制, 浓度为20mg/ml, 共0.1ml。

包装清单:

产品编号	产品名称	包装
S1805	PD 98,059 (MAPKK抑制剂, 20mg/ml)	2mg
—	说明书	1份

保存条件:

-20°C避光保存, 一年有效。

注意事项:

- 本产品对人体有害, 操作时请小心, 并注意有效防护以避免直接接触人体或吸入体内。
- 本PD 98,059在4°C、冰浴等较低温度情况下会凝固而粘在离心管管底、管壁或管盖内, 可以20-25°C水浴温育片刻至全部融解后使用。
- 本产品仅限于专业人员的科学研究用, 不得用于临床诊断或治疗, 不得用于食品或药品, 不得存放于普通住宅内。
- 为了您的安全和健康, 请穿实验服并戴一次性手套操作。

使用说明:

1. PD 98,059常见使用浓度范围为10-100 μ M。具体的最佳工作浓度请参考相关文献, 或根据实验目的, 以及所培养的特定细胞和组织, 通过实验进行摸索和优化。

使用本产品的文献:

1. Tang XD, Zhou X, Zhou KY. Dauricine inhibits insulin-like growth factor-I-induced hypoxia inducible factor 1alpha protein accumulation and vascular endothelial growth factor expression in human breast cancer cells. *Acta Pharmacol Sin.* 2009 May;30(5):605-16.
2. Kui L, Weiwei Z, Ling L, Daikun H, Guoming Z, Linuo Z, Renming H. Ghrelin inhibits apoptosis induced by high glucose and sodium palmitate in adult rat cardiomyocytes through the PI3K-Akt signaling pathway. *Regul Pept.* 2009 Jun 5;155(1-3):62-9.
3. Shi X, Zhou B. The role of Nrf2 and MAPK pathways in PFOS-induced oxidative stress in zebrafish embryos. *Toxicol Sci.* 2010;115(2):391-400.
4. Jin W, Li Q, Lin Y, Lu Y, Li H, Wang L, Hu R, Ma L, Wang J, Pang T. Reversal of Imatinib resistance in BCR-ABL-positive leukemia after inhibition of the Na⁺/H⁺exchanger. *Cancer Lett.* 2011 Sep 1; 308(1):81-90.
5. Gong J, Shen XH, Chen C, Qiu H, Yang RG Down-regulation of HIV-1 infection by inhibition of the MAPK signaling pathway. *Virology.* 2011 Apr;26(2):114-22.
6. Chen S, Liu AR, An FM, Yao WB, Gao XD. Amelioration of neurodegenerative changes in cellular and rat models of diabetes-related Alzheimer's disease by exendin-4. *Age (Dordr).* 2011 Oct;34(5):1211-4.
7. Liu S, Li Y, Lin T, Fan X, Liang Y, Heemann U High dose human insulin and insulin glargine promote T24 bladder cancer cell proliferation viaPI3K-independent activation of Akt. *Diabetes Res Clin Pract.* 2011 Feb;91(2):177-82.
8. Huang Q, Yang Y, Li X, Huang S Transcription suppression of SARI (suppressor of AP-1, regulated by IFN) by BCR-ABL in humanleukemia cells. *Tumour Biol.* 2011 Dec;32(6):1191-7.
9. Fang S, Jin Y, Zheng H, Yan J, Cui Y, Bi H, Jia H, Zhang H, Wang Y, Na L, Gao X, Zhou H High glucose condition upregulated Txnip expression level in rat mesangial cells throughROS/MEK/MAPK pathway. *Mol Cell Biochem.* 2011 Jan;347(1-2):175-82.
10. Li W, Liu Y, Li XX, Yu Y, Wu JJ, Wang Q, Huo H, Wang LM, Yang L. MAPKs are not involved in triptolide-induced cell growth inhibition and apoptosis in prostate cancer cell lines with different p53 status. *Planta Med.* 2011 Jan;77(1):27-31.
11. Feng J, Zhang P, Chen X, He G PI3K and ERK/Nrf2 pathways are involved in oleanolic acid-induced heme oxygenase-1 expressionin rat vascular smooth muscle cells. *J Cell Biochem.* 2011 Jun;112(6):1524-31.
12. Lai T, Wang K, Hou Q, Zhang J, Yuan J, Yuan L, You Z, Xi M. Interleukin 17 induces up-regulation of chemokine and cytokine expression via activation of thenuclear factor κ B and extracellular signal-regulated kinase 1/2 pathways in gynecologic cancer cell lines. *Int J Gynecol Cancer.* 2011 Dec;21(9):1533-9.
13. Sheng X, Li M, Song S, Zhang N, Wang Y, Liang H, Wang W, Ji A Sulfated polysaccharide isolated from the sea cucumber *Stichopus japonicus* promotesneurosphere migration and differentiation via up-regulation of N-Cadherin. *Cell Mol Neurobiol.* 2012 Apr; 32(3):435-42.
14. Xu B, Wang N, Wang X, Tong N, Shao N, Tao J, Li P, Niu X, Feng N, Zhang L, Hua L, Wang Z, Chen M. MiR-146a suppresses tumor growth and progression by targeting EGFR pathway and in a p-ERK-dependent manner in castration-resistant prostate cancer. *Prostate.* 2012 Aug

1;72(11):1171-8.

15. Li ZL, Liu JC, Hu J, Li XQ, Wang SW, Yi DH, Zhao MG. Protective effects of hyperoside against human umbilical vein endothelial cell damage induced by hydrogen peroxide. *J Ethnopharmacol.* 2012 Jan 31;139(2):388-94.
16. Liu L, Wang Y, Fan H, Zhao X, Liu D, Hu Y, Kidd AR 3rd, Bao J, Hou Y. MicroRNA-181a Regulates Local Immune Balance by Inhibiting Proliferation and Immunosuppressive Properties of Mesenchymal Stem Cells. *Stem Cells.* 2012 Aug;30(8):1756-70.
17. Wang LH, Chang GQ, Zhang HJ, Wang J, Lin YN, Jin WN, Li HW, Gao W, Wang RJ, Li QH, Pang TX. Neutrophil gelatinase-associated lipocalin regulates intracellular accumulation of Rh123 in cancer cells. *Genes Cells.* 2012 Mar;17(3):205-17.
18. Song LJ, Liu RJ, Zeng Z, Alper SL, Cui HJ, Lu Y, Zheng L, Yan ZW, Fu GH. Gastrin inhibits a novel, pathological colon cancer signaling pathway involving EGR1, AE2, and P-ERK. *J Mol Med (Berl).* 2012 Jun;90(6):707-18.
19. Hu Y, Bai Y, Chu Z, Wang J, Wang L, Yu M, Lian Z, Hua J. GSK3 inhibitor-BIO regulates proliferation of female germline stem cells from the postnatal mouse ovary. *Cell Prolif.* 2012 Aug;45(4):287-98.
20. Zhong Y, Liu T, Guo Z. Curcumin inhibits ox-LDL-induced MCP-1 expression by suppressing the p38MAPK and NF- κ B pathways in rat vascular smooth muscle cells. *Inflamm Res.* 2012 Jan;61(1):61-7.
21. Tian J, Jiang Y. Insulin upregulates the expression of zinc finger and BTB domain-containing 7A in HepG2 cells. *Mol Med Rep.* 2012 Dec;6(6):1379-84.
22. Liu SY, Liang Y, Lin TX, Su F, Liang WW, Uwe H, Li Y. MEK1 and MEK2 differentially regulate human insulin- and insulin glargine-induced human bladder cancer T24 cell proliferation. *Chin Med J (Engl).* 2012 Dec;125(23):4197-201.
23. Zhong Y, Liu T, Lai W, Tan Y, Tian D, Guo Z. Heme oxygenase-1-mediated reactive oxygen species reduction is involved in the inhibitory effect of curcumin on lipopolysaccharide-induced monocyte chemoattractant protein-1 production in RAW264.7 macrophages. *Mol Med Rep.* 2013 Jan;7(1):242-6.
24. Xu S, Wen H, Jiang H. Urotensin II promotes the proliferation of endothelial progenitor cells through p38 and p44/42 MAPK activation. *Mol Med Rep.* 2012 Jul;6(1):197-200.
25. Hu X, Wang H, Liu J, Fang X, Tao K, Wang Y, Li N, Shi J, Wang Y, Ji P, Cai W, Bai X, Zhu X, Han J, Hu D. The role of ERK and JNK signaling in connective tissue growth factor induced extracellular matrix protein production and scar formation. *Arch Dermatol Res.* 2013 Jul;305(5):433-45.
26. Lv C, Sun W, Sun H, Wei S, Chen R, Wang B, Huang C. Asperolide A, a marine-derived tetranorditerpenoid, induces G2/M arrest in human NCI-H460 lung carcinoma cells, is mediated by p53-p21 stabilization and modulated by Ras/Raf/MEK/ERK signaling pathway. *Mar Drugs.* 2013 Jan 29;11(2):316-31.
27. Yuan Y, Luo Q, Yang L, Song G. Role of FAK-ERK1/2 Signaling Pathway in Proliferation of Rat Bone-marrow Mesenchymal Stem Cells Stimulated by Cyclic Stretching. *Journal of Medical and Biological Engineering.* 2013;33(2):229-38.
28. Liang L, Luo M, Fu Y, Zu Y, Wang W, Gu C, Zhao C, Li C, Efferth T. Cajaninstilbene acid (CSA) exerts cytoprotective effects against oxidative stress through the Nrf2-dependent antioxidant pathway. *Toxicol Lett.* 2013 Jun 7;219(3):254-61.
29. Zhao L, Zhang W, Chen M, Zhang J, Zhang M, Dai K. Aspirin Induces Platelet Apoptosis. *Platelets.* 2013;24(8):637-42.
30. Deng X, Rui W, Zhang F, Ding W. PM2.5 induces Nrf2-mediated defense mechanisms against oxidative stress by activating PIK3/AKT signaling pathway in human lung alveolar epithelial A549 cells. *Cell Biol Toxicol.* 2013 Jun;29(3):143-57.
31. Xuejiao S, Yong X, Ningyu W, Lidan Z, Xuanhong S, Youzhi X, Tinghong Y, Yaojie S, Yongxia Z, Luoting Y. A Novel Benzothiazole Derivative YLT322 Induces Apoptosis via the Mitochondrial Apoptosis Pathway In Vitro with Anti-Tumor Activity in Solid Malignancies. *PLoS One.* 2013 May 30;8(5):e63900.
32. Li ZL, Hu J, Li YL, Xue F, Zhang L, Xie JQ, Liu ZH, Li H, Yi DH, Liu JC, Wang SW. The effect of hyperoside on the functional recovery of the ischemic/reperfused isolated rat heart: potential involvement of the extracellular signal-regulated kinase 1/2 signaling pathway. *Free Radic Biol Med.* 2013 Apr;57:132-40.
33. Zhang K, Meng X, Kong J, Liu FF, Yang JM, Gao F, Zhang Y, Zhang C. Simvastatin increases Prolyl-4-Hydroxylase α 1 expression in atherosclerotic plaque and ox-LDL-stimulated human aortic smooth muscle cells via p38 MAPK and ERK1/2 signaling. *J Mol Cell Cardiol.* 2013 Dec;65:43-50.
34. Feng Y, Niu LL, Wei W, Zhang WY, Li XY, Cao JH, Zhao SH. A feedback circuit between miR-133 and the ERK1/2 pathway involving an exquisite mechanism for regulating myoblast proliferation and differentiation. *Cell Death Dis.* 2013 Nov 28;4:e934.
35. Dong S, Kong J, Kong F, Kong J, Gao J, Ke S, Wang S, Ding X, Sun W, Zheng L. Insufficient radiofrequency ablation promotes epithelial-mesenchymal transition of hepatocellular carcinoma cells through Akt and ERK signaling pathways. *J Transl Med.* 2013 Oct 29;11:273.
36. Wang Y, Huang Y, Guan F, Xiao Y, Deng J, Chen H, Chen X, Li J, Huang H, Shi C. Hypoxia-inducible factor-1 α and MAPK co-regulate activation of hepatic stellate cells upon hypoxia stimulation. *PLoS One.* 2013 Sep 10;8(9):e74051.
37. Li W, Wu Z, Ma Q, Liu J, Xu Q, Han L, Duan W, Lv Y, Wang F, Reindl KM, Wu E. Hyperglycemia regulates TXNIP/TRX/ROS axis via p38 MAPK and ERK pathways in pancreatic cancer. *Curr Cancer Drug Targets.* 2014 May;14(4):348-56.
38. Li L, Han MX, Li S, Xu Y, Wang L. Hypoxia Regulates the Proliferation and Osteogenic Differentiation of Human Periodontal Ligament Cells Under Cyclic Tensile Stress via MAPK Pathways. *J Periodontol.* 2014 Mar;85(3):498-508.
39. Wang P, Zeng Y, Liu T, Zhang C, Yu PW, Hao YX, Luo HX, Liu G. Chloride intracellular channel 1 regulates colon cancer cell migration and invasion through ROS/ERK pathway. *World J Gastroenterol.* 2014 Feb 28;20(8):2071-8.
40. Liu C, Yang J, Fu W, Qi S, Wang C, Quan C, Yang K. Coactivation of the PI3K/Akt and ERK signaling pathways in PCB153-induced NF- κ B activation and caspase inhibition. *Toxicol Appl Pharmacol.* 2014 Jun 15;277(3):270-8.
41. Qian C, Yao J, Wang J, Wang L, Xue M, Zhou T, Liu W, Si J. ERK1/2 inhibition enhances apoptosis induced by JAK2 silencing in human gastric cancer SGC7901 cells. *Mol Cell Biochem.* 2014 Feb;387(1-2):159-70.
42. Qi T, Zhang W, Luan Y, Kong F, Xu D, Cheng G, Wang Y. Proteomic profiling identified multiple short-lived members of the central proteome as the direct targets of the addicted oncogenes in cancer cells. *Mol Cell Proteomics.* 2014 Jan;13(1):49-62.
43. Zheng HL, Yang J, Hou Y, Sun B, Zhang Q, Mou Y, Wand L, Wu C. Oligomer procyanidins (F2) isolated from grape seeds inhibits tumor angiogenesis and cell invasion by targeting HIF-1 α in vitro. *Int J Oncol.* 2015 Feb;46(2):708-20.
44. Xu D, Wang J, Zhou Z, He Z, Zhao Q. Cannabinoid WIN5, 212-2 induces cell cycle arrest and inhibits the proliferation and migration of human BEL7402 hepatocellular carcinoma cells. *Mol Med Rep.* 2015 Dec;12(6):7963-70.
45. Fang J, Wang H, Xi W, Cheng G, Wang S, Su S, Zhang S, Deng Y, Song Z, Xu A, Liu B, Cao J, Wang Z. Downregulation of tNASP inhibits proliferation through regulating cell cycle-related proteins and inactive ERK/MAPK signal pathway in renal cell carcinoma cells. *Tumour Biol.* 2015 Jul;36(7):5209-14.
46. Wang X, Shi L, Deng Y, Qu M, Mao S, Xu L, Xu W, Fang C. Inhibition of leucine aminopeptidase 3 suppresses invasion of ovarian cancer cells through down-regulation of fascin and MMP-2/9. *Eur J Pharmacol.* 2015 Dec 5;768:116-22.
47. Han X, Wang T, Zhang J, Liu X, Li Z, Wang G, Song Q, Pang D, Ouyang H, Tang X. Apolipoprotein CIII regulates lipoprotein-associated phospholipase A2 expression via the MAPK and NF κ B pathways. *Biol Open.* 2015 Apr 2;4(5):661-5.
48. Zhang B, Luo Q, Sun J, Xu B, Ju Y, Yang L, Song G. MGF enhances tenocyte invasion through MMP-2 activity via the FAK-ERK1/2 pathway. *Wound Repair Regen.* 2015 May-Jun;23(3):394-402.
49. Wang J, Xu H, Wang Q, Zhang H, Lin Y, Zhang H, Li Q, Pang T. CIAPIN1 targets Na⁺/H⁺ exchanger 1 to mediate MDA-MB-231 cells' metastasis through regulation of MMPs via ERK1/2 signaling pathway. *Exp Cell Res.* 2015 Apr 10;333(1):60-72.
50. Luo Q, Wu K, Zhang B, Song G. Mechano growth factor E peptide promotes rat bone marrow-derived mesenchymal stem cell migration through CXCR4-ERK1/2. *Growth Factors.* 2015;33(3):210-9.
51. Bao J, Huang B, Zou L, Chen S, Zhang C, Zhang Y, Chen M, Wan JB, Su

- H, Wang Y, He C. Hormetic Effect of Berberine Attenuates the Anticancer Activity of Chemotherapeutic Agents. *PLoS One*. 2015 Sep 30;10(9):e0139298.
52. Xiao M, Tang Y, Chen WW, Wang YL, Yang L, Li X, Song GL, Kuang J. Tubb3 regulation by the Erk and Akt signaling pathways: a mechanism involved in the effect of arginine ADP-ribosyltransferase 1 (Art1) on apoptosis of colon carcinoma CT26 cells. *Tumour Biol*. 2016;37(2):2353-63.
53. Shi M, Ren X, Wang X, Wang H, Liu G, Yuan X, Zheng S, Yu L, Pan S, Song G, Guo Q, Li L, Zhang X, Zhang Z, Ding H, Jiang G. A novel combination of oridonin and valproic acid in enhancement of apoptosis induction of HL-60 leukemia cells. *Int J Oncol*. 2016 Feb;48(2):734-46.
54. Wu X, Yang L, Zheng Z, Li Z, Shi J, Li Y, Han S, Gao J, Tang C, Su L, Hu D. Src promotes cutaneous wound healing by regulating MMP-2 through the ERK pathway. *Int J Mol Med*. 2016 Mar;37(3):639-48.
55. Wan J, Wu W, Chen Y, Kang N, Zhang R. Insufficient radiofrequency ablation promotes the growth of non-small cell lung cancer cells through PI3K/Akt/HIF-1 α signals. *Acta Biochim Biophys Sin (Shanghai)*. 2016 Apr;48(4):371-7.

Version 2016.12.12